

Coloplast Product Guide



Ostomy, Continence, Wound and Skin Care
Effective March, 2021





"I did a lot of research during my early diagnosis and it became clear to me that the Coloplast SenSura® Mio Click would be hugely beneficial. On my first vacation after surgery, the SenSura® Mio Click allowed me to swim in the ocean and participate in activities such as snorkeling and kite surfing. I was confident that my urostomy appliance would not let me down or become a problem while away."

Peter
SenSura® Mio Urostomy User

OSTOMY

Appliances

SenSura® Mio Click

Ostomy Barrier and Pouch - Two Piece System

SenSura® Mio Click Skin Barrier	9
SenSura® Mio Click Drainable Pouch	10
SenSura® Mio Click Drainable Pouch	10
SenSura® Mio Click Closed Pouch	10
SenSura® Mio Click Urostomy Pouch	10
SenSura® Mio Click High Output Pouch	11

SenSura® Mio Flex

Ostomy Barrier and Pouch - Two Piece System

SenSura® Mio Flex Skin Barrier	11
SenSura® Mio Flex Drainable Pouch	12
SenSura® Mio Flex Closed Pouch	12
SenSura® Mio Flex Urostomy Pouch	12
SenSura® Mio Flex High Output Pouch	12

SenSura® Mio

Ostomy Barrier and Pouch - One Piece System

SenSura® Mio Drainable Pouch	12
SenSura® Mio Drainable Pouch	13
SenSura® Mio Drainable Pouch	14
SenSura® Mio Closed Pouch	14
SenSura® Mio Closed Pouch	15
SenSura® Mio One Piece Urostomy	15
SenSura® Mio One Piece High Output Pouch	16
SenSura® Mio One Piece Post-Op Pouch	16

SenSura® Mio Baby

Ostomy Barrier and Pouch - Two Piece System

SenSura® Mio Baby Flex Skin Barrier	16
SenSura® Mio Baby Flex Drainable Pouch	16
SenSura® Mio Baby Adhesive Remover Wipe	16

SenSura® Mio Kids

Ostomy Barrier and Pouch - Two Piece System

SenSura® Mio Kids Flex Skin Barrier	17
SenSura® Mio Kids Flex Drainable Pouch	17
SenSura® Mio Kids Flex Urostomy Pouch	17

Ostomy Barrier and Pouch - One Piece System

SenSura® Mio Kids Drainable Pouch 1 Piece	17
SenSura® Mio Kids Pouch 1 Piece Urostomy	17

SenSura® Click

Ostomy Barrier and Pouch- Two Piece System

SenSura® Click Skin Barrier	18
SenSura® Click Xpro Skin Barrier	18
SenSura® Click Drainable Pouch	19
SenSura® Click Closed Pouch	19

SenSura® Click Urostomy Pouch	19
SenSura® High Output Pouch	19

SenSura® Flex

Ostomy Barrier and Pouch - Two Piece System

SenSura® Flex Skin Barrier	20
SenSura® Flex Xpro Skin Barrier	20
SenSura® Flex Drainable Pouch	20
SenSura® Flex Closed Pouch	20
SenSura® Flex Urostomy Pouch	21

SenSura® - One Piece

Ostomy Barrier and Pouch - One Piece System

SenSura® Drainable Pouch with Easyclose™	21
SenSura® Xpro Drainable Pouch	21
SenSura® Mio Original Closed Pouch	21
SenSura® Closed Pouch	22
SenSura® Urostomy Pouch	22
SenSura® One-Piece High Output Pouch	22
SenSura® Post-Op Pouch	22

Assura®

Ostomy Barrier and Pouch - Two Piece System

Assura® Skin Barrier	23
Assura® Skin Barrier - Extra Extended Wear	24
Assura® 2 pc. High Output Drainable Pouch	24
Assura® 2 pc. High Output Drainable System	25
Assura® Drainable Pouch Clamp Closure	25
Assura® Drainable Pouch with EasiClose™	25
Assura® Closed Pouch	25
Assura® Urostomy Pouch	25

Easiflex®

Ostomy Barrier and Pouch- Two Piece System

Easiflex® Skin Barrier	26
Easiflex® Skin Barrier - Extra Extended Wear	26
Easiflex® Drainable Pouch with EasiClose™	26
Easiflex® Closed Pouch	27
Easiflex® Urostomy Pouch	27

Assura® - One Piece

Ostomy Barrier and Pouch - One Piece System

Assura® Drainable Pouch Clamp Closure	27
Assura® Drainable Pouch with EasiClose™	27
Assura® Drainable Pouch with EasiClose™	28
Assura® Closed Pouch	28
Assura® 1-Piece Post-Op Pouch	28
Assura® Closed Pouch	29
Assura® Urostomy Pouch	29
Assura® Urostomy Pouch - Extra Extended Wear	29

Assura® - Pediatric - Two Piece

Ostomy Barrier and Pouch - Two Piece System - Pediatric

Assura® Pediatric Skin Barrier	30
--------------------------------------	----

Assura® Pediatric Drainable Pouch	30
Assura® Pediatric Closed Pouch	30

EasiFlex® - Pediatric

Ostomy Barrier and Pouch - Two Piece System - Pediatric	
Easiflex® Pediatric Skin Barrier	30
Easiflex® Pediatric Drainable Pouch	30

Assura® - Pediatric - One Piece

Ostomy Barrier and Pouch - One Piece System - Pediatric	
Assura® Pediatric Drainable Pouch	30
Assura® Pediatric Closed Pouch	30
Assura® Pediatric Urostomy Pouch	30

Accessories

Ostomy Accessories

Brava® Protective Sheets.....	31
Brava® Protective Rings	31
Brava® Protective Rings Wide.....	31
Brava® Protective Rings Extra Wide	31
Brava® Mouldable Rings.....	31
Brava® Skin Barrier Spray and Wipe	31
Brava® Adhesive Remover Spray and Wipe.....	31
Brava® Barrier Cream	31
Brava® Ostomy Strip Paste.....	31
Brava® Elastic Barrier Strips.....	32
Brava® Ostomy Powder.....	32
Brava® Mio Belt for All SenSura Mio Products.....	32
Brava® Belt for SenSura Original + Assura Products.....	32
Brava® Lubricating Deodorant.....	32
Brava® Alcohol Free Paste	32
Ostomy Paste (contains alcohol).....	32
Skin Barrier Rings.....	32
Coloplast Prep™ Protective Liquid Skin Barrier.....	32
Coloplast FiltroDor® Pouch Filter	32
Coloplast Pouch Clamp.....	33
Ostomy Scissors	33

Hex-On®

Odour antagonist

Hex-On®	33
---------------	----

Assura® - Ostomy Cap

Coloplast Uro Minicap (100 ml).....	33
Coloplast Stoma Cap	33
Coloplast Stoma Cap - One Piece	33

Assura® Irrigation Products

Assura® Irrigation Set	33
Assura® Irrigation Accessories.....	34
Assura® Irrigation Sleeves.....	34
Coloplast Urostomy Night Bag.....	34

Coloplast Fistula & Wound Management System

Pouch	34
-------------	----

Flexible Lids	34
Accessories	34

CONTINENCE

Hydro Intermittent Catheters

SpeediCath® Flex

Hydrophilic Intermittent Catheter	
SpeediCath® Flex	37

SpeediCath® Standard

Hydrophilic Intermittent Catheter	
SpeediCath® Standard Male Straight	37
SpeediCath® Standard Male Coude	37
SpeediCath® Standard Female Straight	37
SpeediCath® Standard Pediatric Straight	37

SpeediCath® Compact Male

Hydrophilic Intermittent Catheter	
SpeediCath® Compact Male.....	38

SpeediCath® Compact Eve

Hydrophilic Intermittent Catheter	
SpeediCath® Compact Eve.....	38

SpeediCath® Compact Female

Hydrophilic Intermittent Catheter	
SpeediCath® Compact Female.....	38
SpeediCath® Compact Female Plus	38

SpeediCath® Compact Set

Hydrophilic Intermittent Catheter	
SpeediCath® Compact Set Male	39
SpeediCath® Compact Set Female	39

Self-Cath®

Intermittent Catheter

Self-Cath® Female Straight	39
Self-Cath® Pediatric	39
Self-Cath® Extension Tube	39
Self-Cath® Male Straight	40
Self-Cath® Male Soft Straight	40
Self-Cath® Male Straight Curved Packaging	40
Self-Cath® Male Coude Tapered with Guide Stripe	40
Self-Cath® Male Coude Olive with Guide Stripe	40
Self-Cath® Male Closed System Straight	40

External

Conveen® Optima

Male External Catheter - Silicone	
Conveen® Optima Standard	41
Conveen® Optima Short	41

Conveen® Security+

Male External Catheter - SEBS*

Conveen® Security+ Self Sealing	41
Conveen® Security+ Adhesive Liner	41

Clear Advantage® / Freedom Clear™

Male External Catheter - Silicone

Clear Advantage® with Aloe Vera	42
Freedom Clear™	42
Freedom Clear™ "LS" Long Seal	42
Freedom Clear™ "SS" Sport Sheath	42

Conveen®

Male External Catheter - Latex

Conveen® Self-Sealing	42
Conveen® Adhesive Liner	42

Active Cath®

Male External Catheter - Latex

Active Cath®	43
--------------------	----

Freedom Cath®

Male External Catheter - Latex

Freedom Clear®	43
----------------------	----

Conveen® Active

Urine Leg Bag

Conveen® Active Leg Bag, Non Sterile	43
--	----

Conveen® Security+

Urine Leg bag

Conveen® Security+ Leg Bag, Clamp Outlet, Straps, Non Sterile	43
Conveen® Security+ Leg Bag, Clamp Outlet, Sample Port, Straps, Sterile	43
Conveen® Security+ Leg Bag, Lever Outlet, Straps, Sterile	43

Conveen® Standard

Day / Night Urine Bag

Conveen® Standard Day / Night Bag, Clamp Outlet, Non Sterile	44
Conveen® Standard Day / Night Bag, Clamp Outlet, Sample Port, Sterile	44

Conveen® Contour

Urine Leg Bag

Conveen® Contour Leg Bag, Clamp Outlet, Straps, Non Sterile	44
--	----

Conveen® Night Bag

Urine Night Bag

Conveen® Bedside Night Bag, Non-Sterile	44
Conveen® Bedside Night Bag, Sterile	44

Accessories

Conveen® Bag Hanger	44
Conveen® Leg Bag Strap	44
Conveen® Catheter retaining Strap	44

Peristeen

Peristeen® Anal Plug

Peristeen® Anal Plug.....	45
---------------------------	----

Peristeen® Anal Irrigation

Peristeen® Anal Irrigation System	45
Peristeen® Anal Irrigation System - Small	45
Peristeen® Anal Irrigation Accessory unit	45
Peristeen® Anal Irrigation Accessory unit - Small	45
Peristeen® Anal Irrigation Rectal Catheter	45
Peristeen® Anal Irrigation Leg Strap	45
Peristeen® Anal Irrigation Tube	45

WOUND

Biatain® Silicone

Foam Dressing with Silicone Adhesive

Biatain® Silicone	49
-------------------------	----

Biatain® Silicone Lite

Thin Foam Dressing with Silicone Adhesive

Biatain® Silicone Lite.....	49
-----------------------------	----

Biatain® Silicone Ag

Foam Dressing with Silicone Adhesive and Silver

Biatain® Silicone Ag.....	49
---------------------------	----

Biatain®

Foam Dressing

Biatain® Adhesive	50
Biatain® Non-Adhesive	50

Biatain® Soft Hold

Foam Dressing with Gentle Adhesion

Biatain® Soft Hold.....	50
-------------------------	----

Biatain® Ag

Foam Dressing with Silver

Biatain® Ag Adhesive	51
Biatain® Ag Non-Adhesive	51

Biatain® Alginate

Calcium Alginate Dressing

Biatain® Alginate.....	51
------------------------	----

Biatain® Alginate Ag

Calcium Alginate Dressing with Silver

Biatain® Alginate Ag.....	51
---------------------------	----

Biatain® Ibu

Foam Dressing with Ibuprofen

Biatain® Ibu Non-Adhesive	52
Biatain® Ibu Soft-Hold	52

Biatain® Super

Hydrocapillary Pad Dressing

Biatain® Super Adhesive	52
Biatain® Super Non-Adhesive	52

Comfeel® Plus

Hydrocolloid Dressing

Comfeel® Plus	53
---------------------	----

Comfeel® Plus Transparent

Hydrocolloid Dressing

Comfeel® Plus Transparent.....	53
--------------------------------	----

Triad™

Hydrophilic Wound Dressing

Triad™	54
--------------	----

Purilon® Gel

Hydrogel with natural ingredients

Purilon® Gel	54
--------------------	----

Woun'Dres®

Collagen Hydrogel

Woun'Dres® Collagen Hydrogel.....	54
-----------------------------------	----

Sea-Clens®

Wound Cleanser

Sea-Clens® Wound Cleanser	54
---------------------------------	----

SKIN

Sproam®

No-Rinse Spray or Foam Cleanser

Sproam®	57
---------------	----

Bedside-Care® Sensitive Skin

No-Rinse Foam Cleanser

Bedside-Care® Sensitive Skin	57
------------------------------------	----

Bedside-Care® Foam

No-Rinse Foam Cleanser

Bedside-Care®	57
---------------------	----

Bedside-Care®

No-Rinse Spray Cleanser

Bedside-Care®	57
---------------------	----

Bedside-Care® Easi-Cleanse®

No-Rinse Wash Cloth.

Bedside-Care® Easi-Cleanse	58
----------------------------------	----

Gentle Rain® Extra Mild

Rinse Cleanser

Gentle Rain® Extra Mild	58
-------------------------------	----

Atrac-Tain®

Moisturizer

Atrac-Tain®	58
-------------------	----

Sween® 24

Moisturizer

Sween® 24	59
-----------------	----

Sween Cream®

Moisturizer

Sween Cream®	59
--------------------	----

Sween® Lotion Unscented

Moisturizer

Sween® Lotion Unscented	59
-------------------------------	----

Critic-Aid® Clear

Skin Protectant

Critic-Aid® Clear	60
-------------------------	----

Baza® Protect II

Skin Protectant

Baza® Protect II	60
------------------------	----

Baza Cleanse and Protect®

Skin Protectant

Baza Cleanse and Protect®	60
---------------------------------	----

InterDry®

Moisture-wicking textile with antimicrobial silver complex

InterDry®	61
-----------------	----

Isagel®

Hand Sanitizer Gel

Isagel®	61
---------------	----

Sween® Body Powder

All-body powder

Sween® Body Powder	61
--------------------------	----

Accessories

Coloplast® Wall Dispenser for Bag	62
---	----

Coloplast® Wall Dispenser Drip Tray.....	62
--	----

Embrace™ Dispenser Bracket.....	62
---------------------------------	----


Ostomy

SenSura® Mio Click

Ostomy Barrier and Pouch - Two Piece System

Click mechanical coupling with a lock ring for security. Provides a tight, secure fit to individual body shapes with optimal discretion. Elastic adhesive, neutral grey textile and full-circle filter. All 2pc Mio Click Flanges have 4 belt attachment points.

SenSura® Mio Click Skin Barrier					
Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	40mm (1 9/16") ●	10-35mm (3/8 - 1 3/8")	10502	5
		50mm (2") ●	10-45mm (3/8 - 1 3/4")	10512	5
		60mm (2 3/8") ●	10-55mm (3/8 - 2 1/8")	10522	5
		70mm (2 3/4") ●	10-65mm (3/8 - 2 1/2")	10531	5
Pre-Cut	Non-Convex	40mm (1 9/16") ●	20mm (3/4")	10503	5
			25mm (1")	10504	5
		50mm (2") ●	30mm (1 1/8")	10513	5
			35mm (1 3/8")	10515	5
Cut-to-Fit	Convex Light	40mm (1 9/16") ●	10-20mm (3/8 - 3/4")	16901	5
		50mm (2") ●	15-30mm (5/8 - 1 1/8")	16911	5
		60mm (2 3/8") ●	15-40mm (5/8 - 1 9/16")	16921	5
		70mm (2 3/4") ●	15-53mm (5/8 - 2 1/16")	16931	5
Pre-Cut	Convex Light	40mm (1 9/16") ●	18mm (3/4")	16903	5
		50mm (2") ●	25mm (1")	16914	5
			28mm (1 1/8")	16915	5
		60mm (2 3/8") ●	31mm (1 1/4")	16923	5
			35mm (1 3/8")	16924	5
Cut-to-Fit	Convex Deep	40mm (1 9/16") ●	10-20mm (3/8 - 3/4")	16941	5
		50mm (2") ●	15-30mm (5/8 - 1 1/8")	16951	5
		60mm (2 3/8") ●	15-40mm (5/8 - 1 9/16")	16961	5
		70mm (2 3/4") ●	15-53mm (5/8 - 2 1/16")	16971	5
Pre-Cut	Convex Deep	40mm (1 9/16") ●	18mm (3/4")	16942	5
		50mm (2") ●	25mm (1")	16954	5
			28mm (1 1/8")	16955	5
Cut-to-Fit	Concave	40mm (1 9/16") ●	10-35mm (3/8 - 1 3/8")	18500	5
		50mm (2") ●	10-45mm (3/8 - 1 3/4")	18510	5
		60mm (2 3/8") ●	10-55mm (3/8 - 2 1/8")	18520	5
		70mm (2 3/4") ●	10-65mm (3/8 - 2 1/2")	18530	5
Pre-Cut	Concave	50mm (2") ●	25-45mm (1" - 1 3/4")	18511	5
		50mm (2") ●	30-45mm (1 1/8" - 1 3/4")	18512	5
		50mm (2") ●	35-45mm (1 3/8" - 1 3/4")	18513	5
	Concave	60mm (2 3/8") ●	30-55mm (1 1/8" - 2 1/8")	18521	5
		60mm (2 3/8") ●	35-55mm (1 3/8" - 2 1/8")	18522	5
		60mm (2 3/8") ●	40-55mm (1 9/16" - 2 1/8")	18523	5
		60mm (2 3/8") ●	45-55mm (1 3/4" - 2 1/8")	18524	5
	Concave	70mm (2 3/4") ●	45-65mm (1 3/4" - 2 1/2")	18531	5
		70mm (2 3/4") ●	50-65mm (2" - 2 1/2")	18532	5



SenSura® Mio Click Drainable Pouch with Easiclese™ WIDE-Outlet

Size	Colour	Coupling Size		Product Code	Units
Midi 26cm (10 1/4")	Neutral grey with inspection window	40mm (1 9/16")		11461	10
		50mm (2")		11471	10
		40mm (1 9/16")		11462	10
		50mm (2")		11472	10
Maxi 28cm (11")	Neutral grey with inspection window	60mm (2 3/8")		11482	10
		70mm (2 3/4")		11486	10
		40mm (1 9/16")	 Narrow outlet with ileo filter	11432	10
		50mm (2")	 Narrow outlet with ileo filter	11443	10



SenSura® Mio Click Drainable Pouch with Easiclese™ WIDE-Outlet

Size	Colour	Coupling Size		Product Code	Units
Maxi 28cm (11")	Transparent	40mm (1 9/16")		11463	10
		50mm (2")		11473	10
		50mm (2")	 No filter	11474	10
		60mm (2 3/8")		11483	10
		60mm (2 3/8")	 No filter	11484	10
		70mm (2 3/4")		11487	10



SenSura® Mio Click Closed Pouch

Size	Colour	Coupling Size	Product Code	Units
Midi 18cm (7")	Neutral grey with inspection window	40mm (1 9/16")	11401	30
		50mm (2")	11412	30
		60mm (2 3/8")	11421	30
Maxi 22cm (8 1/2")	Neutral grey with inspection window	40mm (1 9/16")	11402	30
		50mm (2")	11413	30
		60mm (2 3/8")	11422	30
		70mm (2 3/4")	11425	30



SenSura® Mio Click Urostomy Pouch

Size	Colour	Coupling Size	Product Code	Units
Midi 24cm (9 1/2")	Neutral grey	50mm (2")	11496	10
		40mm (1 9/16")	11491	10
		50mm (2")	11497	10
Maxi 26cm (10 1/4")	Neutral grey with inspection window	60mm (2 3/8")	11499	10
		40mm (1 9/16")	11492	10
	Transparent	50mm (2")	11498	10
		60mm (2 3/8")	11500	10



SenSura® Mio Click High Output Pouch

Outlet	Colour	Coupling Size	Volume	Product Code	Units
Wide-lumen Hard Outlet 32cm (12 11/16")	Neutral grey with inspection window	40mm (1 9/16")	1160 mL	18600	10
		50mm (2")	1120 mL	18601	10
		60mm (2 3/8")	1080 mL	18602	10
		70mm (2 3/4")	1040 mL	18603	10
	Transparent	40mm (1 9/16")	1160 mL	18610	10
		50mm (2")	1120 mL	18611	10
		60mm (2 3/8")	1080 mL	18612	10
		70mm (2 3/4")	1040 mL	18613	10
Mini Pouch Wide-lumen Hard Outlet 24cm (9 1/2")	Neutral grey with inspection window	40mm (1 9/16")	650 mL	18620	10
		50mm (2")	581 mL	18621	10
		60mm (2 3/8")	546 mL	18622	10
Soft Outlet Universal Connector 32cm (12 11/16")	Neutral grey with inspection window	50mm (2")	1120 mL	18630	10
		60mm (2 3/8")	1080 mL	18631	10
		70mm (2 3/4")	1040 mL	18632	10
	Transparent	50mm (2")	1120 mL	18640	10
		60mm (2 3/8")	1080 mL	18641	10
		70mm (2 3/4")	1040 mL	18642	10



Pouches with
Wide Lumen
Hard Outlet
will connect
to Coloplast
drainage bag
14010

SenSura® Mio Flex

Ostomy Barrier and Pouch - Two Piece System

Flex adhesive coupling allows flexibility and freedom of movement. Provides a tight, secure fit to individual body shapes with optimal discretion. Elastic adhesive, neutral grey textile and full-circle filter. Only Mio Flex Convex Flanges have 4 belt attachment points.

SenSura® Mio Flex Skin Barrier

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	35mm (1 3/8")	10-33mm (3/8 - 1 1/4")	10551	10
		50mm (2")	10-48mm (3/8 - 2")	10561	10
		70mm (2 3/4")	10-68mm (3/8 - 2 3/4")	10571	10
Pre-Cut	Non-Convex	50mm (2")	25mm (1")	10562	10
			30mm (1 1/4")	10563	10
			35mm (1 3/8")	10564	10
Cut-to-Fit	Convex Light	35mm (1 3/8")	10-20mm (3/8 - 3/4")	16471	5
		50mm (2")	15-30mm (5/8 - 1 1/8")	16481	5
			15-40mm (5/8 - 1 9/16")	16483	5
		70mm (2 3/4")	15-53mm (5/8 - 2 1/16")	16491	5
Pre-Cut	Convex Light	35mm (1 3/8")	18mm (3/4")	16472	5
		50mm (2")	25mm (1")	16485	5
			28mm (1 1/8")	16486	5
			31mm (1 1/4")	16487	5
			35mm (1 3/8")	16488	5



SenSura® Mio Flex Drainable Pouch with Easyclose™ WIDE-Outlet

Size	Colour	Coupling Size		Product Code	Units
Midi 26cm (10 1/4")	Neutral grey with inspection window	50mm (2")	●	12271	10
		35mm (1 3/8")	●	12232	10
		35mm (1 3/8")	●	12262	10
	Neutral grey with inspection window	50mm (2")	●	12243	10
		50mm (2")	●	12272	10
		70mm (2 3/4")	●	12282	10
Maxi 28cm (11")	Transparent	50mm (2")	●	12274	10
		70mm (2 3/4")	●	12284	10



SenSura® Mio Flex Closed Pouch

Size	Colour	Coupling Size		Product Code	Units
Midi 18cm (7")	Neutral grey with inspection window	35mm (1 3/8")	●	12201	30
		50mm (2")	●	12212	30
		35mm (1 3/8")	●	12202	30
Maxi 22cm (8 1/2")	Neutral grey with inspection window	50mm (2")	●	12213	30
		70mm (2 3/4")	●	12222	30



SenSura® Mio Flex Urostomy Pouch

Size	Colour	Coupling Size		Product Code	Units
Maxi 26cm (10 1/4")	Neutral grey	35mm (1 3/8")	●	12291	10
		50mm (2")	●	12297	10
	Transparent	35mm (1 3/8")	●	12292	10
		50mm (2")	●	12298	10



SenSura® Mio Flex High Output Pouch

Outlet	Colour	Coupling Size	Volume	Product Code	Units
	Neutral grey with inspection window	50mm (2")	● 980 mL	18650	10
		70mm (2 3/4")	● 906 mL	18651	10
Wide-lumen Hard Outlet 32cm (12 11/16")		50mm (2")	● 980 mL	18655	10
	Transparent	70mm (2 3/4")	● 906 mL	18656	10



Pouches with Wide Lumen Hard Outlet will connect to Coloplast drainage bag 14010

SenSura® Mio

Ostomy Barrier and Pouch - One Piece System

Elastic adhesive, neutral grey textile and full-circle filter. Provides a tight, secure fit to individual body shapes with optimal discretion. Only 1pc Convex systems have 4 belt attachment points.

SenSura® Mio Drainable Pouch with Easyclose™ WIDE-Outlet

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 26cm (10 1/4")	Cut-to-Fit	Non-Convex	Neutral grey with inspection window	10-45mm (3/8 - 1 3/4")	10461	10
				10-55mm (3/8 - 2 1/8")	10451	10
	Cut-to-Fit	Non-Convex	Neutral grey with inspection window	Narrow outlet with ileo filter	10471	10
				10-55mm (3/8 - 2 1/8")	10481	10
			Transparent	10-55mm (3/8 - 2 1/8")	10489 No filter	10
Maxi 28cm (11")				25mm (1")	10472	10
				30mm (1 1/8")	10474	10
	Pre-Cut	Non-Convex	Neutral grey with inspection window	35mm (1 3/8")	10476	10
				40mm (1 9/16")	10479	10
				45mm (1 3/4")	10477	10
				50mm (2")	10478	10



SenSura® Mio Drainable Pouch with Easyclose™ WIDE-Outlet

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 28cm (11")	Pre-Cut	Non-Convex	Transparent	25mm (1")	10482	10
				30mm (1 1/8")	10483	10
				35mm (1 3/8")	10484	10
				40mm (1 9/16")	10487	10
Maxi 28cm (11")	Cut-to-Fit	Convex Soft	Neutral grey with inspection window	15-33mm (5/8 - 1 1/4")	16705	10
				10-50mm (3/8 - 2")	16706	10
			Transparent	15-33mm (5/8 - 1 1/4")	16715	10
				10-50mm (3/8 - 2")	16716	10
	Cut-to-Fit	Convex Light	Neutral grey with inspection window	10-23mm (3/8 - 7/8")	16725	10
				15-33mm (5/8 - 1 1/4")	16726	10
				10-43mm (3/8 - 1 3/4")	16733	10
				10-23mm (3/8 - 7/8")	16737	10
			Transparent	15-33mm (5/8 - 1 1/4")	16738	10
				10-43mm (3/8 - 1 3/4")	16739	10
				10-43mm (3/8 - 1 3/4")	16748	No filter 10
				25mm (1")	16731	10
	Pre-Cut	Convex Light	Neutral grey with inspection window	28mm (1 1/8")	16732	10
				31mm (1 1/4")	16727	10
				35mm (1 3/8")	16734	10
				10-23mm (3/8 - 7/8")	16755	10
	Cut-to-Fit	Convex Deep	Neutral grey with inspection window	15-33mm (5/8 - 1 1/4")	16756	10
				10-43mm (3/8 - 1 3/4")	16768	10
			Transparent	10-43mm (3/8 - 1 3/4")	16767	10
				10-43mm (3/8 - 1 3/4")	16778	No filter 10
	Pre-Cut	Convex Deep	Neutral grey with inspection window	21mm (7/8")	16760	10
				25mm (1")	16761	10
				28mm (1 1/8")	16762	10
				31mm (1 1/4")	16757	10
			Transparent	35mm (1 3/8")	16763	10
				25mm (1")	16771	10
				28mm (1 1/8")	16772	10
				31mm (1 1/4")	16773	10
Midi 26cm (10 1/4")	Cut-to-fit	Concave	Neutral grey with inspection window	10-40mm (3/8 - 1 9/16")	18300	10
				10-50mm (3/8 - 2")	18310	10
Maxi 28cm (11")	Pre-Cut	Concave	Neutral grey with inspection window	25-40mm (1 - 1 9/16")	18301	10
				30-40mm (1 1/8 - 1 9/16")	18302	10
	Cut-to-fit	Concave	Transparent	10-40mm (3/8 - 1 9/16")	18305	10
				10-50mm (3/8 - 2")	18312	10
Maxi 28cm (11")	Pre-Cut	Concave	Neutral grey with inspection window	25-50mm (1 - 2")	18313	10
				30-50mm (1 1/8 - 2")	18314	10
				35-50mm (1 3/8 - 2")	18315	10
	Cut-to-fit	Concave	Transparent	10-50mm (3/8 - 2")	No filter 18320	10
				10-50mm (3/8 - 2")	18321	10
				25-50mm (1 - 2")	18322	10
				30-50mm (1 1/8 - 2")	18323	10
	Pre-Cut	Concave	Transparent	35-50mm (1 3/8 - 2")	18324	10



SenSura® Mio Drainable Pouch with Easyclose™ WIDE-Outlet

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 28cm (11")	Cut-to-fit	Concave	Neutral grey with inspection window	10-60mm ($\frac{2}{8}$ - $2\frac{3}{8}$ "	18330	10
	Pre-Cut	Concave		40-60mm ($1\frac{9}{16}$ - $2\frac{3}{8}$ "	18331	10
				45-60mm ($1\frac{3}{4}$ - $2\frac{3}{8}$ "	18332	10
				50-60mm ($2''$ - $2\frac{3}{8}$ "	18333	10
	Cut-to-fit	Concave	Transparent	10-60mm ($\frac{3}{8}$ - $2\frac{3}{8}$ " <i>No filter</i>	18340	10
				10-60mm ($\frac{2}{8}$ - $2\frac{3}{8}$ "	18341	10
	Pre-Cut	Concave		40-60mm ($1\frac{9}{16}$ - $2\frac{3}{8}$ "	18342	10
				45-60mm ($1\frac{3}{4}$ - $2\frac{3}{8}$ "	18343	10
				50-60mm ($2''$ - $2\frac{3}{8}$ "	18344	10



SenSura® Mio Closed Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Mini 14cm	Cut-to-fit	Non-Convex	Neutral grey with inspection window	15-45mm ($\frac{5}{8}$ - $1\frac{3}{4}$ ")	10870	30
Midi 18cm (7")	Cut-to-fit	Non-Convex	Neutral grey with inspection window	15-45mm ($\frac{5}{8}$ - $1\frac{3}{4}$ ")	10871	30
	Pre-Cut	Non-Convex	Neutral grey with inspection window	25mm (1")	10872	30
				30mm ($1\frac{1}{8}$ ")	10873	30
				35mm ($1\frac{3}{8}$ ")	10874	30
Maxi 22cm ($8\frac{1}{2}$ ")	Cut-to-Fit	Non-Convex	Neutral grey with inspection window	15-45mm ($\frac{5}{8}$ - $1\frac{3}{4}$ ")	10881	30
	Pre-Cut	Non-Convex	Neutral grey with inspection window	40-55mm ($1\frac{9}{16}$ - $2\frac{1}{8}$ ")	10885	30
				30mm ($1\frac{1}{8}$ ")	10883	30
				35mm ($1\frac{3}{8}$ ")	10884	30
Midi 18cm (7")	Cut-to-Fit	Convex Soft	Neutral grey with inspection window	15-33mm ($\frac{5}{8}$ - $1\frac{1}{4}$ ")	16301	10
	Cut-to-Fit	Convex Light	Neutral grey with inspection window	15-33mm ($\frac{5}{8}$ - $1\frac{1}{4}$ ")	16322	10
	Pre-Cut	Convex Light	Neutral grey with inspection window	31mm ($1\frac{1}{4}$ ")	16323	10
	Cut-to-Fit	Convex Soft	Neutral grey with inspection window	15-33mm ($\frac{5}{8}$ - $1\frac{1}{4}$ ")	16305	10
			Neutral grey with inspection window	10-50mm ($\frac{3}{8}$ - 2 ")	16306	10
			Transparent	10-50mm ($\frac{3}{8}$ - 2 ")	16315	10
			Transparent	25mm (1")	16310	10
	Pre-Cut	Convex Soft	Neutral grey with inspection window	28mm ($1\frac{1}{8}$ ")	16311	10
				31mm ($1\frac{1}{4}$ ")	16312	10
				35mm ($1\frac{3}{8}$ ")	16313	10
				10-23mm ($\frac{3}{8}$ - $\frac{7}{8}$ ")	16325	10
	Cut-to-Fit	Convex Light	Neutral grey with inspection window	15-33mm ($\frac{5}{8}$ - $1\frac{1}{4}$ ")	16326	10
			Neutral grey with inspection window	10-43mm ($\frac{3}{8}$ - $1\frac{3}{4}$ ")	16332	10
			Transparent	10-43mm ($\frac{3}{8}$ - $1\frac{3}{4}$ ")	16335	10
			Transparent	25mm (1")	16330	10
	Pre-Cut	Convex Light	Neutral grey with inspection window	28mm ($1\frac{1}{8}$ ")	16331	10
				31mm ($1\frac{1}{4}$ ")	16327	10
				35mm ($1\frac{3}{8}$ ")	16333	10
				10-23mm ($\frac{3}{8}$ - $\frac{7}{8}$ ")	16345	10
Maxi 22cm ($8\frac{1}{2}$ ")	Cut-to-Fit	Convex Deep	Neutral grey with inspection window	15-33mm ($\frac{5}{8}$ - $1\frac{1}{4}$ ")	16346	10
			Neutral grey with inspection window	10-43mm ($\frac{3}{8}$ - $1\frac{3}{4}$ ")	16352	10
			Transparent	25mm (1")	16350	10
			Transparent	28mm ($1\frac{1}{8}$ ")	16351	10
	Pre-Cut	Convex Deep	Neutral grey with inspection window	31mm ($1\frac{1}{4}$ ")	16347	10
				35mm ($1\frac{3}{8}$ ")	16353	10
				35mm ($1\frac{3}{8}$ ")	16353	10



SenSura® Mio Closed Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 18cm (7")	Cut-to-Fit	Concave	Neutral grey with inspection window	10-40mm ($\frac{3}{8}$ - 1 $\frac{9}{16}$ "	18151	10
				10-50mm ($\frac{3}{8}$ - 2")	18160	10
	Pre-Cut	Concave		25-50mm (1 - 2")	18161	10
				30-50mm (1 $\frac{1}{8}$ " - 2")	18162	10
				35-50mm (1 $\frac{3}{8}$ " - 2")	18163	10
Maxi 22cm (8 $\frac{1}{2}$ ")	Cut-to-Fit	Concave	Neutral grey with inspection window	10-40mm ($\frac{3}{8}$ - 1 $\frac{9}{16}$ "	18152	10
Maxi 22cm (8 $\frac{1}{2}$ ")	Cut-to-Fit	Concave	Neutral grey with inspection window	10-50mm ($\frac{3}{8}$ - 2")	18170	10
	Pre-Cut	Concave		25-50mm (1 - 2")	18171	10
				30-50mm (1 $\frac{1}{8}$ " - 2")	18172	10
				35-50mm (1 $\frac{3}{8}$ " - 2")	18173	10
	Maxi 22cm (8 $\frac{1}{2}$ ")	Cut-to-Fit		Concave	Neutral grey with inspection window	10-60mm (3/8 - 2 $\frac{3}{8}$ ")
Pre-Cut		Concave	40-60mm (1 $\frac{9}{16}$ - 2 $\frac{3}{8}$ ")	18181		10
			45-60mm (1 $\frac{3}{4}$ - 2 $\frac{3}{8}$ ")	18182		10
		Cut-to-Fit	Concave	Transparent	10-60mm ($\frac{3}{8}$ - 2 $\frac{3}{8}$ ")	18185



SenSura® Mio One Piece Urostomy

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 26cm (10 $\frac{1}{4}$ ")	Cut-to-Fit	Non-Convex	Neutral grey	10-45mm ($\frac{3}{8}$ - $1\frac{3}{4}$ ")	10582	10
			Transparent	10-45mm ($\frac{3}{8}$ - $1\frac{3}{4}$ ")	10586	10
				20mm ($\frac{3}{4}$ ")	10583	10
	Pre-Cut	Non-Convex	Neutral grey	25mm (1")	10584	10
				30mm ($1\frac{1}{8}$ ")	10585	10
Maxi 26cm (10 $\frac{1}{4}$ ")	Cut-to-Fit	Convex Soft	Neutral grey	15-33mm ($\frac{5}{8}$ - $1\frac{1}{4}$ ")	16805	10
			Transparent	10-50mm ($\frac{3}{8}$ - 2")	16810	10
		Convex Light	Neutral grey	10-23mm ($\frac{3}{8}$ - $\frac{7}{8}$ ")	16825	10
				15-33mm ($\frac{5}{8}$ - $1\frac{1}{4}$ ")	16826	10
		Convex Light	Transparent	10-23mm ($\frac{3}{8}$ - $\frac{7}{8}$ ")	16835	10
				15-33mm ($\frac{5}{8}$ - $1\frac{1}{4}$ ")	16836	10
	Pre-Cut	Convex Light	Neutral grey	10-43mm ($\frac{3}{8}$ - $1\frac{3}{4}$ ")	16837	10
				31mm ($1\frac{1}{4}$ ")	16827	10
				21mm ($\frac{7}{8}$ ")	16840	10
		Convex Deep	Neutral grey	25mm (1")	16841	10
				28mm ($1\frac{1}{8}$ ")	16842	10
				10-23mm ($\frac{3}{8}$ - $\frac{7}{8}$ ")	16855	10
	Cut-to-Fit	Convex Deep	Neutral grey	15-33mm ($\frac{5}{8}$ - $1\frac{1}{4}$ ")	16856	10
				10-23mm ($\frac{3}{8}$ - $\frac{7}{8}$ ")	16865	10
				15-33mm ($\frac{5}{8}$ - $1\frac{1}{4}$ ")	16866	10
		Convex Deep	Transparent	10-43mm ($\frac{3}{8}$ - $1\frac{3}{4}$ ")	16867	10
				31mm ($1\frac{1}{4}$ ")	16857	10
				21mm ($\frac{7}{8}$ ")	16870	10
	Pre-Cut	Convex Deep	Transparent	25mm (1")	16871	10
				28mm ($1\frac{1}{8}$ ")	16872	10
				10-40mm ($\frac{3}{8}$ - $1\frac{9}{16}$ ")	18401	10
Maxi 26cm (10 $\frac{1}{4}$ ")	Pre-Cut	Concave	Neutral grey	20-40mm ($\frac{3}{4}$ " - $1\frac{9}{16}$ ")	18402	10
				25-40mm (1 " - $1\frac{9}{16}$ ")	18403	10
				30-40mm ($1\frac{1}{8}$ " - $1\frac{9}{16}$ ")	18404	10
	Cut-to-Fit	Concave	Transparent	10-40mm ($\frac{3}{8}$ - $1\frac{9}{16}$ ")	18405	10



SenSura® Mio One Piece High Output Pouch

Cut/Outlet	Convexity	Colour	Volume	Stoma Size	Product Code	Units
Cut-to-fit/ Soft Outlet Universal Connector 28cm (11 3/16")	Non-Convex	Neutral grey with inspection window	630 mL	10-70mm (3/8 - 2 3/4")	18660	10
			725 mL	10-100mm (3/8 - 4")	18661	10
		Transparent	630 mL	10-70mm (3/8 - 2 3/4")	18665	10
			725 mL	10-100mm (3/8 - 4")	18666	10
	Convex Soft	Neutral grey with inspection window	630 mL	10-50mm (3/8 - 2")	18670	10
		Transparent	630 mL	10-50mm (3/8 - 2")	18671	10
	Convex Light	Neutral grey with inspection window	630 mL	10-43mm (3/8 - 1 3/4")	18675	10
		Transparent	630 mL	10-43mm (3/8 - 1 3/4")	18676	10



Does not connect to Coloplast drainage bag 14010

SenSura® Mio One Piece Post-Op Pouch

Size	Cut	Convexity	Volume	Stoma Size	Product Code	Units
28cm (10 15/16")	Cut-to-Fit	Transparent	680 mL	10-70mm (3/8 - 2 3/4")	18680 Sterile	5
29cm (11 1/2")	Cut-to-Fit	Transparent	725 mL	10-100mm (3/8 - 4")	18681 Sterile	5
28cm (10 15/16")	Cut-to-Fit	Transparent	680 mL	10-70mm (3/8 - 2 3/4")	18690	5
29cm (11 1/2")	Cut-to-Fit	Transparent	725 mL	10-100mm (3/8 - 4")	18691	5



SenSura® Mio Baby

Ostomy Barrier and Pouch - Two Piece System

SenSura® Mio Baby Flex Skin Barrier

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	35mm (1 3/8")	0-35mm	18705	5



SenSura® Mio Baby Flex Drainable Pouch

Size	Colour	Coupling Size	Product Code	Units
20cm (7 3/4")	Transparent	35mm (1 3/8")	18700	10



SenSura® Mio Baby Adhesive Remover Wipe

	Product Code	Units
Single use	12014	30



SenSura[®] Mio Kids

Ostomy Barrier and Pouch - Two Piece System

SenSura® Mio Kids Flex Skin Barrier						
Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units	
Cut-to-Fit	Non-Convex	27 mm	0-27 mm	18735	5	
		35 mm	0-35 mm	18747	5	
SenSura® Mio Kids Flex Drainable Pouch						
Cut	Colour	Coupling Size		Product Code	Units	
Mini 16cm (6 1⁄8")	Neutral grey	27 mm		18720	10	
	Transparent	27 mm		18721	10	
Midi 17cm (6 3⁄4")	Neutral grey	27 mm		18722	10	
	Transparent	27 mm		18723	10	
	Neutral grey	35 mm		18724	10	
	Transparent	35 mm		18725	10	
SenSura® Mio Kids Flex Urostomy Pouch						
Cut	Colour	Coupling Size	Stoma Size	Product Code	Units	
Midi 17cm (6 3⁄4")	Neutral grey	27 mm		18730	10	
		35 mm		18731	10	

Ostomy Barrier and Pouch - One Piece System

SenSura® Mio Kids Drainable Pouch 1 Piece							
Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units	
Midi 17cm (6 3/4")	Cut-to-fit	Non-Convex	Neutral grey	0-35mm	18710	10	
			Transparent	0-35mm	18711	10	
SenSura® Mio Kids Pouch 1 Piece Urostomy							
Cut		Convexity	Colour	Stoma Size	Product Code	Units	
Cut-to-fit 19cm (7 1/4")		Non-Convex	Transparent	10-35mm	18712	10	
			Neutral grey	10-35mm	18713	10	

SenSura® Click

Ostomy Barrier and Pouch- Two Piece System

Click mechanical coupling with a lock ring for security. Provides a tight, secure fit to individual body shapes. All 2pc SenSura Click Flanges have 2 belt attachment points.

SenSura® Click Skin Barrier					
Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	40mm (1 9/16")	10-35mm (3/8 - 1 3/8")	10011	5
		50mm (2")	10-45mm (3/8 - 1 3/4")	10021	5
		60mm (2 3/8")	10-55mm (3/8 - 2 1/8")	10031	5
		70mm (2 3/4")	10-65mm (3/8 - 2 1/2")	10041	5
	Convex Light	40mm (1 9/16")	15-23mm (5/8 - 7/8")	11011	5
		50mm (2")	15-33mm (5/8 - 1 1/4")	11021	5
		60mm (2 3/8")	15-43mm (5/8 - 1 3/4")	11031	5
		70mm (2 3/4")	15-53mm (5/8 - 2 1/16")	11041	5
Pre-Cut	Non-Convex	40mm (1 9/16")	20mm (3/4")	10012	5
			25mm (1")	10013	5
		50mm (2")	30mm (1 1/8")	10022	5
			35mm (1 3/8")	10023	5
		60mm (2 3/8")	40mm (1 9/16")	10032	5
		40mm (1 9/16")	21mm (7/8")	11012	5
	Convex Light	50mm (2")	25mm (1")	11022	5
			28mm (1 1/8")	11023	5
		60mm (2 3/8")	31mm (1 1/4")	11024	5
			35mm (1 3/8")	11032	5
		40mm (1 9/16")	20mm (3/4")	10016	5
			25mm (1")	10017	5
Cut-to-Fit	Non-Convex	40mm (1 9/16")	10-35mm (3/8 - 1 3/8")	10015	5
		50mm (2")	10-45mm (3/8 - 1 3/4")	10025	5
		60mm (2 3/8")	10-55mm (3/8 - 2 1/8")	10035	5
		70mm (2 3/4")	10-65mm (3/8 - 2 1/2")	10045	5
	Convex Light	40mm (1 9/16")	15-23mm (5/8 - 7/8")	11015	5
		50mm (2")	15-33mm (5/8 - 1 1/4")	11025	5
		60mm (2 3/8")	15-43mm (5/8 - 1 3/4")	11035	5
		70mm (2 3/4")	15-53mm (5/8 - 2 1/16")	11045	5
	Non-Convex	40mm (1 9/16")	20mm (3/4")	10016	5
			25mm (1")	10017	5
		50mm (2")	30mm (1 1/8")	10026	5
			35mm (1 3/8")	10027	5
		60mm (2 3/8")	40mm (1 9/16")	10036	5
		40mm (1 9/16")	21mm (7/8")	11016	5
Pre-Cut	Convex Light	50mm (2")	25mm (1")	11026	5
			28mm (1 1/8")	11027	5
		60mm (2 3/8")	31mm (1 1/4")	11028	5
			35mm (1 3/8")	11036	5
		40mm (1 9/16")	20mm (3/4")	10016	5
			25mm (1")	10017	5



SenSura® Click Drainable Pouch with Easyclose™ WIDE-Outlet

Pouch Size	Pouch Colour	Coupling Size	Product Code	Units
Midi 27cm (10 1/2")	Opaque	40mm (1 9/16") ●	11114	10
		50mm (2") ●	11115	10
Maxi 30cm (12")	Opaque	40mm (1 9/16") ●	11124	10
		50mm (2") ●	11125	10
		50mm (2") ● No filter	11135	10
		60mm (2 3/8") ●	11126	10
		60mm (2 3/8") ● No filter	11136	10
		70mm -2 3/4" ●	11127	10
	Transparent	40mm (1 9/16") ●	11194	10
		50mm (2") ●	11195	10
		50mm (2") ● No filter	11185	10
		60mm (2 3/8") ●	11196	10
		60mm (2 3/8") ● No filter	11186	10
		70mm -2 3/4" ●	11197	10



SenSura® Click Closed Pouch

Pouch Size	Pouch Colour	Coupling Size	Product Code	Units
Midi 18cm (7")	Opaque	40mm (1 9/16") ●	10154	30
		50mm (2") ●	10155	30
		60mm (2 3/8") ●	10156	30
Maxi 22cm (8 1/2")	Opaque	40mm (1 9/16") ●	10164	30
		50mm (2") ●	10165	30
		60mm (2 3/8") ●	10166	30
		70mm (2 3/4") ●	10167	30



SenSura® Click Urostomy Pouch

Pouch Size	Pouch Colour	Coupling Size	Product Code	Units
Midi 24cm (9 1/4")	Opaque	40mm (1 9/16") ●	11841	10
		50mm (2") ●	11842	10
	Transparent	40mm (1 9/16") ●	11851	10
		50mm (2") ●	11852	10
Maxi 26cm (10 1/4")	Opaque	40mm (1 9/16") ●	11844	10
		50mm (2") ●	11845	10
		60mm (2 3/8") ●	11846	10
	Transparent	40mm (1 9/16") ●	11854	10
		50mm (2") ●	11855	10
		60mm (2 3/8") ●	11856	10



SenSura® High Output Pouch

Colour	Coupling Size	Volume	Product Code	Units
Opaque	40mm (1 9/16") ●	1160 mL	19040	10
	50mm (2") ●	1120 mL	19041	10
	60mm (2 3/8") ●	1080 mL	19042	10
Transparent	40mm (1 9/16") ●	1160 mL	19050	10
	50mm (2") ●	1120 mL	19051	10
	60mm (2 3/8") ●	1080 mL	19052	10



Will connect to
Coloplast Drainage
Bag 14010

SenSura® Flex

Ostomy Barrier and Pouch - Two Piece System

Low profile adhesive coupling system that provide flexibility and freedom of movement. Double layer adhesive for secure, safe, long lasting wear. Drainable pouch with velcro closure system. Only SenSura Flex Convex Flanges have 2 belt attachment points.

SenSura® Flex Skin Barrier						
Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units	
Cut-to-Fit	Non-Convex	35mm (1 3⁄8")	10-33mm (3⁄8- 1 1⁄4")	10101	10	
		50mm (2")	10-48mm (3⁄8- 2")	10102	10	
		70mm (2 3⁄4")	10-68mm (3⁄8- 2 3⁄4")	10103	10	
		90mm - 3 1⁄2"	10-88mm (3⁄8- 3 1⁄2")	10108	5	
	Convex light	35mm (1 3⁄8")	15-23mm (5⁄8- 7⁄8")	11301	5	
		50mm (2")	15-33mm (5⁄8- 1 1⁄4")	11302	5	
		70mm (2 3⁄4")	15-43mm (5⁄8- 1 3⁄4")	11303	5	
Pre-Cut	Non-Convex	35mm (1 3⁄8")	21mm (7⁄8")	10111	10	
			28mm (1 1⁄8")	10113	10	
			25mm (1")	10114	10	
		50mm (2")	31mm (1 1⁄4")	10116	10	
	Convex light		35mm (1 3⁄8")	10117	10	
		35mm (1 3⁄8")	21mm (7⁄8")	11312	5	
			25mm (1")	11313	5	
		50mm (2")	28mm (1 1⁄8")	11314	5	
			31mm (1 1⁄4")	11315	5	
			35mm (1 3⁄8")	11316	5	

SenSura® Flex Xpro Skin Barrier						
Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units	
Cut-to-Fit	Non-Convex	35mm (1 3⁄8")	10-33mm (3⁄8- 1 1⁄4")	10104	10	
		50mm (2")	10-48mm (3⁄8- 2")	10105	10	
		70mm (2 3⁄4")	10-68mm (3⁄8- 2 3⁄4")	10106	10	
	Convex light	35mm (1 3⁄8")	15-23mm (5⁄8- 7⁄8")	11304	5	
		50mm (2")	15-33mm (5⁄8- 1 1⁄4")	11305	5	
			15-43mm (5⁄8- 1 3⁄4")	11306	5	
		70mm (2 3⁄4")	15-56mm (5⁄8- 2 1⁄8")	11308	5	

SenSura® Flex Drainable Pouch with Easyclose™ WIDE-Outlet						
Size	Colour	Coupling Size		Product Code	Units	
Midi 27cm (10 1⁄2")	Opaque	50mm (2")		11505	20	
Maxi 30cm (12")	Opaque	35mm (1 3⁄8")		11514	20	
		50mm (2")		11515	20	
		50mm (2")	No filter	11615	20	
		70mm (2 3⁄4")		11518	20	
	Transparent	50mm (2")	No filter	11612	20	
		90mm (3 1⁄2")	No filter	11690	20	

SenSura® Flex Closed Pouch						
Size	Colour	Coupling Size		Product Code	Units	
Midi 18cm (7")	Opaque	35mm (1 3⁄8")		10904	30	
		50mm (2")		10905	30	
Maxi 22cm (8 1⁄2")	Opaque	35mm (1 3⁄8")		10914	30	
		50mm (2")		10915	30	
		70mm (2 3⁄4")		10918	30	

SenSura® Flex Urostomy Pouch

Size	Colour	Coupling Size	Product Code	Units
Midi 24mm (9 1/2")	Opaque	35mm (1 3/8") ●	11821	20
		50mm (2") ●	11822	20
Maxi 26cm (10 1/4")	Opaque	35mm (1 3/8") ●	11824	20
		50mm (2") ●	11825	20
	Transparent	35mm (1 3/8") ●	11827	20
		50mm (2") ●	11828	20



SenSura®

Ostomy Barrier and Pouch - One Piece System

Easiclese™ system with Velcro® to ensure security. Double layer adhesive provides secure adhesion, flexibility and ease of removal. Only SenSura 1pc Convex Systems have 2 belt attachment points.

SenSura® Drainable Pouch with Easiclese™ WIDE-Outlet

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 30cm (12")	Cut-to-Fit	Non-Convex	Opaque	10-76mm (3⁄8-3")	15531	10
			Transparent	10-76mm (3⁄8-3")	15521	10
				10-76mm (3⁄8-3")	15821	10
		Convex Light	Opaque	15-33mm (5⁄8 - 1 1⁄4")	15607	10
				15-43mm (5⁄8 - 1 3⁄4")	15608	10
			Transparent	15-33mm (5⁄8 - 1 1⁄4")	15605	10
				15-43mm (5⁄8 - 1 3⁄4")	15606	10
		Pre-Cut	Non-Convex	Opaque	25mm (1")	15532
	30mm (1 1⁄8")				15533	10
	35mm (1 3⁄8")				15534	10
	Convex Light		Opaque	25mm (1")	15633	10
				28mm (1 1⁄8")	15634	10
				31mm (1 1⁄4")	15635	10
	Non-Convex		Transparent	25mm (1")	15522	10
				30mm (1 1⁄8")	15523	10
		35mm (1 3⁄8")		15524	10	



SenSura® Xpro Drainable Pouch with Easiclese™ WIDE-Outlet

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 30cm (12")	Cut-to-Fit	Non-Convex	Transparent	10-76mm (3/8- 3")	15981	20
		Convex Light	Transparent	15-33mm (5/8 - 1 1/4")	15695	10
				15-43mm (5/8 - 1 3/4")	15696	10



SenSura® Mio Original Closed Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Mini 14cm (5 1/2")	Cut-to-Fit	Non-Convex	Opaque with inspection window	15-45mm (5/8 - 1 3/4")	10811	30
Midi 18cm (7")	Cut-to-Fit	Non-Convex	Opaque with inspection window	15-45mm (5/8 - 1 3/4")	10821	30
	Pre-Cut	Non-Convex	Opaque with inspection window	25mm (1")	10822	30
				30mm (1 1/8")	10823	30
				35mm (1 3/8")	10824	30
Maxi 22cm (8 1/2")	Cut-to-Fit	Non-Convex	Opaque with inspection window	15-45mm (5/8 - 1 3/4")	10831	30
				40-55mm (1 9/16 - 2 1/8")	10835	30
	Pre-Cut	Non-Convex	Opaque with inspection window	30mm (1 1/8")	10833	30
				35mm (1 3/8")	10834	30



SenSura® Mio Original Closed Pouches will be discontinued in 2021 and will be replaced with SenSura® Mio Closed Pouches (pg.14)

SenSura® Closed Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 18cm (7")	Cut-to-Fit	Non-Convex	Opaque	10-66mm ($\frac{3}{8}$ - 2 $\frac{5}{8}$ ")	15450	30
			Transparent	10-66mm ($\frac{3}{8}$ - 2 $\frac{5}{8}$ ")	15441	30
		Convex Light	Opaque	15-33mm ($\frac{5}{8}$ - 1 $\frac{1}{4}$ ")	15403	10
				15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	15404	10
	Pre-Cut	Convex Light	Opaque	25mm (1")	15413	10
				28mm (1 $\frac{1}{8}$ ")	15414	10
				31mm (1 $\frac{1}{4}$ ")	15415	10
Maxi 22cm (8 $\frac{1}{2}$ ")	Cut-to-Fit	Non-Convex	Opaque	10-76mm ($\frac{3}{8}$ - 3")	15480	30
			Transparent	10-76mm ($\frac{3}{8}$ - 3")	15470	30
		Convex Light	Opaque	15-33mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	15407	10
				15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	15408	10
			Transparent	15-33mm ($\frac{5}{8}$ - 1 $\frac{1}{4}$ ")	15405	10
				15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	15406	10
	Pre-Cut	Non-Convex	Opaque	25mm (1")	15482	30
				30mm (1 $\frac{1}{8}$ ")	15483	30
				35mm (1 $\frac{3}{8}$ ")	15484	30
			Transparent	25mm (1")	15472	30
				30mm (1 $\frac{1}{8}$ ")	15473	30
				35mm (1 $\frac{3}{8}$ ")	15474	30



SenSura® Urostomy Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 24cm (9 $\frac{1}{2}$ ")	Cut-to-Fit	Non-Convex	Opaque	10-66mm ($\frac{3}{8}$ - 2 $\frac{5}{8}$ ")	11801	10
			Transparent	10-66mm ($\frac{3}{8}$ - 2 $\frac{5}{8}$ ")	11803	10
		Non-Convex	Opaque	10-76mm ($\frac{3}{8}$ - 3")	11802	10
			Transparent	10-76mm ($\frac{3}{8}$ - 3")	11804	10
Maxi 26cm (10 $\frac{1}{4}$ ")	Cut-to-Fit	Non-Convex	Opaque	15-33mm ($\frac{5}{8}$ - 1 $\frac{1}{4}$ ")	11811	10
				15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	11812	10
		Convex Light	Transparent	15-33mm ($\frac{5}{8}$ - 1 $\frac{1}{4}$ ")	11814	10
				15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	11815	10



SenSura® One-Piece High Output Pouch

Cut	Sterile	Stoma Size	Volume	Product Code	Units
Cut-to-Fit	Sterile	10-76mm ($\frac{3}{8}$ -3")	630 mL	19001	5
	Non-Sterile	10-76mm ($\frac{3}{8}$ -3")	630 mL	19003	6



SenSura® Post-Op Pouch


Cut	Sterile	Stoma Size	Window Size	Volume	Product Code	Units
Cut-to-Fit	Sterile	10-76mm ($\frac{3}{8}$ -3")	70mm	680 mL	19010	5
		10-115mm ($\frac{3}{8}$ - 4 $\frac{1}{2}$ ")	100mm	725 mL	19011	5
	Non-Sterile	10-76mm ($\frac{3}{8}$ -3")	70mm	680 mL	19020	6
		10-115mm ($\frac{3}{8}$ - 4 $\frac{1}{2}$ ")	100mm	725 mL	19021	6
		10-76mm ($\frac{3}{8}$ -3")	No window	680 mL	19030	6



Ostomy Barrier and Pouch - Two Piece System


Click mechanical coupling with a lock ring for security. Skin barrier with spiral technology for secure, long lasting wear. All 2pc Assura Click Flanges have 2 belt attachment points.

Assura® Skin Barrier					
Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
		50mm (2")	● 10-45mm ($\frac{3}{8}$ - $1\frac{3}{4}$ ")	2882	5
		60mm ($2\frac{3}{8}$ ")	● 10-55mm ($\frac{3}{8}$ - $2\frac{1}{8}$ ")	2883	5
	Convex Light	40mm ($1\frac{9}{16}$ ")	● 15-23mm ($\frac{5}{8}$ - $\frac{7}{8}$ ")	14261	5
		50mm (2")	● 15-33mm ($\frac{5}{8}$ - $1\frac{1}{4}$ ")	14262	5
		60mm ($2\frac{3}{8}$ ")	● 15-43mm ($\frac{5}{8}$ - $1\frac{3}{4}$ ")	14263	5
	Convex Deep	50mm (2")	● 15-33mm ($\frac{5}{8}$ - $1\frac{1}{4}$ ")	12716	5
		60mm ($2\frac{3}{8}$ ")	● 15-43mm ($\frac{5}{8}$ - $1\frac{3}{4}$ ")	12719	5
	Non-Convex	40mm ($1\frac{9}{16}$ ")	● 20mm ($\frac{3}{4}$ ")	12841	5
			25mm (1")	12842	5
		50mm (2")	● 30mm ($1\frac{1}{8}$ ")	12844	5
			35mm ($1\frac{3}{8}$ ")	12845	5
		60mm ($2\frac{3}{8}$ ")	● 45mm ($1\frac{3}{4}$ ")	12847	5
		60mm ($2\frac{3}{8}$ ")	● 45mm ($1\frac{3}{4}$ ")	12848	5
		40mm ($1\frac{9}{16}$ ")	● 18mm ($\frac{3}{4}$ ")	14271	5
			21mm ($\frac{7}{8}$ ")	14272	5
		50mm (2")	● 25mm (1")	14273	5
			28mm ($1\frac{1}{8}$ ")	14274	5
Pre-Cut	Convex Light		31mm ($1\frac{1}{4}$ ")	14275	5
			35mm ($1\frac{3}{8}$ ")	14276	5
		60mm ($2\frac{3}{8}$ ")	● 38mm ($1\frac{1}{2}$ ")	14277	5
	Convex Deep		41mm ($1\frac{5}{8}$ ")	14278	5
			15mm ($\frac{5}{8}$ ")	12701	5
		40mm ($1\frac{9}{16}$ ")	● 18mm ($\frac{3}{4}$ ")	12702	5
			21mm ($\frac{7}{8}$ ")	12703	5
			25mm (1")	12704	5
		50mm (2")	● 28mm ($1\frac{1}{8}$ ")	12705	5
			31mm ($1\frac{1}{4}$ ")	12706	5
			35mm ($1\frac{3}{8}$ ")	12707	5
		60mm ($2\frac{3}{8}$ ")	● 38mm ($1\frac{1}{2}$ ")	12708	5
			41mm ($1\frac{5}{8}$ ")	12709	5



Assura® Skin Barrier - Extra Extended Wear

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	40mm (1 9/16")	10-35mm (3/8 - 1 3/8")	2831	5
		50mm (2")	10-45mm (3/8 - 1 3/4")	2832	5
		60mm (2 3/8")	10-55mm (3/8 - 2 1/8")	2833	5
	Convex Light	40mm (1 9/16")	15-23mm (5/8 - 7/8")	14281	5
		50mm (2")	15-33mm (5/8 - 1 1/4")	14282	5
		60mm (2 3/8")	15-43mm (5/8 - 1 3/4")	14283	5
	Convex Deep	40mm - (1 9/16")	15-23mm (5/8 - 7/8")	14243	5
		50mm (2")	15-33mm (5/8 - 1 1/4")	14246	5
		60mm (2 3/8")	15-43mm (5/8 - 1 3/4")	14249	5
Pre-Cut	Non-Convex	40mm (1 9/16")	20mm (3/4")	14251	5
			25mm (1")	14252	5
		50mm (2")	30mm (1 1/8")	14254	5
			35mm (1 3/8")	14255	5
		60mm (2 3/8")	45mm - 1 3/4"	14258	5
	Convex Light	40mm (1 9/16")	18mm (3/4")	14291	5
			21mm (7/8")	14292	5
		50mm (2")	25mm (1")	14293	5
			28mm (1 1/8")	14294	5
			31mm (1 1/4")	14295	5
		60mm (2 3/8")	35mm (1 3/8")	14296	5
	Convex Deep	40mm (1 9/16")	18mm (3/4")	14232	5
			21mm (7/8")	14233	5
		50mm (2")	25mm (1")	14234	5
			28mm (1 1/8")	14235	5
			31mm (1 1/4")	14236	5
			35mm (1 3/8")	14237	5
		60mm (2 3/8")	38mm (1 1/2")	14238	5
			41mm (1 5/8")	14239	5



Assura® 2 pc. High Output Drainable Pouch

Size	Colour	Coupling Size	Product Code	Units
Maxi 32cm (12 1/2")	Opaque	50mm (2")	2846	10
		60mm (2 3/8")	2847	10



Assura® 2 pc. High Output Drainable System, Ileostomy Pouch 2L

Size	Colour	Coupling Size	Product Code	Units
2 L	Opaque	50mm (2") ●	2836	5
		60mm (2 3/8") ●	2837	5



Assura® Drainable Pouch Clamp Closure

Size	Colour	Coupling Size	Product Code	Units
Maxi 30cm (12 1/2")	Transparent	50mm (2") ●	2792	10
		60mm (2 3/8") ●	2793	10



Assura® Drainable Pouch with EasiClose™

Size	Colour	Coupling Size	Product Code	Units
Mini 23cm (9")	Opaque	40mm (1 9/16") ●	13924	10
		50mm (2") ●	13925	10
Midi 26cm (10 3/4")	Opaque with Graphic	40mm (1 9/16") ●	13964	10
		50mm (2") ●	13965	10
		60mm (2 3/8") ●	13966	10
Maxi 28cm (11")	Opaque	40mm (1 9/16") ●	13984	10
		40mm (1 9/16") ● No filter	14494	10
		50mm (2") ●	13985	10
		60mm (2 3/8") ●	13986	10
	Transparent	60mm (2 3/8") ● No filter	14498	10
		40mm (1 9/16") ●	13974	10
		50mm (2") ●	13975	10
		60mm (2 3/8") ●	13976	10



Assura® Closed Pouch

Size	Colour	Coupling Size	Product Code	Units
Midi 18cm (7")	Opaque	40mm (1 9/16") ●	12354	30
		50mm (2") ●	12355	30
		60mm (2 3/8") ●	12356	30
	Transparent	40mm (1 9/16") ●	12344	30
		60mm (2 3/8") ●	12346	30
		40mm (1 9/16") ●	12384	30
Maxi 22cm (8 1/2")	Opaque	50mm (2") ●	12385	30
		60mm (2 3/8") ●	12386	30
		40mm (1 9/16") ●	12374	30
	Transparent	50mm (2") ●	12375	30
		60mm (2 3/8") ●	12376	30
		40mm (1 9/16") ●	12374	30



Assura® Urostomy Pouch

Size	Colour	Coupling Size	Product Code	Units
Maxi 25cm (10")	Opaque	40mm (1 9/16") ●	14224	10
		50mm (2") ●	14225	10
		60mm (2 3/8") ●	14226	10
	Transparent	40mm (1 9/16") ●	14227	10
		50mm (2") ●	14228	10
		60mm (2 3/8") ●	14229	10



Ostomy Barrier and Pouch- Two Piece System

Soft and flexible adhesive coupling allows for extra comfort. Easy to apply the pouch to the skin barrier and easy to remove. All Easiflex Convex Flanges have 2 belt attachment points.

Easiflex® Skin Barrier					
Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	40mm (1 9/16")	10-35mm (3/8 - 1 3/8")	14301	10
		50mm (2")	10-50mm (3/8 - 2")	14302	10
		70mm (2 3/4")	10-68mm (3/8 - 2 3/8")	14303	10
		90mm - (3 1/2")	15-88mm (5/8 - 3 1/2")	14309	5
	Convex Light	40mm (1 9/16")	15-23mm (5/8 - 7/8")	14401	5
		50mm (2")	15-33mm (5/8 - 1 1/4")	14402	5
		50mm (2")	15-43mm (5/8 - 1 3/4")	14403	5
		70mm (2 3/4")	15-61mm (5/8 - 2 1/8")	14406	5
	Pre-Cut	40mm (1 9/16")	21mm (7/8")	14642	5
		50mm (2")	25mm (1")	14643	5
		50mm (2")	28mm (1 1/8")	14644	5



Easiflex® Skin Barrier - Extra Extended Wear					
Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	40mm (1 9/16")	10-35mm (3/8 - 1 3/8")	14304	10
		50mm (2")	10-48mm (3/8 - 2")	14305	10
		70mm (2 3/4")	10-68mm (3/8 - 2 3/8")	14306	10
	Convex Light	40mm (1 9/16")	15-23mm (5/8 - 7/8")	14601	5
		50mm (2")	15-33mm (5/8 - 1 1/4")	14602	5
		50mm (2")	15-43mm (5/8 - 1 3/4")	14603	5
		70mm (2 3/4")	15-61mm (5/8 - 2 1/8")	14408	5

Easiflex® Drainable Pouch with EasiClose™					
Size	Colour	Coupling Size		Product Code	Units
Midi 26cm (10 3/4")	Opaque	40mm (1 9/16")		14346	20
		50mm (2")		14347	20
	Transparent	50mm (2")		14342	20
Maxi 28cm (11")	Opaque	40mm (1 9/16")		14361	20
		40mm (1 9/16")	 No filter	14524	20
		50mm (2")		14362	20
		50mm (2")	 No filter	14525	20
		70mm (2 3/4")		14349	20
		40mm (1 9/16")		14356	20
	Transparent	40mm (1 9/16")	 No filter	14520	20
		50mm (2")		14357	20
		50mm (2")	 No filter	14521	20
		70mm (2 3/4")		14344	20
		70mm (2 3/4")	 No filter	14528	20
		90mm - (3 1/2")	 No filter	14527	20
XL 30cm (12")	Opaque	70mm (2 3/4")		14364	20
	Transparent	70mm (2 3/4")		14359	20



Easiflex® Closed Pouch

Size	Colour	Coupling Size	Product Code	Units
Midi 18cm (7")	Opaque	40mm (1 9/16")	14316	30
		50mm (2")	14317	30
Maxi 21cm (8 1/4")	Opaque	40mm (1 9/16")	14331	30
		50mm (2")	14332	30
		70mm (2 3/4")	14319	30
		90mm - (3 1/2")	14355	30
XL 26cm (10 1/4")	Transparent	70mm (2 3/4")	14329	30



Easiflex® Urostomy Pouch

Size	Colour	Coupling Size	Product Code	Units
Maxi 25cm (10")	Opaque	40mm (1 9/16")	14553	20
		50mm (2")	14555	20
	Transparent	40mm (1 9/16")	14552	20
		50mm (2")	14554	20



Assura®

Ostomy Barrier and Pouch - One Piece System

Skin Barrier with spiral technology for secure, long lasting wear. Different convexity and pouch shape. All Assura 1pc Convex Systems have 2 belt attachment points.

Assura® Drainable Pouch Clamp Closure

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 30cm (12")	Cut-to-Fit	Non-Convex	Opaque	10-70mm (3/8 - 2 3/4")	2520 No filter	10
			Transparent	10-70mm (3/8 - 2 3/4")	2490 No filter	10
		Convex Deep	Opaque	15-43mm (5/8 - 1 3/4")	12533 No filter	10
			Transparent	15-43mm (5/8 - 1 3/4")	13706 No filter	10



Assura® Drainable Pouch with EasiClose™

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 26cm (10 1/4")	Cut-to-Fit	Non-Convex	Opaque	10-55mm (3/8 - 2 1/8")	13840	10
		Convex Light	Opaque	15-33mm (5/8 - 1 1/4")	14404	10
				15-43mm (5/8 - 1 3/4")	14405	10
		Convex Deep	Opaque	15-43mm (5/8 - 1 3/4")	14196	10
	Pre-Cut	Non-Convex	Opaque	25mm (1")	13844	10
				30mm (1 1/8")	13845	10
Maxi 28cm (11")	Cut-to-Fit	Non-Convex	Opaque	35mm (1 3/8")	13846	10
				10-70mm (3/8 - 2 3/4")	13870	10
		Convex Light	Transparent	10-70mm (3/8 - 2 3/4")	13860	10
			Opaque	15-33mm (5/8 - 1 1/4")	14414	10
				15-43mm (5/8 - 1 3/4")	14415	10
			Transparent	15-33mm (5/8 - 1 1/4")	14411	10
				15-43mm (5/8 - 1 3/4")	14412	10



Assura® Drainable Pouch with EasiClose™

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 28cm (11")	Cut-to-Fit	Convex Deep	Opaque	15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	14106	10
			Transparent	15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	14103	10
	Pre-Cut	Non-Convex	Opaque	25mm (1")	13874	10
				30mm (1 $\frac{1}{8}$ ")	13875	10
				35mm (1 $\frac{3}{8}$ ")	13876	10
				25mm (1")	14422	10
		Convex Light	Opaque	28mm (1 $\frac{1}{8}$ ")	14423	10
				31mm (1 $\frac{1}{4}$ ")	14424	10
				21mm ($\frac{7}{8}$ ")	14173	10
				25mm (1")	14174	10
		Convex Deep	Opaque	28mm (1 $\frac{1}{8}$ ")	14175	10
				31mm (1 $\frac{1}{4}$ ")	14176	10
				21mm ($\frac{7}{8}$ ")	14163	10
			Transparent	25mm (1")	14164	10
				28mm (1 $\frac{1}{8}$ ")	14165	10
				31mm (1 $\frac{1}{4}$ ")	14166	10



Assura® Closed Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Mini 12cm (4 $\frac{3}{4}$ ")	Cut-to-Fit	Non-Convex	Opaque	20-55mm ($\frac{3}{8}$ - 2 $\frac{1}{8}$ ")	12110	30
				20-55mm ($\frac{3}{8}$ - 2 $\frac{1}{8}$ ")	12140	30
	Cut-to-Fit	Non-Convex	Transparent	20-55mm ($\frac{3}{8}$ - 2 $\frac{1}{8}$ ")	12130	30
Midi 18cm (7")	Pre-Cut	Non-Convex	Opaque	25mm (1")	12144	30
				30mm (1 $\frac{1}{8}$ ")	12145	30
				35mm (1 $\frac{3}{8}$ ")	12146	30
				40mm (1 $\frac{9}{16}$ ")	12147	30




Assura® 1-Piece Post-Op Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Mini 31cm (12 $\frac{1}{4}$ ")	Cut-to-Fit	Non-Convex	Transparent	13-70mm ($\frac{1}{2}$ - 2 $\frac{3}{4}$ ")	12805	5
				13-70mm ($\frac{1}{2}$ - 2 $\frac{3}{4}$ ")	12800	5
Midi 33cm (13")	Cut-to-Fit	Non-Convex	Transparent	13-100mm ($\frac{1}{2}$ - 4")	12802	5
				13-70mm ($\frac{1}{2}$ - 2 $\frac{3}{4}$ ")	12808 Sterile	5
				13-100mm ($\frac{1}{2}$ - 4")	12810 Sterile	5



Assura® Closed Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 18cm (7")	Pre-Cut	Convex Light	Opaque	25mm (1")	14452	10
				28mm (1 1/8")	14453	10
				31mm (1 1/4")	14454	10
Maxi 22cm (8 3/4")	Cut-to-Fit	Non-Convex	Opaque	20-70mm (3/4 - 2 3/4")	12170	30
		Convex Light	Opaque	15-43mm (5/8 - 1 3/4")	14445	10
			Transparent	15-33mm (5/8 - 1 1/4")	14441	10
	Pre-Cut	Non-Convex	Opaque	25mm (1")	12174	30
				30mm (1 1/8")	12175	30
				35mm (1 3/8")	12176	30
				40mm (1 9/16")	12177	30



Assura® Urostomy Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 27cm (10 1/2")	Cut-to-Fit	Non-Convex	Opaque	10-55mm (3/8 - 2 1/8")	5580	10
			Transparent	10-55mm (3/8 - 2 1/8")	5585	10
		Convex Deep	Transparent	15-33mm (5/8 - 1 1/4")	12595	10
				15-43mm (5/8 - 1 3/4")	12596	10
	Pre-Cut	Convex Deep	Transparent	15mm (5/8")	12991	10
				18mm (3/4")	12992	10
				21mm (7/8")	12993	10
				25mm (1")	12994	10
				28mm (1 1/8")	12995	10
				31mm (1 1/4")	12996	10
				35mm (1 3/8")	12997	10



Assura® Urostomy Pouch - Extra Extended Wear

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 25cm (10")	Cut-to-Fit	Non-Convex	Opaque	10-55mm (3/8 - 2 1/8")	14201	10
			Transparent	10-55mm (3/8 - 2 1/8")	14202	10
		Convex Light	Opaque	15-43mm (5/8 - 1 3/4")	14701	10
			Transparent	15-43mm (5/8 - 1 3/4")	14702	10
Maxi 27cm (10 1/2")	Cut-to-Fit	Non-Convex	Opaque	10-55mm (3/8 - 2 1/8")	12471	10
			Transparent	10-55mm (3/8 - 2 1/8")	12472	10
		Convex Deep	Transparent	15-33mm (5/8 - 1 1/4")	12432	10
				15-43mm (5/8 - 1 3/4")	12433	10



Ostomy Barrier and Pouch - Two Piece System - Pediatric

Click mechanical coupling.

Assura® Pediatric Skin Barrier					
Cut	Convexity	Stoma Size	Product Code	Units	
Cut-to-Fit	Non-Convex	10-35mm ($\frac{3}{8}$ - $1\frac{3}{8}$ ")	2181	5	

Assura® Pediatric Drainable Pouch				
Size	Colour	Product Code	Units	
13cm (5")	Opaque	2151	10	

Assura® Pediatric Closed Pouch				
Size	Colour	Product Code	Units	
13cm (5")	Transparent	2166	30	

EasiFlex®

Ostomy Barrier and Pouch - Two Piece System - Pediatric

Flex adhesive coupling.

Easiflex® Pediatric Skin Barrier						
Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units	
Cut-to-Fit	Non-Convex	20mm ($\frac{3}{4}$ ")	0-15mm ($0\frac{5}{8}$ ")	14307	10	
		30mm ($1\frac{1}{8}$ ")	0-25mm (0-1")	14308	10	

Easiflex® Pediatric Drainable Pouch					
Size	Colour	Coupling Size	Product Code	Units	
105 mL	Teddy Print	20mm ($\frac{3}{4}$ ")	14691	10	
	Transparent	20mm ($\frac{3}{4}$ ")	14681	10	
150 mL	Teddy Print	30mm ($1\frac{1}{8}$ ")	14692	10	
	Transparent	30mm ($1\frac{1}{8}$ ")	14682	10	

Assura®

Ostomy Barrier and Pouch - One Piece System - Pediatric

Assura® Pediatric Drainable Pouch							
Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units	
15cm (6")	Cut-to-Fit	Non-Convex	Transparent	10-35mm ($\frac{3}{8}$ - $1\frac{3}{8}$ ")	2115	10	

Assura® Pediatric Closed Pouch							
Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units	
13cm (5")	Cut-to-Fit	Non-Convex	Transparent	10-35mm ($\frac{3}{8}$ - $1\frac{3}{8}$ ")	2125	30	

Assura® Pediatric Urostomy Pouch							
Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units	
15cm (6")	Cut-to-Fit	Non-Convex	Transparent	10-35mm ($\frac{3}{8}$ - $1\frac{3}{8}$ ")	8009	10	

Ostomy Accessories

Brava® Protective Sheets

Size (cms)	Size (inches)	Product Code	Units
10 x 10 cm	4 x 4"	3210	10
15 x 15 cm	6 x 6"	3215	5
20 x 20 cm	8 x 8"	3220	5



Brava® Protective Rings

Size (mm) Internal diameter	Size (mm) External diameter	Size (mm) thickness	Product Code	Units
18 mm	48 mm	2.5 mm	12035	10
18 mm	48 mm	4.2 mm	12045	10
27 mm	57 mm	2.5 mm	12037	10
27 mm	57 mm	4.2 mm	12047	10
34 mm	64 mm	2.5 mm	12039	10
34 mm	64 mm	4.2 mm	12049	10



Brava® Protective Rings Wide

Size (mm) Internal diameter	Size (mm) External diameter	Size (mm) thickness	Product Code	Units
18 mm	57 mm	2.5 mm	12036	10
18 mm	64 mm	2.5 mm	12038	10
18 mm	57 mm	4.2 mm	12046	10
18 mm	64 mm	4.2 mm	12048	10



Brava® Protective Rings Extra Wide

Size (mm) Internal diameter	Size (mm) External diameter	Size (mm) thickness	Product Code	Units
18 mm	76 mm	2.5 mm	12032	10
27 mm	76 mm	2.5 mm	12033	10



Brava® Protective Ring Convex

Size (mm) Internal diameter	Size (mm) External diameter	Size (mm) thickness	Product Code	Units
20 mm	60 mm	8.0 mm	12090	10
25 mm	65 mm	8.0 mm	12091	10
30 mm	70 mm	8.0 mm	12092	10
35 mm	75 mm	8.0 mm	12093	10
40 mm	80 mm	8.0 mm	12094	10
45 mm	85 mm	8.0 mm	12095	10



Brava® Mouldable Rings

Size (mm) Internal diameter	Size (mm) thickness	Product Code	Units
18 mm	2.0 mm	12030	10
18 mm	4.2 mm	12042	10



Brava® Skin Barrier Spray and Wipe

Size	Product Code	Units
50 mL Spray	12020	1
Single unit wipe	12021	30



Brava® Adhesive Remover Spray and Wipe

Size	Product Code	Units
50 mL Spray	12010	1
Single unit wipe	12011	30



Brava® Barrier Cream

Size	Product Code	Units
60g Tube	12000	1



Brava® Ostomy Strip Paste			
Size	Product Code	Units	
57g / 2oz	2655	10	
Brava® Elastic Barrier Strips			
Size	Product Code	Units	
Brava Elastic Barrier Strip Original (a.)	12070	20	
Brava Elastic Barrier Strip Straight (b.)	12074	20	
Brava Elastic Barrier Strip Belt (c.)	12072	20	
Brava Elastic Barrier Strip XL (d.)	12076	20	
Brava® Ostomy Powder			
Size	Product Code	Units	
25g / 1oz bottle	1907	16	
Brava® Mio Belt for All SenSura Mio Products			
New design for easier attachment and improved security			
Size	Product Code	Units	
Standard - 100 cm - 40"	0423	15	
XXL - 155cm - 61"	0424	15	
Brava® Belt for SenSura Original + Assura Products			
Size	Product Code	Units	
Standard - 100 cm - 40"	0421	15	
XL - 135cm - 54"	0422	15	
Brava® Lubricating Deodorant			
Size	Product Code	Units	
7.5 mL sachet	12060	20	
240 mL bottle	12061	1	
Brava® Alcohol Free Paste			
Size	Product Code	Units	
60 g tube	12050	1	
Ostomy Paste (contains alcohol)			
Size	Product Code	Units	
60 g tube	2650	12	
Skin Barrier Rings			
Pre-cut Size	Product Code	Units	
20mm (3/4")	2320	30	
25mm (1")	2325	30	
30mm (1 1/8")	2330	30	
40mm (1 9/16")	2340	30	
50mm (2")	2350	30	
Coloplast Prep™ Protective Liquid Skin Barrier (contains alcohol)			
Size	Product Code	Units	
Single Application Wipe	2041	12 boxes/54	
59 mL - 2oz. dab bottle	0925	12	

Coloplast Filtrodor® Pouch Filter

	Product Code	Units
Charcoal-activated deodorizing filter with self-adhering foam center, adheres to front of Assura & SenSura Pouches	0509	5 pouches (10 units per pouch)



Coloplast Pouch Clamp

	Product Code	Units
Simple closure clamp for Assura Drainable Pouches	9500	20



Ostomy Scissors

	Product Code	Units
Stainless steel medical grade ostomy scissors with curved blades and rounded ends - right handed	9505	30
Stainless steel medical grade ostomy scissors with curved blades and rounded ends - left handed	9506	25



Hex-On®

Odour antagonist

Eliminates odour with a light, mild scent.

Hex-On®

Size (ounces)	Size (mL)	Product Code	Units
2 oz.	59 mL	7583	12



Ostomy Cap

Designed for discreet situations. 2 piece compatible with SenSura Click and Assura 2 piece skin barriers.

**Note Uro Minicap and Coloplast Stoma cap do not fit SenSura Mio products*

Coloplast Uro Minicap (100 ml)

Size	Colour	Coupling Size	Product Code	Units
Mini	Opaque	40mm (1 9/16") 	2807	30
		50mm (2") 	2808	30




Coloplast Stoma Cap

Size	Colour	Coupling Size	Product Code	Units
Mini	Opaque	40mm (1 9/16") 	2801	30
		50mm (2") 	2802	30
		60mm (2 3/8") 	2803	30



Coloplast Stoma Cap - One Piece

Size	Colour	Cut to Fit - Stoma size	Product Code	Units
Mini	Opaque	20-55mm (3/4 - 2 1/8")	2501	30



Assura® Irrigation Products

Assura® Irrigation Set

Description	Product Code	Units
Includes: - 1 Water Bag/Regulator - 1 Cone	1500	1



Assura® Irrigation Accessories

Size	Product Code	Units
Stoma Cone - Soft, flexible stoma cone with tubing	1110	4



Assura® Irrigation Sleeves

Size	Product Code	Units
Transparent, disposable irrigation sleeve with adhesive back. - 1-pc sleeve	1003	10



The irrigation sleeves below must be used with their corresponding colour coded flange or faceplate. Consult this chart to order the correct transparent sleeves.

Description	Coupling Size	Use with Assura and SenSura Barrier Product Codes (will not work with SenSura Mio)	Product Code	Units
Irrigation sleeves	40mm (1 9/16") ●	Will couple to any green 40 mm barriers within the Assura and SenSura lines	12834	5
	50mm (2") ●	Will couple to any red 50 mm barriers within the Assura and SenSura lines	12835	5
	60mm (2 3/8") ●	Will couple to any blue 60 mm barriers within the Assura and SenSura lines	12836	5



Coloplast Urostomy Night Bag

Description	Product Code	Units
Includes a fixed connector that fits on to the urostomy pouch. The bag holds up to 2 litres and has a 120cm long flexible tube with anti-reflux valve. Not Compatible with External Catheters or Leg Bags	21365	10



Coloplast Fistula & Wound Management System

Pouch

Size	Cutting Surface	Product Code	Units
Mini	104-159mm (4 1/8 - 6 1/4")	14050	6
Midi	156-228mm (6 3/16 - 9")	14060	3
Maxi	208-297mm (8 1/4 - 11 3/4")	14070	



Flexible Lids

Size	For	Product Code	Units
Mini	Mini Fistula & Wound Management System	14021	10
Midi	Midi Fistula & Wound Management System	14031	
Maxi	Maxi Fistula & Wound Management System	14041	



Accessories

Size	For	Product Code	Units
Bed Drainage Bag	All Fistula & Wound Management systems (Capacity 2000 mL)	14010	10
Drain Port	All Fistula & Wound Management systems	14015	12



Continence

SpeediCath® Flex

Hydrophilic Intermittent Catheter

A soft catheter for men with a dry-sleeve and a flexible tip with all the benefits of the SpeediCath® catheter. Offered with a loop system that fits into a discreet packaging which is 100% aluminum free. Latex, phthalate and PVC free. Sterile.

SpeediCath® Flex				
Length (cm)	Length (inches)	Size	Product Code	Units
44 cm	17"	10 FR ●	28920	30
		12 FR ○	28922	30
		14 FR ●	28924	30
		16 FR ●	28916	30



SpeediCath® Standard

Hydrophilic Intermittent Catheter

Instantly ready to use catheter. Featuring a smooth and uniform hydrophilic coating and polished eyelets to ensure a convenient and simple catheterization. Latex, phthalate and PVC free. Sterile.

SpeediCath® Standard Male Straight				
Length (cm)	Length (inches)	Size	Product Code	Units
40 cm	16"	8 FR ●	27408	30
		10 FR ●	27410	30
		12 FR ○	27412	30
		14 FR ●	27414	30
		16 FR ●	27416	30
		18 FR ●	27418	30

SpeediCath® Standard Male Coude				
Length (cm)	Length (inches)	Size	Product Code	Units
40 cm	16"	10 FR ●	27490	30
		12 FR ○	27492	30
		14 FR ●	27494	30
		16 FR ●	27496	30

SpeediCath® Standard Female Straight				
Length (cm)	Length (inches)	Size	Product Code	Units
20 cm	8"	8 FR ●	27508	30
		10 FR ●	27510	30
		12 FR ○	27512	30
		14 FR ●	27514	30
		16 FR ●	27516	30

SpeediCath® Standard Pediatric Straight				
Length (cm)	Length (inches)	Size	Product Code	Units
20 cm	8"	6 FR ●	27706	30
		8 FR ●	27708	30
		10 FR ●	27710	30



SpeediCath® Compact Male

Hydrophilic Intermittent Catheter

Developed to offer discretion and practicality with all the benefits of the SpeediCath® catheter and more. Telescopic design makes the packaging half the size of a standard male catheter. Latex, phthalate and PVC free. Sterile.

SpeediCath® Compact Male

Length (cm)	Length (inches)	Size	Product Code	Units
33 cm	13"	12/18 FR (one size)	28692	30



SpeediCath® Compact Eve

Hydrophilic Intermittent Catheter

Stylish and discreet compact catheter for women. Designed for modern life, delivering a one-step easy opening and secure reclosing. Triangular shape ensures easy handling and grip. Latex, phthalate and PVC free. Sterile.

SpeediCath® Compact Eve

Length (cm)	Length (inches)	Size	Product Code	Units
8.3 cm	3.3"	10 FR ●	28110	30
		12 FR ○	28112	30
		14 FR ●	28114	30



SpeediCath® Compact Female

Hydrophilic Intermittent Catheter

Compact and discreet solution that is instantly ready to use. The catheter can easily be kept in a pocket or handbag. Latex, phthalate and PVC free. Sterile.

SpeediCath® Compact Female

Length (cm)	Length (inches)	Size	Product Code	Units
7 cm	2.75"	8 FR ●	28578	30
		10 FR ●	28580	30
		12 FR ○	28582	30
		14 FR ●	28584	30



SpeediCath® Compact Female Plus

Length (cm)	Length (inches)	Size	Product Code	Units
9 cm	3.5"	10 FR ●	28810	30
		12 FR ○	28812	30
		14 FR ●	28814	30



SpeediCath® Compact Set

Hydrophilic Intermittent Catheter

Compact and discreet catheter and bag solution that is instantly ready to use. Can easily be kept in a pocket or handbag. Latex, phthalate and PVC free. Sterile.

SpeediCath® Compact Set Male					
Length (cm)	Length (inches)	Size		Product Code	Units
33 cm	13"	12/18 FR	(one size)	28422	20



SpeediCath® Compact Set Female					
Length (cm)	Length (inches)	Size		Product Code	Units
9 cm	3.5"	10 FR	●	28520	20
		12 FR	○	28522	20
		14 FR	●	28524	20



Self-Cath®

Intermittent Catheter

Siliconized surface for smooth insertion. Made of flexible medical grade PVC with smooth fire polished eyelets for patient comfort and ease of use. Latex free. Sterile.

Self-Cath® Female Straight						
Length (cm)	Length (inches)	Size	Reference Code		Product Code	Units
15 cm	6"	8 FR	●	504600	208	30
		10 FR	●	504610	210	30
		12 FR	○	504620	212	30
		14 FR	●	504630	214	30
		14 FR	●	501960	without connector	240



Self-Cath® Pediatric							
Length (cm)	Length (inches)	Size	Reference Code		Product Code	Units	
25 cm	10"	5 FR		504400	305	30	
		6 FR		504410	306	30	
		8 FR		504420	without connector	308	30
		10 FR		504430	without connector	310	30



Self-Cath® Extension Tube						
Length (cm)	Length (inches)	Size		Reference Code	Product Code	Units
61 cm	24"	n/a		507850	475	12



Self-Cath® Male Straight

Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
40 cm	16"	8 FR	● 504500	408	30
		10 FR	● 504510	410	30
		12 FR	○ 504520	412	30
		14 FR	● 504530	414	30
		16 FR	● 504540	416	30
		18 FR	● 504550	418	30



Self-Cath® Male Soft Straight

Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
40 cm	16"	12 FR	○ 501810	112	30
		14 FR	● 501820	114	30
		16 FR	● 501830	116	30

Self-Cath® Male Straight Curved Packaging

Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
40 cm	16"	12 FR	○ 502200	460	50
		14 FR	● 502210	450	50

Self-Cath® Male Coude Tapered with Guide Stripe

Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
40 cm	16"	8 FR	● 504650	608	30
		10 FR	● 504660	610	30
		12 FR	○ 504670	612	30
		14 FR	● 504680	614	30



Self-Cath® Male Coude Olive with Guide Stripe

Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
40 cm	16"	8 FR	● 504710	808	30
		10 FR	● 504720	810	30
		12 FR	○ 504730	812	30
		14 FR	● 504740	814	30
		16 FR	● 504750	816	30
		18 FR	● 504760	818	30



Self-Cath® Male Closed System Straight

Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
40 cm	16"	14 FR	● 502730	1114	50
		16 FR	● 502740	1116	50



Conveen® Optima

Male External Catheter - Silicone

Discreet, easy to use, reliable all day self-adhering external catheter. Double-grip strip ensures easy and smooth roll-out for a secure application. Push ring with sure grip ribs provides secure connection to the collecting bag. Anti-kink bellow stop kinks ensuring a free flow. Compact easy-to-open packaging. Latex free.

Conveen® Optima Standard

Length (cm)	Length (inches)	Size	Product Code	Units
7.8 cm	3"	25 mm	22025	30
		28 mm	22028	30
		30 mm	22030	30
		35 mm	22035	30
		40 mm	22040	30

Conveen® Optima Short

Length (cm)	Length (inches)	Size	Product Code	Units
4.6 cm	2"	21 mm	22121	30
		25 mm	22125	30
		30 mm	22130	30
		35 mm	22135	30



Conveen® Security+

Male External Catheter - SEBS*

Easy to use, reliable all day self-sealing external catheter. Applicator grip and single-grip strip ensures easy and smooth roll-out for a secure application. Push ring with sure grip ribs provides secure connection to the collecting bag. Anti-kink bellow stop kinks ensuring a free flow. Latex free.

Conveen® Security+ Self Sealing

Length (cm)	Length (inches)	Size	Product Code	Units
8.2 cm	3"	21 mm	5221	35
		25 mm	5225	35
		25 mm	5225H	100
		30 mm	5230	35
		30 mm	5230H	100
		35 mm	5235	35
		35 mm	5235H	100



Conveen® Security+ Adhesive Liner

Length (cm)	Length (inches)	Size	Product Code	Units
8.2 cm	3"	30 mm	5030	35
		30 mm	5030H	100
		35 mm	5035	35
		35 mm	5035H	100



*Styrene Ethylene Butylene Styrene

Clear Advantage® / Freedom Clear™

Male External Catheter - Silicone

Self-adhering male external catheter. Latex-free.

Clear Advantage® with Aloe Vera

Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
7.9 cm	3"	23 mm	504311	C6100	100
		28 mm	504331	6200	100
		31 mm	504341	6300	100
		35 mm	504351	6400	100

Freedom Clear™

Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
7.9 cm	3"	23 mm	505711	5100	100
		28 mm	505731	C5200	100
		31 mm	505741	5300	100
		35 mm	505751	5400	100

Freedom Clear™ "LS" Long Seal

Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
8.9 cm	3.5"	28 mm	505421	5290	100
		31 mm	505431	5390	100
		35 mm	505341	5490	100
		40 mm	505451	5590	100

Freedom Clear™ "SS" Sport Sheath

Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
5.5 cm	2"	23 mm	505311	5110	100
		28 mm	505321	C5210	100
		31 mm	505331	C5310	100
		35 mm	505441	C5410	100



Conveen®

Male External Catheter - Latex

Conveen® Self-Sealing

Length (cm)	Length (inches)	Size	Product Code	Units
10.5 cm	4"	20 mm	5212	35
		25 mm	5200	35
		30 mm	5205	35
		35 mm	5210	35
		40 mm	5215	35



Conveen® Adhesive Liner

Length (cm)	Length (inches)	Size	Product Code	Units
10.5 cm	4"	30 mm	5130	35
		35 mm	5135	35



Active Cath®

Male External Catheter - Latex

Self-sealing male external catheter.

Active Cath®					
Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
7.9 cm	3"	23 mm	506220	8100	100
		28 mm	506200	8300	100
		31 mm	506240	8305	100
		35 mm	506260	8500	100



Freedom Cath®

Male External Catheter - Latex

Self-sealing male external catheter.

Freedom Clear®					
Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
7.9 cm	3"	23 mm	506110	8000	100
		28 mm	506130	8200	100
		31 mm	506140	8205	100
		35 mm	506150	8400	100




Conveen® Active

Urine Leg Bag

Developed for an active lifestyle. Integrated in a comfortable elastic band that wraps around the thigh for ease of use, mobility and discretion. Hide-away outlet tap concealed in elastic band when not in use. Soft fabric band for increased wearing comfort and a no-slip silicone feature inside straps for security. Latex and PVC free.

Conveen® Active Leg Bag, Non Sterile				
Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
9 oz	250 mL	3"	25501	10



Conveen® Security+

Urine Leg bag

Secure and reliable urine leg bag with anti-kink tubing to reduce the risk of accidental leakage. Soft backing for increased wearing comfort (except codes 21034 and 21054). Latex and PVC free.

Conveen® Security+ Leg Bag, Clamp Outlet, Straps, Non Sterile				
Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
17 oz	500 mL	20"	5161	10
26 oz	750 mL	20"	5167	10

Conveen® Security+ Leg Bag, Clamp Outlet, Sample Port, Straps, Sterile				
Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
17 oz	500 mL	20"	5163	10

Conveen® Security+ Leg Bag, Lever Outlet, Straps, Sterile				
Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
17 oz	500 mL	20"	21034	10
34 oz	1000 mL	20"	21054	10



Conveen® Standard

Day / Night Urine Bag

Secure and reliable urine bag for bedside drainage or large capacity leg bag. Clamp outlet makes it easy to verify proper closure and reduce the risk of accidental leakage. Standard tubing. Leg bag straps sold separately. Soft backing for increased wearing comfort. Latex free.

Conveen® Standard Day / Night Bag, Clamp Outlet, Non Sterile

Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
51 oz	1500 mL	36"	5062	10

Conveen® Standard Day / Night Bag, Clamp Outlet, Sample Port, Sterile

Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
51 oz	1500 mL	36"	5063	10



Conveen® Contour

Urine Leg Bag

Secure and reliable urine leg bag. Enhanced chambered design for added discretion. Soft backing for increased wearing comfort. Latex free.

Conveen® Contour Leg Bag, Clamp Outlet, Straps, Non Sterile

Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
21 oz	600 mL	18"	5170	10

28 oz	800 mL	18"	5174	10
-------	--------	-----	------	----



Conveen® Night Bag

Urine Night Bag

Secure and reliable urine night bag with anti-kink tubing to prevent accidental leakage. Large capacity allows the user to feel secure all night for a full night's sleep. Latex and PVC free. Not Compatible with Urostomy Pouches.

Conveen® Bedside Night Bag, Non-Sterile

Capacity (oz)	Capacity (L)	Product Code	Units
68 oz	2L	21346	30

Conveen® Bedside Night Bag, Sterile

Capacity (oz)	Capacity (L)	Product Code	Units
68 oz	2L	21356	10



Accessories

Conveen® Bag Hanger

Length (cm)	Length (inches)	Product Code	Units
18 cm	7"	5070	1



Conveen® Leg Bag Strap

Length (cm)	Length (inches)	Product Code	Units
63 cm	25"	5050	10



Conveen® Catheter retaining Strap

Length (cm)	Length (inches)	Product Code	Units
90 cm	35"	383004	2



Peristeen® Anal Plug

Simple, safe and discreet aid for fecal incontinence. Inserted like a suppository. It can be used for up to 12 hours, preventing the uncontrolled loss of solid stool. Each box comes with a tube of lubricant. Latex free.

Peristeen® Anal Plug			
Length (cm)	Length (inches)	Product Code	Units
3.7 cm]	1.5"	1450	20
4.5 cm	1.75"	1451	20



Peristeen® Anal Irrigation

Transanal irrigation system for fecal incontinence and constipation. Minimizes the likelihood of involuntary bowel leakage and/or constipation.

Peristeen® Anal Irrigation System		
Items	Description	Product Code
- 1 control unit		
- 2 sterile rectal catheters		
- 1 bag (with lid)	** 90 irrigations = 6 months	29121
- 2 straps		
- 1 toiletries bag		



Peristeen® Anal Irrigation System - Small		
Items	Description	Product Code
- 1 control unit		
- 2 sterile rectal catheters - small		
- 1 bag (with lid)	** 90 irrigations = 6 months	29126
- 2 straps		
- 1 toiletries bag		

Peristeen® Anal Irrigation Accessory unit		
Items	Description	Product Code
- 15 sterile rectal catheters	**15 irrigations = 1 month	29122
- 1 bag (no lid)		



Peristeen® Anal Irrigation Accessory unit - Small		
Items	Description	Product Code
- 15 sterile rectal catheters - small	**15 irrigations = 1 month	29127
- 1 bag (no lid)		



Peristeen® Anal Irrigation Rectal Catheter	
Items	Product Code
- 10 sterile rectal catheters	29123



Peristeen® Anal Irrigation Leg Strap	
Items	Product Code
- 1 pair	29124



Peristeen® Anal Irrigation Tube	
Items	Product Code
- 2 tubes with blue connectors	29125



** if used every other day

Wound

Biatain® Silicone

Foam Dressing with Silicone Adhesive

Soft, flexible and absorbent foam dressing with gentle silicone adhesive for use on a broad range of exuding wounds, making it the obvious choice for moist wound healing of commonly occurring wounds. For more information please visit www.biatainsiliconeca.com.



reddot award 2017
winner



Biatain® Silicone

Size (cm)	Size (inches)	Product Code	Units
7.5 x 7.5 cm	3 x 3"	33434	10
10 x 10 cm	4 x 4"	33435	10
12.5 x 12.5 cm	5 x 5"	33436	10
15 x 15 cm	6 x 6"	33437	5
17.5 x 17.5 cm	7 x 7"	33438	5
10 x 20 cm	3.9 x 7.9"	33400	5
10 x 30 cm	3.9 x 11.8"	33401	5
18 x 18 cm	7.1 x 7.1" Heel	33406	5
25 x 25 cm	9.8 x 9.8" Large Sacral	39000*	5
15 x 19 cm	5.9 x 7.5" Small Sacral	39001*	5
14 x 19.5 cm	5.5 x 7.7" Multi-shape	39010*	5



* This product has the Pressure Injury Prevention indication.

Biatain® Silicone Lite

Thin Foam Dressing with Silicone Adhesive

Thin and ultra flexible foam dressing with gentle silicone adhesive for use on a wide range of low to moderate exuding wounds.



Biatain® Silicone Lite

Size (cm)	Size (inches)	Product Code	Units
5 x 5 cm	2 x 2"	33452	5
5 x 12.5 cm	2 x 5"	33453	5
7.5 x 7.5 cm	3 x 3"	33444	10
10 x 10 cm	4 x 4"	33445	10
12.5 x 12.5 cm	5 x 5"	33446	10



Biatain® Silicone Ag

Foam Dressing with Silicone Adhesive and Silver

Flexible and absorbent foam dressing with gentle silicone adhesive and a patented silver complex for use on a broad range of exuding infected wounds or wounds at risk of infection. For more information please visit biatainsiliconeca.com.



Biatain® Silicone Ag

Size (cm)	Size (inches)	Product Code	Units
7.5 x 7.5 cm	3 x 3"	39636	5
10 x 10 cm	4 x 4"	39637	5
12.5 x 12.5 cm	5 x 5"	39638	5
15 x 15 cm	6 x 6"	39639	5
17.5 x 17.5 cm	7 x 7"	39640	5
10 x 20 cm	3.9 x 7.9"	39644	5
10 x 30 cm	3.9 x 11.8"	39645	5
15 x 19 cm	5.9 x 3.5" Small Sacral	39650	5
25 x 25 cm	9.8 x 9.8" Large Sacral	39651	5
18 x 18 cm	7.1 x 7.1" Heel	39652	5



Biatain®

Foam Dressing

Soft and conformable polyurethane foam dressing that effectively absorbs and retains wound exudate. This ensures a balance that is optimal for management of exuding wounds. Available in adhesive and non-adhesive. The adhesive has a hydrocolloid border – for a secure adhesion. This product is also used in high friction and shear areas such as the coccyx.



Biatain® Adhesive

Size (cm)	Size (inches)		Product Code	Units
7.5 x 7.5 cm	3 x 3"		3462	10
10 x 10 cm	4 x 4"		3430	10
12.5 x 12.5 cm	5 x 5"		3420	10
17 x 17 cm	7 x 7"	Small Sacral	3483	5
19 x 20 cm	7.5 x 8"	Heel	3488	5
23 x 23 cm	9 x 9"	Large Sacral	3485	5



Biatain® Non-Adhesive

Size (cm)	Size (inches)	Product Code	Units
5 x 7 cm	2 x 2.8"	6105	10
10 x 10 cm	4 x 4"	3410	10
10 x 20 cm	4 x 8"	3412	5
15 x 15 cm	6 x 6"	3413	5
20 x 20 cm	8 x 8"	3416	5



Biatain® Soft Hold

Foam Dressing with Gentle Adhesion

Soft and conformable polyurethane foam dressing with a gentle adherent layer covering less than 50% of the foam surface. This product is used under compression. This allows both hands to be free during dressing application and removal.



Biatain® Soft Hold

Size (cm)	Size (inches)	Product Code	Units
10 x 10 cm	4 x 4"	3470	5
10 x 20 cm	4 x 8"	3472	5
15 x 15 cm	6 x 6"	3475	5



Biatain® Ag

Foam Dressing with Silver

Soft and conformable polyurethane foam dressing with a patented silver complex for use on infected wounds or wounds at risk of infection. Available in adhesive and non-adhesive. The adhesive has a hydrocolloid border – for a secure adhesion. This product is also used in high friction and shear areas such as the coccyx.



Biatain® Ag Adhesive

Size (cm)	Size (inches)		Product Code	Units
7.5 x 7.5 cm	3 x 3"		9631	5
12.5 x 12.5 cm	5 x 5"		9632	5
18 x 18 cm	7 x 7"		9635	5
23 x 23 cm	9 x 9"	Sacral	9641	5
19 x 20 cm	7.5 x 8"	Heel	9643	5



Biatain® Ag Non-Adhesive

Size (cm)	Size (inches)		Product Code	Units
5 x 7 cm	2 x 2.8"		5105	5
10 x 10 cm	4 x 4"		9622	5
10 x 20 cm	4 x 8"		9623	5
15 x 15 cm	6 x 6"		9625	5
20 x 20 cm	8 x 8"		9626	5



Biatain® Alginate

Calcium Alginate Dressing

Highly absorbent alginate dressing for use on moderate to heavily exuding wounds of any size, shape and depth.

Biatain® Alginate

Size (cm)	Size (inches)		Product Code	Units
5 x 5 cm	2 x 2"		3705	30
10 x 10 cm	4 x 4"		3710	10
15 x 15 cm	6 x 6"		3715	10
3 x 44 cm	1 x 17.5"	Ribbon	3740	6



Biatain® Alginate Ag

Calcium Alginate Dressing with Silver

Highly absorbent alginate dressing with silver for use on moderate to heavily exuding infected wounds or wounds at risk of infection of any size, shape and depth.

Biatain® Alginate Ag

Size (cm)	Size (inches)		Product Code	Units
5 x 5 cm	2 x 2"		3755	30
10 x 10 cm	4 x 4"		3760	10
15 x 15 cm	6 x 6"		3765	10
3 x 44 cm	1 x 17.5"	Ribbon	3780	10



Biatain® Ibu

Foam Dressing with Ibuprofen

Soft and flexible absorbent foam dressing containing Ibuprofen for use on painful wounds. Available in non-adhesive and soft hold.



Biatain® Ibu Non-Adhesive				
Size (cm)	Size (inches)	Product Code	Units	
10 x 10 cm	4 x 4"	4110	5	
10 x 20 cm	4 x 8"	4112	5	
15 x 15 cm	6 x 6"	4115	5	
20 x 20 cm	8 x 8"	4120	5	
Biatain® Ibu Soft-Hold				
Size (cm)	Size (inches)	Product Code	Units	
10 x 10 cm	4 x 4"	4140	5	
10 x 20 cm	4 x 8"	4142	5	

Biatain® Super

Hydrocapillary Pad Dressing

Highly absorbent hydrocapillary pad dressing for use on heavily exuding wounds. Available in non-adhesive and adhesive.


Biatain® Super Adhesive				
Size (cm)	Size (inches)	Product Code	Units	
10 x 10 cm	4 x 4"	4610	10	
12.5 x 12.5 cm	5 x 5"	4612	10	
12 x 20 cm	4.5 x 8"	4625	10	
15 x 15 cm	6 x 6"	4615	10	
20 x 20 cm	8 x 8"	4620	10	
Biatain® Super Non-Adhesive				
Size (cm)	Size (inches)	Product Code	Units	
10 x 10 cm	4 x 4"	4630	10	
12.5 x 12.5 cm	5 x 5"	4632	10	
12 x 20 cm	4.5 x 8"	4645	10	
15 x 15 cm	6 x 6"	4635	10	
20 x 20 cm	8 x 8"	4639	10	

Comfeel® Plus

Hydrocolloid Dressing

Hydrocolloid dressing for use on low to moderate exuding wounds. Provides faster wound healing by sealing and protecting the wounds.

Comfeel® Plus				
Size (cm)	Size (inches)	Product Code	Units	
4 x 6 cm	1.5 x 2.5"	33146	30	
10 x 10 cm	4 x 4"	33110	10	
15 x 15 cm	6 x 6"	33115	5	
20 x 20 cm	8 x 8"	33120	5	
18 x 20 cm	7 x 8" Sacral	33285	5	
9 x 11 cm	23 x 28" Contour	33283	5	




Comfeel® Plus Transparent

Hydrocolloid Dressing

Thin and flexible hydrocolloid dressing for use on dry or low exuding wounds. Provides faster wound healing by sealing and protecting the wound. Transparency allows for visual monitoring of the wound.

Comfeel® Plus Transparent				
Size (cm)	Size (inches)	Product Code	Units	
5 x 7 cm	2 x 2.75"	33530	10	
10 x 10 cm	4 x 4"	33533	10	
9 x 14 cm	3.5 x 5.5"	33536	10	
15 x 15 cm	6 x 6"	33539	5	
9 x 25 cm	3.5 x 10"	33537	5	
5 x 15 cm	2 x 6"	33547	10	
5 x 25 cm	2 x 10"	33548	10	



Hydrophilic Wound Dressing

Zinc-oxide based hydrophilic sterile coating. Used for autolytic debridement. The ideal alternative for difficult-to-dress areas. Sterile.

Triad™			
Size (ounces)	Size (grams)	Product Code	Units
2.5 oz.	71 g tube	1964	12
6 oz.	170 g tube	1967	12



Purilon® Gel

Hydrogel with natural ingredients

Effective and gentle debridement.

Purilon® Gel			
Size (ounces)	Size (grams)	Product Code	Units
28 oz.	8 g tube	3906	10
5 oz.	15 g tube	3900	10



Woun'Dres®

Collagen Hydrogel

Rehydrates dry wounds and eschar.

Woun'Dres® Collagen Hydrogel			
Size (ounces)	Size (grams)	Product Code	Units
1 oz.	28 g tube	1166	36
3 oz.	84 g tube	7690	12



Sea-Clens®

Wound Cleanser

Saline based solution with a surfactant for cleansing and irrigating acute and chronic wounds.

Sea-Clens® Wound Cleanser			
Size (ounces)	Size (mL)	Product Code	Units
6 oz.	178 mL	1063	12
12 oz.	355 mL	1061	12



Skin

Sproam®



No-Rinse Spray or Foam Cleanser

Spray or foam body wash, shampoo and incontinence cleanser.

Sproam®		
Size	Product Code	Units
175 ml	0897	12
350 ml	0901	12



Bedside-Care® Sensitive Skin



No-Rinse Foam Cleanser

Unscented foaming body wash, shampoo and incontinence cleanser.

Bedside-Care® Sensitive Skin		
Size	Product Code	Units
125 mL	7300	36
240 mL	7301	12



Bedside-Care® Foam



No-Rinse Foam Cleanser

Scented foaming body wash, shampoo and incontinence cleanser.

Bedside-Care®		
Size	Product Code	Units
240 mL	67145	12



Bedside-Care®



No-Rinse Spray Cleanser

Spray body wash, shampoo and incontinence cleanser.

Bedside-Care®		
Size	Product Code	Units
120 mL Scented	61768	36
240 mL Scented	61762	12
240 mL Unscented	61761	12



No-Rinse Cleanser



Rinse Cleanser



Moisturize



Protect



CHG Compatible



Neonate Safe

Bedside-Care® Easi-Cleanse®



No-Rinse Wash Cloth.

Body wash, shampoo, and incontinence cleanser in a self-sudsing disposable wash cloth.

Bedside-Care® Easi-Cleanse		
Size	Product Code	Units
5 per package	7056	100 pkgs/case



Gentle Rain® Extra Mild



Rinse Cleanser

Body wash and shampoo.

Gentle Rain® Extra Mild		
Size	Product Code	Units
5 mL sachet	7231	300
50 mL	7230	36
600 mL (accessories available, see pg 62)	7233	12
1000 mL (accessories available, see pg 62)	7234	12



Atrac-Tain®



Moisturizer

Skin moisturizer for severely dry, scaly and fissured skin with 4% Alphahydroxyacid (AHA).

Atrac-Tain®		
Size	Product Code	Units
2 g sachet	1843	300
60 mL	1802	12
140 mL	1814	12



Sween® 24



Moisturizer

Once a day skin moisturizer for itchy, irritated, dry and flaky skin.

Sween® 24		
Size	Product Code	Units
4 g sachet	7090	300
90 g	7094	12



DIN 02250322

Sween Cream®



Moisturizer

Soothing body cream for irritated dry skin.

Sween Cream®		
Size	Product Code	Units
90 mL	7067	12
190 mL	7068	12



Sween® Lotion Unscented



Moisturizer

Everyday moisturizer. For all body use.

Sween® Lotion Unscented		
Size	Product Code	Units
60 mL	7222	36
120 mL	7223	36
620 mL	7224	12



No-Rinse Cleanser



Rinse Cleanser



Moisturize



Protect



CHG Compatible



Neonate Safe

Critic-Aid® Clear



Skin Protectant

Protects skin and treats irritation due to incontinence.

Critic-Aid® Clear		
Size	Product Code	Units
4 g sachet	7565	300
71 g	7573	12



DIN 02277778

Baza® Protect II



Skin Protectant

Protects skin and treats irritation due to incontinence. Helps seal out wetness & relieves sore, irritated skin.

Baza® Protect II		
Size	Product Code	Units
60 g	7086	12
140 g	7087	12



DIN 02315793

Baza Cleanse and Protect®



Skin Protectant

All-in-one cleanser, moisturizer and a breathable skin protectant for sore, irritated, chafed skin due to incontinence.

Baza Cleanse and Protect®		
Size	Product Code	Units
4 g sachet (fragrance free)	7710	300
240 mL (fragrance free)	7712	12
240 mL (with fragrance)	7725	12



DIN 02243966



InterDry®

Moisture-wicking textile with antimicrobial silver complex

Effectively manage complications associated with skin folds. Appropriate for a variety of skin-to-skin and skin-to-device application sites. For more information please visit www.interdry.ca.

InterDry®			
Size (cm)	Size (Inches)	Product Code	Units
25 x 46 cm	10 x 18"	7915	40
25 x 92 cm	10 x 36"	7912	10
25 x 366 cm	10 x 144" textile on a roll	7910	10



Isagel®



Hand Sanitizer Gel

Sanitizing hand gel with added moisturizer.

Isagel®			
Size	Product Code	Units	
115 mL	1644	36	
621 mL (Accessories available. Please see page 62.)	1645	12	
800 mL (Accessories available. Please see page 62.)	7041	12	



Sween® Body Powder

All-body powder

Provides excellent absorption to keep patient dry.

Sween® Body Powder			
Size	Product Code	Units	
90 g	0506	36	
227 g	0505	36	



No-Rinse Cleanser



Rinse Cleanser



Moisturize



Protect



CHG Compatible



Neonate Safe

Accessories

Coloplast® Wall Dispenser for Bag

Description	Product Code	Units
Compatible with 800 mL (Isagel®) and 1000 mL (Gentle Rain® Extra Mild) sized bags (7234 and 7041)*	7251	12



Coloplast® Wall Dispenser Drip Tray

Description	Product Code	Units
To be used with Wall Dispenser to prevent spillage (7251)*	7262	60



Embrace™ Dispenser Bracket

Description	Product Code	Units
Compatible with 600 mL (Gentle Rain® Extra Mild) and 621 mL (Isagel®) sized pump bottles (1645, 7224, 7233)*	0089	12



*Product not included.

Brief Statements

Biatain® Family of Dressings

Biatain® dressings are single-use only. See the device manual for detailed information regarding complete instructions, contraindications, warnings, precautions and potential adverse events. For further questions, call Coloplast Canada Corporation at 1-888-880-8605 and/or consult the company website at www.coloplast.ca.

InterDry®

InterDry® is a skin protectant for the management of skin folds and other skin-to-skin contact areas. It reduces microbial colonization in the fabric. This single use fabric provides moisture transportation to keep skin dry while the antimicrobial in the fabric reduces odor. The safe use of InterDry fabric during pregnancy, lactation and on children has not been demonstrated. Although many patients benefit from the use of this device, results may vary. See the device labeling for detailed information including the indication, contraindications, warnings, precautions, and potential complications/adverse events. For further questions, call Coloplast Canada Corporation at 1-888-880-8605 and/or consult the company website at www.coloplast.ca.

Peristeen Anal Plug®

The Peristeen Anal Plug is intended as a barrier to faeces in the rectum. The anal plug can be used by people who are faecal incontinent or suffer from involuntary passing of air (flatus). It functions as an effective barrier to faeces in the rectum.

Contraindication: Peristeen Anal Plug must not be used, if you suffer from: hemorrhoids (3rd and 4th degree), and anal stenosis

Caution: This product should not be used without assessment by your health care professional. Contact your health care professional before you use Peristeen Anal Plug for the first time, as professional guidance in the use of the plug is vital. The anal plug should not be used if you are suffering from diarrhoea. If placed in the mouth, the Peristeen Anal Plug will expand and may cause choking. Remove the anal plug when you feel the urge to defecate. If you have no sensation in the rectal area and do not feel the urge to defecate, changing plug at regular intervals is recommended.

The Peristeen Anal Plug must be removed no later than 12 hours after insertion. If you experience difficulties removing the Peristeen Anal Plug, you must contact your health care professional. If you repeatedly experience discomfort when using the Peristeen Anal Plug, stop using it and contact your health care professional. Keep out of reach of children. See the device manual for detailed information regarding instructions, contraindications, warnings, precautions, and potential complications/adverse events or, call Coloplast Canada Corporation at 1-866-293-6349 and/or consult the company website at www.coloplast.ca.

Peristeen®

The Peristeen® Anal Irrigation System is intended to instill water into the colon through a rectal catheter, which incorporates an inflatable balloon, inserted into the rectum to promote evacuation of the contents of the lower colon. The Peristeen Anal Irrigation System is indicated for use by children (3 years - <12 years old), adolescent (12 years - < 18 years old), transitional adolescent (18 - <21 years old) and adult patients with neurogenic bowel dysfunction who suffer from fecal incontinence, chronic constipation, and/or time-consuming bowel management procedures.

Contraindications: Peristeen Anal Irrigation must not be used in the following situations: known obstruction of the large bowel due to strictures or tumors, acute inflammatory bowel disease, diverticulitis, complex diverticular disease, abdominal or anal surgery within the last 3 months, spinal shock phase, and in patients who are pregnant and have not used the system before (if the patient is pregnant and has never used anal irrigation before, they should not start the irrigation procedure during pregnancy).

Warning: Anal irrigation procedure should always be carried out with care. Bowel perforation is an extremely rare, but serious and potentially lethal complication to anal irrigation and will require immediate admission to a hospital, often requiring surgery.

See the device manual for detailed information regarding the procedure, contraindications, warnings, precautions, and potential complications/adverse events. For further information, call Coloplast Canada Corporation at 1-866-293-6349 and/or consult the company website at www.coloplast.ca.

SpeediCath®

SpeediCath® catheters are indicated for use by patients with chronic urine retention and patients with a post void residual volume (PVR) due to neurogenic and non-neurogenic voiding dysfunction. The catheter is inserted into the urethra to reach the bladder allowing urine to drain. There is a separate SpeediCath device intended for either males or females. This is a single-use only product. Patients are to contact their healthcare professional before performing self-catheterization for the first time and follow the instructions carefully. See the device manual for further information or, call Coloplast Canada Corporation at 1-866-293-6349 and/or consult the company website at www.coloplast.ca.

Self-Cath®

Self-Cath® Catheters are intended for use in male, female, and pediatric patients (neonates, infants, children, adolescents, and transitional adolescents) requiring bladder drainage as determined by their physician. The catheter is indicated for those individuals unable to promote a natural urine flow or for those individuals who have a significant volume of residual urine following a natural bladder-voiding episode. Single-use only. See the device manual for detailed information regarding the procedure, contraindications, warnings, precautions, and potential complications/adverse events. For further information, call Coloplast Canada Corporation at 1-866-293-6349 and/or consult the company website at www.coloplast.ca.

Triad™

Triad™ hydrophilic wound dressing is indicated for the local management of pressure and venous stasis ulcers, dermal lesions/injuries, superficial wounds, scrapes, first-and second-degree burns, partial-and full-thickness wounds. Suitable for use on neonates.

Contraindications: Triad™ is contraindicated for third degree burns or infected wounds. See device label for complete information including warnings and precautions. For further questions, call Coloplast Canada Corporation at 1-888-880-8605 and/or consult the company website at www.coloplast.ca.

Product Code	Page
0089	62
112	40
114	40
116	40
208	39
210	39
212	39
214	39
240	39
305	39
306	39
308	39
310	39
408	40
410	40
412	40
414	40
416	40
418	40
0421	32
0422	32
0423	32
0424	32
450	40
460	40
475	39
0505	61
0506	61
0509	32
608	40
610	40
612	40
614	40
808	40
810	40
812	40
814	40
816	40
818	40
0897	57
0901	57
0925	32
1003	34
1061	54
1063	54
1110	34
1114	40
1116	40
1166	54
1450	45
1451	45
1500	33
1644	61
1645	61
1802	58
1814	58

Product Code	Page
1843	58
1907	32
1964	54
1967	54
2041	32
2115	30
2125	30
2151	30
2166	30
2181	30
2320	32
2325	32
2330	32
2340	32
2350	32
2490	27
2501	33
2520	27
2650	32
2655	31
2792	25
2793	25
2801	33
2802	33
2803	33
2807	33
2808	33
2831	24
2832	24
2833	24
2836	25
2837	25
2846	24
2847	24
2882	23
2883	23
3210	31
3215	31
3220	31
3410	50
3412	50
3413	50
3416	50
3420	50
3430	50
3462	50
3470	50
3472	50
3475	50
3483	50
3485	50
3488	50
3705	51
3710	51
3715	51
3740	51

Product Code	Page
3755	51
3760	51
3765	51
3780	51
3900	54
3906	54
4110	52
4112	52
4115	52
4120	52
4140	52
4142	52
4610	52
4612	52
4615	52
4620	52
4625	52
4630	52
4632	52
4635	52
4639	52
4645	52
5030	41
5030H	41
5035	41
5035H	41
5050	44
5062	44
5063	44
5070	44
5100	42
5105	51
5110	42
5130	42
5135	42
5161	43
5163	43
5167	43
5170	44
5174	44
5200	42
5205	42
5210	42
5212	42
5215	42
5221	41
5225	41
5225H	41
5230	41
5230H	41
5235	41
5235H	41
5290	42
5300	42
5390	42
5400	42

Product Code	Page
5490	42
5580	29
5585	29
5590	42
6105	50
6200	42
6300	42
6400	42
7041	61
7056	58
7067	59
7068	59
7086	60
7087	60
7090	59
7094	59
7222	59
7223	59
7224	59
7230	58
7231	58
7233	58
7234	58
7251	62
7262	62
7300	57
7301	57
7565	60
7573	60
7583	33
7690	54
7710	60
7712	60
7725	60
7910	61
7912	61
7915	61
8000	43
8009	30
8100	43
8200	43
8205	43
8300	43
8305	43
8400	43
8500	43
9500	33
9505	33
9506	33
9622	51
9623	51
9625	51
9626	51
9631	51
9632	51
9635	51

Product Code	Page
9641	51
9643	51
10011	18
10012	18
10013	18
10015	18
10016	18
10017	18
10021	18
10022	18
10023	18
10025	18
10026	18
10027	18
10031	18
10032	18
10035	18
10036	18
10041	18
10045	18
10101	20
10102	20
10103	20
10104	20
10105	20
10106	20
10108	20
10111	20
10113	20
10114	20
10116	20
10117	20
10154	19
10155	19
10156	19
10164	19
10165	19
10166	19
10167	19
10451	12
10461	12
10471	12
10472	12
10474	12
10476	12
10477	12
10478	12
10479	12
10481	12
10482	13
10483	13
10484	13
10487	13
10489	12
10502	9
10503	9

Product Code	Page
10504	9
10512	9
10513	9
10515	9
10522	9
10531	9
10551	11
10561	11
10562	11
10563	11
10564	11
10571	11
10582	15
10583	15
10584	15
10585	15
10586	15
10811	21
10821	21
10822	21
10823	21
10824	21
10831	21
10833	21
10834	21
10835	21
10870	14
10871	14
10872	14
10873	14
10874	14
10881	14
10883	14
10884	14
10885	14
10904	20
10905	20
10914	20
10915	20
10918	20
11011	18
11012	18
11015	18
11016	18
11021	18
11022	18
11023	18
11024	18
11025	18
11026	18
11027	18
11028	18
11031	18
11032	18
11035	18
11036	18

Product Code	Page
11041	18
11045	18
11114	19
11115	19
11124	19
11125	19
11126	19
11127	19
11135	19
11136	19
11185	19
11186	19
11194	19
11195	19
11196	19
11197	19
11301	20
11302	20
11303	20
11304	20
11305	20
11306	20
11307	20
11308	20
11312	20
11313	20
11314	20
11315	20
11316	20
11401	10
11402	10
11412	10
11413	10
11421	10
11422	10
11425	10
11432	10
11443	10
11461	10
11462	10
11463	10
11471	10
11472	10
11473	10
11474	10
11482	10
11483	10
11484	10
11486	10
11487	10
11491	10
11492	10
11496	10
11497	10
11498	10
11499	10

Product Code	Page
11500	10
11505	20
11514	20
11515	20
11518	20
11612	20
11615	20
11690	20
11801	22
11802	22
11803	22
11804	22
11811	22
11812	22
11814	22
11815	22
11821	21
11822	21
11824	21
11825	21
11827	21
11828	21
11841	19
11842	19
11844	19
11845	19
11846	19
11851	19
11852	19
11854	19
11855	19
11856	19
12000	31
12010	31
12011	31
12014	16
12020	31
12021	31
12030	31
12032	31
12033	31
12035	31
12036	31
12037	31
12038	31
12039	31
12042	31
12045	31
12046	31
12047	31
12048	31
12049	31
12050	32
12060	32
12061	32
12070	32

Product Code	Page
12072	32
12074	32
12076	32
12110	28
12130	28
12140	28
12144	28
12145	28
12146	28
12147	28
12170	29
12174	29
12175	29
12176	29
12177	29
12201	12
12202	12
12212	12
12213	12
12222	12
12232	12
12243	12
12262	12
12271	12
12272	12
12274	12
12282	12
12284	12
12291	12
12292	12
12297	12
12298	12
12344	25
12346	25
12354	25
12355	25
12356	25
12374	25
12375	25
12376	25
12384	25
12385	25
12386	25
12432	29
12433	29
12471	29
12472	29
12533	27
12595	29
12596	29
12701	23
12702	23
12703	23
12704	23
12705	23
12706	23

Product Code	Page
12707	23
12708	23
12709	23
12716	23
12719	23
12800	28
12802	28
12805	28
12808	28
12810	28
12834	34
12835	34
12836	34
12841	23
12842	23
12844	23
12845	23
12847	23
12848	23
12991	29
12992	29
12993	29
12994	29
12995	29
12996	29
12997	29
13706	27
13840	27
13844	27
13845	27
13846	27
13860	27
13870	27
13874	28
13875	28
13876	28
13924	25
13925	25
13964	25
13965	25
13966	25
13974	25
13975	25
13976	25
13984	25
13985	25
13986	25
14010	34
14015	34
14021	34
14031	34
14041	34
14050	34
14060	34
14070	34
14103	28

Product Code	Page
14106	28
14163	28
14164	28
14165	28
14166	28
14173	28
14174	28
14175	28
14176	28
14196	27
14201	29
14202	29
14224	25
14225	25
14226	25
14227	25
14228	25
14229	25
14232	24
14233	24
14234	24
14235	24
14236	24
14237	24
14238	24
14239	24
14243	24
14246	24
14249	24
14251	24
14252	24
14254	24
14255	24
14258	24
14261	23
14262	23
14263	23
14271	23
14272	23
14273	23
14274	23
14275	23
14276	23
14277	23
14278	23
14281	24
14282	24
14283	24
14291	24
14292	24
14293	24
14294	24
14295	24
14296	24
14301	26
14302	26

Product Code	Page
14303	26
14304	26
14305	26
14306	26
14307	30
14308	30
14309	26
14316	27
14317	27
14319	27
14329	27
14331	27
14332	27
14342	26
14344	26
14346	26
14347	26
14349	26
14355	27
14356	26
14357	26
14359	26
14361	26
14362	26
14364	26
14401	26
14402	26
14403	26
14404	27
14405	27
14406	26
14408	26
14411	27
14412	27
14414	27
14415	27
14422	28
14423	28
14424	28
14441	29
14445	29
14452	29
14453	29
14454	29
14494	25
14498	25
14520	26
14521	26
14524	26
14525	26
14527	26
14528	26
14552	27
14553	27
14554	27
14555	27

Product Code	Page
14601	26
14602	26
14603	26
14642	26
14643	26
14644	26
14681	30
14682	30
14691	30
14692	30
14701	29
14702	29
15403	22
15404	22
15405	22
15406	22
15407	22
15408	22
15413	22
15414	22
15415	22
15441	22
15450	22
15470	22
15472	22
15473	22
15474	22
15480	22
15482	22
15483	22
15484	22
15521	21
15522	21
15523	21
15524	21
15531	21
15532	21
15533	21
15534	21
15605	21
15606	21
15607	21
15608	21
15633	21
15634	21
15635	21
15695	21
15696	21
15821	21
15981	21
16301	14
16305	14
16306	14
16310	14
16311	14
16312	14

Product Code	Page
16313	14
16315	14
16322	14
16323	14
16325	14
16326	14
16327	14
16330	14
16331	14
16332	14
16333	14
16335	14
16345	14
16346	14
16347	14
16350	14
16351	14
16352	14
16353	14
16471	11
16472	11
16481	11
16483	11
16485	11
16486	11
16487	11
16488	11
16491	11
16705	13
16706	13
16715	13
16716	13
16725	13
16726	13
16727	13
16731	13
16732	13
16733	13
16734	13
16737	13
16738	13
16739	13
16748	13
16755	13
16756	13
16757	13
16760	13
16761	13
16762	13
16763	13
16767	13
16768	13
16771	13
16772	13
16773	13
16778	13

Product Code	Page
16805	15
16810	15
16825	15
16826	15
16827	15
16835	15
16836	15
16837	15
16840	15
16841	15
16842	15
16855	15
16856	15
16857	15
16865	15
16866	15
16867	15
16870	15
16871	15
16872	15
16901	9
16903	9
16911	9
16914	9
16915	9
16921	9
16923	9
16924	9
16931	9
16941	9
16942	9
16951	9
16954	9
16955	9
16961	9
16971	9
18151	15
18152	15
18160	15
18161	15
18162	15
18163	15
18170	15
18171	15
18172	15
18173	15
18180	15
18181	15
18182	15
18185	15
18300	13
18301	13
18302	13
18305	13
18310	13
18312	13

Product Code	Page	Product Code	Page	Product Code	Page	Product Code	Page	Product Code	Page
18313	13	18651	12	22040	41	33146	53	505711	42
18314	13	18655	12	22121	41	33283	53	505731	42
18315	13	18656	12	22125	41	33285	53	505741	42
18320	13	18660	16	22130	41	33400	49	505751	42
18321	13	18661	16	22135	41	33401	49	506110	43
18322	13	18665	16	25501	43	33406	49	506130	43
18323	13	18666	16	27408	37	33434	49	506140	43
18324	13	18670	16	27410	37	33435	49	506150	43
18330	14	18671	16	27412	37	33436	49	506200	43
18331	14	18675	16	27414	37	33437	49	506220	43
18332	14	18676	16	27416	37	33438	49	506240	43
18333	14	18680	16	27418	37	33444	49	506260	43
18340	14	18681	16	27490	37	33445	49	C5200	42
18341	14	18690	16	27492	37	33446	49	C5210	42
18342	14	18691	16	27494	37	33452	49	C5310	42
18343	14	18700	16	27496	37	33453	49	C5410	42
18344	14	18705	16	27508	37	33530	53	C6100	42
18401	15	18710	17	27510	37	33533	53		
18402	15	18711	17	27512	37	33536	53		
18403	15	18712	17	27514	37	33537	53		
18404	15	18713	17	27516	37	33539	53		
18405	15	18720	17	27706	37	33547	53		
18500	9	18721	17	27708	37	33548	53		
18510	9	18722	17	27710	37	39000*	49		
18511	9	18723	17	28110	38	39001*	49		
18512	9	18724	17	28112	38	39010*	49		
18513	9	18725	17	28114	38	39636	49		
18520	9	18730	17	28422	39	39637	49		
18521	9	18731	17	28520	39	39638	49		
18522	9	18735	17	28522	39	39639	49		
18523	9	18747	17	28524	39	39640	49		
18524	9	19001	22	28578	38	39644	49		
18530	9	19003	22	28580	38	39645	49		
18531	9	19010	22	28582	38	39650	49		
18532	9	19011	22	28584	38	39651	49		
18600	11	19020	22	28692	38	39652	49		
18601	11	19021	22	28810	38	61761	57		
18602	11	19030	22	28812	38	61762	57		
18603	11	19040	19	28814	38	61768	57		
18610	11	19041	19	28916	37	67145	57		
18611	11	19042	19	28920	37	383004	44		
18612	11	19050	19	28922	37	504311	42		
18613	11	19051	19	28924	37	504331	42		
18620	11	19052	19	29121	45	504341	42		
18621	11	21034	43	29122	45	504351	42		
18622	11	21054	43	29123	45	505311	42		
18630	11	21346	44	29124	45	505321	42		
18631	11	21356	44	29125	45	505331	42		
18632	11	21365	34	29126	45	505341	42		
18640	11	22025	41	29127	45	505421	42		
18641	11	22028	41	33110	53	505431	42		
18642	11	22030	41	33115	53	505441	42		
18650	12	22035	41	33120	53	505451	42		

Coloplast® Care

Support for your ostomy or continence care needs



Website and emails

A dedicated website and regular emails with articles and reliable advice

Self-Assessment Tools

Helping users keep track of their ostomy or continence health through regular self-assessment

Phone support

Advisors ready to answer questions and provide support by phone

Product guidance

Access to product samples on request

Enrolling in Care is easy!



Online at www.coloplast.ca/careform



By phone at 1-866-293-6349



Request Coloplast Care enrollment forms from your Coloplast Territory Manager

Velcro® is a registered trademark of Velcro Industries B.V.

Coloplast Canada Corporation, 2401 Bristol Circle, Suite A205, Oakville, ON
www.coloplast.ca The Coloplast logo is a registered trademark of Coloplast A/S.
©2021-03 All rights reserved Coloplast A/S

