

Coloplast Product Guide



Ostomy, Continence, Wound and Skin Care
Effective June, 2022





"I was so excited to try this technique because it would give me the independence I desperately needed. It was so easy I could hardly believe it. As long as I have my standard SpeediCath® or my SpeediCath® Compact Set (which is a catheter and bag all in one) I am good to go. Both of these products are small, require no lubrication, and are easy to use!"

Brittney
Alberta, Canada
Speedicath® User

OSTOMY

Appliances

SenSura® Mio Click

Ostomy Barrier and Pouch - Two Piece System

SenSura® Mio Click Skin Barrier	9
SenSura® Mio Click Drainable Pouch	10
SenSura® Mio Click Drainable Pouch	10
SenSura® Mio Click Closed Pouch	10
SenSura® Mio Click Urostomy Pouch	11
SenSura® Mio Click High Output Pouch	11

SenSura® Mio Flex

Ostomy Barrier and Pouch - Two Piece System

SenSura® Mio Flex Skin Barrier	12
SenSura® Mio Flex Drainable Pouch	12
SenSura® Mio Flex Closed Pouch	13
SenSura® Mio Flex Urostomy Pouch	13
SenSura® Mio Flex High Output Pouch	13

SenSura® Mio

Ostomy Barrier and Pouch - One Piece System

SenSura® Mio Drainable Pouch	13
SenSura® Mio Drainable Pouch	14
SenSura® Mio Drainable Pouch	15
SenSura® Mio Closed Pouch	16
SenSura® Mio Closed Pouch	17
SenSura® Mio One Piece Urostomy	18
SenSura® Mio One Piece High Output Pouch	18
SenSura® Mio One Piece Post-Op Pouch	18

SenSura® Mio Baby

Ostomy Barrier and Pouch - Two Piece System

SenSura® Mio Baby Flex Skin Barrier	19
SenSura® Mio Baby Flex Drainable Pouch	19
SenSura® Mio Baby Adhesive Remover Wipe	19

SenSura® Mio Kids

Ostomy Barrier and Pouch - Two Piece System

SenSura® Mio Kids Flex Skin Barrier	19
SenSura® Mio Kids Flex Drainable Pouch	19
SenSura® Mio Kids Flex Urostomy Pouch	19

Ostomy Barrier and Pouch - One Piece System

SenSura® Mio Kids Drainable Pouch 1 Piece	20
SenSura® Mio Kids Pouch 1 Piece Urostomy	20

SenSura® Click

Ostomy Barrier and Pouch- Two Piece System

SenSura® Click Skin Barrier	20
SenSura® Click Xpro Skin Barrier	21
SenSura® Click Drainable Pouch	21
SenSura® Click Closed Pouch	21
SenSura® Click Urostomy Pouch	22
SenSura® High Output Pouch	22

SenSura® Flex

Ostomy Barrier and Pouch - Two Piece System

SenSura® Flex Skin Barrier	22
SenSura® Flex Xpro Skin Barrier	23
SenSura® Flex Drainable Pouch	23
SenSura® Flex Closed Pouch	23
SenSura® Flex Urostomy Pouch	23

SenSura® - One Piece

Ostomy Barrier and Pouch - One Piece System

SenSura® Drainable Pouch with Easicle™	24
SenSura® Xpro Drainable Pouch	24
SenSura® Closed Pouch	25
SenSura® Urostomy Pouch	25
SenSura® One-Piece High Output Pouch	25
SenSura® Post-Op Pouch	25

Assura®

Ostomy Barrier and Pouch - Two Piece System

Assura® Skin Barrier	26
Assura® Skin Barrier - Extra Extended Wear	27
Assura® 2 pc. High Output Drainable Pouch	27
Assura® 2 pc. High Output Drainable System	27
Assura® Drainable Pouch Clamp Closure	28
Assura® Drainable Pouch with EasiClose™	28
Assura® Closed Pouch	28
Assura® Urostomy Pouch	28

Easiflex®

Ostomy Barrier and Pouch- Two Piece System

Easiflex® Skin Barrier	29
Easiflex® Skin Barrier - Extra Extended Wear	29
Easiflex® Drainable Pouch with EasiClose™	29
Easiflex® Closed Pouch	30
Easiflex® Urostomy Pouch	30

Assura® - One Piece

Ostomy Barrier and Pouch - One Piece System

Assura® Drainable Pouch Clamp Closure	30
Assura® Drainable Pouch with EasiClose™	30
Assura® Drainable Pouch with EasiClose™	31
Assura® Closed Pouch	31
Assura® 1-Piece Post-Op Pouch	31
Assura® Closed Pouch	32
Assura® Urostomy Pouch	32
Assura® Urostomy Pouch - Extra Extended Wear	32

Assura® - Pediatric - Two Piece

Ostomy Barrier and Pouch - Two Piece System - Pediatric

Assura® Pediatric Skin Barrier	32
Assura® Pediatric Drainable Pouch	33
Assura® Pediatric Closed Pouch	33

EasiFlex® - Pediatric

Ostomy Barrier and Pouch - Two Piece System - Pediatric	
Easiflex® Pediatric Skin Barrier	33
Easiflex® Pediatric Drainable Pouch	33

Assura® - Pediatric - One Piece

Ostomy Barrier and Pouch - One Piece System - Pediatric	
Assura® Pediatric Drainable Pouch	33
Assura® Pediatric Closed Pouch	33
Assura® Pediatric Urostomy Pouch	33

Accessories

Ostomy Accessories

Brava® Protective Sheets.....	33
Brava® Protective Rings	34
Brava® Protective Rings Wide.....	34
Brava® Protective Rings Extra Wide.....	34
Brava® Protective Ring Convex.....	34
Brava® Mouldable Rings.....	34
Brava® Skin Barrier Spray and Wipe	34
Brava® Adhesive Remover Spray and Wipe.....	34
Brava® Barrier Cream	34
Brava® Ostomy Strip Paste.....	34
Brava® Elastic Barrier Strips.....	35
Brava® Ostomy Powder.....	35
Brava® Mio Belt for All SenSura Mio Products.....	35
Brava® Belt for SenSura Original + Assura Products.....	35
Brava® Lubricating Deodorant.....	35
Brava® Alcohol Free Paste	35
Ostomy Paste (contains alcohol).....	35
Skin Barrier Rings.....	35
Coloplast Prep™ Protective Liquid Skin Barrier	35
Coloplast Filtrodor® Pouch Filter	35
Coloplast Pouch Clamp.....	36
Ostomy Scissors	36

Hex-On®

Odour antagonist	
Hex-On®	36

Assura® - Ostomy Cap

Coloplast Uro Minicap (100 ml).....	36
Coloplast Stoma Cap	36
Coloplast Stoma Cap - One Piece	36

Assura® Irrigation Products

Assura® Irrigation Set	37
Assura® Irrigation Accessories.....	37
Assura® Irrigation Sleeves.....	37
Coloplast Urostomy Night Bag.....	37

Coloplast Fistula & Wound Management System

Pouch	38
Flexible Lids	38
Accessories	38

CONTINENCE

Hydro Intermittent Catheters

SpeediCath® Flex

Hydrophilic Intermittent Catheter	
SpeediCath® Flex	41

SpeediCath® Standard

Hydrophilic Intermittent Catheter	
SpeediCath® Standard Male Straight	41
SpeediCath® Standard Male Coude	41
SpeediCath® Standard Female Straight	41
SpeediCath® Standard Pediatric Straight	41

SpeediCath® Compact Male

Hydrophilic Intermittent Catheter	
SpeediCath® Compact Male.....	42

SpeediCath® Compact Eve

Hydrophilic Intermittent Catheter	
SpeediCath® Compact Eve.....	42

SpeediCath® Compact Female

Hydrophilic Intermittent Catheter	
SpeediCath® Compact Female.....	42
SpeediCath® Compact Female Plus	42

SpeediCath® Compact Set

Hydrophilic Intermittent Catheter	
SpeediCath® Compact Set Male	43
SpeediCath® Compact Set Female	43

Self-Cath®

Intermittent Catheter

Self-Cath® Female Straight	43
Self-Cath® Pediatric	43
Self-Cath® Extension Tube	43
Self-Cath® Male Straight	44
Self-Cath® Male Soft Straight	44
Self-Cath® Male Straight Curved Packaging	44
Self-Cath® Male Coude Tapered with Guide Stripe	44
Self-Cath® Male Coude Olive with Guide Stripe	44
Self-Cath® Male Closed System Straight	44

External

Conveen® Optima

Male External Catheter - Silicone

Conveen® Optima Standard	45
Conveen® Optima Short	45

Conveen® Security+

Male External Catheter - SEBS*

Conveen® Security+ Self Sealing	45
Conveen® Security+ Adhesive Liner	45

Clear Advantage® / Freedom Clear™

Male External Catheter - Silicone

Clear Advantage® with Aloe Vera	46
Freedom Clear™	46
Freedom Clear™ "LS" Long Seal	46
Freedom Clear™ "SS" Sport Sheath	46

Conveen®

Male External Catheter - Latex

Conveen® Self-Sealing	46
Conveen® Adhesive Liner	46

Active Cath®

Male External Catheter - Latex

Active Cath®	47
--------------------	----

Freedom Cath®

Male External Catheter - Latex

Freedom Clear®	47
----------------------	----

Conveen® Active

Urine Leg Bag

Conveen® Active Leg Bag, Non Sterile	47
--	----

Conveen® Security+

Urine Leg bag

Conveen® Security+ Leg Bag, Clamp Outlet, Straps, Non Sterile	47
Conveen® Security+ Leg Bag, Clamp Outlet, Sample Port, Straps, Sterile	47
Conveen® Security+ Leg Bag, Lever Outlet, Straps, Sterile	47

Conveen® Standard

Day / Night Urine Bag

Conveen® Standard Day / Night Bag, Clamp Outlet, Non Sterile	48
Conveen® Standard Day / Night Bag, Clamp Outlet, Sample Port, Sterile	48

Conveen® Contour

Urine Leg Bag

Conveen® Contour Leg Bag, Clamp Outlet, Straps, Non Sterile	48
---	----

Conveen® Night Bag

Urine Night Bag

Conveen® Bedside Night Bag, Non-Sterile	48
Conveen® Bedside Night Bag, Sterile	48

Accessories

Conveen® Bag Hanger	48
Conveen® Leg Bag Strap	48
Conveen® Catheter retaining Strap	48

Peristeen

Peristeen® Anal Plug

Peristeen® Anal Plug.....	49
---------------------------	----

Peristeen Plus® Anal Irrigation

Peristeen Plus® Anal Irrigation System	49
Peristeen Plus® Anal Irrigation System - Small	49
Peristeen Plus® Anal Irrigation Accessory unit	49
Peristeen Plus® Anal Irrigation Accessory unit - Small	49
Peristeen Plus® Anal Irrigation Rectal Catheter	49
Peristeen Plus® Anal Irrigation Leg Strap	49
Peristeen Plus® Anal Irrigation Tube	49

WOUND

Biatain® Silicone

Foam Dressing with Silicone Adhesive

Biatain® Silicone	53
-------------------------	----

Biatain® Silicone Non-Border

Biatain® Silicone Non-Border.....	53
-----------------------------------	----

Biatain® Silicone Lite

Thin Foam Dressing with Silicone Adhesive

Biatain® Silicone Lite.....	53
-----------------------------	----

Biatain® Silicone Ag

Foam Dressing with Silicone Adhesive and Silver

Biatain® Silicone Ag.....	54
---------------------------	----

Biatain®

Foam Dressing

Biatain® Adhesive	54
Biatain® Non-Adhesive	54

Biatain® Soft Hold

Foam Dressing with Gentle Adhesion

Biatain® Soft Hold.....	55
-------------------------	----

Biatain® Ag

Foam Dressing with Silver

Biatain® Ag Adhesive	55
Biatain® Ag Non-Adhesive	55

Biatain® Fiber

Biatain® Fiber	55
----------------------	----

Biatain® Alginate

Calcium Alginate Dressing

Biatain® Alginate.....	56
------------------------	----

Biatain® Alginate Ag

Calcium Alginate Dressing with Silver

Biatain® Alginate Ag.....	56
---------------------------	----

Biatain® Ibu

Foam Dressing with Ibuprofen

Biatain® Ibu Non-Adhesive	56
Biatain® Ibu Soft-Hold	56

Biatain® Super

Hydrocapillary Pad Dressing

Biatain® Super Adhesive 57

Biatain® Super Non-Adhesive 57

Comfeel® Plus

Hydrocolloid Dressing

Comfeel® Plus 57

Comfeel® Plus Transparent

Hydrocolloid Dressing

Comfeel® Plus Transparent..... 57

Triad™

Hydrophilic Wound Dressing

Triad™ 58

Purilon® Gel

Hydrogel with natural ingredients

Purilon® Gel 58

Woun'Dres®

Collagen Hydrogel

Woun'Dres® Collagen Hydrogel 58

Sea-Clens®

Wound Cleanser

Sea-Clens® Wound Cleanser 58

SKIN

Sproam®

No-Rinse Spray or Foam Cleanser

Sproam® 61

Bedside-Care® Sensitive Skin

No-Rinse Foam Cleanser

Bedside-Care® Sensitive Skin 61

Bedside-Care® Foam

No-Rinse Foam Cleanser

Bedside-Care® 61

Bedside-Care®

No-Rinse Spray Cleanser

Bedside-Care® 61

EasiCleanse® Bath

No-Rinse Wash Cloth

EasiCleanse® Bath 62

Gentle Rain® Extra Mild

Rinse Cleanser

Gentle Rain® Extra Mild 62

Gentle Rain® Extra Mild Fragrance Free

Rinse Cleanser

Gentle Rain® Extra Mild Fragrance Free 62

Atrac-Tain®

Moisturizer

Atrac-Tain® 62

Sween® 24

Moisturizer

Sween® 24 63

Sween Cream®

Moisturizer

Sween Cream® 63

Sween® Lotion Unscented

Moisturizer

Sween® Lotion Unscented 63

Critic-Aid® Clear

Skin Protectant

Critic-Aid® Clear 64

Baza® Protect II

Skin Protectant

Baza® Protect II 64

Baza Cleanse and Protect®

Skin Protectant

Baza Cleanse and Protect® 64

InterDry®

Moisture-wicking textile with antimicrobial silver complex

InterDry® 65

Isagel®

Hand Sanitizer Gel

Isagel® 65

Sween® Body Powder

All-body powder

Sween® Body Powder 65

Accessories

Coloplast® Wall Dispenser for Bag 66

Coloplast® Wall Dispenser Drip Tray 66

Embrace™ Dispenser Bracket 66


Ostomy

SenSura® Mio Click

Ostomy Barrier and Pouch - Two Piece System

Click mechanical coupling with a lock ring for security. Provides a tight, secure fit to individual body shapes with optimal discretion. Elastic adhesive, neutral grey textile and full-circle filter. All 2pc Mio Click Flanges have 4 belt attachment points.

SenSura® Mio Click Skin Barrier					
Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	40mm (1 9/16")	10-35mm (3/8 - 1 3/8")	10502	5
		50mm (2")	10-45mm (3/8 - 1 3/4")	10512	5
		60mm (2 3/8")	10-55mm (3/8 - 2 1/8")	10522	5
		70mm (2 3/4")	10-65mm (3/8 - 2 1/2")	10531	5
Pre-Cut	Non-Convex	40mm (1 9/16")	20mm (3/4")	10503	5
			22mm (7/8")	10506 <i>New!</i>	5
			25mm (1")	10504	5
			28mm (1 1/8")	10505 <i>New!</i>	5
		50mm (2")	30mm (1 1/8")	10513	5
			32mm (1 1/4")	10514 <i>New!</i>	5
			35mm (1 3/8")	10515	5
			38mm (1 1/2")	10516 <i>New!</i>	5
		60mm (2 3/8")	40mm (1 9/16")	10523 <i>New!</i>	5
			45mm (1 3/4")	10524 <i>New!</i>	5
			50mm (2")	10525 <i>New!</i>	5
Cut-to-Fit	Convex Light	40mm (1 9/16")	10-20mm (3/8 - 3/4")	16901	5
		50mm (2")	15-30mm (5/8 - 1 1/8")	16911	5
		60mm (2 3/8")	15-40mm (5/8 - 1 9/16")	16921	5
		70mm (2 3/4")	15-53mm (5/8 - 2 1/16")	16931	5
Pre-Cut	Convex Light	40mm (1 9/16")	18mm (3/4")	16903	5
			21mm (13/16")	16913 <i>New!</i>	5
			25mm (1")	16914	5
		50mm (2")	28mm (1 1/8")	16915	5
			31mm (1 1/4")	16923	5
			35mm (1 3/8")	16924	5
		70mm (2 3/4")	38mm (1 1/2")	16933 <i>New!</i>	5
			41mm (1 5/8")	16934 <i>New!</i>	5
Cut-to-Fit	Convex Deep	40mm (1 9/16")	10-20mm (3/8 - 3/4")	16941	5
		50mm (2")	15-30mm (5/8 - 1 1/8")	16951	5
		60mm (2 3/8")	15-40mm (5/8 - 1 9/16")	16961	5
		70mm (2 3/4")	15-53mm (5/8 - 2 1/16")	16971	5
Pre-Cut	Convex Deep	40mm (1 9/16")	18mm (3/4")	16942	5
			21mm (13/16")	16953 <i>New!</i>	5
			25mm (1")	16954	5
		50mm (2")	28mm (1 1/8")	16955	5
			31mm (1 1/4")	16963 <i>New!</i>	5
			35mm (1 3/8")	16964 <i>New!</i>	5
		70mm (2 3/4")	38mm (1 1/2")	16973 <i>New!</i>	5
			41mm (1 5/8")	16974 <i>New!</i>	5



SenSura® Mio Click Skin Barrier

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Concave	40mm (1 9/16") ●	10-35mm (3/8 - 1 3/8")	18500	5
		50mm (2") ●	10-45mm (3/8 - 1 3/4")	18510	5
		60mm (2 3/8") ●	10-55mm (3/8 - 2 1/8")	18520	5
		70mm (2 3/4") ●	10-65mm (3/8 - 2 1/2")	18530	5
Pre-Cut	Concave	40mm (1 9/16") ●	20mm (3/4")	18501 <i>New!</i>	5
			25mm (1")	18511	5
			30mm (1 1/8")	18512	5
		50mm (2") ●	35mm (1 3/8")	18513	5
			30mm (1 1/8")	18521	5
			35mm (1 3/8")	18522	5
		60mm (2 3/8") ●	40mm (1 9/16")	18523	5
			45mm (1 3/4")	18524	5
			45mm (1 3/4")	18531	5
		70mm (2 3/4") ●	50mm (2")	18532	5



SenSura® Mio Click Drainable Pouch with Easiclese™ WIDE-Outlet

Size	Colour	Coupling Size		Product Code	Units
Midi 26cm (10 1/4")	Neutral grey with inspection window	40mm (1 9/16")		11461	10
		50mm (2")		11471	10
Maxi 28cm (11")	Neutral grey with inspection window	40mm (1 9/16")		11462	10
		50mm (2")		11472	10
		60mm (2 3/8")		11482	10
		70mm (2 3/4")		11486	10
		40mm (1 9/16")	 Narrow outlet with Ileo filter	11432	10
		50mm (2")	 Narrow outlet with Ileo filter	11443	10
Maxi 28cm (11")	Transparent	40mm (1 9/16")		11463	10
		50mm (2")		11473	10
		50mm (2")	 No filter	11474	10
		60mm (2 3/8")		11483	10
		60mm (2 3/8")	 No filter	11484	10
		70mm (2 3/4")		11487	10



SenSura® Mio Click Closed Pouch

Size	Colour	Coupling Size	Product Code	Units
Midi 18cm (7")	Neutral grey with inspection window	40mm (1 9/16") ●	11401	30
		50mm (2") ●	11412	30
		60mm (2 3/8") ●	11421	30
Maxi 22cm (8 1/2")	Neutral grey with inspection window	40mm (1 9/16") ●	11402	30
		50mm (2") ●	11413	30
		60mm (2 3/8") ●	11422	30
		70mm (2 3/4") ●	11425	30



SenSura® Mio Click Urostomy Pouch

Size	Colour	Coupling Size		Product Code	Units
Midi 24cm (9 1/2")	Neutral grey	50mm (2")	●	11496	10
		40mm (1 9/16")	●	11491	10
	Neutral grey	50mm (2")	●	11497	10
	Neutral grey with inspection window	60mm (2 3/8")	●	11499	10
Maxi 26cm (10 1/4")		40mm (1 9/16")	●	11492	10
	Transparent	50mm (2")	●	11498	10
		60mm (2 3/8")	●	11500	10



SenSura® Mio Click High Output Pouch

Outlet	Colour	Coupling Size	Volume	Product Code	Units
		40mm (1 9/16")	● 1160 mL	18600	10
	Neutral grey with inspection window	50mm (2")	● 1120 mL	18601	10
		60mm (2 3/8")	● 1080 mL	18602	10
Wide-lumen Hard Outlet 32cm (12 11/16")		70mm (2 3/4")	● 1040 mL	18603	10
		40mm (1 9/16")	● 1160 mL	18610	10
	Transparent	50mm (2")	● 1120 mL	18611	10
		60mm (2 3/8")	● 1080 mL	18612	10
		70mm (2 3/4")	● 1040 mL	18613	10
Mini Pouch Wide-lumen Hard Outlet 24cm (9 1/2")	Neutral grey with inspection window	40mm (1 9/16")	● 650 mL	18620	10
		50mm (2")	● 581 mL	18621	10
		60mm (2 3/8")	● 546 mL	18622	10
		50mm (2")	● 1120 mL	18630	10
	Neutral grey with inspection window	60mm (2 3/8")	● 1080 mL	18631	10
		70mm (2 3/4")	● 1040 mL	18632	10
Soft Outlet Universal Connector 32cm (12 11/16")		50mm (2")	● 1120 mL	18640	10
	Transparent	60mm (2 3/8")	● 1080 mL	18641	10
		70mm (2 3/4")	● 1040 mL	18642	10



Pouches with Wide Lumen Hard Outlet will connect to Coloplast drainage bag 14010

SenSura[®] Mio Flex

Ostomy Barrier and Pouch - Two Piece System

Flex adhesive coupling allows flexibility and freedom of movement. Provides a tight, secure fit to individual body shapes with optimal discretion. Elastic adhesive, neutral grey textile and full-circle filter. Only Mio Flex Convex Flanges have 4 belt attachment points.

SenSura® Mio Flex Skin Barrier						
Cut	Convexity	Coupling Size		Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	35mm (1 3⁄8")	●	10-33mm (3⁄8 - 1 1⁄4")	10551	10
		50mm (2")	●	10-48mm (3⁄8 - 2")	10561	10
		70mm (2 3⁄4")	●	10-68mm (3⁄8 - 2 3⁄4")	10571	10
		90mm (3 1⁄2")	●	10-88mm (3⁄8" - 3 1⁄2")	18390 <i>New!</i>	10
Pre-Cut	Non-Convex			25mm (1")	10562	10
				28mm (1 1⁄8")	10566 <i>New!</i>	10
				30mm (1 1⁄4")	10563	10
		50mm (2")	●	32mm (1 1⁄4")	10567 <i>New!</i>	10
				35mm (1 3⁄8")	10564	10
				38mm (1 1⁄2")	10568 <i>New!</i>	10
				40mm (1 9⁄16")	10565 <i>New!</i>	10
		70mm (2 3⁄4")	●	45mm (1 3⁄4")	10572 <i>New!</i>	5
				50mm (2")	10573 <i>New!</i>	5
		Cut-to-Fit	Convex Light	35mm (1 3⁄8")	●	10-20mm (3⁄8 - 3⁄4")
50mm (2")	●			15-30mm (5⁄8 - 1 1⁄8")	16481	5
				15-40mm (5⁄8 - 1 9⁄16")	16483	5
70mm (2 3⁄4")	●			15-53mm (5⁄8 - 2 1⁄16")	16491	5
Pre-Cut	Convex Light	35mm (1 3⁄8")	●	18mm (3⁄4")	16472	5
				21mm (13⁄16")	16484 <i>New!</i>	5
				25mm (1")	16485	5
		50mm (2")	●	28mm (1 1⁄8")	16486	5
				31mm (1 1⁄4")	16487	5
				35mm (1 3⁄8")	16488	5
		70mm (2 3⁄4")	●	38mm (1 1⁄2")	16493 <i>New!</i>	5
				41mm (1 5⁄8")	16494 <i>New!</i>	5



SenSura® Mio Flex Drainable Pouch with Easyclose™ WIDE-Outlet						
Size	Colour	Coupling Size			Product Code	Units
Midi 26cm (10 1⁄4")	Neutral grey with inspection window	50mm (2")	●		12271	10
		35mm (1 3⁄8")	●	Narrow outlet with ileo filter	12232	10
		35mm (1 3⁄8")	●		12262	10
Maxi 28cm (11")	Neutral grey with inspection window	50mm (2")	●	Narrow outlet with ileo filter	12243	10
		50mm (2")	●		12272	10
		70mm (2 3⁄4")	●		12282	10
	Transparent	50mm (2")	●	No filter	12274	10
		70mm (2 3⁄4")	●	No filter	12284	10
		90mm (3 1⁄2")	●		18386 <i>New!</i>	10





SenSura® Mio Flex Closed Pouch

Size	Colour	Coupling Size	Product Code	Units
Midi 18cm (7")	Neutral grey with inspection window	35mm (1 3/8") ●	12201	30
		50mm (2") ●	12212	30
Maxi 22cm (8 1/2")	Neutral grey with inspection window	35mm (1 3/8") ●	12202	30
		50mm (2") ●	12213	30
		70mm (2 3/4") ●	12222	30
		90mm (3 1/2") ●	18389 <i>New!</i>	10



SenSura® Mio Flex Urostomy Pouch

Size	Colour	Coupling Size	Product Code	Units
Maxi 26cm (10 1/4")	Neutral grey	35mm (1 3/8") ●	12291	10
		50mm (2") ●	12297	10
	Transparent	35mm (1 3/8") ●	12292	10
		50mm (2") ●	12298	10



SenSura® Mio Flex High Output Pouch

Outlet	Colour	Coupling Size	Volume	Product Code	Units
Wide-lumen Hard Outlet 32cm (12 11/16")	Neutral grey with inspection window	50mm (2") ●	980 mL	18650	10
		70mm (2 3/4") ●	906 mL	18651	10
	Transparent	50mm (2") ●	980 mL	18655	10
		70mm (2 3/4") ●	906 mL	18656	10
Soft outlet 32 cm (12 1/2")	Transparent	50mm (2") ●	1080 mL	18657 <i>New!</i>	10
	Neutral grey	50mm (2") ●	1080 mL	18652 <i>New!</i>	10



Pouches with Wide Lumen Hard Outlet will connect to Coloplast drainage bag 14010

SenSura® Mio

Ostomy Barrier and Pouch - One Piece System

Elastic adhesive, neutral grey textile and full-circle filter. Provides a tight, secure fit to individual body shapes with optimal discretion. Only 1pc Convex systems have 4 belt attachment points.

SenSura® Mio Drainable Pouch with Easyclose™ WIDE-Outlet

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 24cm (9 1/2")	Cut-to-Fit	Soft Convex	Neutral grey with inspection window	10-33mm (3/8 - 1 5/16")	16801 <i>New!</i>	10
Midi 26cm (10 1/4")	Cut-to-Fit	Non-Convex	Neutral grey with inspection window	10-45mm (3/8 - 1 3/4")	10461	10
				10-55mm (3/8 - 2 1/8") <i>Narrow outlet with ileo filter</i>	10451	10
	Cut-to-Fit	Non-Convex	Neutral grey with inspection window	10-55mm (3/8 - 2 1/8")	10471	10
				10-55mm (3/8 - 2 1/8")	10481	10
			Transparent	10-55mm (3/8 - 2 1/8")	10489 <i>No filter</i>	10
				25mm (1")	10472	10
	Pre-Cut	Non-Convex	Neutral grey with inspection window	28mm (1 1/8")	10473 <i>New!</i>	10
				30mm (1 1/8")	10474	10
				32mm (1 1/4")	10475 <i>New!</i>	10
				35mm (1 3/8")	10476	10
				38mm (1 1/2")	10480 <i>New!</i>	10
				40mm (1 9/16")	10479	10
				45mm (1 3/4")	10477	10
				50mm (2")	10478	10
Maxi 28cm (11")						



SenSura® Mio Drainable Pouch with Easyclose™ WIDE-Outlet

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 28cm (11")	Pre-Cut	Non-Convex	Transparent	25mm (1")	10482	10
				30mm (1 1/8")	10483	10
				35mm (1 3/8")	10484	10
				40mm (1 9/16")	10487	10
Maxi 28cm (11")	Cut-to-Fit	Convex Soft	Neutral grey with inspection window	15-33mm (5/8 - 1 1/4")	16705	10
				10-50mm (3/8 - 2")	16706	10
			Transparent	15-33mm (5/8 - 1 1/4")	16715	10
				10-50mm (3/8 - 2")	16716	10
	Pre-Cut	Convex Soft	Neutral grey with inspection window	21mm (7/8")	13671 <i>New!</i>	10
				25 mm (1")	13672 <i>New!</i>	10
				28 mm (1 1/8")	13673 <i>New!</i>	10
				31mm (1 1/4")	13674 <i>New!</i>	10
				35 mm (1 3/8")	13675 <i>New!</i>	10
				38 mm (1 1/2")	13676 <i>New!</i>	10
	Cut-to-Fit	Convex Light	Neutral grey with inspection window	10-23mm (3/8 - 7/8")	16725	10
				15-33mm (5/8 - 1 1/4")	16726	10
				10-43mm (3/8 - 1 3/4")	16733	10
			Transparent	10-23mm (3/8 - 7/8")	16737	10
				15-33mm (5/8 - 1 1/4")	16738	10
				10-43mm (3/8 - 1 3/4")	16739	10
	Pre-Cut	Convex Light	Neutral grey with inspection window	10-43mm (3/8 - 1 3/4")	16748 <i>No filter</i>	10
				21mm (13/16")	16730 <i>New!</i>	10
				25mm (1")	16731	10
				28mm (1 1/8")	16732	10
				31mm (1 1/4")	16727	10
				35mm (1 3/8")	16734	10
				38mm (1 1/2")	16735 <i>New!</i>	10
				10-23mm (3/8 - 7/8")	16755	10
				15-33mm (5/8 - 1 1/4")	16756	10
				10-43mm (3/8 - 1 3/4")	16768	10
	Cut-to-Fit	Convex Deep	Neutral grey with inspection window	10-43mm (3/8 - 1 3/4")	16767	10
				10-43mm (3/8 - 1 3/4")	16778 <i>No filter</i>	10
			Transparent	21mm (7/8")	16760	10
				25mm (1")	16761	10
	Pre-Cut	Convex Deep	Neutral grey with inspection window	28mm (1 1/8")	16762	10
				31mm (1 1/4")	16757	10
				35mm (1 3/8")	16763	10
			Transparent	38mm (1 1/2")	16764 <i>New!</i>	10
				21mm (13/16")	16770 <i>New!</i>	10
				25mm (1")	16771	10
				28mm (1 1/8")	16772	10
				31mm (1 1/4")	16773	10



SenSura® Mio Drainable Pouch with Easyclose™ WIDE-Outlet

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 26cm (10 1/4")	Cut-to-fit	Concave	Neutral grey with inspection window	10-40mm (3/8 - 1 9/16")	18300	10
				10-50mm (3/8 - 2")	18310	10
				25mm (1")	18301	10
Maxi 28cm (11")	Pre-Cut	Concave	Neutral grey with inspection window	30mm (1 1/8")	18302	10
				10-40mm (3/8 - 1 9/16")	18305	10
				10-50mm (3/8 - 2")	18312	10
Maxi 28cm (11")	Cut-to-fit	Concave	Transparent	25mm (1")	18313	10
				30mm (1 1/8")	18314	10
				35mm (1 3/8")	18315	10
	Pre-Cut	Concave	Neutral grey with inspection window	10-50mm (3/8 - 2") <i>No filter</i>	18320	10
				10-50mm (3/8 - 2")	18321	10
				25mm (1")	18322	10
	Cut-to-fit	Concave	Transparent	30mm (1 1/8")	18323	10
				35mm (1 3/8")	18324	10
				10-60mm (3/8 - 2 3/8")	18330	10
Maxi 28cm (11")	Pre-Cut	Concave	Neutral grey with inspection window	40mm (1 9/16")	18331	10
				45mm (1 3/4")	18332	10
				50mm (2")	18333	10
	Cut-to-fit	Concave	Transparent	10-60mm (3/8 - 2 3/8") <i>No filter</i>	18340	10
				10-60mm (3/8 - 2 3/8")	18341	10
				40mm (1 9/16")	18342	10
	Pre-Cut	Concave	Transparent	45mm (1 3/4")	18343	10
				50mm (2")	18344	10



SenSura® Mio Closed Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Mini 14cm	Cut-to-fit	Non-Convex	Neutral grey with inspection window	15-45mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	10870	30
Midi 18cm (7")	Cut-to-fit	Non-Convex	Neutral grey with inspection window	15-45mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	10871	30
	Pre-Cut	Non-Convex	Neutral grey with inspection window	25mm (1")	10872	30
				30mm (1 $\frac{1}{8}$ ")	10873	30
				35mm (1 $\frac{3}{8}$ ")	10874	30
Maxi 22cm (8 $\frac{1}{2}$ ")	Cut-to-Fit	Non-Convex	Neutral grey with inspection window	15-45mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	10881	30
	Pre-Cut	Non-Convex	Neutral grey with inspection window	40-55mm (1 $\frac{9}{16}$ - 2 $\frac{1}{8}$ ")	10885	30
				30mm (1 $\frac{1}{8}$ ")	10883	30
				35mm (1 $\frac{3}{8}$ ")	10884	30
Midi 18cm (7")	Cut-to-Fit	 Soft	Neutral grey with inspection window	15-33mm ($\frac{5}{8}$ - 1 $\frac{1}{4}$ ")	16301	10
	Cut-to-Fit	 Light	Neutral grey with inspection window	15-33mm ($\frac{5}{8}$ - 1 $\frac{1}{4}$ ")	16322	10
	Pre-Cut	 Light	Neutral grey with inspection window	31mm (1 $\frac{1}{4}$ ")	16323	10
Maxi 22cm (8 $\frac{1}{2}$ ")	Cut-to-Fit	 Soft	Neutral grey with inspection window	15-33mm ($\frac{5}{8}$ - 1 $\frac{1}{4}$ ")	16305	10
				10-50mm ($\frac{3}{8}$ - 2")	16306	10
			Transparent	10-50mm ($\frac{3}{8}$ - 2")	16315	10
	Pre-Cut	 Soft	Neutral grey with inspection window	25mm (1")	16310	10
				28mm (1 $\frac{1}{8}$ ")	16311	10
				31mm (1 $\frac{1}{4}$ ")	16312	10
35mm (1 $\frac{3}{8}$ ")				16313	10	
Maxi 28cm (11")	Pre-Cut	 Soft	Neutral grey with inspection window	38 mm (1 $\frac{1}{2}$ ")	13684 <i>New!</i>	10
				41 mm (1 $\frac{5}{8}$ ")	13685 <i>New!</i>	10
Maxi 22cm (8 $\frac{1}{2}$ ")	Cut-to-Fit	 Light	Neutral grey with inspection window	10-23mm ($\frac{3}{8}$ - $\frac{7}{8}$ ")	16325	10
				15-33mm ($\frac{5}{8}$ - 1 $\frac{1}{4}$ ")	16326	10
				10-43mm ($\frac{3}{8}$ - 1 $\frac{3}{4}$ ")	16332	10
				Transparent	10-43mm ($\frac{3}{8}$ - 1 $\frac{3}{4}$ ")	16335
	Pre-Cut	 Light	Neutral grey with inspection window	25mm (1")	16330	10
				28mm (1 $\frac{1}{8}$ ")	16331	10
				31mm (1 $\frac{1}{4}$ ")	16327	10
				35mm (1 $\frac{3}{8}$ ")	16333	10
Maxi 28cm (11")	Pre-Cut	 Light	Neutral grey with inspection window	41 mm (1 $\frac{5}{8}$ ")	13683 <i>New!</i>	10
Maxi 22cm (8 $\frac{1}{2}$ ")	Cut-to-Fit	 Deep	Neutral grey with inspection window	10-23mm ($\frac{3}{8}$ - $\frac{7}{8}$ ")	16345	10
				15-33mm ($\frac{5}{8}$ - 1 $\frac{1}{4}$ ")	16346	10
				10-43mm ($\frac{3}{8}$ - 1 $\frac{3}{4}$ ")	16352	10
	Pre-Cut	 Deep	Neutral grey with inspection window	25mm (1")	16350	10
				28mm (1 $\frac{1}{8}$ ")	16351	10
31mm (1 $\frac{1}{4}$ ")				16347	10	
				35mm (1 $\frac{3}{8}$ ")	16353	10



SenSura® Mio Closed Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 18cm (7")	Cut-to-Fit	Concave	Neutral grey with inspection window	10-40mm ($\frac{3}{8}$ - 1 $\frac{9}{16}$ "	18151	10
				10-50mm ($\frac{3}{8}$ - 2")	18160	10
	Pre-Cut	Concave		25mm (1")	18161	10
				30mm (1 $\frac{1}{8}$ "	18162	10
				35mm (1 $\frac{3}{8}$ "	18163	10
Maxi 22cm (8 $\frac{1}{2}$ ")	Cut-to-Fit	Concave	Neutral grey with inspection window	10-40mm ($\frac{3}{8}$ - 1 $\frac{9}{16}$ "	18152	10
Maxi 22cm (8 $\frac{1}{2}$ ")	Cut-to-Fit	Concave	Neutral grey with inspection window	10-50mm ($\frac{3}{8}$ - 2")	18170	10
				25mm (1")	18171	10
	Pre-Cut	Concave		30mm (1 $\frac{1}{8}$ "	18172	10
				35mm (1 $\frac{3}{8}$ "	18173	10
				Maxi 22cm (8 $\frac{1}{2}$ ")	Cut-to-Fit	Concave
Pre-Cut	Concave	40mm (1 $\frac{9}{16}$ "	18181		10	
		45mm (1 $\frac{3}{4}$ "	18182		10	
	Cut-to-Fit	Concave	Transparent		10-60mm ($\frac{3}{8}$ - 2 $\frac{3}{8}$ "	18185



SenSura® Mio One Piece Urostomy

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 26cm (10 3/4")	Cut-to-Fit	Non-Convex	Neutral grey	10-45mm (3/8 - 1 3/4")	10582	10
			Transparent	10-45mm (3/8 - 1 3/4")	10586	10
	Pre-Cut	Non-Convex	Neutral grey	20mm (3/4")	10583	10
				25mm (1")	10584	10
				30mm (1 1/8")	10585	10
Maxi 28cm (11")	Cut-to-Fit	 Convex Soft	Neutral grey	10-50 mm (3/8" - 2")	13678 <i>New!</i>	10
Maxi 26cm (10 3/4")	Cut-to-Fit	 Convex Soft	Neutral grey	15-33mm (5/8 - 1 1/4")	16805	10
			Transparent	10-50mm (3/8 - 2")	16810	10
Maxi 28cm (11")	Cut-to-Fit	 Convex Soft	Transparent	10-33 mm (3/8" - 1 5/16")	13679 <i>New!</i>	10
Maxi 26cm (10 3/4")	Cut-to-Fit	 Convex Light	Neutral grey	10-23mm (3/8 - 7/8")	16825	10
				15-33mm (5/8 - 1 1/4")	16826	10
	Cut-to-Fit	 Convex Light	Transparent	10-23mm (3/8 - 7/8")	16835	10
				15-33mm (5/8 - 1 1/4")	16836	10
				10-43mm (3/8 - 1 3/4")	16837	10
Maxi 28cm (11")	Pre-Cut	 Convex Light	Neutral grey	25 mm (1")	13680 <i>New!</i>	10
	Pre-Cut	 Convex Light	Neutral grey	31mm (1 1/4")	16827	10
			21mm (7/8")	16840	10	
			25mm (1")	16841	10	
			28mm (1 1/8")	16842	10	
	Maxi 26cm (10 3/4")	Cut-to-Fit	 Convex Deep	Neutral grey	10-23mm (3/8 - 7/8")	16855
15-33mm (5/8 - 1 1/4")				16856	10	
10-23mm (3/8 - 7/8")				16865	10	
15-33mm (5/8 - 1 1/4")				16866	10	
10-43mm (3/8 - 1 3/4")				16867	10	
Pre-Cut		 Convex Deep	Neutral grey	31mm (1 1/4")	16857	10
			21mm (7/8")	16870	10	
			Transparent	25mm (1")	16871	10
			28mm (1 1/8")	16872	10	



SenSura® Mio One Piece Urostomy

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 26cm (10 1/4")	Cut-to-Fit	Concave	Neutral grey	10-40mm (3/8- 1 9/16")	18401	10
				20mm (3/4")	18402	10
	Pre-Cut	Concave	Neutral grey	25mm (1")	18403	10
				30mm (1 1/8")	18404	10
	Cut-to-Fit	Concave	Transparent	10-40mm (3/8- 1 9/16")	18405	10



SenSura® Mio One Piece High Output Pouch

Cut/Outlet	Convexity	Colour	Volume	Stoma Size	Product Code	Units
	Non-Convex	Neutral grey with inspection window	630 mL	10-70mm (3/8- 2 3/4")	18660	10
			725 mL	10-100mm (3/8- 4")	18661	10
		Transparent	630 mL	10-70mm (3/8- 2 3/4")	18665	10
			725 mL	10-100mm (3/8- 4")	18666	10
Cut-to-fit/ Soft Outlet Universal Connector 28cm (11 3/16")	Convex Soft	Neutral grey with inspection window	630 mL	10-50mm (3/8- 2")	18670	10
		Transparent	630 mL	10-50mm (3/8- 2")	18671	10
	Convex Light	Neutral grey with inspection window	630 mL	10-43mm (3/8- 1 3/4")	18675	10
		Transparent	630 mL	10-43mm (3/8- 1 3/4")	18676	10



Does not connect to
Coloplast drainage
bag 14010

SenSura® Mio One Piece Post-Op Pouch

Size	Cut	Convexity	Volume	Stoma Size	Product Code	Units
28cm (10 15/18")	Cut-to-Fit	Transparent	680 mL	10-70mm (3/8- 2 3/4")	18680	Sterile 5
29cm (11 1/2")	Cut-to-Fit	Transparent	725 mL	10-100mm (3/8- 4")	18681	Sterile 5
28cm (10 15/18")	Cut-to-Fit	Transparent	680 mL	10-70mm (3/8- 2 3/4")	18690	5
29cm (11 1/2")	Cut-to-Fit	Transparent	725 mL	10-100mm (3/8- 4")	18691	5



SenSura[®] Mio Baby

Ostomy Barrier and Pouch - Two Piece System

SenSura[®] Mio Baby Flex Skin Barrier

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	35mm (1 3/8")	0-35mm	18705	5



SenSura[®] Mio Baby Flex Drainable Pouch

Size	Colour	Coupling Size	Product Code	Units
20cm (7 3/4")	Transparent	35mm (1 3/8")	18700	10



SenSura[®] Mio Baby Adhesive Remover

	Product Code	Units
Single use	12014	30



SenSura[®] Mio Kids

Ostomy Barrier and Pouch - Two Piece System

SenSura[®] Mio Kids Flex Skin Barrier

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	27 mm	0-27 mm	18735	5
		35 mm	0-35 mm	18747	5



SenSura[®] Mio Kids Flex Drainable Pouch

Cut	Colour	Coupling Size	Product Code	Units
Mini 16cm (6 1/8")	Neutral grey	27 mm	18720	10
	Transparent	27 mm	18721	10
Midi 17cm (6 3/4")	Neutral grey	27 mm	18722	10
	Transparent	27 mm	18723	10
	Neutral grey	35 mm	18724	10
	Transparent	35 mm	18725	10



SenSura[®] Mio Kids Flex Urostomy Pouch

Cut	Colour	Coupling Size	Stoma Size	Product Code	Units
Midi 17cm (6 3/4")	Neutral grey	27 mm		18730	10
		35 mm		18731	10



Ostomy Barrier and Pouch - One Piece System

SenSura® Mio Kids Drainable Pouch 1 Piece

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 17cm (6 3/4")	Cut-to-fit	Non-Convex	Neutral grey	0-35mm	18710	10
			Transparent	0-35mm	18711	10



SenSura® Mio Kids Pouch 1 Piece Urostomy

Cut	Convexity	Colour	Stoma Size	Product Code	Units
Cut-to-fit 19cm (7 1/4")	Non-Convex	Transparent	10-35mm	18712	10
		Neutral grey	10-35mm	18713	10



SenSura® Click

Ostomy Barrier and Pouch- Two Piece System

Click mechanical coupling with a lock ring for security. Provides a tight, secure fit to individual body shapes. All 2pc SenSura Click Flanges have 2 belt attachment points.

SenSura® Click Skin Barrier

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	40mm (1 9/16")	10-35mm (3/8 - 1 3/8")	10011	5
		50mm (2")	10-45mm (3/8 - 1 3/4")	10021	5
		60mm (2 3/8")	10-55mm (3/8 - 2 1/8")	10031	5
		70mm (2 3/4")	10-65mm (3/8 - 2 1/2")	10041	5
	Convex Light	40mm (1 9/16")	15-23mm (5/8 - 7/8")	11011	5
		50mm (2")	15-33mm (5/8 - 1 1/4")	11021	5
		60mm (2 3/8")	15-43mm (5/8 - 1 3/4")	11031	5
		70mm (2 3/4")	15-53mm (5/8 - 2 1/16")	11041	5
Pre-Cut	Non-Convex	40mm (1 9/16")	20mm (3/4")	10012	5
			25mm (1")	10013	5
		50mm (2")	30mm (1 1/8")	10022	5
			35mm (1 3/8")	10023	5
		60mm (2 3/8")	40mm (1 9/16")	10032	5
			45mm (1 1/2")	10033	5
	Convex Light	40mm (1 9/16")	21mm (7/8")	11012	5
			25mm (1")	11022	5
		50mm (2")	28mm (1 1/8")	11023	5
			31mm (1 1/4")	11024	5
		60mm (2 3/8")	35mm (1 3/8")	11032	5
			40mm (1 9/16")	11033	5



SenSura® Click Xpro Skin Barrier

(Engineered to offer extra protection for liquid ileostomies and urostomies)

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	40mm (1 9/16")	10-35mm (3/8 - 1 3/8")	10015	5
		50mm (2")	10-45mm (3/8 - 1 3/4")	10025	5
		60mm (2 3/8")	10-55mm (3/8 - 2 1/8")	10035	5
		70mm (2 3/4")	10-65mm (3/8 - 2 1/2")	10045	5
	Convex Light	40mm (1 9/16")	15-23mm (5/8 - 7/8")	11015	5
		50mm (2")	15-33mm (5/8 - 1 1/4")	11025	5
		60mm (2 3/8")	15-43mm (5/8 - 1 3/4")	11035	5
Pre-Cut	Non-Convex	70mm (2 3/4")	15-53mm (5/8 - 2 1/16")	11045	5
		40mm (1 9/16")	20mm (3/4")	10016	5
			25mm (1")	10017	5
		50mm (2")	30mm (1 1/8")	10026	5
	Convex Light		35mm (1 3/8")	10027	5
		60mm (2 3/8")	40mm (1 9/16")	10036	5
		40mm (1 9/16")	21mm (7/8")	11016	5
			25mm (1")	11026	5
		50mm (2")	28mm (1 1/8")	11027	5
			31mm (1 1/4")	11028	5
		60mm (2 3/8")	35mm (1 3/8")	11036	5



SenSura® Click Drainable Pouch with Easyclose™ WIDE-Outlet

Pouch Size	Pouch Colour	Coupling Size		Product Code	Units
Midi 27cm (10 1/2")	Opaque	40mm (1 9/16")		11114	10
		50mm (2")		11115	10
	Opaque	40mm (1 9/16")		11124	10
		50mm (2")		11125	10
50mm (2")		 No filter	11135	10	
60mm (2 3/8")			11126	10	
60mm (2 3/8")		 No filter	11136	10	
70mm -2 3/4"			11127	10	
Maxi 30cm (12")	Transparent	40mm (1 9/16")		11194	10
		50mm (2")		11195	10
		50mm (2")	 No filter	11185	10
		60mm (2 3/8")		11196	10
		60mm (2 3/8")	 No filter	11186	10
		70mm -2 3/4"		11197	10



SenSura® Click Closed Pouch

Pouch Size	Pouch Colour	Coupling Size	Product Code	Units
Midi 18cm (7")	Opaque	40mm (1 9/16")	10154	30
		50mm (2")	10155	30
		60mm (2 3/8")	10156	30
Maxi 22cm (8 1/2")	Opaque	40mm (1 9/16")	10164	30
		50mm (2")	10165	30
		60mm (2 3/8")	10166	30
		70mm (2 3/4")	10167	30



SenSura® Click Urostomy Pouch

Pouch Size	Pouch Colour	Coupling Size	Product Code	Units
Midi 24cm (9 1/4")	Opaque	40mm (1 9/16") ●	11841	10
		50mm (2") ●	11842	10
	Transparent	40mm (1 9/16") ●	11851	10
		50mm (2") ●	11852	10
Maxi 26cm (10 1/4")	Opaque	40mm (1 9/16") ●	11844	10
		50mm (2") ●	11845	10
		60mm (2 3/8") ●	11846	10
	Transparent	40mm (1 9/16") ●	11854	10
		50mm (2") ●	11855	10
		60mm (2 3/8") ●	11856	10



SenSura® High Output Pouch

Colour	Coupling Size	Volume	Product Code	Units
Opaque	40mm (1 9/16") ●	1160 mL	19040	10
	50mm (2") ●	1120 mL	19041	10
	60mm (2 3/8") ●	1080 mL	19042	10
Transparent	40mm (1 9/16") ●	1160 mL	19050	10
	50mm (2") ●	1120 mL	19051	10
	60mm (2 3/8") ●	1080 mL	19052	10



Will connect to
Coloplast Drainage
Bag 14010

SenSura® Flex

Ostomy Barrier and Pouch - Two Piece System

Low profile adhesive coupling system that provide flexibility and freedom of movement. Double layer adhesive for secure, safe, long lasting wear. Drainable pouch with velcro closure system. Only SenSura Flex Convex Flanges have 2 belt attachment points.

SenSura® Flex Skin Barrier

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	35mm (1 3/8") ●	10-33mm (3/8 - 1 1/4")	10101	10
		50mm (2") ●	10-48mm (3/8 - 2")	10102	10
		70mm (2 3/4") ●	10-68mm (3/8 - 2 3/4")	10103	10
		90mm - 3 1/2" ●	10-88mm (3/8 - 3 1/2")	10108	5
Pre-Cut	Convex light	35mm (1 3/8") ●	15-23mm (5/8 - 7/8")	11301	5
		50mm (2") ●	15-33mm (5/8 - 1 1/4")	11302	5
			15-43mm (5/8 - 1 3/4")	11303	5
		70mm (2 3/4") ●	15-56mm (5/8 - 2 1/8")	11307	5
	Non-Convex	35mm (1 3/8") ●	21mm (7/8")	10111	10
			28mm (1 1/8")	10113	10
			25mm (1")	10114	10
		50mm (2") ●	31mm (1 1/4")	10116	10
Pre-Cut	Convex light		35mm (1 3/8")	10117	10
		35mm (1 3/8") ●	21mm (7/8")	11312	5
			25mm (1")	11313	5
			28mm (1 1/8")	11314	5
			31mm (1 1/4")	11315	5
			35mm (1 3/8")	11316	5



SenSura® Flex Xpro Skin Barrier

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	35mm (1 3/8")	10-33mm (3/8 - 1 1/4")	10104	10
		50mm (2")	10-48mm (3/8 - 2")	10105	10
		70mm (2 3/4")	10-68mm (3/8 - 2 3/4")	10106	10
	Convex light	35mm (1 3/8")	15-23mm (5/8 - 7/8")	11304	5
		50mm (2")	15-33mm (5/8 - 1 1/4")	11305	5
		70mm (2 3/4")	15-43mm (5/8 - 1 3/4")	11306	5
		70mm (2 3/4")	15-56mm (5/8 - 2 1/8")	11308	5



SenSura® Flex Drainable Pouch with Easyclose™ WIDE-Outlet

Size	Colour	Coupling Size			Product Code	Units
Midi 27cm (10 1/2")	Opaque	50mm (2")			11505	20
Maxi 30cm (12")	Opaque	35mm (1 3/8")			11514	20
		50mm (2")			11515	20
		50mm (2")		No filter	11615	20
		70mm (2 3/4")			11518	20
	Transparent	50mm (2")		No filter	11612	20
		90mm (3 1/2")		No filter	11690	20



SenSura® Flex Closed Pouch

Size	Colour	Coupling Size	Product Code	Units
Midi 18cm (7")	Opaque	35mm (1 3/8")	10904	30
		50mm (2")	10905	30
Maxi 22cm (8 1/2")	Opaque	35mm (1 3/8")	10914	30
		50mm (2")	10915	30
		70mm (2 3/4")	10918	30



SenSura® Flex Urostomy Pouch

Size	Colour	Coupling Size	Product Code	Units
Midi 24cm (9 1/2")	Opaque	35mm (1 3/8")	11821	20
		50mm (2")	11822	20
Maxi 26cm (10 1/4")	Opaque	35mm (1 3/8")	11824	20
		50mm (2")	11825	20
	Transparent	35mm (1 3/8")	11827	20
		50mm (2")	11828	20



Ostomy Barrier and Pouch - One Piece System

Easiclose™ system with Velcro® to ensure security. Double layer adhesive provides secure adhesion, flexibility and ease of removal. Only SenSura 1pc Convex Systems have 2 belt attachment points.

SenSura® Drainable Pouch with Easiclose™ WIDE-Outlet						
Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 30cm (12")	Cut-to-Fit	Non-Convex	Opaque	10-76mm ($\frac{3}{8}$ -3")	15531	10
			Transparent	10-76mm ($\frac{3}{8}$ -3")	15521	10
				10-76mm ($\frac{3}{8}$ -3")	15821	10
		Convex Light	Opaque	15-33mm ($\frac{5}{8}$ -1 $\frac{1}{4}$ ")	15607	10
				15-43mm ($\frac{5}{8}$ -1 $\frac{3}{4}$ ")	15608	10
			Transparent	15-33mm ($\frac{5}{8}$ -1 $\frac{1}{4}$ ")	15605	10
				15-43mm ($\frac{5}{8}$ -1 $\frac{3}{4}$ ")	15606	10
	Pre-Cut	Non-Convex	Opaque	25mm (1")	15532	10
				30mm (1 $\frac{1}{8}$ ")	15533	10
				35mm (1 $\frac{3}{8}$ ")	15534	10
		Convex Light	Opaque	25mm (1")	15633	10
				28mm (1 $\frac{1}{8}$ ")	15634	10
				31mm (1 $\frac{1}{4}$ ")	15635	10
				25mm (1")	15522	10
		Non-Convex	Transparent	30mm (1 $\frac{1}{8}$ ")	15523	10
				35mm (1 $\frac{3}{8}$ ")	15524	10



SenSura® Xpro Drainable Pouch with Easiclose™ WIDE-Outlet

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 30cm (12")	Cut-to-Fit	Non-Convex	Transparent	10-76mm ($\frac{3}{8}$ - 3")	15981	20
		Convex Light	Transparent	15-33mm ($\frac{5}{8}$ - 1 $\frac{1}{4}$ ")	15695	10
				15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	15696	10

SenSura® Closed Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 18cm (7")	Cut-to-Fit	Non-Convex	Opaque	10-66mm ($\frac{3}{8}$ - 2 $\frac{5}{8}$ ")	15450	30
			Transparent	10-66mm ($\frac{3}{8}$ - 2 $\frac{5}{8}$ ")	15441	30
		Convex Light	Opaque	15-33mm ($\frac{5}{8}$ - 1 $\frac{1}{4}$ ")	15403	10
				15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	15404	10
	Pre-Cut	Convex Light	Opaque	25mm (1")	15413	10
				28mm (1 $\frac{1}{8}$ ")	15414	10
				31mm (1 $\frac{1}{4}$ ")	15415	10
	Cut-to-Fit	Non-Convex	Opaque	10-76mm ($\frac{3}{8}$ - 3")	15480	30
				10-76mm ($\frac{3}{8}$ - 3")	15470	30
		Convex Light	Opaque	15-33mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	15407	10
				15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	15408	10
			Transparent	15-33mm ($\frac{5}{8}$ - 1 $\frac{1}{4}$ ")	15405	10
				15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	15406	10
Maxi 22cm (8 $\frac{1}{2}$ ")	Pre-Cut	Non-Convex	Opaque	25mm (1")	15482	30
				30mm (1 $\frac{1}{8}$ ")	15483	30
				35mm (1 $\frac{3}{8}$ ")	15484	30
				25mm (1")	15472	30
			Transparent	30mm (1 $\frac{1}{8}$ ")	15473	30
				35mm (1 $\frac{3}{8}$ ")	15474	30
				25mm (1")	15472	30
				30mm (1 $\frac{1}{8}$ ")	15473	30



SenSura® Urostomy Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 24cm (9 $\frac{1}{2}$ ")	Cut-to-Fit	Non-Convex	Opaque	10-66mm ($\frac{3}{8}$ - 2 $\frac{5}{8}$ ")	11801	10
			Transparent	10-66mm ($\frac{3}{8}$ - 2 $\frac{5}{8}$ ")	11803	10
		Non-Convex	Opaque	10-76mm ($\frac{3}{8}$ - 3")	11802	10
			Transparent	10-76mm ($\frac{3}{8}$ - 3")	11804	10
Maxi 26cm (10 $\frac{1}{4}$ ")	Cut-to-Fit	Convex Light	Opaque	15-33mm ($\frac{5}{8}$ - 1 $\frac{1}{4}$ ")	11811	10
				15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	11812	10
			Transparent	15-33mm ($\frac{5}{8}$ - 1 $\frac{1}{4}$ ")	11814	10
				15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	11815	10



SenSura® One-Piece High Output Pouch

Cut	Sterile	Stoma Size	Volume	Product Code	Units
Cut-to-Fit	Sterile	10-76mm ($\frac{3}{8}$ -3")	630 mL	19001	5
	Non-Sterile	10-76mm ($\frac{3}{8}$ -3")	630 mL	19003	6



SenSura® Post-Op Pouch


Cut	Sterile	Stoma Size	Window Size	Volume	Product Code	Units
Cut-to-Fit	Sterile	10-76mm ($\frac{3}{8}$ -3")	70mm	680 mL	19010	5
		10-115mm ($\frac{3}{8}$ - 4 $\frac{1}{2}$ ")	100mm	725 mL	19011	5
	Non-Sterile	10-76mm ($\frac{3}{8}$ -3")	70mm	680 mL	19020	6
		10-115mm ($\frac{3}{8}$ - 4 $\frac{1}{2}$ ")	100mm	725 mL	19021	6
		10-76mm ($\frac{3}{8}$ -3")	No window	680 mL	19030	6



Ostomy Barrier and Pouch - Two Piece System


Click mechanical coupling with a lock ring for security. Skin barrier with spiral technology for secure, long lasting wear. All 2pc Assura Click Flanges have 2 belt attachment points.

Assura® Skin Barrier					
Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Pre-Cut		50mm (2")	● 10-45mm ($\frac{3}{8}$ - $1\frac{3}{4}$ ")	2882	5
		60mm ($2\frac{3}{8}$ ")	● 10-55mm ($\frac{3}{8}$ - $2\frac{1}{8}$ ")	2883	5
	Convex Light	40mm ($1\frac{9}{16}$ ")	● 15-23mm ($\frac{5}{8}$ - $\frac{7}{8}$ ")	14261	5
		50mm (2")	● 15-33mm ($\frac{5}{8}$ - $1\frac{1}{4}$ ")	14262	5
		60mm ($2\frac{3}{8}$ ")	● 15-43mm ($\frac{5}{8}$ - $1\frac{3}{4}$ ")	14263	5
	Convex Deep	50mm (2")	● 15-33mm ($\frac{5}{8}$ - $1\frac{1}{4}$ ")	12716	5
		60mm ($2\frac{3}{8}$ ")	● 15-43mm ($\frac{5}{8}$ - $1\frac{3}{4}$ ")	12719	5
	Non-Convex	40mm ($1\frac{9}{16}$ ")	● 20mm ($\frac{3}{4}$ ")	12841	5
			● 25mm (1")	12842	5
		50mm (2")	● 30mm ($1\frac{1}{8}$ ")	12844	5
			● 35mm ($1\frac{3}{8}$ ")	12845	5
		60mm ($2\frac{3}{8}$ ")	● 40mm ($1\frac{5}{8}$ ")	12847	5
		60mm ($2\frac{3}{8}$ ")	● 45mm ($1\frac{3}{4}$ ")	12848	5
		40mm ($1\frac{9}{16}$ ")	● 18mm ($\frac{3}{4}$ ")	14271	5
			● 21mm ($\frac{7}{8}$ ")	14272	5
		50mm (2")	● 25mm (1")	14273	5
			● 28mm ($1\frac{1}{8}$ ")	14274	5
			● 31mm ($1\frac{1}{4}$ ")	14275	5
			● 35mm ($1\frac{3}{8}$ ")	14276	5
	Convex Light	60mm ($2\frac{3}{8}$ ")	● 38mm ($1\frac{1}{2}$ ")	14277	5
			● 41mm ($1\frac{5}{8}$ ")	14278	5
		40mm ($1\frac{9}{16}$ ")	● 15mm ($\frac{5}{8}$ ")	12701	5
			● 18mm ($\frac{3}{4}$ ")	12702	5
			● 21mm ($\frac{7}{8}$ ")	12703	5
		50mm (2")	● 25mm (1")	12704	5
			● 28mm ($1\frac{1}{8}$ ")	12705	5
			● 31mm ($1\frac{1}{4}$ ")	12706	5
			● 35mm ($1\frac{3}{8}$ ")	12707	5
		60mm ($2\frac{3}{8}$ ")	● 38mm ($1\frac{1}{2}$ ")	12708	5
			● 41mm ($1\frac{5}{8}$ ")	12709	5



Assura® Skin Barrier - Extra Extended Wear

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	40mm (1 9/16")	10-35mm (3/8 - 1 3/8")	2831	5
		50mm (2")	10-45mm (3/8 - 1 3/4")	2832	5
		60mm (2 3/8")	10-55mm (3/8 - 2 1/8")	2833	5
	Convex Light	40mm (1 9/16")	15-23mm (5/8 - 7/8")	14281	5
		50mm (2")	15-33mm (5/8 - 1 1/4")	14282	5
		60mm (2 3/8")	15-43mm (5/8 - 1 3/4")	14283	5
	Convex Deep	40mm - (1 9/16")	15-23mm (5/8 - 7/8")	14243	5
		50mm (2")	15-33mm (5/8 - 1 1/4")	14246	5
		60mm (2 3/8")	15-43mm (5/8 - 1 3/4")	14249	5
Pre-Cut	Non-Convex	40mm (1 9/16")	20mm (3/4")	14251	5
			25mm (1")	14252	5
		50mm (2")	30mm (1 1/8")	14254	5
			35mm (1 3/8")	14255	5
		60mm (2 3/8")	45mm - 1 3/4"	14258	5
	Convex Light	40mm (1 9/16")	18mm (3/4")	14291	5
			21mm (7/8")	14292	5
		50mm (2")	25mm (1")	14293	5
			28mm (1 1/8")	14294	5
		60mm (2 3/8")	31mm (1 1/4")	14295	5
			35mm (1 3/8")	14296	5
	Convex Deep	40mm (1 9/16")	18mm (3/4")	14232	5
			21mm (7/8")	14233	5
		50mm (2")	25mm (1")	14234	5
			28mm (1 1/8")	14235	5
		60mm (2 3/8")	31mm (1 1/4")	14236	5
			35mm (1 3/8")	14237	5
		60mm (2 3/8")	38mm (1 1/2")	14238	5
			41mm (1 5/8")	14239	5



Assura® 2 pc. High Output Drainable Pouch

Size	Colour	Coupling Size	Product Code	Units
Maxi 32cm (12 1/2")	Opaque	50mm (2")	2846	10
		60mm (2 3/8")	2847	10



Assura® 2 pc. High Output Drainable System, Ileo Night Pouch 2L

Size	Colour	Coupling Size	Product Code	Units
2 L	Opaque	50mm (2")	2836	5
		60mm (2 3/8")	2837	5



Assura® Drainable Pouch Clamp Closure

Size	Colour	Coupling Size	Product Code	Units
Maxi 30cm (12 1/2")	Transparent	50mm (2") ●	2792	10
		60mm (2 3/8") ●	2793	10

Assura® Drainable Pouch with EasiClose™

Size	Colour	Coupling Size	Product Code	Units
Mini 23cm (9")	Opaque	40mm (1 9/16") ●	13924	10
		50mm (2") ●	13925	10
Midi 26cm (10 1/4")	Opaque with Graphic	40mm (1 9/16") ●	13964	10
		50mm (2") ●	13965	10
		60mm (2 3/8") ●	13966	10
Maxi 28cm (11")	Opaque	40mm (1 9/16") ●	13984	10
		40mm (1 9/16") ● No filter	14494	10
		50mm (2") ●	13985	10
		60mm (2 3/8") ●	13986	10
		60mm (2 3/8") ● No filter	14498	10
		40mm (1 9/16") ●	13974	10
	Transparent	50mm (2") ●	13975	10
		60mm (2 3/8") ●	13976	10



Assura® Closed Pouch

Size	Colour	Coupling Size	Product Code	Units
Midi 18cm (7")	Opaque	40mm (1 9/16") ●	12354	30
		50mm (2") ●	12355	30
		60mm (2 3/8") ●	12356	30
	Transparent	40mm (1 9/16") ●	12344	30
		60mm (2 3/8") ●	12346	30
		60mm (2 3/8") ●	12346	30
Maxi 22cm (8 1/2")	Opaque	40mm (1 9/16") ●	12384	30
		50mm (2") ●	12385	30
		60mm (2 3/8") ●	12386	30
	Transparent	40mm (1 9/16") ●	12374	30
		50mm (2") ●	12375	30
		60mm (2 3/8") ●	12376	30



Assura® Urostomy Pouch

Size	Colour	Coupling Size	Product Code	Units
Maxi 25cm (10")	Opaque	40mm (1 9/16") ●	14224	10
		50mm (2") ●	14225	10
		60mm (2 3/8") ●	14226	10
	Transparent	40mm (1 9/16") ●	14227	10
		50mm (2") ●	14228	10
		60mm (2 3/8") ●	14229	10



Please note: The adult Easiflex® line of ostomy appliances will be discontinued in 2022 and will be replaced with the SenSura® Mio Flex line (page 12). Please call Coloplast Care 1-866-293-6349 for any questions.

Ostomy Barrier and Pouch- Two Piece System

Soft and flexible adhesive coupling allows for extra comfort. Easy to apply the pouch to the skin barrier and easy to remove. All Easiflex Convex Flanges have 2 belt attachment points.

Easiflex® Skin Barrier					
Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	40mm (1 9/16")	10-35mm (3/8 - 1 3/8")	14301	10
		50mm (2")	10-50mm (3/8 - 2")	14302	10
		70mm (2 3/4")	10-68mm (3/8 - 2 3/8")	14303	10
		90mm - (3 1/2")	15-88mm (5/8 - 3 1/2")	14309	5
	Convex Light	40mm (1 9/16")	15-23mm (5/8 - 7/8")	14401	5
		50mm (2")	15-33mm (5/8 - 1 1/4")	14402	5
		50mm (2")	15-43mm (5/8 - 1 3/4")	14403	5
		70mm (2 3/4")	15-61mm (5/8 - 2 1/8")	14406	5
	Pre-Cut	40mm (1 9/16")	21mm (7/8")	14642	5
		50mm (2")	25mm (1")	14643	5
		50mm (2")	28mm (1 1/8")	14644	5



Easiflex® Skin Barrier - Extra Extended Wear					
Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	40mm (1 9/16")	10-35mm (3/8 - 1 3/8")	14304	10
		50mm (2")	10-48mm (3/8 - 2")	14305	10
		70mm (2 3/4")	10-68mm (3/8 - 2 3/8")	14306	10
	Convex Light	40mm (1 9/16")	15-23mm (5/8 - 7/8")	14601	5
		50mm (2")	15-33mm (5/8 - 1 1/4")	14602	5
		50mm (2")	15-43mm (5/8 - 1 3/4")	14603	5
		70mm (2 3/4")	15-61mm (5/8 - 2 1/8")	14408	5

Easiflex® Drainable Pouch with EasiClose™					
Size	Colour	Coupling Size		Product Code	Units
Midi 26cm (10 1/4")	Opaque	40mm (1 9/16")		14346	20
		50mm (2")		14347	20
	Transparent	50mm (2")		14342	20
		40mm (1 9/16")		14361	20
Maxi 28cm (11")	Opaque	40mm (1 9/16")	No filter	14524	20
		50mm (2")		14362	20
		50mm (2")	No filter	14525	20
		70mm (2 3/4")		14349	20
		40mm (1 9/16")		14356	20
		40mm (1 9/16")	No filter	14520	20
	Transparent	50mm (2")		14357	20
		50mm (2")	No filter	14521	20
		70mm (2 3/4")		14344	20
		70mm (2 3/4")	No filter	14528	20
		90mm - (3 1/2")	No filter	14527	20
		70mm (2 3/4")		14364	20
XL 30cm (12")	Transparent	70mm (2 3/4")		14359	20



Easiflex® Closed Pouch

Size	Colour	Coupling Size	Product Code	Units
Midi 18cm (7")	Opaque	40mm (1 9/16")	14316	30
		50mm (2")	14317	30
Maxi 21cm (8 1/4")	Opaque	40mm (1 9/16")	14331	30
		50mm (2")	14332	30
		70mm (2 3/4")	14319	30
		90mm - (3 1/2")	14355	30
XL 26cm (10 1/4")	Transparent	70mm (2 3/4")	14329	30



Easiflex® Urostomy Pouch

Size	Colour	Coupling Size	Product Code	Units
Maxi 25cm (10")	Opaque	40mm (1 9/16")	14553	20
		50mm (2")	14555	20
	Transparent	40mm (1 9/16")	14552	20
		50mm (2")	14554	20



Assura®

Ostomy Barrier and Pouch - One Piece System

Skin Barrier with spiral technology for secure, long lasting wear. Different convexity and pouch shape. All Assura 1pc Convex Systems have 2 belt attachment points.

Assura® Drainable Pouch Clamp Closure

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 30cm (12")	Cut-to-Fit	Non-Convex	Opaque	10-70mm (3/8 - 2 3/4")	2520 No filter	10
			Transparent	10-70mm (3/8 - 2 3/4")	2490 No filter	10
		Convex Deep	Opaque	15-43mm (5/8 - 1 3/4")	12533 No filter	10
			Transparent	15-43mm (5/8 - 1 3/4")	13706 No filter	10



Assura® Drainable Pouch with EasiClose™

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 26cm (10 1/4")	Cut-to-Fit	Non-Convex	Opaque	10-55mm (3/8 - 2 1/8")	13840	10
		Convex Light	Opaque	15-33mm (5/8 - 1 1/4")	14404	10
				15-43mm (5/8 - 1 3/4")	14405	10
		Convex Deep	Opaque	15-43mm (5/8 - 1 3/4")	14196	10
	Pre-Cut	Non-Convex	Opaque	25mm (1")	13844	10
				30mm (1 1/8")	13845	10
Maxi 28cm (11")	Cut-to-Fit	Non-Convex	Opaque	10-70mm (3/8 - 2 3/4")	13870	10
			Transparent	10-70mm (3/8 - 2 3/4")	13860	10
		Convex Light	Opaque	15-33mm (5/8 - 1 1/4")	14414	10
				15-43mm (5/8 - 1 3/4")	14415	10
			Transparent	15-33mm (5/8 - 1 1/4")	14411	10
				15-43mm (5/8 - 1 3/4")	14412	10



Assura® Drainable Pouch with EasiClose™

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 28cm (11")	Cut-to-Fit	Convex Deep	Opaque	15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	14106	10
			Transparent	15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	14103	10
	Pre-Cut	Non-Convex	Opaque	25mm (1")	13874	10
				30mm (1 $\frac{1}{8}$ ")	13875	10
				35mm (1 $\frac{3}{8}$ ")	13876	10
				25mm (1")	14422	10
		Convex Light	Opaque	28mm (1 $\frac{1}{8}$ ")	14423	10
				31mm (1 $\frac{1}{4}$ ")	14424	10
				21mm ($\frac{7}{8}$ ")	14173	10
				25mm (1")	14174	10
		Convex Deep	Opaque	28mm (1 $\frac{1}{8}$ ")	14175	10
				31mm (1 $\frac{1}{4}$ ")	14176	10
				21mm ($\frac{7}{8}$ ")	14163	10
			Transparent	25mm (1")	14164	10
				28mm (1 $\frac{1}{8}$ ")	14165	10
				31mm (1 $\frac{1}{4}$ ")	14166	10



Assura® Closed Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Mini 12cm (4 $\frac{3}{4}$ ")	Cut-to-Fit	Non-Convex	Opaque	20-55mm ($\frac{3}{8}$ - 2 $\frac{1}{8}$ ")	12110	30
			Transparent	20-55mm ($\frac{3}{8}$ - 2 $\frac{1}{8}$ ")	12130	30
Midi 18cm (7")	Pre-Cut	Non-Convex	Opaque	25mm (1")	12144	30
				30mm (1 $\frac{1}{8}$ ")	12145	30
				35mm (1 $\frac{3}{8}$ ")	12146	30
				40mm (1 $\frac{9}{16}$ ")	12147	30




Assura® 1-Piece Post-Op Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Mini 31cm (12 $\frac{1}{4}$ ")	Cut-to-Fit	Non-Convex	Transparent	13-70mm ($\frac{1}{2}$ - 2 $\frac{3}{4}$ ")	12805	5
				13-70mm ($\frac{1}{2}$ - 2 $\frac{3}{4}$ ")	12800	5
Midi 33cm (13")	Cut-to-Fit	Non-Convex	Transparent	13-100mm ($\frac{1}{2}$ - 4")	12802	5
				13-70mm ($\frac{1}{2}$ - 2 $\frac{3}{4}$ ")	12808 Sterile	5
				13-100mm ($\frac{1}{2}$ - 4")	12810 Sterile	5



Assura® Closed Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 18cm (7")	Pre-Cut	Convex Light	Opaque	25mm (1")	14452	10
				28mm (1 1/8")	14453	10
				31mm (1 1/4")	14454	10
Maxi 22cm (8 3/4")	Cut-to-Fit	Non-Convex	Opaque	20-70mm (3/4 - 2 3/4")	12170	30
		Convex Light	Opaque	15-43mm (5/8 - 1 3/4")	14445	10
			Transparent	15-33mm (5/8 - 1 1/4")	14441	10
	Pre-Cut	Non-Convex	Opaque	25mm (1")	12174	30
				30mm (1 1/8")	12175	30
				35mm (1 3/8")	12176	30
				40mm (1 9/16")	12177	30



Assura® Urostomy Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 27cm (10 1/2")	Cut-to-Fit	Non-Convex	Opaque	10-55mm (3/8 - 2 1/8")	5580	10
			Transparent	10-55mm (3/8 - 2 1/8")	5585	10
		Convex Deep	Transparent	15-33mm (5/8 - 1 1/4")	12595	10
				15-43mm (5/8 - 1 3/4")	12596	10
	Pre-Cut	Convex Deep	Transparent	15mm (5/8")	12991	10
				18mm (3/4")	12992	10
				21mm (7/8")	12993	10
				25mm (1")	12994	10
				28mm (1 1/8")	12995	10
				31mm (1 1/4")	12996	10
				35mm (1 3/8")	12997	10



Assura® Urostomy Pouch - Extra Extended Wear

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 25cm (10")	Cut-to-Fit	Non-Convex	Opaque	10-55mm (3/8 - 2 1/8")	14201	10
			Transparent	10-55mm (3/8 - 2 1/8")	14202	10
		Convex Light	Opaque	15-43mm (5/8 - 1 3/4")	14701	10
			Transparent	15-43mm (5/8 - 1 3/4")	14702	10
Maxi 27cm (10 1/2")	Cut-to-Fit	Non-Convex	Opaque	10-55mm (3/8 - 2 1/8")	12471	10
			Transparent	10-55mm (3/8 - 2 1/8")	12472	10
		Convex Deep	Transparent	15-33mm (5/8 - 1 1/4")	12432	10
				15-43mm (5/8 - 1 3/4")	12433	10



Assura®

Ostomy Barrier and Pouch - Two Piece System - Pediatric

Click mechanical coupling.

Assura® Pediatric Skin Barrier

Cut	Convexity	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	10-35mm (3/8 - 1 3/8")	2181	5



Assura® Pediatric Drainable Pouch

Size	Colour	Product Code	Units
13cm (5")	Opaque	2151	10

Assura® Pediatric Closed Pouch

Size	Colour	Product Code	Units
13cm (5")	Transparent	2166	30



EasiFlex®

Ostomy Barrier and Pouch - Two Piece System - Pediatric

Flex adhesive coupling.

Easiflex® Pediatric Skin Barrier

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	20mm ($\frac{3}{4}$ "	0-15mm ($0-\frac{5}{8}$ "	14307	10
		30mm ($1\frac{1}{8}$ "	0-25mm (0-1"	14308	10



Easiflex® Pediatric Drainable Pouch

Size	Colour	Coupling Size	Product Code	Units
105 mL	Teddy Print	20mm ($\frac{3}{4}$ "	14691	10
	Transparent	20mm ($\frac{3}{4}$ "	14681	10
150 mL	Teddy Print	30mm ($1\frac{1}{8}$ "	14692	10
	Transparent	30mm ($1\frac{1}{8}$ "	14682	10



Assura®

Ostomy Barrier and Pouch - One Piece System - Pediatric

Assura® Pediatric Drainable Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
15cm (6")	Cut-to-Fit	Non-Convex	Transparent	10-35mm ($\frac{3}{8}$ - $1\frac{3}{8}$ "	2115	10

Assura® Pediatric Closed Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
13cm (5")	Cut-to-Fit	Non-Convex	Transparent	10-35mm ($\frac{3}{8}$ - $1\frac{3}{8}$ "	2125	30

Assura® Pediatric Urostomy Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
15cm (6")	Cut-to-Fit	Non-Convex	Transparent	10-35mm ($\frac{3}{8}$ - $1\frac{3}{8}$ "	8009	10



Ostomy Accessories

Brava® Protective Sheets

Size (cms)	Size (inches)	Product Code	Units
10 x 10 cm	4 x 4"	3210	10
15 x 15 cm	6 x 6"	3215	5
20 x 20 cm	8 x 8"	3220	5



Brava® Protective Rings					
Size (mm) Internal diameter	Size (mm) External diameter	Size (mm) thickness	Product Code	Units	
18 mm	48 mm	2.5 mm	12035	10	
18 mm	48 mm	4.2 mm	12045	10	
27 mm	57 mm	2.5 mm	12037	10	
27 mm	57 mm	4.2 mm	12047	10	
34 mm	64 mm	2.5 mm	12039	10	
34 mm	64 mm	4.2 mm	12049	10	
Brava® Protective Rings Wide					
Size (mm) Internal diameter	Size (mm) External diameter	Size (mm) thickness	Product Code	Units	
18 mm	57 mm	2.5 mm	12036	10	
18 mm	64 mm	2.5 mm	12038	10	
18 mm	57 mm	4.2 mm	12046	10	
18 mm	64 mm	4.2 mm	12048	10	
Brava® Protective Rings Extra Wide					
Size (mm) Internal diameter	Size (mm) External diameter	Size (mm) thickness	Product Code	Units	
18 mm	76 mm	2.5 mm	12032	10	
27 mm	76 mm	2.5 mm	12033	10	
Brava® Protective Ring Convex					
Size (mm) Internal diameter	Size (mm) External diameter	Size (mm) thickness	Product Code	Units	
20 mm	60 mm	8.0 mm	12090	10	
25 mm	65 mm	8.0 mm	12091	10	
30 mm	70 mm	8.0 mm	12092	10	
35 mm	75 mm	8.0 mm	12093	10	
40 mm	80 mm	8.0 mm	12094	10	
45 mm	85 mm	8.0 mm	12095	10	
Brava® Mouldable Rings					
Size (mm) Internal diameter	Size (mm) thickness		Product Code	Units	
18 mm	2.0 mm		12030	10	
18 mm	4.2 mm		12042	10	
Brava® Skin Barrier Spray and Wipe					
Size	Product Code		Units		
50 mL Spray	12020		1		
Single unit wipe	12021		30		
Brava® Adhesive Remover Spray and Wipe					
Size	Product Code		Units		
50 mL Spray	12010		1		
Single unit wipe	12011		30		
Brava® Barrier Cream					
Size	Product Code		Units		
60g Tube	12000		1		
Brava® Ostomy Strip Paste					
Size	Product Code		Units		
57g / 2oz	2655		10		

Brava® Elastic Barrier Strips

Size	Product Code	Units
Brava Elastic Barrier Strip Original (a.)	12070	20
Brava Elastic Barrier Strip Straight (b.)	12074	20
Brava Elastic Barrier Strip Belt (c.)	12072	20
Brava Elastic Barrier Strip XL (d.)	12076	20



Brava® Ostomy Powder

Size	Product Code	Units
25g / 1oz bottle	1907	16



Brava® Mio Belt for All SenSura Mio Products

New design for easier attachment and improved security

Size	Product Code	Units
Standard - 100 cm - 40"	0423	15
XXL - 155cm - 61"	0424	15



Brava® Belt for SenSura Original + Assura Products

Size	Product Code	Units
Standard - 100 cm - 40"	0421	15
XL - 135cm - 54"	0422	15



Brava® Lubricating Deodorant

Size	Product Code	Units
7.5 mL sachet	12060	20
240 mL bottle	12061	1



Brava® Alcohol Free Paste

Size	Product Code	Units
60 g tube	12050	1



Ostomy Paste (contains alcohol)

Size	Product Code	Units
60 g tube	2650	12



Skin Barrier Rings

Pre-cut Size	Product Code	Units
20mm ($\frac{3}{4}$ ")	2320	30
25mm (1")	2325	30
30mm (1 $\frac{1}{8}$ ")	2330	30
40mm (1 $\frac{9}{16}$ ")	2340	30
50mm (2")	2350	30



Coloplast Prep™ Protective Liquid Skin Barrier (contains alcohol)

Size	Product Code	Units
Single Application Wipe	2041	12 boxes/54
59 mL - 2oz. dab bottle	0925	12



Coloplast Filtrodor® Pouch Filter

	Product Code	Units
Charcoal-activated deodorizing filter with self-adhering foam center, adheres to front of Assura & SenSura Pouches	0509	5 pouches (10 units per pouch)



Coloplast Pouch Clamp

Product Code Units

Simple closure clamp for Assura Drainable Pouches

9500 20



Ostomy Scissors

Product Code Units

Stainless steel medical grade ostomy scissors with curved blades and rounded ends - right handed

9505 30



Hex-On®

Odour antagonist

Eliminates odour with a light, mild scent.

Hex-On®

Size (ounces) Size (mL) Product Code Units

2 oz. 59 mL 7583 12



Ostomy Cap

Designed for discreet situations. 2 piece compatible with SenSura Click and Assura 2 piece skin barriers.

**Note Uro Minicap and Coloplast Stoma cap do not fit SenSura Mio products*

Coloplast Uro Minicap (100 ml)

Size	Colour	Coupling Size	Product Code	Units
Mini	Opaque	40mm (1 9/16") 	2807	30
		50mm (2") 	2808	30




Coloplast Stoma Cap

Size	Colour	Coupling Size	Product Code	Units
Mini	Opaque	40mm (1 9/16") 	2801	30
		50mm (2") 	2802	30
		60mm (2 3/8") 	2803	30



Coloplast Stoma Cap - One Piece

Size	Colour	Cut to Fit - Stoma size	Product Code	Units
Mini	Opaque	20-55mm (3/4 - 2 1/8")	2501	30



Assura® Irrigation Products

Assura® Irrigation Set

Description	Product Code	Units
Includes: - 1 Water Bag/Regulator - 1 Cone	1500	1



Assura® Irrigation Accessories

Size	Product Code	Units
Stoma Cone - Soft, flexible stoma cone with tubing	1110	4



Assura® Irrigation Sleeves

Size	Product Code	Units
Transparent, disposable irrigation sleeve with adhesive back. - 1-pc sleeve	1003	10



The irrigation sleeves below must be used with their corresponding colour coded flange or faceplate. Consult this chart to order the correct transparent sleeves.

Description	Coupling Size	Use with Assura and SenSura Barrier Product Codes (will not work with SenSura Mio)	Product Code	Units
Irrigation sleeves	40mm (1 9/16") ●	Will couple to any green 40 mm barriers within the Assura and SenSura lines	12834	5
	50mm (2") ●	Will couple to any red 50 mm barriers within the Assura and SenSura lines	12835	5
	60mm (2 3/8") ●	Will couple to any blue 60 mm barriers within the Assura and SenSura lines	12836	5



Coloplast Urostomy Night Bag

Description	Product Code	Units
Includes a fixed connector that fits on to the urostomy pouch. The bag holds up to 2 litres and has a 120cm long flexible tube with anti-reflux valve. Not Compatible with External Catheters or Leg Bags	21365	10



Coloplast Fistula & Wound Management System

Pouch				
Size	Cutting Surface	Product Code	Units	
Mini	104-159mm (4 1/8 - 6 1/4")	14050	6	
Midi	156-228mm (6 3/16 - 9")	14060	3	
Maxi	208-297mm (8 1/4 - 11 3/4")	14070		
Flexible Lids				
Size	For	Product Code	Units	
Mini	Mini Fistula & Wound Management System	14021	10	
Midi	Midi Fistula & Wound Management System	14031		
Maxi	Maxi Fistula & Wound Management System	14041		
Accessories				
Size	For	Product Code	Units	
Bed Drainage Bag	All Fistula & Wound Management systems (Capacity 2000 mL)	14010	10	
Drain Port	All Fistula & Wound Management systems	14015	12	

Continence

SpeediCath® Flex

Hydrophilic Intermittent Catheter

A soft catheter for men with a dry-sleeve and a flexible tip with all the benefits of the SpeediCath® catheter. Offered with a loop system that fits into a discreet packaging which is 100% aluminum free. Latex, phthalate and PVC free. Sterile.

SpeediCath® Flex				
Length (cm)	Length (inches)	Size	Product Code	Units
44 cm	17"	10 FR ●	28920	30
		12 FR ○	28922	30
		14 FR ●	28924	30
		16 FR ●	28916	30



SpeediCath® Standard

Hydrophilic Intermittent Catheter

Instantly ready to use catheter. Featuring a smooth and uniform hydrophilic coating and polished eyelets to ensure a convenient and simple catheterization. Latex, phthalate and PVC free. Sterile.

SpeediCath® Standard Male Straight				
Length (cm)	Length (inches)	Size	Product Code	Units
40 cm	16"	8 FR ●	27408	30
		10 FR ●	27410	30
		12 FR ○	27412	30
		14 FR ●	27414	30
		16 FR ●	27416	30
		18 FR ●	27418	30



SpeediCath® Standard Male Coude				
Length (cm)	Length (inches)	Size	Product Code	Units
40 cm	16"	10 FR ●	27490	30
		12 FR ○	27492	30
		14 FR ●	27494	30
		16 FR ●	27496	30



SpeediCath® Standard Female Straight				
Length (cm)	Length (inches)	Size	Product Code	Units
20 cm	8"	8 FR ●	27508	30
		10 FR ●	27510	30
		12 FR ○	27512	30
		14 FR ●	27514	30
		16 FR ●	27516	30



SpeediCath® Standard Pediatric Straight				
Length (cm)	Length (inches)	Size	Product Code	Units
20 cm	8"	6 FR ●	27706	30
		8 FR ●	27708	30
		10 FR ●	27710	30



SpeediCath® Compact Male

Hydrophilic Intermittent Catheter

Developed to offer discretion and practicality with all the benefits of the SpeediCath® catheter and more. Telescopic design makes the packaging half the size of a standard male catheter. Latex, phthalate and PVC free. Sterile.

SpeediCath® Compact Male					
Length (cm)	Length (inches)	Size	Product Code	Units	
33 cm	13"	12/18 FR (one size)	28692	30	



SpeediCath® Compact Eve

Hydrophilic Intermittent Catheter

Stylish and discreet compact catheter for women. Designed for modern life, delivering a one-step easy opening and secure reclosing. Triangular shape ensures easy handling and grip. Latex, phthalate and PVC free. Sterile.

SpeediCath® Compact Eve					
Length (cm)	Length (inches)	Size	Product Code	Units	
8.3 cm	3.3"	10 FR ●	28110	30	
		12 FR ○	28112	30	
		14 FR ●	28114	30	



SpeediCath® Compact Female

Hydrophilic Intermittent Catheter

Compact and discreet solution that is instantly ready to use. The catheter can easily be kept in a pocket or handbag. Latex, phthalate and PVC free. Sterile.

SpeediCath® Compact Female					
Length (cm)	Length (inches)	Size	Product Code	Units	
7 cm	2.75"	8 FR ●	28578	30	
		10 FR ●	28580	30	
		12 FR ○	28582	30	
		14 FR ●	28584	30	



SpeediCath® Compact Female Plus					
Length (cm)	Length (inches)	Size	Product Code	Units	
9 cm	3.5"	10 FR ●	28810	30	
		12 FR ○	28812	30	
		14 FR ●	28814	30	



SpeediCath® Compact Set

Hydrophilic Intermittent Catheter

Compact and discreet catheter and bag solution that is instantly ready to use. Can easily be kept in a pocket or handbag. Latex, phthalate and PVC free. Sterile.

SpeediCath® Compact Set Male

Length (cm)	Length (inches)	Size		Product Code	Units
33 cm	13"	12/18 FR	(one size)	28422	20



SpeediCath® Compact Set Female

Length (cm)	Length (inches)	Size		Product Code	Units
9 cm	3.5"	10 FR	●	28520	20
		12 FR	○	28522	20
		14 FR	●	28524	20



Self-Cath®

Intermittent Catheter

Siliconized surface for smooth insertion. Made of flexible medical grade PVC with smooth fire polished eyelets for patient comfort and ease of use. Latex free. Sterile.

Self-Cath® Female Straight

Length (cm)	Length (inches)	Size		Reference Code	Product Code	Units
15 cm	6"	8 FR	●	504600	208	30
		10 FR	●	504610	210	30
		12 FR	○	504620	212	30
		14 FR	●	504630	214	30
		14 FR	●	501960	without connector 240	30



Self-Cath® Pediatric

Length (cm)	Length (inches)	Size		Reference Code	Product Code	Units
25 cm	10"	5 FR	○	504400	305	30
		6 FR	●	504410	306	30
		8 FR	●	504420	without connector 308	30
		10 FR	●	504430	without connector 310	30



Self-Cath® Extension Tube

Length (cm)	Length (inches)	Size		Reference Code	Product Code	Units
61 cm	24"	n/a		507850	475	12



Self-Cath® Male Straight

Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
40 cm	16"	8 FR	● 504500	408	30
		10 FR	● 504510	410	30
		12 FR	○ 504520	412	30
		14 FR	● 504530	414	30
		16 FR	● 504540	416	30
		18 FR	● 504550	418	30



Self-Cath® Male Soft Straight

Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
40 cm	16"	12 FR	○ 501810	112	30
		14 FR	● 501820	114	30
		16 FR	● 501830	116	30

Self-Cath® Male Straight Curved Packaging

Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
40 cm	16"	12 FR	○ 502200	460	50
		14 FR	● 502210	450	50

Self-Cath® Male Coude Tapered with Guide Stripe

Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
40 cm	16"	8 FR	● 504650	608	30
		10 FR	● 504660	610	30
		12 FR	○ 504670	612	30
		14 FR	● 504680	614	30



Self-Cath® Male Coude Olive with Guide Stripe

Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
40 cm	16"	8 FR	● 504710	808	30
		10 FR	● 504720	810	30
		12 FR	○ 504730	812	30
		14 FR	● 504740	814	30
		16 FR	● 504750	816	30
		18 FR	● 504760	818	30



Self-Cath® Male Closed System Straight

Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
40 cm	16"	14 FR	● 502730	1114	50
		16 FR	● 502740	1116	50



Conveen® Optima

Male External Catheter - Silicone

Discreet, easy to use, reliable all day self-adhering external catheter. Double-grip strip ensures easy and smooth roll-out for a secure application. Push ring with sure grip ribs provides secure connection to the collecting bag. Anti-kink bellow stop kinks ensuring a free flow. Compact easy-to-open packaging. Latex free.

Conveen® Optima Standard

Length (cm)	Length (inches)	Size	Product Code	Units
7.8 cm	3"	25 mm	22025	30
		28 mm	22028	30
		30 mm	22030	30
		35 mm	22035	30
		40 mm	22040	30

Conveen® Optima Short

Length (cm)	Length (inches)	Size	Product Code	Units
4.6 cm	2"	21 mm	22121	30
		25 mm	22125	30
		30 mm	22130	30
		35 mm	22135	30



Conveen® Security+

Male External Catheter - SEBS*

Easy to use, reliable all day self-sealing external catheter. Applicator grip and single-grip strip ensures easy and smooth roll-out for a secure application. Push ring with sure grip ribs provides secure connection to the collecting bag. Anti-kink bellow stop kinks ensuring a free flow. Latex free.

Conveen® Security+ Self Sealing

Length (cm)	Length (inches)	Size	Product Code	Units
8.2 cm	3"	21 mm	5221	35
		25 mm	5225	35
		25 mm	5225H	100
		30 mm	5230	35
		30 mm	5230H	100
		35 mm	5235	35
		35 mm	5235H	100



Conveen® Security+ Adhesive Liner

Length (cm)	Length (inches)	Size	Product Code	Units
8.2 cm	3"	30 mm	5030	35
		30 mm	5030H	100
		35 mm	5035	35
		35 mm	5035H	100



*Styrene Ethylene Butylene Styrene

Clear Advantage® / Freedom Clear™

Male External Catheter - Silicone

Self-adhering male external catheter. Latex-free.

Clear Advantage® with Aloe Vera						
Length (cm)	Length (inches)	Size		Reference Code	Product Code	Units
7.9 cm	3"	23 mm		504311	C6100	100
		28 mm		504331	6200	100
		31 mm		504341	6300	100
		35 mm		504351	6400	100
Freedom Clear™						
Length (cm)	Length (inches)	Size		Reference Code	Product Code	Units
7.9 cm	3"	23 mm		505711	5100	100
		28 mm		505731	C5200	100
		31 mm		505741	5300	100
		35 mm		505751	5400	100
Freedom Clear™ "LS" Long Seal						
Length (cm)	Length (inches)	Size		Reference Code	Product Code	Units
8.9 cm	3.5"	28 mm		505421	5290	100
		31 mm		505431	5390	100
		35 mm		505341	5490	100
		40 mm		505451	5590	100
Freedom Clear™ "SS" Sport Sheath						
Length (cm)	Length (inches)	Size		Reference Code	Product Code	Units
5.5 cm	2"	23 mm		505311	5110	100
		28 mm		505321	C5210	100
		31 mm		505331	C5310	100
		35 mm		505441	C5410	100



Conveen®

Male External Catheter - Latex

Conveen® Self-Sealing					
Length (cm)	Length (inches)	Size		Product Code	Units
10.5 cm	4"	20 mm		5212	35
		25 mm		5200	35
		30 mm		5205	35
		35 mm		5210	35
		40 mm		5215	35
Conveen® Adhesive Liner					
Length (cm)	Length (inches)	Size		Product Code	Units
10.5 cm	4"	30 mm		5130	35
		35 mm		5135	35

Active Cath[®]

Male External Catheter - Latex

Self-sealing male external catheter.

Active Cath [®]						
Length (cm)	Length (inches)	Size		Reference Code	Product Code	Units
7.9 cm	3"	23 mm		506220	8100	100
		28 mm		506200	8300	100
		31 mm		506240	8305	100
		35 mm		506260	8500	100



Freedom Cath[®]

Male External Catheter - Latex

Self-sealing male external catheter.

Freedom Clear [®]						
Length (cm)	Length (inches)	Size		Reference Code	Product Code	Units
7.9 cm	3"	23 mm		506110	8000	100
		28 mm		506130	8200	100
		31 mm		506140	8205	100
		35 mm		506150	8400	100




Conveen[®] Active

Urine Leg Bag

Developed for an active lifestyle. Integrated in a comfortable elastic band that wraps around the thigh for ease of use, mobility and discretion. Hide-away outlet tap concealed in elastic band when not in use. Soft fabric band for increased wearing comfort and a no-slip silicone feature inside straps for security. Latex and PVC free.

Conveen [®] Active Leg Bag, Non Sterile				
Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
9 oz	250 mL	3"	25501	10



Conveen[®] Security+

Urine Leg bag

Secure and reliable urine leg bag with anti-kink tubing to reduce the risk of accidental leakage. Soft backing for increased wearing comfort (except codes 21034 and 21054). Latex and PVC free.

Conveen [®] Security+ Leg Bag, Clamp Outlet, Straps, Non Sterile				
Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
17 oz	500 mL	20"	5161	10
26 oz	750 mL	20"	5167	10

Conveen [®] Security+ Leg Bag, Clamp Outlet, Sample Port, Straps, Sterile				
Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
17 oz	500 mL	20"	5163	10

Conveen [®] Security+ Leg Bag, Lever Outlet, Straps, Sterile				
Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
17 oz	500 mL	20"	21034	10
34 oz	1000 mL	20"	21054	10



Conveen® Standard

Day / Night Urine Bag

Secure and reliable urine bag for bedside drainage or large capacity leg bag. Clamp outlet makes it easy to verify proper closure and reduce the risk of accidental leakage. Standard tubing. Leg bag straps sold separately. Soft backing for increased wearing comfort. Latex free.

Conveen® Standard Day / Night Bag, Clamp Outlet, Non Sterile

Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
51 oz	1500 mL	36"	5062	10

Conveen® Standard Day / Night Bag, Clamp Outlet, Sample Port, Sterile

Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
51 oz	1500 mL	36"	5063	10



Conveen® Contour

Urine Leg Bag

Secure and reliable urine leg bag. Enhanced chambered design for added discretion. Soft backing for increased wearing comfort. Latex free.

Conveen® Contour Leg Bag, Clamp Outlet, Straps, Non Sterile

Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
21 oz	600 mL	18"	5170	10

28 oz	800 mL	18"	5174	10
-------	--------	-----	------	----



Conveen® Night Bag

Urine Night Bag

Secure and reliable urine night bag with anti-kink tubing to prevent accidental leakage. Large capacity allows the user to feel secure all night for a full night's sleep. Latex and PVC free. Not Compatible with Urostomy Pouches.

Conveen® Bedside Night Bag, Non-Sterile

Capacity (oz)	Capacity (L)	Product Code	Units
68 oz	2L	21346	30

Conveen® Bedside Night Bag, Sterile

Capacity (oz)	Capacity (L)	Product Code	Units
68 oz	2L	21356	10



Accessories

Conveen® Bag Hanger

Length (cm)	Length (inches)	Product Code	Units
18 cm	7"	5070	1



Conveen® Leg Bag Strap

Length (cm)	Length (inches)	Product Code	Units
63 cm	25"	5050	10



Conveen® Catheter retaining Strap

Length (cm)	Length (inches)	Product Code	Units
90 cm	35"	383004	2



Peristeen® Anal Plug

Simple, safe and discreet aid for fecal incontinence. Inserted like a suppository. It can be used for up to 12 hours, preventing the uncontrolled loss of solid stool. Each box comes with a tube of lubricant. Latex free.

Peristeen® Anal Plug			
Length (cm)	Length (inches)	Product Code	Units
3.7 cm]	1.5"	1450	20
4.5 cm	1.75"	1451	20



Peristeen Plus® Anal Irrigation

Transanal irrigation system for fecal incontinence and constipation. Minimizes the likelihood of involuntary bowel leakage and/or constipation.

Peristeen Plus® Anal Irrigation System		
Items	Description	Product Code
- 1 control unit		
- 2 sterile rectal catheters		
- 1 bag (with lid)	** 90 irrigations = 6 months	29140
- 2 straps		
- 1 toiletries bag		



Peristeen Plus® Anal Irrigation System - Small		
Items	Description	Product Code
- 1 control unit		
- 2 sterile rectal catheters - small		
- 1 bag (with lid)	** 90 irrigations = 6 months	29147
- 2 straps		
- 1 toiletries bag		

Peristeen Plus® Anal Irrigation Accessory unit		
Items	Description	Product Code
- 15 sterile rectal catheters	**15 irrigations = 1 month	29142
- 1 bag (no lid)		



Peristeen Plus® Anal Irrigation Accessory unit - Small		
Items	Description	Product Code
- 15 sterile rectal catheters - small	**15 irrigations = 1 month	29149
- 1 bag (no lid)		

Peristeen Plus® Anal Irrigation Rectal Catheter		
Items	Product Code	
- 10 sterile rectal catheters	29143	



Peristeen Plus® Anal Irrigation Leg Strap		
Items	Product Code	
- 1 pair	29145	



Peristeen Plus® Anal Irrigation Tube		
Items	Product Code	
- 2 tubes with blue connectors	29146	



** if used every other day

Wound

Biatain® Silicone

Foam Dressing with Silicone Adhesive

Soft, flexible and absorbent foam dressing with gentle silicone adhesive for use on a broad range of exuding wounds, making it the obvious choice for moist wound healing of commonly occurring wounds. For more information please visit www.biatainsiliconeca.com.



reddot award 2017
winner



Biatain® Silicone			
Size (cm)	Size (inches)	Product Code	Units
7.5 x 7.5 cm	3 x 3"	33434	10
10 x 10 cm	4 x 4"	33435	10
12.5 x 12.5 cm	5 x 5"	33436	10
15 x 15 cm	6 x 6"	33437	5
17.5 x 17.5 cm	7 x 7"	33438	5
10 x 20 cm	3.9 x 7.9"	33400	5
10 x 30 cm	3.9 x 11.8"	33401	5
18 x 18 cm	7.1 x 7.1" Heel	33406	5
25 x 25 cm	9.8 x 9.8" Large Sacral	39000*	5
15 x 19 cm	5.9 x 7.5" Small Sacral	39001*	5
14 x 19.5 cm	5.5 x 7.7" Multi-shape	39010*	5



* This product has the Pressure Injury Prevention indication.

Biatain® Silicone Non-Border

Designed for high absorption and retention under compression to minimise the risk of leakage and maceration of fragile periwound skin. Can also be cut to fit areas that are difficult to dress, e.g. on the heel or in between toes and fingers. For high absorption and retention under compression. 3DFit Technology minimizes risk of leakage and maceration.

Biatain® Silicone			
Size (cm)	Size (inches)	Product Code	Units
5 x 7.5 cm	2 x 3"	39020	10
7.5 x 7.5 cm	3 x 3"	39021	10
10 x 10 cm	4 x 4"	39022	10
12.5 x 12.5 cm	5 x 5"	39023	10
10 x 20 cm	4 x 8"	39024	10
15 x 15 cm	6 x 6"	39025	10
20 x 20 cm	8 x 8"	39026	10



Biatain® Silicone Lite

Thin Foam Dressing with Silicone Adhesive

Thin and ultra flexible foam dressing with gentle silicone adhesive for use on a wide range of low to moderate exuding wounds.



Biatain® Silicone Lite			
Size (cm)	Size (inches)	Product Code	Units
5 x 5 cm	2 x 2"	33452	5
5 x 12.5 cm	2 x 5"	33453	5
7.5 x 7.5 cm	3 x 3"	33444	10
10 x 10 cm	4 x 4"	33445	10
12.5 x 12.5 cm	5 x 5"	33446	10



Biatain® Silicone Ag

Foam Dressing with Silicone Adhesive and Silver

Flexible and absorbent foam dressing with gentle silicone adhesive and a patented silver complex for use on a broad range of exuding infected wounds or wounds at risk of infection.



Biatain® Silicone Ag				
Size (cm)	Size (inches)	Product Code	Units	
7.5 x 7.5 cm	3 x 3"	39636	5	     
10 x 10 cm	4 x 4"	39637	5	
12.5 x 12.5 cm	5 x 5"	39638	5	
15 x 15 cm	6 x 6"	39639	5	
17.5 x 17.5 cm	7 x 7"	39640	5	
10 x 20 cm	3.9 x 7.9"	39644	5	
10 x 30 cm	3.9 x 11.8"	39645	5	
15 x 19 cm	5.9 x 3.5" Small Sacral	39650	5	
25 x 25 cm	9.8 x 9.8" Large Sacral	39651	5	
18 x 18 cm	7.1 x 7.1" Heel	39652	5	

Biatain®

Foam Dressing

Soft and conformable polyurethane foam dressing that effectively absorbs and retains wound exudate. This ensures a balance that is optimal for management of exuding wounds. Available in adhesive and non-adhesive. The adhesive has a hydrocolloid border – for a secure adhesion. This product is also used in high friction and shear areas such as the coccyx.



Biatain® Adhesive				
Size (cm)	Size (inches)		Product Code	Units
7.5 x 7.5 cm	3 x 3"		3462	10
10 x 10 cm	4 x 4"		3430	10
12.5 x 12.5 cm	5 x 5"		3420	10
17 x 17 cm	6.7 x 6.7"	Small Sacral	3483	5
18 x 18 cm	7 x 7"		3423	5
19 x 20 cm	7.5 x 8"	Heel	3488	5
23 x 23 cm	9 x 9"	Large Sacral	3485	5





Biatain® Non-Adhesive				
Size (cm)	Size (inches)		Product Code	Units
5 x 7 cm	2 x 2.8"		6105	10
10 x 10 cm	4 x 4"		3410	10
10 x 20 cm	4 x 8"		3412	5
15 x 15 cm	6 x 6"		3413	5
20 x 20 cm	8 x 8"		3416	5



Biatain® Soft Hold

Foam Dressing with Gentle Adhesion

Soft and conformable polyurethane foam dressing with a gentle adherent layer covering less than 50% of the foam surface. This product is used under compression. This allows both hands to be free during dressing application and removal.



Biatain® Soft Hold			
Size (cm)	Size (inches)	Product Code	Units
10 x 10 cm	4 x 4"	3470	5
10 x 20 cm	4 x 8"	3472	5
15 x 15 cm	6 x 6"	3475	5



Biatain® Ag

Foam Dressing with Silver

Soft and conformable polyurethane foam dressing with a patented silver complex for use on infected wounds or wounds at risk of infection. Available in adhesive and non-adhesive. The adhesive has a hydrocolloid border – for a secure adhesion. This product is also used in high friction and shear areas such as the coccyx.



Biatain® Ag Adhesive			
Size (cm)	Size (inches)	Product Code	Units
7.5 x 7.5 cm	3 x 3"	9631	5
12.5 x 12.5 cm	5 x 5"	9632	5
18 x 18 cm	7 x 7"	9635	5
23 x 23 cm	9 x 9"	Sacral 9641	5
19 x 20 cm	7.5 x 8"	Heel 9643	5



Biatain® Ag Non-Adhesive			
Size (cm)	Size (inches)	Product Code	Units
5 x 7 cm	2 x 2.8"	5105	5
10 x 10 cm	4 x 4"	9622	5
10 x 20 cm	4 x 8"	9623	5
15 x 15 cm	6 x 6"	9625	5
20 x 20 cm	8 x 8"	9626	5




Biatain® Fiber

Suitable for highly exuding wounds and locks in exudate and bacteria. Effectively absorbs and retains exudate and bacteria, thereby reducing exudate pooling. HexaLock Technology is the new gelling fiber for all exuding wounds.



reddot award 2017
winner

Biatain® Fiber			
Size (cm)	Size (inches)	Product Code	Units
5 x 5 cm	2 x 2"	33501	10
10 x 10 cm	4 x 4"	33503	10
15 x 15 cm	6 x 6"	33506	10
19 x 25 cm	7.5 x 9.8"	33507	10
2.5 x 46 cm	1.5 x 18"	33508	10



Biatain® Alginate

Calcium Alginate Dressing

Highly absorbent alginate dressing for use on moderate to heavily exuding wounds of any size, shape and depth.

Biatain® Alginate					
Size (cm)	Size (inches)		Product Code	Units	
5 x 5 cm	2 x 2"		3705	30	
10 x 10 cm	4 x 4"		3710	10	
15 x 15 cm	6 x 6"		3715	10	
3 x 44 cm	1 x 17.5"	Ribbon	3740	6	



Biatain® Alginate Ag

Calcium Alginate Dressing with Silver

Highly absorbent alginate dressing with silver for use on moderate to heavily exuding infected wounds or wounds at risk of infection of any size, shape and depth.

Biatain® Alginate Ag					
Size (cm)	Size (inches)		Product Code	Units	
5 x 5 cm	2 x 2"		3755	30	
10 x 10 cm	4 x 4"		3760	10	
15 x 15 cm	6 x 6"		3765	10	
3 x 44 cm	1 x 17.5"	Ribbon	3780	10	



Biatain® Ibu

Foam Dressing with Ibuprofen

Soft and flexible absorbent foam dressing containing Ibuprofen for use on painful wounds. Available in non-adhesive and soft hold.



Biatain® Ibu Non-Adhesive					
Size (cm)	Size (inches)		Product Code	Units	
10 x 10 cm	4 x 4"		4110	5	
10 x 20 cm	4 x 8"		4112	5	
15 x 15 cm	6 x 6"		4115	5	
20 x 20 cm	8 x 8"		4120	5	

Biatain® Ibu Soft-Hold					
Size (cm)	Size (inches)		Product Code	Units	
10 x 10 cm	4 x 4"		4140	5	
10 x 20 cm	4 x 8"		4142	5	



Biatain® Super

Hydrocapillary Pad Dressing

Highly absorbent hydrocapillary pad dressing for use on heavily exuding wounds. Available in non-adhesive and adhesive.

Biatain® Super Adhesive				
Size (cm)	Size (inches)	Product Code	Units	
10 x 10 cm	4 x 4"	4610	10	
12.5 x 12.5 cm	5 x 5"	4612	10	
12 x 20 cm	4.5 x 8"	4625	10	
15 x 15 cm	6 x 6"	4615	10	
20 x 20 cm	8 x 8"	4620	10	
Biatain® Super Non-Adhesive				
Size (cm)	Size (inches)	Product Code	Units	
10 x 10 cm	4 x 4"	4630	10	
12.5 x 12.5 cm	5 x 5"	4632	10	
12 x 20 cm	4.5 x 8"	4645	10	
15 x 15 cm	6 x 6"	4635	10	
20 x 20 cm	8 x 8"	4639	10	

Comfeel® Plus

Hydrocolloid Dressing

Hydrocolloid dressing for use on low to moderate exuding wounds. Provides faster wound healing by sealing and protecting the wounds.

Comfeel® Plus				
Size (cm)	Size (inches)	Product Code	Units	
4 x 6 cm	1.5 x 2.5"	33146	30	
10 x 10 cm	4 x 4"	33110	10	
15 x 15 cm	6 x 6"	33115	5	
20 x 20 cm	8 x 8"	33120	5	
18 x 20 cm	7 x 8"	Sacral	33285	5
9 x 11 cm	3 x 4"	Contour	33283	5

Comfeel® Plus Transparent

Hydrocolloid Dressing

Thin and flexible hydrocolloid dressing for use on dry or low exuding wounds. Provides faster wound healing by sealing and protecting the wound. Transparency allows for visual monitoring of the wound.

Comfeel® Plus Transparent				
Size (cm)	Size (inches)	Product Code	Units	
5 x 7 cm	2 x 2.75"	33530	10	
10 x 10 cm	4 x 4"	33533	10	
9 x 14 cm	3.5 x 5.5"	33536	10	
15 x 15 cm	6 x 6"	33539	5	
9 x 25 cm	3.5 x 10"	33537	5	
5 x 15 cm	2 x 6"	33547	10	
5 x 25 cm	2 x 10"	33548	10	

Hydrophilic Wound Dressing

Zinc-oxide based hydrophilic sterile coating. Used for autolytic debridement. The ideal alternative for difficult-to-dress areas. Sterile. For more information please visit <https://triad.coloplat.ca/>.

Triad™			
Size (ounces)	Size (grams)	Product Code	Units
2.5 oz.	71 g tube	1964	12
6 oz.	170 g tube	1967	12



Purilon® Gel

Hydrogel with natural ingredients

Effective and gentle debridement.

Purilon® Gel			
Size (ounces)	Size (grams)	Product Code	Units
28 oz.	8 g tube	3906	10
5 oz.	15 g tube	3900	10



Woun'Dres®

Collagen Hydrogel

Rehydrates dry wounds and eschar.

Woun'Dres® Collagen Hydrogel			
Size (ounces)	Size (grams)	Product Code	Units
1 oz.	28 g tube	1166	36
3 oz.	84 g tube	7690	12



Sea-Clens®

Wound Cleanser

Saline based solution with a surfactant for cleansing and irrigating acute and chronic wounds.

Sea-Clens® Wound Cleanser			
Size (ounces)	Size (mL)	Product Code	Units
6 oz.	178 mL	1063	12
12 oz.	355 mL	1061	12



Skin

Sproam®



No-Rinse Spray or Foam Cleanser

Spray or foam body wash, shampoo and incontinence cleanser.

Sproam®		
Size	Product Code	Units
175 ml	0897	12
350 ml	0901	12



Bedside-Care® Sensitive Skin



No-Rinse Foam Cleanser

Unscented foaming body wash, shampoo and incontinence cleanser.

Bedside-Care® Sensitive Skin		
Size	Product Code	Units
125 mL	7300	36
240 mL	7301	12



Bedside-Care® Foam



No-Rinse Foam Cleanser

Scented foaming body wash, shampoo and incontinence cleanser.

Bedside-Care®		
Size	Product Code	Units
240 mL	67145	12



Bedside-Care®



No-Rinse Spray Cleanser

Spray body wash, shampoo and incontinence cleanser.

Bedside-Care®		
Size	Product Code	Units
120 mL Scented	61768	36
240 mL Scented	61762	12
240 mL Unscented	61761	12



No-Rinse Cleanser



Rinse Cleanser



Moisturize



Protect



CHG Compatible



Neonate Safe

Easi-Cleanse® Bath



No-Rinse Wash Cloth.

Body wash, shampoo, and incontinence cleanser in a self-sudsing disposable wash cloth.

Easi-Cleanse		
Size	Product Code	Units
5 per package	7056	100 pkgs/case



Gentle Rain® Extra Mild



Rinse Cleanser

Body wash and shampoo.

Gentle Rain® Extra Mild		
Size	Product Code	Units
5 mL sachet	7231	300
50 mL	7230	36
600 mL (accessories available, see pg 62)	7233	12
1000 mL (accessories available, see pg 62)	7234	12



Gentle Rain® Extra Mild Fragrance Free



Rinse Cleanser

Body wash and shampoo.

Gentle Rain® Extra Mild Fragrance Free		
Size	Product Code	Units
5 mL sachet	7073	300
50 mL	7074	36
600 mL (accessories available, see pg 62)	7075	12



Atrac-Tain®



Moisturizer

Skin moisturizer for severely dry, scaly and fissured skin with 4% Alphahydroxyacid (AHA).

Atrac-Tain®		
Size	Product Code	Units
2 g sachet	1843	300
60 mL	1802	12
140 mL	1814	12



Sween® 24



Moisturizer

Once a day skin moisturizer for itchy, irritated, dry and flaky skin.

Sween® 24		
Size	Product Code	Units
4 g sachet	7090	300
90 g	7094	12



DIN 02250322



DIN 02250322

Sween Cream®



Moisturizer

Soothing body cream for irritated dry skin.

Sween Cream®		
Size	Product Code	Units
90 mL	7067	12
190 mL	7068	12





Sween® Lotion Unscented



Moisturizer

Everyday moisturizer. For all body use.

Sween® Lotion Unscented		
Size	Product Code	Units
60 mL	7222	36
120 mL	7223	36
620 mL	7224	12





No-Rinse Cleanser



Rinse Cleanser



Moisturize



Protect



CHG Compatible



Neonate Safe

Critic-Aid® Clear



Skin Protectant

Protects skin and treats irritation due to incontinence.

Critic-Aid® Clear		
Size	Product Code	Units
4 g sachet	7565	300
71 g	7573	12



DIN 02277778

Baza® Protect II



Skin Protectant

Protects skin and treats irritation due to incontinence. Helps seal out wetness & relieves sore, irritated skin.

Baza® Protect II		
Size	Product Code	Units
60 g	7086	12
140 g	7087	12



DIN 02315793

Baza Cleanse and Protect®



Skin Protectant

All-in-one cleanser, moisturizer and a breathable skin protectant for sore, irritated, chafed skin due to incontinence.

Baza Cleanse and Protect®		
Size	Product Code	Units
4 g sachet (fragrance free)	7710	300
240 mL (fragrance free)	7712	12
240 mL (with fragrance)	7725	12



DIN 02243966



InterDry®

Moisture-wicking textile with antimicrobial silver complex

Effectively manage complications associated with skin folds. Appropriate for a variety of skin-to-skin and skin-to-device application sites. For more information please visit www.interdry.ca.

InterDry®			
Size (cm)	Size (Inches)	Product Code	Units
25 x 46 cm	10 x 18"	7915	40
25 x 92 cm	10 x 36"	7912	10
25 x 366 cm	10 x 144" textile on a roll	7910	10



Isagel®



Hand Sanitizer Gel

Sanitizing hand gel with added moisturizer.

Isagel®		
Size	Product Code	Units
115 mL	1644	36
621 mL (Accessories available. Please see page 62.)	1645	12



Sween® Body Powder

All-body powder

Provides excellent absorption to keep patient dry.

Sween® Body Powder		
Size	Product Code	Units
90 g	0506	36
227 g	0505	36



No-Rinse Cleanser



Rinse Cleanser



Moisturize



Protect



CHG Compatible



Neonate Safe

Accessories

Coloplast® Wall Dispenser for Bag

Description	Product Code	Units
Compatible with 800 mL (Isagel®) and 1000 mL (Gentle Rain® Extra Mild) sized bags (7234 and 7041)*	7251	12



Coloplast® Wall Dispenser Drip Tray

Description	Product Code	Units
To be used with Wall Dispenser to prevent spillage (7251)*	7262	60



Embrace™ Dispenser Bracket

Description	Product Code	Units
Compatible with 600 mL (Gentle Rain® Extra Mild) and 621 mL (Isagel®) sized pump bottles (1645, 7224, 7233)*	0089	12



*Product not included.

Brief Statements

Biatain® Family of Dressings

Biatain® dressings are single-use only. See the device manual for detailed information regarding complete instructions, contraindications, warnings, precautions and potential adverse events. For further questions, call Coloplast Canada Corporation at 1-888-880-8605 and/or consult the company website at www.coloplast.ca.

InterDry®

InterDry® is a skin protectant for the management of skin folds and other skin-to-device contact areas. It reduces microbial colonization in the fabric. This single use fabric provides moisture transportation to keep skin dry while the antimicrobial in the fabric reduces odor. The safe use of InterDry fabric during pregnancy, lactation and on children has not been demonstrated. Although many patients benefit from the use of this device, results may vary. See the device labeling for detailed information including the indication, contraindications, warnings, precautions, and potential complications/adverse events. For further questions, call Coloplast Canada Corporation at 1-888-880-8605 and/or consult the company website at www.coloplast.ca.

Peristeen Anal Plug®

The Peristeen Anal Plug is intended as a barrier to faeces in the rectum. The anal plug can be used by people who are faecal incontinent or suffer from involuntary passing of air (flatus). It functions as an effective barrier to faeces in the rectum.

Contraindication: Peristeen Anal Plug must not be used, if you suffer from: hemorrhoids (3rd and 4th degree), and anal stenosis

Caution: This product should not be used without assessment by your health care professional. Contact your health care professional before you use Peristeen Anal Plug for the first time, as professional guidance in the use of the plug is vital. The anal plug should not be used if you are suffering from diarrhoea. If placed in the mouth, the Peristeen Anal Plug will expand and may cause choking. Remove the anal plug when you feel the urge to defecate. If you have no sensation in the rectal area and do not feel the urge to defecate, changing plug at regular intervals is recommended.

The Peristeen Anal Plug must be removed no later than 12 hours after insertion. If you experience difficulties removing the Peristeen Anal Plug, you must contact your health care professional. If you repeatedly experience discomfort when using the Peristeen Anal Plug, stop using it and contact your health care professional. Keep out of reach of children. See the device manual for detailed information regarding instructions, contraindications, warnings, precautions, and potential complications/adverse events or, call Coloplast Canada Corporation at 1-866-293-6349 and/or consult the company website at www.coloplast.ca.

Peristeen®

The Peristeen® Anal Irrigation System is intended to instill water into the colon through a rectal catheter, which incorporates an inflatable balloon, inserted into the rectum to promote evacuation of the contents of the lower colon. The Peristeen Anal Irrigation System is indicated for use by children (3 years - <12 years old), adolescent (12 years - < 18 years old), transitional adolescent (18 - <21 years old) and adult patients with neurogenic bowel dysfunction who suffer from fecal incontinence, chronic constipation, and/or time-consuming bowel management procedures.

Contraindications: Peristeen Anal Irrigation must not be used in the following situations: known obstruction of the large bowel due to strictures or tumors, acute inflammatory bowel disease, diverticulitis, complex diverticular disease, abdominal or anal surgery within the last 3 months, spinal shock phase, and in patients who are pregnant and have not used the system before (if the patient is pregnant and has never used anal irrigation before, they should not start the irrigation procedure during pregnancy).

Warning: Anal irrigation procedure should always be carried out with care. Bowel perforation is an extremely rare, but serious and potentially lethal complication to anal irrigation and will require immediate admission to a hospital, often requiring surgery.

See the device manual for detailed information regarding the procedure, contraindications, warnings, precautions, and potential complications/adverse events. For further information, call Coloplast Canada Corporation at 1-866-293-6349 and/or consult the company website at www.coloplast.ca.

SpeediCath®

SpeediCath® catheters are indicated for use by patients with chronic urine retention and patients with a post void residual volume (PVR) due to neurogenic and non-neurogenic voiding dysfunction. The catheter is inserted into the urethra to reach the bladder allowing urine to drain. There is a separate SpeediCath device intended for either males or females. This is a single-use only product. Patients are to contact their healthcare professional before performing self-catheterization for the first time and follow the instructions carefully. See the device manual for further information or, call Coloplast Canada Corporation at 1-866-293-6349 and/or consult the company website at www.coloplast.ca.

Self-Cath®

Self-Cath® Catheters are intended for use in male, female, and pediatric patients (neonates, infants, children, adolescents, and transitional adolescents) requiring bladder drainage as determined by their physician. The catheter is indicated for those individuals unable to promote a natural urine flow or for those individuals who have a significant volume of residual urine following a natural bladder-voiding episode. Single-use only. See the device manual for detailed information regarding the procedure, contraindications, warnings, precautions, and potential complications/adverse events. For further information, call Coloplast Canada Corporation at 1-866-293-6349 and/or consult the company website at www.coloplast.ca.

Triad™

Triad™ hydrophilic wound dressing is indicated for the local management of pressure and venous stasis ulcers, dermal lesions/injuries, superficial wounds, scrapes, first-and second-degree burns, partial-and full-thickness wounds. Suitable for use on neonates.

Contraindications: Triad™ is contraindicated for third degree burns or infected wounds. See device label for complete information including warnings and precautions. For further questions, call Coloplast Canada Corporation at 1-888-880-8605 and/or consult the company website at www.coloplast.ca.

Product Code	Page
0089	66
112	44
114	44
116	44
208	43
210	43
212	43
214	43
240	43
305	43
306	43
308	43
310	43
408	44
410	44
412	44
414	44
416	44
418	44
0421	35
0422	35
0423	35
0424	35
450	44
460	44
475	43
0505	65
0506	65
0509	35
608	44
610	44
612	44
614	44
808	44
810	44
812	44
814	44
816	44
818	44
0897	61
0901	61
0925	35
1003	37
1061	58
1063	58
1110	37
1114	44
1116	44
1166	58
1450	49
1451	49
1500	37
1644	65
1645	65
1802	62
1814	62

Product Code	Page
1843	62
1907	35
1964	58
1967	58
2041	35
2115	33
2125	33
2151	33
2166	33
2181	32
2320	35
2325	35
2330	35
2340	35
2350	35
2490	30
2501	36
2520	30
2650	35
2655	34
2792	28
2793	28
2801	36
2802	36
2803	36
2807	36
2808	36
2831	27
2832	27
2833	27
2836	27
2837	27
2846	27
2847	27
2882	26
2883	26
3210	33
3215	33
3220	33
3410	54
3412	54
3413	54
3416	54
3420	54
3423	54
3430	54
3462	54
3470	55
3472	55
3475	55
3483	54
3485	54
3488	54
3705	56
3710	56
3715	56

Product Code	Page
3740	56
3755	56
3760	56
3765	56
3780	56
3900	58
3906	58
4110	56
4112	56
4115	56
4120	56
4140	56
4142	56
4610	57
4612	57
4615	57
4620	57
4625	57
4630	57
4632	57
4635	57
4639	57
4645	57
5030	45
5030H	45
5035	45
5035H	45
5050	48
5062	48
5063	48
5070	48
5100	46
5105	55
5110	46
5130	46
5135	46
5161	47
5163	47
5167	47
5170	48
5174	48
5200	46
5205	46
5210	46
5212	46
5215	46
5221	45
5225	45
5225H	45
5230	45
5230H	45
5235	45
5235H	45
5290	46
5300	46
5390	46

Product Code	Page
5400	46
5490	46
5580	32
5585	32
5590	46
6105	54
6200	46
6300	46
6400	46
7056	62
7067	63
7068	63
7073	62
7074	62
7075	62
7086	64
7087	64
7090	63
7094	63
7222	63
7223	63
7224	63
7230	62
7231	62
7233	62
7234	62
7251	66
7262	66
7300	61
7301	61
7565	64
7573	64
7583	36
7690	58
7710	64
7712	64
7725	64
7910	65
7912	65
7915	65
8000	47
8009	33
8100	47
8200	47
8205	47
8300	47
8305	47
8400	47
8500	47
9500	36
9505	36
9622	55
9623	55
9625	55
9626	55
9631	55

Product Code	Page
9632	55
9635	55
9641	55
9643	55
10011	20
10012	20
10013	20
10015	21
10016	21
10017	21
10021	20
10022	20
10023	20
10025	21
10026	21
10027	21
10031	20
10032	20
10035	21
10036	21
10041	20
10045	21
10101	22
10102	22
10103	22
10104	23
10105	23
10106	23
10108	22
10111	22
10113	22
10114	22
10116	22
10117	22
10154	21
10155	21
10156	21
10164	21
10165	21
10166	21
10167	21
10451	13
10461	13
10471	13
10472	13
10473	13
10474	13
10475	13
10476	13
10477	13
10478	13
10479	13
10480	13
10481	13
10482	14
10483	14

Product Code	Page
10484	14
10487	14
10489	13
10502	9
10503	9
10504	9
10505	9
10506	9
10512	9
10513	9
10514	9
10515	9
10516	9
10522	9
10523	9
10524	9
10525	9
10531	9
10551	12
10561	12
10562	12
10563	12
10564	12
10565	12
10566	12
10567	12
10568	12
10571	12
10572	12
10573	12
10582	17
10583	17
10584	17
10585	17
10586	17
10870	16
10871	16
10872	16
10873	16
10874	16
10881	16
10883	16
10884	16
10885	16
10904	23
10905	23
10914	23
10915	23
10918	23
11011	20
11012	20
11015	21
11016	21
11021	20
11022	20
11023	20

Product Code	Page
11024	20
11025	21
11026	21
11027	21
11028	21
11031	20
11032	20
11035	21
11036	21
11041	20
11045	21
11114	21
11115	21
11124	21
11125	21
11126	21
11127	21
11135	21
11136	21
11185	21
11186	21
11194	21
11195	21
11196	21
11197	21
11301	22
11302	22
11303	22
11304	23
11305	23
11306	23
11307	22
11308	23
11312	22
11313	22
11314	22
11315	22
11316	22
11505	23
11514	23
11515	23
11518	23
11612	23
11615	23
11690	23
11801	25
11802	25
11803	25
11804	25
11811	25
11812	25
11814	25
11815	25
11821	23
11822	23
11824	23

Product Code	Page
11825	23
11827	23
11828	23
11841	22
11842	22
11844	22
11845	22
11846	22
11851	22
11852	22
11854	22
11855	22
11856	22
12000	34
12010	34
12011	34
12014	19
12020	34
12021	34
12030	34
12032	34
12033	34
12035	34
12036	34
12037	34
12038	34
12039	34
12042	34
12045	34
12046	34
12047	34
12048	34
12049	34
12050	35
12060	35
12061	35
12070	35
12072	35
12074	35
12076	35
12090	34
12091	34
12092	34
12093	34
12094	34
12095	34
12110	31
12130	31
12140	31
12144	31
12145	31
12146	31
12147	31
12170	32
12174	32
12175	32

Product Code	Page
12176	32
12177	32
12201	13
12202	13
12212	13
12213	13
12222	13
12232	12
12243	12
12262	12
12271	12
12272	12
12274	12
12282	12
12284	12
12291	13
12292	13
12297	13
12298	13
12344	28
12346	28
12354	28
12355	28
12356	28
12374	28
12375	28
12376	28
12384	28
12385	28
12386	28
12432	32
12433	32
12471	32
12472	32
12533	30
12595	32
12596	32
12701	26
12702	26
12703	26
12704	26
12705	26
12706	26
12707	26
12708	26
12709	26
12716	26
12719	26
12800	31
12802	31
12805	31
12808	31
12810	31
12834	37
12835	37
12836	37

Product Code	Page
12841	26
12842	26
12844	26
12845	26
12847	26
12848	26
12991	32
12992	32
12993	32
12994	32
12995	32
12996	32
12997	32
13671	14
13672	14
13673	14
13674	14
13675	14
13676	14
13678	17
13679	17
13680	17
13706	30
13840	30
13844	30
13845	30
13846	30
13860	30
13870	30
13874	31
13875	31
13876	31
13924	28
13925	28
13964	28
13965	28
13966	28
13974	28
13975	28
13976	28
13984	28
13985	28
13986	28
14010	38
14015	38
14021	38
14031	38
14041	38
14050	38
14060	38
14070	38
14103	31
14106	31
14163	31
14164	31
14165	31

Product Code	Page
14166	31
14173	31
14174	31
14175	31
14176	31
14196	30
14201	32
14202	32
14224	28
14225	28
14226	28
14227	28
14228	28
14229	28
14232	27
14233	27
14234	27
14235	27
14236	27
14237	27
14238	27
14239	27
14243	27
14246	27
14249	27
14251	27
14252	27
14254	27
14255	27
14258	27
14261	26
14262	26
14263	26
14271	26
14272	26
14273	26
14274	26
14275	26
14276	26
14277	26
14278	26
14281	27
14282	27
14283	27
14291	27
14292	27
14293	27
14294	27
14295	27
14296	27
14301	29
14302	29
14303	29
14304	29
14305	29
14306	29

Product Code	Page
14307	33
14308	33
14309	29
14316	30
14317	30
14319	30
14329	30
14331	30
14332	30
14342	29
14344	29
14346	29
14347	29
14349	29
14355	30
14356	29
14357	29
14359	29
14361	29
14362	29
14364	29
14401	29
14402	29
14403	29
14404	30
14405	30
14406	29
14408	29
14411	30
14412	30
14414	30
14415	30
14422	31
14423	31
14424	31
14441	32
14445	32
14452	32
14453	32
14454	32
14494	28
14498	28
14520	29
14521	29
14524	29
14525	29
14527	29
14528	29
14552	30
14553	30
14554	30
14555	30
14601	29
14602	29
14603	29
14642	29

Product Code	Page
14643	29
14644	29
14681	33
14682	33
14691	33
14692	33
14701	32
14702	32
15403	25
15404	25
15405	25
15406	25
15407	25
15408	25
15413	25
15414	25
15415	25
15441	25
15450	25
15470	25
15472	25
15473	25
15474	25
15480	25
15482	25
15483	25
15484	25
15521	24
15522	24
15523	24
15524	24
15531	24
15532	24
15533	24
15534	24
15605	24
15606	24
15607	24
15608	24
15633	24
15634	24
15635	24
15695	24
15696	24
15821	24
15981	24
16301	16
16305	16
16306	16
16310	16
16311	16
16312	16
16313	16
16315	16
16322	16
16323	16

Product Code	Page
16325	16
16326	16
16327	16
16330	16
16331	16
16332	16
16333	16
16335	16
16345	16
16346	16
16347	16
16350	16
16351	16
16352	16
16353	16
16471	12
16472	12
16481	12
16483	12
16484	12
16485	12
16486	12
16487	12
16488	12
16491	12
16493	12
16494	12
16705	14
16706	14
16715	14
16716	14
16725	14
16726	14
16727	14
16730	14
16731	14
16732	14
16733	14
16734	14
16735	14
16737	14
16738	14
16739	14
16748	14
16755	14
16756	14
16757	14
16760	14
16761	14
16762	14
16763	14
16764	14
16767	14
16768	14
16770	14
16771	14

Product Code	Page
16772	14
16773	14
16778	14
16801	13
16805	17
16810	17
16825	17
16826	17
16827	17
16835	17
16836	17
16837	17
16840	17
16841	17
16842	17
16855	17
16856	17
16857	17
16865	17
16866	17
16867	17
16870	17
16871	17
16872	17
16901	9
16903	9
16911	9
16913	9
16914	9
16915	9
16921	9
16923	9
16924	9
16931	9
16933	9
16934	9
16941	9
16942	9
16951	9
16953	9
16954	9
16955	9
16961	9
16963	9
16964	9
16971	9
16973	9
16974	9
18151	17
18152	17
18160	17
18161	17
18162	17
18163	17
18170	17
18171	17

Product Code	Page
18172	17
18173	17
18180	17
18181	17
18182	17
18185	17
18300	15
18301	15
18302	15
18305	15
18310	15
18312	15
18313	15
18314	15
18315	15
18320	15
18321	15
18322	15
18323	15
18324	15
18330	15
18331	15
18332	15
18333	15
18340	15
18341	15
18342	15
18343	15
18344	15
18386	12
18389	13
18390	12
18401	18
18402	18
18403	18
18404	18
18405	18
18500	10
18501	10
18510	10
18511	10
18512	10
18513	10
18520	10
18521	10
18522	10
18523	10
18524	10
18530	10
18531	10
18532	10
18650	13
18651	13

Product Code	Page
18652	13
18655	13
18656	13
18657	13
18660	18
18661	18
18665	18
18666	18
18670	18
18671	18
18675	18
18676	18
18680	18
18681	18
18690	18
18691	18
18700	19
18705	19
18710	20
18711	20
18712	20
18713	20
18720	19
18721	19
18722	19
18723	19
18724	19
18725	19
18730	19
18731	19
18735	19
18747	19
19001	25
19003	25
19010	25
19011	25
19020	25
19021	25
19030	25
19040	22
19041	22
19042	22
19050	22
19051	22
19052	22
21034	47
21054	47
21346	48
21356	48
21365	37
22025	45
22028	45
22030	45


Product Code	Page
22035	45
22040	45
22121	45
22125	45
22130	45
22135	45
25501	47
27408	41
27410	41
27412	41
27414	41
27416	41
27418	41
27490	41
27492	41
27494	41
27496	41
27508	41
27510	41
27512	41
27514	41
27516	41
27706	41
27708	41
27710	41
28110	42
28112	42
28114	42
28422	43
28520	43
28522	43
28524	43
28578	42
28580	42
28582	42
28584	42
28692	42
28810	42
28812	42
28814	42
28916	41
28920	41
28922	41
28924	41
29140	49
29142	49
29143	49
29145	49
29146	49
29147	49
29149	49
33110	57
33115	57

Product Code	Page
33120	57
33146	57
33283	57
33285	57
33400	53
33401	53
33406	53
33434	53
33435	53
33436	53
33437	53
33438	53
33444	53
33445	53
33446	53
33452	53
33453	53
33501	55
33503	55
33506	55
33507	55
33508	55
33530	57
33533	57
33536	57
33537	57
33539	57
33547	57
33548	57
39000*	53
39001*	53
39010*	53
39020	53
39021	53
39022	53
39023	53
39024	53
39025	53
39026	53
39636	54
39637	54
39638	54
39639	54
39640	54
39644	54
39645	54
39650	54
39651	54
39652	54
61761	61
61762	61
61768	61
67145	61

Product Code	Page
383004	48
504311	46
504331	46
504341	46
504351	46
505311	46
505321	46
505331	46
505341	46
505421	46
505431	46
505441	46
505451	46
505711	46
505731	46
505741	46
505751	46
506110	47
506130	47
506140	47
506150	47
506200	47
506220	47
506240	47
506260	47
C5200	46
C5210	46
C5310	46
C5410	46
C6100	46

Coloplast® Care

Support for your ostomy or continence care needs



Website and emails

A dedicated website and regular emails with articles and reliable advice

Self-Assessment Tools

Helping users keep track of their ostomy or continence health through regular self-assessment

Phone support

Advisors ready to answer questions and provide support by phone

Product guidance

Access to product samples on request

Enrolling in Care is easy!



Online at www.coloplast.ca/careform



By phone at 1-866-293-6349



Request Coloplast Care enrollment forms from your Coloplast Territory Manager

Velcro® is a registered trademark of Velcro Industries B.V.

[Ostomy Care](#) / [Continence Care](#) / [Wound & Skin Care](#) / [Interventional Urology](#)

Coloplast Corp. Oakville, ON L6H 5S9 / 1-888-880-8605
www.coloplast.ca

The Coloplast logo is a registered trademark of Coloplast A/S. © 2022 Coloplast Corp. All rights reserved.

