

SpeediCath®

Setting the standard





We didn't invent the catheter. We just constantly reinvent it.

When it comes to the design of any SpeediCath catheter, our starting point is always users' needs. The commitment to understanding the daily lives of the people who use our products helped us discover the important connection between design and product satisfaction when it comes to improving compliance¹. We know that

when people have a product they are confident in using, and which suits their lifestyle, they are more likely to use it. This approach has helped us set the standard of catheter design and performance for over 15 years. And while we'll never pretend to understand all your patients' needs, we'll never stop trying to meet them.

SpeediCath® Our range,



For men

a SpeediCath Flex

The first dry-sleeve catheter with a flexible tip

b SpeediCath Compact Male

The first compact male catheter

c SpeediCath Compact Set

The first compact all-in-one catheter and bag solution

d SpeediCath Standard

The proven and reliable standard catheter

your choice



f SpeediCath Compact Female Plus SpeediCath Compact Female The most compact female catheter

g SpeediCath Compact Set
The first compact all-in-one catheter and bag solution

h SpeediCath Standard The proven and reliable catheter

Why SpeediCath® catheters are the **most preferred** in Europe²

SpeediCath catheters have one thing in common – they are all designed with users' needs in mind. This focused dedication has helped make SpeediCath the number one catheter brand in Europe².

A commitment to users' needs

No two catheter users have precisely the same needs. Factors such as hand dexterity and mobility, and concerns over hygiene and discretion can all play critical roles in patient acceptance and compliance. We believe that by providing a range of products designed to meet the different needs of your patients, we give them the best chance of succeeding with intermittent catheterization.


A commitment to quality and safety

Regardless of which SpeediCath you choose, our commitment to quality and safety ensures high performance. All catheters have a unique and uniform hydrophilic coating that also covers the polished eyelets. This ensures a smooth catheterization during both


insertion and withdrawal.^{3,4,5} All SpeediCath catheters are PVC- and phthalate-free, ensuring patient health is not compromised and limiting the environmental impact.

A commitment to compliance

At Coloplast, we know that when a user is satisfied with their product and it fits into their daily lives, it greatly enhances the likelihood of compliance. That's why when we talk about design, we don't mean just the way something looks, but more the way it works. For us, good design is reflected in a user's seamless experience with one of our catheters. If a product is easy to use, the user quickly becomes confident using it and trusts that it is working optimally. That's our goal with every catheter.



SpeediCath® A history of *innovation*





SpeediCath® Flex

The only dry-sleeve catheter with a flexible tip

To help men get off to a good start with intermittent catheterization, SpeediCath Flex is designed to make it easier to teach and perform a clean technique. From opening to disposal, SpeediCath Flex has a range of features that make every step of catheterization easy – the first time and every time.

The benefits of SpeediCath Flex

- Flexible tip for easy and gentle navigation
- Bacteria proof closed-loop system protecting the catheter
- Dry-sleeve for clean and hygienic catheterization
- Double re-closure for discreet and hygienic disposal





SpeediCath® Compact

The first compact catheter for men

The revolutionary and award-winning design of SpeediCath Compact Male has shown to improve the ability to be compliant⁶. It provides discretion without compromising on performance, safety or hygiene. Thanks to the patented telescopic design, the catheter is straightforward to navigate and control during insertion.

The benefits of SpeediCath Compact Male

- Improves ability to be compliant⁶
- Improves quality of life related to intermittent catheterization⁷
- No-touch insertion^{8, 9, 10}
- Re-closable after use for discreet and hygienic disposal



SpeediCath® Compact Set

The first compact all-in-one catheter and bag solution for men

The simple, intuitive and award winning design of SpeediCath Compact Set, provides users with more flexibility and freedom in terms of when and where to catheterize. The compact and discreet design marks a paradigm shift in catheter set solutions.

The benefits of SpeediCath Compact Set

- · Compact catheter with an integrated, sterile bag
- Improves quality of life related to intermittent catheterization⁷
- No-touch insertion^{8, 9, 10}
- Re-closable after use for discreet and hygienic disposal



red<mark>dot</mark> design award winner 2013

SpeediCath® Standard

The proven and reliable standard catheter

SpeediCath Standard is the first instantly ready to use hydrophilic-coated catheter and has been setting the standards for over 15 years. It is the proven and reliable choice.

The benefits of SpeediCath Standard

- Available in all sizes including paediatric
- Available both with Straight and Coude tip
- Easy opening with ring-pull feature
- Adhesive dot for stable placement





SpeediCath® Compact Eve

The first triangular catheter design

By listening to the emotional needs of women, SpeediCath Compact Eve is taking female catheters to a new level, by combining function with an esthetic design. This positively influences acceptance of intermittent catheterization for women.

The benefits of SpeediCath Compact Eve

- Increases acceptance of intermittent catheterization¹¹
- Improves quality of life related to intermittent catheterization⁷
- One-step easy opening and secure reclosing
- Triangular shape for easy handling and better grip
- · Remains sterile when connecting a urine bag



SpeediCath® Compact

The most compact catheter for women

No bigger than a lipstick, SpeediCath Compact Female was a revolution in catheter design when it was launched. Still one of the most popular catheters in the speediCath range, it continues to set the standards for design and performance.

The benefits of SpeediCath Compact Female

- Improves quality of life related to Intermittent catheterization⁷
- Compact, discreet design
- Available with extra length and in paediatric sizes



SpeediCath® Compact Set

The first compact all-in-one catheter and bag solution for women

The simple, intuitive and award winning design of SpeediCath Compact Set, provides users with more flexibility and freedom in terms of when and where to catheterize. The compact and discreet design marks a paradigm shift in catheter set solutions.

The benefits of SpeediCath Compact Set

- Compact catheter with an integrated, sterile bag
- Improves quality of life related to intermittent catheterization⁷
- Re-closable after use for discreet and hygienic disposal









SpeediCath® Standard

The proven and reliable standard catheter

SpeediCath Standard is the first instantly ready to use hydrophilic-coated catheter and has been setting the standards for over 15 years. It is the proven and reliable choice.

The benefits of SpeediCath Standard

- Available in all sizes including paediatric
- · Easy opening with ring-pull feature
- Adhesive dot for stable placement



Product overview

SpeediCath® assortment - Men

Туре		FR	Length (inches/cms)	Product code	Units
SpeediCath Flex	Male	10	17/44	28920	30
	Male	12	17/44	28922	30
	Male	14	17/44	28924	30
	Male	16	17/44	28916	30
SpeediCath® Standard	Male	8	16/40	27408	30
	Male	10	16/40	27410	30
	Male	12	16/40	27412	30
	Male	14	16/40	27414	30
	Male	16	16/40	27416	30
	Male	18	16/40	27418	30
	Male Coude	10	16/40	27490	30
	Male Coude	12	16/40	27492	30
	Male Coude	14	16/40	27494	30
	Male Coude	16	16/40	27496	30
	Paediatric	6	8/20	27706	30
	Paediatric	8	8/20	27708	30
	Paediatric	10	8/20	27710	30
SpeediCath® Compact	Male	12/18	12/30	28692	30
SpeediCath® Compact Set	Male	12/18	12/30	28422	20

SpeediCath® assortment - Women

Туре		FR	Length (inches/cms)	Product code	Units
SpeediCath® Standard	Female	8	8/20	27508	30
	Female	10	8/20	27510	30
	Female	12	8/20	27512	30
	Female	14	8/20	27514	30
	Female	16	8/20	27516	30
	Paediatric	6	8/20	27706	30
	Paediatric	8	8/20	27708	30
	Paediatric	10	8/20	27710	30
SpeediCath® Compact Eve	Female	10	3.3/8.3	28110	30
	Female	12	3.3/8.3	28112	30
	Female	14	3.3/8.3	28114	30
SpeediCath® Compact Set	Female	10	3.5/9	28520	20
	Female	12	3.5/9	28522	20
	Female	14	3.5/9	28524	20
			3.5/5		

Brief statement: SpeediCath catheters are indicated for use by patients with chronic urine retention and patients with a post void residual volume (PVR) due to neurogenic and non-neurogenic voiding dysfunction. The catheter is inserted into the urethra to reach the bladder allowing urine to drain. This is a single use only product. Patients are to contact their healthcare professional before performing self-catheterization for the first time and follow the instructions carefully. See the device manual for further information or, call Coloplast Canada Corporation at 1-877-820-7008 and/or consult the company website at www.coloplast.ca.

- 1. Barbosa et al. 2012
- Coloplast sales data, GERS, IMS, Assobiomedica, Nefemed, PCA, 2015/2016 Cardenas DD et al. PM&R Vol. 3 408-417, May 2011 De Ridder DJ et al. Eur Urol 2005; 48:991-995 Chartier-Kastler E. and Denys P. DOI 10.1002/nau.20929 SCCM Evaluation (n=482), November 2016, Data on file

- Chartier- Kastler E. et al. J Urol. 2013 Sep;190(3):942-7
 Vahr S et al. European Association of Urology Nurse 2013
 Chartier- Kastler E. Spinalcord 2011, 49(7), 844-850
 Bagi P et al. Urol Int 2011; 86:179-184
 SCCE Evaluation (n=161), September 2015, Data on file

The Coloplast story began back in 1954. Our company reflects the passion, ambition and commitment of the people who made it happen. Elise Sørensen, a nurse, invented the first disposable, self-adhesive ostomy pouch, because she was determined to help her sister out of isolation. Aage Louis-Hansen added his engineering ability and entrepreneurial drive, and through strong commitment and resilience he founded Coloplast.

Still today, people with intimate healthcare needs often live in isolation.

We are fighting to change that.

