

Coloplast Product Guide



Leanne,
2024 Paralympic
Paratriathlete for
Team Canada.
SpeediCath® Compact
Female user.

Ostomy, Continence, Advanced Wound Dressings
Effective June, 2024





Leanne,
2024 Paralympic
Paratriathlete for
Team Canada.
SpeediCath® Compact
Female user.

Leanne is 31 years old from England, and moved to Canada when she was 10. She is now a proud Canadian living in Winnipeg, Manitoba. In 2018, a cycling accident left her completely paralyzed from the waist down. Following her spinal cord injury, she fell in love with paratriathlon and went on to complete her first race just 8 months after the accident. When a fellow para-athlete introduced her to SpeediCath® catheters, she was so happy to see that they were designed with accessibility in mind. Leanne's favourite product is the SpeediCath® Compact Female catheter because it comes pre-lubricated so she can use it right out to the package. After years of training, Leanne officially qualified for the Paris 2024 Summer Paralympic Games and will be representing Team Canada this September.

Leanne
Winnipeg, Canada
SpeediCath® User

OSTOMY

Appliances

SenSura® Mio Click

Ostomy Barrier and Pouch - Two Piece System

SenSura® Mio Click Skin Barrier	9
SenSura® Mio Click Drainable Pouch	10
SenSura® Mio Click Drainable Pouch	10
SenSura® Mio Click Closed Pouch	10
SenSura® Mio Click Urostomy Pouch	11
SenSura® Mio Click High Output Pouch	11

SenSura® Mio Flex

Ostomy Barrier and Pouch - Two Piece System

SenSura® Mio Flex Skin Barrier	14
SenSura® Mio Flex Drainable Pouch	15
SenSura® Mio Flex Closed Pouch	15
SenSura® Mio Flex Urostomy Pouch	15
SenSura® Mio Flex High Output Pouch	15

SenSura® Mio

Ostomy Barrier and Pouch - One Piece System

SenSura® Mio Drainable Pouch	16
SenSura® Mio 1 pc, Closed	16
SenSura® Mio Drainable Pouch	17
SenSura® Mio Closed Pouch	19
SenSura® Mio One Piece Urostomy	20
SenSura® Mio One Piece High Output Pouch	21
SenSura® Mio One Piece Post-Op Pouch	21

SenSura® Mio Baby

Ostomy Barrier and Pouch - Two Piece System

SenSura® Mio Baby Flex Skin Barrier	22
SenSura® Mio Baby Flex Drainable Pouch	22
Brava® Baby Adhesive Remover Wipe	21

SenSura® Mio Kids

Ostomy Barrier and Pouch - Two Piece System

SenSura® Mio Kids Flex Skin Barrier	22
SenSura® Mio Kids Flex Drainable Pouch	22
SenSura® Mio Kids Flex Urostomy Pouch	22

Ostomy Barrier and Pouch - One Piece System

SenSura® Mio Kids Drainable Pouch 1 Piece	23
SenSura® Mio Kids Pouch 1 Piece Urostomy	23

SenSura® Click

Ostomy Barrier and Pouch- Two Piece System

SenSura® Click Skin Barrier	23
SenSura® Click Xpro Skin Barrier	24
SenSura® Click Drainable Pouch	24
SenSura® Click Closed Pouch	24
SenSura® Click Urostomy Pouch	25
SenSura® High Output Pouch	25

SenSura® Flex

Ostomy Barrier and Pouch - Two Piece System

SenSura® Flex Skin Barrier	25
SenSura® Flex Xpro Skin Barrier	26
SenSura® Flex Drainable Pouch	26
SenSura® Flex Closed Pouch	26
SenSura® Flex Urostomy Pouch	26

SenSura® - One Piece

Ostomy Barrier and Pouch - One Piece System

SenSura® Drainable Pouch with Easyclose™	27
SenSura® Xpro Drainable Pouch	27
SenSura® Closed Pouch	28
SenSura® Urostomy Pouch	28
SenSura® One-Piece High Output Pouch	28
SenSura® Post-Op Pouch	28

Assura®

Ostomy Barrier and Pouch - Two Piece System

Assura® Skin Barrier	29
Assura® Skin Barrier - Extra Extended Wear	30
Assura® 2 pc. High Output Drainable Pouch	30
Assura® 2 pc. High Output Drainable System	30
Assura® Drainable Pouch Clamp Closure	31
Assura® Drainable Pouch with EasiClose™	31
Assura® Closed Pouch	31
Assura® Urostomy Pouch	31

Assura® - One Piece

Ostomy Barrier and Pouch - One Piece System

Assura® Drainable Pouch Clamp Closure	32
Assura® Drainable Pouch with EasiClose™	32
Assura® Drainable Pouch with EasiClose™	32
Assura® Closed Pouch	33
Assura® 1-Piece Post-Op Pouch	33
Assura® Closed Pouch	33
Assura® Urostomy Pouch	34
Assura® Urostomy Pouch - Extra Extended Wear	34
Coloplast Urostomy Night Bag	34

Assura® - Pediatric - Two Piece

Ostomy Barrier and Pouch - Two Piece System - Pediatric

Assura® Pediatric Skin Barrier	34
Assura® Pediatric Drainable Pouch	34
Assura® Pediatric Closed Pouch	34

EasiFlex® - Pediatric

Ostomy Barrier and Pouch - Two Piece System - Pediatric

Easiflex® Pediatric Skin Barrier	35
Easiflex® Pediatric Drainable Pouch	35

Assura® - Pediatric - One Piece

Ostomy Barrier and Pouch - One Piece System - Pediatric	
Assura® Pediatric Drainable Pouch	35
Assura® Pediatric Closed Pouch	35
Assura® Pediatric Urostomy Pouch	35

Accessories

Ostomy Accessories

Brava® Protective Sheets	38
Brava® Protective Rings	38
Brava® Protective Rings Wide	38
Brava® Protective Rings Extra Wide	38
Brava® Protective Ring Convex	38
Brava® Mouldable Rings	38
Brava® Skin Barrier Spray and Wipe	38
Brava® Adhesive Remover Spray and Wipe	39
Brava® Barrier Cream	39
Brava® Ostomy Strip Paste	39
Brava® Elastic Barrier Strips	39
Brava® Ostomy Powder	39
Brava® Mio Belt for All SenSura Mio Products	39
Brava® Belt for SenSura Original + Assura Products	39
Brava® Lubricating Deodorant	39
Brava® Alcohol Free Paste	39
Ostomy Paste (contains alcohol)	39
Skin Barrier Rings	40
Coloplast Prep™ Protective Liquid Skin Barrier	40
Coloplast Filtrodor® Pouch Filter	40
Coloplast Pouch Clamp	40
Ostomy Scissors	40

Hex-On®

Odour antagonist	
Hex-On®	40

Assura® - Ostomy Cap

Coloplast Uro Minicap (100 ml)	41
Coloplast Stoma Cap	41
Coloplast Stoma Cap - One Piece	41

Mio Irrigation Products

Mio Irrigation One Piece Sleeve	41
---------------------------------------	----

Assura® Irrigation Products

Assura® Irrigation Set	41
Assura® Irrigation Accessories	41
Assura® Irrigation Sleeves	41

Coloplast Fistula & Wound Management System

Pouch	42
Flexible Lids	42
Accessories	42

CONTINENCE

Hydro Intermittent Catheters

SpeediCath® Flex

Hydrophilic Intermittent Catheter	
SpeediCath® Flex	47

SpeediCath® Standard

Hydrophilic Intermittent Catheter	
SpeediCath® Standard Male Straight	47
SpeediCath® Standard Male Coude	47
SpeediCath® Standard Female Straight	47
SpeediCath® Standard Pediatric Straight	47

SpeediCath® Compact Male

Hydrophilic Intermittent Catheter	
SpeediCath® Compact Male	48

SpeediCath® Compact Eve

Hydrophilic Intermittent Catheter	
SpeediCath® Compact Eve	48

SpeediCath® Compact Female

Hydrophilic Intermittent Catheter	
SpeediCath® Compact Female	48
SpeediCath® Compact Female Plus	48

SpeediCath® Compact Set

Hydrophilic Intermittent Catheter	
SpeediCath® Compact Set Male	49
SpeediCath® Compact Set Female	49

Self-Cath®

Intermittent Catheter	
Self-Cath® Female Straight	49
Self-Cath® Pediatric	49
Self-Cath® Extension Tube	49
Self-Cath® Male Straight	50
Self-Cath® Male Soft Straight	50
Self-Cath® Male Straight Curved Packaging	50
Self-Cath® Male Coude Tapered with Guide Stripe	50
Self-Cath® Male Coude Olive with Guide Stripe	50
Self-Cath® Male Closed System Straight	50

External

Conveen® Optima

Male External Catheter - Silicone	
Conveen® Optima Standard	53
Conveen® Optima Short	53

Conveen® Security+

Male External Catheter - SEBS*	
Conveen® Security+ Self Sealing	53
Conveen® Security+ Adhesive Liner	53

Conveen®

Male External Catheter - Latex

Conveen® Self-Sealing	54
Conveen® Adhesive Liner	54

Active Cath®

Male External Catheter - Latex

Active Cath®	54
--------------------	----

Freedom Cath®

Male External Catheter - Latex

Freedom Clear®	54
----------------------	----

Conveen® Active

Urine Leg Bag

Conveen® Active Leg Bag, Non Sterile	55
--	----

Conveen® Security+

Urine Leg bag

Conveen® Security+ Leg Bag, Clamp Outlet, Non Sterile	55
Conveen® Security+ Leg Bag, Clamp Outlet, Sample Port, Sterile	55
Conveen® Security+ Leg Bag, Lever Outlet, Sterile	55

Conveen® Standard

Day / Night Urine Bag

Conveen® Standard Day / Night Bag, Clamp Outlet, Non Sterile	55
Conveen® Standard Day / Night Bag, Clamp Outlet, Sample Port, Sterile	55

Conveen® Contour

Urine Leg Bag

Conveen® Contour Leg Bag, Clamp Outlet, Straps, Non Sterile	56
--	----

Conveen® Night Bag

Urine Night Bag

Conveen® Bedside Night Bag, Non-Sterile	56
Conveen® Bedside Night Bag, Sterile	56

Accessories

Conveen® Bag Hanger	56
Conveen® Leg Bag Strap	56
Conveen® Catheter retaining Strap	56

Peristeen Plus® Anal Irrigation

Peristeen Plus® Anal Irrigation System	57
Peristeen Plus® Anal Irrigation System - Small	57
Peristeen Plus® Anal Irrigation Accessory unit	57
Peristeen Plus® Anal Irrigation Accessory unit - Small	57
Peristeen Plus® Anal Irrigation Rectal Catheter	57
Peristeen Plus® Anal Irrigation Leg Strap	57
Peristeen Plus® Anal Irrigation Tube	57

ADVANCED WOUND DRESSINGS

Biatain® Silicone

Foam Dressing with Silicone Adhesive

Biatain® Silicone	61
-------------------------	----

Biatain® Silicone Non-Border

Biatain® Silicone Non-Border	61
------------------------------------	----

Biatain® Silicone Lite

Thin Foam Dressing with Silicone Adhesive

Biatain® Silicone Lite	61
------------------------------	----

Biatain® Silicone Ag

Foam Dressing with Silicone Adhesive and Silver

Biatain® Silicone Ag	62
----------------------------	----

Biatain®

Foam Dressing

Biatain® Adhesive	62
Biatain® Non-Adhesive	62

Biatain® Soft Hold

Foam Dressing with Gentle Adhesion

Biatain® Soft Hold	63
--------------------------	----

Biatain® Ag

Foam Dressing with Silver

Biatain® Ag Non-Adhesive	63
--------------------------------	----

Biatain® Fiber

Biatain® Fiber	63
----------------------	----

Biatain® Alginate

Calcium Alginate Dressing

Biatain® Alginate	64
-------------------------	----

Biatain® Alginate Ag

Calcium Alginate Dressing with Silver

Biatain® Alginate Ag	64
----------------------------	----

Biatain® Ibu

Foam Dressing with Ibuprofen

Biatain® Ibu Non-Adhesive	64
Biatain® Ibu Soft-Hold	64

Biatain® Super

Hydrocapillary Pad Dressing

Biatain® Super Adhesive	64
Biatain® Super Non-Adhesive	64

Comfeel® Plus

Hydrocolloid Dressing

Comfeel® Plus	65
---------------------	----

Comfeel® Plus Transparent

Hydrocolloid Dressing

Comfeel® Plus Transparent.....65

Triad®

Hydrophilic Wound Dressing

Triad®65

Purilon® Gel

Hydrogel with natural ingredients

Purilon® Gel66

InterDry®

Moisture-wicking textile with antimicrobial silver complex

InterDry®66

Ostomy

Is your patient getting the right fit?

A poor fit to the peristomal body profile is the main cause of leakage and skin complications.^{1,2}



BodyFit Technology® is an elastic adhesive technology unique to Sensura® Mio and is the foundation for achieving the right fit.




**Secure body
contact for
protection**



**Flexible fit
during
movement**



1. Rolstad, B. S. & Erwin-Toth P. L. Peristomal Skin Complications: Prevention and Management. Ostomy Wound Manage. 2004;50(9):68-77.
2. Colwell, J. C. et al. International Consensus Results. Development of Practice Guidelines for Assessment of Peristomal Body and Stoma Profiles, Patient Engagement, and Patient Follow-up. Journal of Wound, Ostomy and Continence Nursing: November/December 2019 - Volume 46 - Issue 6 - p 497-504.



Sensura® Mio

SenSura® Mio Click

Ostomy Barrier and Pouch - Two Piece System

Click mechanical coupling with a lock ring for security.

- All 2pc Mio Click Flanges have 4 belt attachment points.
- Not made with natural rubber latex.

SenSura® Mio Click Skin Barrier					
Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	40mm (1 9/16") ●	10-35mm (3/8 - 1 3/8")	10502	5
		50mm (2") ●	10-45mm (3/8 - 1 3/4")	10512	5
		60mm (2 3/8") ●	10-55mm (3/8 - 2 1/8")	10522	5
		70mm (2 3/4") ●	10-65mm (3/8 - 2 1/2")	10531	5
Pre-Cut	Non-Convex	40mm (1 9/16") ●	20mm (3/4")	10503	5
			22mm (7/8")	10506	5
			25mm (1")	10504	5
			28mm (1 1/8")	10505	5
		50mm (2") ●	30mm (1 1/8")	10513	5
			32mm (1 1/4")	10514	5
			35mm (1 3/8")	10515	5
			38mm (1 1/2")	10516	5
		60mm (2 3/8") ●	40mm (1 9/16")	10523	5
			45mm (1 3/4")	10524	5
			50mm (2")	10525	5
Cut-to-Fit	Convex Light	40mm (1 9/16") ●	10-20mm (3/8 - 3/4")	16901	5
		50mm (2") ●	15-30mm (5/8 - 1 1/8")	16911	5
		60mm (2 3/8") ●	15-40mm (5/8 - 1 9/16")	16921	5
		70mm (2 3/4") ●	15-53mm (5/8 - 2 1/16")	16931	5
Pre-Cut	Convex Light	40mm (1 9/16") ●	18mm (3/4")	16903	5
			21mm (13/16")	16913	5
			25mm (1")	16914	5
		50mm (2") ●	28mm (1 1/8")	16915	5
			31mm (1 1/4")	16923	5
			35mm (1 3/8")	16924	5
		70mm (2 3/4") ●	38mm (1 1/2")	16933	5
			41mm (1 5/8")	16934	5
Cut-to-Fit	Convex Deep	40mm (1 9/16") ●	10-20mm (3/8 - 3/4")	16941	5
		50mm (2") ●	15-30mm (5/8 - 1 1/8")	16951	5
		60mm (2 3/8") ●	15-40mm (5/8 - 1 9/16")	16961	5
		70mm (2 3/4") ●	15-53mm (5/8 - 2 1/16")	16971	5
Pre-Cut	Convex Deep	40mm (1 9/16") ●	18mm (3/4")	16942	5
			21mm (13/16")	16953	5
			25mm (1")	16954	5
		50mm (2") ●	28mm (1 1/8")	16955	5
			31mm (1 1/4")	16963	5
			35mm (1 3/8")	16964	5
		70mm (2 3/4") ●	38mm (1 1/2")	16973	5
			41mm (1 5/8")	16974	5



SenSura® Mio Click Skin Barrier

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Concave	40mm (1 9/16")	10-35mm (3/8 - 1 3/8")	18500	5
		50mm (2")	10-45mm (3/8 - 1 3/4")	18510	5
		60mm (2 3/8")	10-55mm (3/8 - 2 1/8")	18520	5
		70mm (2 3/4")	10-65mm (3/8 - 2 1/2")	18530	5
Pre-Cut	Concave	40mm (1 9/16")	20mm (3/4")	18501	5
			25mm (1")	18511	5
		50mm (2")	30mm (1 1/8")	18512	5
			35mm (1 3/8")	18513	5
		60mm (2 3/8")	30mm (1 1/8")	18521	5
			35mm (1 3/8")	18522	5
			40mm (1 9/16")	18523	5
			45mm (1 3/4")	18524	5
		70mm (2 3/4")	45mm (1 3/4")	18531	5
			50mm (2")	18532	5



SenSura® Mio Click Drainable Pouch with Easyclose™ WIDE-Outlet

Size	Colour	Coupling Size	Product Code	Units
Midi 26cm (10 1/4") 485 mL	Neutral grey with inspection window	40mm (1 9/16")	11461	10
		50mm (2")	11471	10
Maxi 28cm (11") 655 mL	Neutral grey with inspection window	40mm (1 9/16")	11462	10
		50mm (2")	11472	10
		60mm (2 3/8")	11482	10
		70mm (2 3/4")	11486	10
		40mm (1 9/16")	Narrow outlet with Ileo filter 11432	10
		50mm (2")	Narrow outlet with Ileo filter 11443	10
Maxi 28cm (11") 655 mL	Transparent	40mm (1 9/16")	11463	10
		50mm (2")	11473	10
		50mm (2")	No filter 11474	10
		60mm (2 3/8")	11483	10
		60mm (2 3/8")	No filter 11484	10
		70mm (2 3/4")	11487	10



SenSura® Mio Click Closed Pouch

Size	Colour	Coupling Size	Product Code	Units
Midi 18cm (7") 175 mL	Neutral grey with inspection window	40mm (1 9/16")	11401	30
		50mm (2")	11412	30
		60mm (2 3/8")	11421	30
Maxi 22cm (8 1/2") 570 mL	Neutral grey with inspection window	40mm (1 9/16")	11402	30
		50mm (2")	11413	30
		60mm (2 3/8")	11422	30
		70mm (2 3/4")	11425	30



SenSura® Mio Click Urostomy Pouch

Size	Colour	Coupling Size		Product Code	Units
Midi 24cm (9 1/2") 330 mL	Neutral grey	50mm (2")	●	11496	10
		40mm (1 9/16")	●	11491	10
		50mm (2")	●	11497	10
Maxi 26cm (10 1/4") 480 mL	Neutral grey	60mm (2 3/8")	●	11499	10
		40mm (1 9/16")	●	11492	10
	Transparent	50mm (2")	●	11498	10
		60mm (2 3/8")	●	11500	10



SenSura® Mio Click High Output Pouch

Outlet	Colour	Coupling Size	Volume	Product Code	Units
Wide-lumen Hard Outlet 32cm (12 11/16") 1080 mL	Neutral grey with inspection window	40mm (1 9/16")	● 1160 mL	18600	10
		50mm (2")	● 1120 mL	18601	10
		60mm (2 3/8")	● 1080 mL	18602	10
	Transparent	70mm (2 3/4")	● 1040 mL	18603	10
		40mm (1 9/16")	● 1160 mL	18610	10
		50mm (2")	● 1120 mL	18611	10
		60mm (2 3/8")	● 1080 mL	18612	10
		70mm (2 3/4")	● 1040 mL	18613	10
		70mm (2 3/4")	● 1040 mL	18613	10
Mini Pouch Wide-lumen Hard Outlet 24cm (9 1/2") 650 mL	Neutral grey with inspection window	40mm (1 9/16")	● 650 mL	18620	10
		50mm (2")	● 581 mL	18621	10
		60mm (2 3/8")	● 546 mL	18622	10
Soft Outlet Universal Connector 32cm (12 11/16") 1080 mL	Neutral grey with inspection window	50mm (2")	● 1120 mL	18630	10
		60mm (2 3/8")	● 1080 mL	18631	10
		70mm (2 3/4")	● 1040 mL	18632	10
	Transparent	50mm (2")	● 1120 mL	18640	10
		60mm (2 3/8")	● 1080 mL	18641	10
		70mm (2 3/4")	● 1040 mL	18642	10



Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast fistula drain bag #14010

Soft outlet provides effective drainage; connects to a standard drain bag or Coloplast fistula drain bag #14010

The SenSura® Mio range with BodyFit Technology®, makes it easier to provide an appliance with the right fit to body



The SenSura® Mio range with BodyFit Technology® makes it easier to provide an appliance with the right fit to body




BodyFit Technology® is an elastic adhesive technology that is unique to SenSura® Mio. BodyFit Technology® provides secure body contact for protection and a flexible fit during movement.



SenSura® Mio

A flat SenSura® Mio adhesive will usually be adequate for a person with a stoma opening positioned above skin surface. In case of dips, creases and folds in the skin surface, a protective seal can be added to help create an even surface to ensure a secure seal.



SenSura® Mio Convex

For a better fit to an inward peristomal body profile, consider a SenSura® Mio Convex solution:

- Convex Soft, for stomas above skin level that need support due to dips and creases.
- Convex Light, for slightly deep-seated areas, and for stomas with openings in level with the skin, to help the stoma protrude.
- Convex Deep, for when the stoma opening is in or below skin level to make the stoma protrude.

Soft Light Deep



SenSura® Mio Concave

For people with an outward body profile, SenSura® Mio Concave with a curved fit can help obtain a secure seal without creases and folds in the adhesive.

Three more reasons why SenSura® Mio is a great choice for your patients

Discreet pouch colour


The all-textile material has a neutral colour designed to be discreet and unnoticeable, even under white clothing.

A filter you can trust

A unique, full-circle pre-filter provides greater protection against clogging.

Proven to reduce incidents of ballooning² by

61%



A choice of couplings

Click coupling has a lock ring designed for reassurance and a feeling of security. An audible 'click' confirms when the lock ring is securely locking the bag to the baseplate.

Flex coupling is flexible enough to adapt to body contours and comfortable during wear, with a very low profile.



SenSura® Mio is available in many different variants for people with colostomy, ileostomy and urostomy



2. Coloplast, Clinically, 2010 (PAP 001128) for people with ileostomy, urostomy. Statistical Case

SenSura[®] Mio Flex

Ostomy Barrier and Pouch - Two Piece System

Flex adhesive coupling allows flexibility and freedom of movement.

- Only Mio Flex Convex Flanges have 4 belt attachment points.
- Not made with natural rubber latex.

SenSura [®] Mio Flex Skin Barrier					
Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	35mm (1 3/8")	10-33mm (3/8 - 1 1/4")	10551	10
		50mm (2")	10-48mm (3/8 - 2")	10561	10
		70mm (2 3/4")	10-68mm (3/8 - 2 3/4")	10571	10
		90mm (3 1/2")	10-88mm (3/8 - 3 1/2")	18390	10
Pre-Cut	Non-Convex		25mm (1")	10562	10
			28mm (1 1/8")	10566	10
			30mm (1 1/4")	10563	10
		50mm (2")	32mm (1 1/4")	10567	10
			35mm (1 3/8")	10564	10
			38mm (1 1/2")	10568	10
			40mm (1 9/16")	10565	10
		70mm (2 3/4")	45mm (1 3/4")	10572	5
			50mm (2")	10573	5
Cut To Fit	Convex Soft	35mm (1 3/8")	10-20 mm	16880 <i>New</i>	5
Pre-Cut	Convex Soft	35mm (1 3/8")	18mm	16881 <i>New</i>	5
Cut To Fit	Convex Soft	50mm (2")	15-30 mm	16882 <i>New</i>	5
Pre-Cut	Convex Soft		21 mm	16883 <i>New</i>	5
		50mm (2")	25 mm	16884 <i>New</i>	5
			28 mm	16885 <i>New</i>	5
Cut To Fit	Convex Soft	50mm (2")	15-40 mm	16886 <i>New</i>	5
Pre-Cut	Convex Soft		31 mm	16887 <i>New</i>	5
		50mm (2")	35 mm	16888 <i>New</i>	5
Cut To Fit	Convex Soft	70mm (2 3/4")	15-53 mm	16889 <i>New</i>	5
Pre-Cut	Convex Soft		38 mm	16890 <i>New</i>	5
		70mm (2 3/4")	41 mm	16891 <i>New</i>	5
Cut-to-Fit	Convex Light	35mm (1 3/8")	10-20mm (3/8 - 3/4")	16471	5
			15-30mm (5/8 - 1 1/8")	16481	5
		50mm (2")	15-40mm (5/8 - 1 9/16")	16483	5
		70mm (2 3/4")	15-53mm (5/8 - 2 1/16")	16491	5
Pre-Cut	Convex Light	35mm (1 3/8")	18mm (3/4")	16472	5
			21mm (13/16")	16484	5
			25mm (1")	16485	5
		50mm (2")	28mm (1 1/8")	16486	5
			31mm (1 1/4")	16487	5
			35mm (1 3/8")	16488	5
			38mm (1 1/2")	16493	5
		70mm (2 3/4")	41mm (1 5/8")	16494	5



SenSura® Mio Flex Drainable Pouch with Easyclose™ WIDE-Outlet

Size	Colour	Coupling Size		Product Code	Units
Midi 26cm (10 1/4") 485 mL	Neutral grey with inspection window	50mm (2")	●	12271	10
		35mm (1 3/8")	●	Narrow outlet with Ileo filter 12232	10
		35mm (1 3/8")	●	12262	10
	Neutral grey with inspection window	50mm (2")	●	Narrow outlet with Ileo filter 12243	10
		50mm (2")	●	12272	10
		70mm (2 3/4")	●	12282	10
Maxi 28cm (11") 655 mL		50mm (2")	●	No filter 12274	10
	Transparent	70mm (2 3/4")	●	No filter 12284	10
		90mm (3 1/2")	●	18389	10



SenSura® Mio Flex Closed Pouch

Size	Colour	Coupling Size		Product Code	Units
Midi 18cm (7") 275 mL	Neutral grey with inspection window	35mm (1 3/8")	●	12201	30
		50mm (2")	●	12212	30
		35mm (1 3/8")	●	12202	30
Maxi 22cm (8 1/2") 570 mL	Neutral grey with inspection window	50mm (2")	●	12213	30
		70mm (2 3/4")	●	12222	30
		90mm (3 1/2")	●	18386	10



SenSura® Mio Flex Urostomy Pouch

Size	Colour	Coupling Size		Product Code	Units
Maxi 26cm (10 1/4") 400mL	Neutral grey	35mm (1 3/8")	●	12291	10
		50mm (2")	●	12297	10
	Transparent	35mm (1 3/8")	●	12292	10
		50mm (2")	●	12298	10



SenSura® Mio Flex High Output Pouch

Outlet	Colour	Coupling Size	Volume	Product Code	Units
	Neutral grey with inspection window	50mm (2")	● 980 mL	18650	10
		70mm (2 3/4")	● 906 mL	18651	10
Wide-lumen Hard Outlet 32cm (12 11/16") 1080 mL	Transparent	50mm (2")	● 980 mL	18655	10
		70mm (2 3/4")	● 906 mL	18656	10
Soft outlet 32 cm (12 11/16") 1080 mL	Neutral grey with inspection window	50mm (2")	● 1080 mL	18652	10
	Transparent	50mm (2")	● 1080 mL	18657	10



Pouches with Wide Lumen Hard Outlet will connect to Coloplast drainage bag 14010

SenSura[®] Mio

Ostomy Barrier and Pouch - One Piece System

Elastic adhesive, neutral grey textile and full-circle filter.

- Only 1pc Convex systems have 4 belt attachment points.
- Not made with natural rubber latex.

SenSura[®] Mio Drainable Pouch with Easyclose™ WIDE-Outlet

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 26cm (10 1/4") 485 mL	Cut-to-Fit	Non-Convex	Neutral grey with inspection window	10-45mm (3/8 - 1 3/4")	10461	10
				10-55mm (3/8 - 2 1/8") <i>Narrow outlet with Ileo filter</i>	10451	10
	Cut-to-Fit	Non-Convex	Neutral grey with inspection window	10-55mm (3/8 - 2 1/8")	10471	10
				10-55mm (3/8 - 2 1/8")	10481	10
				10-55mm (3/8 - 2 1/8")	10489 <i>No filter</i>	10
Maxi 28cm (11") 655 mL	Pre-Cut	Non-Convex	Neutral grey with inspection window	25mm (1")	10472	10
				28mm (1 1/8")	10473	10
				30mm (1 1/8")	10474	10
				32mm (1 1/4")	10475	10
				35mm (1 3/8")	10476	10
				38mm (1 1/2")	10480	10
				40mm (1 9/16")	10479	10
				45mm (1 3/4")	10477	10
				50mm (2")	10478	10



SenSura[®] Mio 1 pc, Closed, Max Cuts with Circle Filter

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi	Cut To Fit	Non Convex	Opaque	10-70 mm	18371 <i>New</i>	30
				10-100 mm	18373 <i>New</i>	30



SenSura[®] Mio 1 pc, Open, Max Cuts with Circle Filter, Wide-Outlet

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi	Cut To Fit	Non Convex	Transparent	10-70 mm	18382 <i>New</i>	10
				10-100 mm	18383 <i>New</i>	10



SenSura® Mio Drainable Pouch with Easyclose™ WIDE-Outlet

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units		
Maxi 28cm (11") 655 mL	Pre-Cut	Non-Convex	Transparent	25mm (1")	10482	10		
				30mm (1 1⁄8")	10483	10		
				35mm (1 3⁄8")	10484	10		
				40mm (1 9⁄16")	10487	10		
Maxi 28cm (11") 655 mL	Cut-to-Fit		Neutral grey with inspection window	15-33mm (5⁄8 - 1 1⁄4")	16705	10		
				10-50mm (3⁄8 - 2")	16706	10		
			Transparent	15-33mm (5⁄8 - 1 1⁄4")	16715	10		
				10-50mm (3⁄8 - 2")	16716	10		
			Pre-Cut		Neutral grey with inspection window	21mm (7⁄8")	13671	10
						25 mm (1")	13672	10
	28 mm (1 1⁄8")	13673				10		
	31mm (1 1⁄4")	13674				10		
	35 mm (1 3⁄8")	13675				10		
	38 mm (1 1⁄2")	13676				10		
	Cut-to-Fit		Neutral grey with inspection window	10-23mm (3⁄8 - 7⁄8")	16725	10		
				15-33mm (5⁄8 - 1 1⁄4")	16726	10		
				10-43mm (3⁄8 - 1 3⁄4")	16733	10		
			Transparent	10-23mm (3⁄8 - 7⁄8")	16737	10		
				15-33mm (5⁄8 - 1 1⁄4")	16738	10		
				10-43mm (3⁄8 - 1 3⁄4")	16739	10		
				10-43mm (3⁄8 - 1 3⁄4")	16748	No filter 10		
			Pre-Cut		Neutral grey with inspection window	21mm (13⁄16")	16730	10
						25mm (1")	16731	10
						28mm (1 1⁄8")	16732	10
	31mm (1 1⁄4")	16727				10		
	35mm (1 3⁄8")	16734				10		
	38mm (1 1⁄2")	16735				10		
	Cut-to-Fit		Neutral grey with inspection window	10-23mm (3⁄8 - 7⁄8")	16755	10		
				15-33mm (5⁄8 - 1 1⁄4")	16756	10		
				10-43mm (3⁄8 - 1 3⁄4")	16768	10		
			Transparent	10-43mm (3⁄8 - 1 3⁄4")	16767	10		
				10-43mm (3⁄8 - 1 3⁄4")	16778	No filter 10		
Pre-Cut		Neutral grey with inspection window	21mm (7⁄8")	16760	10			
			25mm (1")	16761	10			
			28mm (1 1⁄8")	16762	10			
			31mm (1 1⁄4")	16757	10			
			35mm (1 3⁄8")	16763	10			
			38mm (1 1⁄2")	16764	10			
		Transparent	21mm (13⁄16")	16770	10			
			25mm (1")	16771	10			
			28mm (1 1⁄8")	16772	10			
			31mm (1 1⁄4")	16773	10			



SenSura® Mio Drainable Pouch with Easiclese™ WIDE-Outlet

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 26cm (10 1/4") 485 mL	Cut-to-fit	Concave 	Neutral grey with inspection window	10-40mm (3/8 - 1 9/16")	18300	10
				10-50mm (3/8 - 2")	18310	10
			Neutral grey with inspection window	10-50mm (3/8 - 2")	18312	10
				10-60mm (3/8 - 2 3/8")	18330	10
				10-40mm (3/8 - 1 9/16")	18305	10
		Concave 	Transparent	10-50mm (3/8 - 2")	18321	10
				10-60mm (3/8 - 2 3/8")	18340 <i>No filter</i>	10
			Transparent	10-60mm (3/8 - 2 3/8")	18341	10
				25-40mm (1 - 1 9/16")	18301	10
				30-40mm (1 1/8 - 1 9/16")	18302	10
Maxi 28cm (11") 655 mL	Pre-Cut	Concave 	Neutral grey with inspection window	25-50mm (1 - 2")	18313	10
				30-50mm (1 1/8 - 2")	18314	10
				35-50mm (1 3/8 - 2")	18315	10
				40-60mm (1 9/16 - 2 3/8")	18331	10
				45-60mm (1 3/4 - 2 3/8")	18332	10
				50-60mm (2" - 2 3/8")	18333	10
			Transparent	25-50mm (1 - 2")	18322	10
				30-50mm (1 1/8 - 2")	18323	10
				35-50mm (1 3/8 - 2")	18324	10
				40-60mm (1 9/16 - 2 3/8")	18342	10
		Concave 	Transparent	45-60mm (1 3/4 - 2 3/8")	18343	10
				50-60mm (2" - 2 3/8")	18344	10



SenSura® Mio Closed Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Mini 14cm (5 1/2") 110 mL	Cut-to-fit	Non-Convex	Neutral grey with inspection window	15-45mm (5/8 - 1 3/4")	10870	30
	Cut-to-fit	Non-Convex	Neutral grey with inspection window	15-45mm (5/8 - 1 3/4")	10871	30
	Pre-Cut	Non-Convex	Neutral grey with inspection window	25mm (1")	10872	30
				30mm (1 1/8")	10873	30
35mm (1 3/8")				10874	30	
Midi 18cm (7") 340 mL	Cut-to-Fit	Convex Soft	Neutral grey with inspection window	15-33mm (5/8 - 1 1/4")	16301	10
		Convex Light	Neutral grey with inspection window	15-33mm (5/8 - 1 1/4")	16322	10
	Pre-Cut	Convex Light	Neutral grey with inspection window	31mm (1 1/4")	16323	10
		Cut-to-Fit	Non-Convex	Neutral grey with inspection window	15-45mm (5/8 - 1 3/4")	10881
Pre-Cut		Non-Convex	Neutral grey with inspection window	40-55mm (1 9/16 - 2 1/8")	10885	30
				30mm (1 1/8")	10883	30
				35mm (1 3/8")	10884	30
				15-33mm (5/8 - 1 1/4")	16305	10
Cut-to-Fit		Convex Soft	Neutral grey with inspection window	10-50mm (3/8 - 2")	16306	10
			Transparent	10-50mm (3/8 - 2")	16315	10
			25mm (1")	16310	10	
Pre-Cut		Convex Soft	Neutral grey with inspection window	28mm (1 1/8")	16311	10
				31mm (1 1/4")	16312	10
				35mm (1 3/8")	16313	10
				38 mm (1 1/2")	13684	10
				41 mm (1 5/8")	13685	10
				10-23mm (3/8 - 7/8")	16325	10
				15-33mm (5/8 - 1 1/4")	16326	10
Cut-to-Fit		Convex Light	10-43mm (3/8 - 1 3/4")	16332	10	
			Transparent	10-43mm (3/8 - 1 3/4")	16335	10
			25mm (1")	16330	10	
			28mm (1 1/8")	16331	10	
Pre-Cut		Convex Light	Neutral grey with inspection window	31mm (1 1/4")	16327	10
	35mm (1 3/8")			16333	10	
	41 mm (1 5/8")			13683	10	
	10-23mm (3/8 - 7/8")			16345	10	
Cut-to-Fit	Convex Deep	Neutral grey with inspection window	15-33mm (5/8 - 1 1/4")	16346	10	
			10-43mm (3/8 - 1 3/4")	16352	10	
			25mm (1")	16350	10	
Pre-Cut	Convex Deep	Neutral grey with inspection window	28mm (1 1/8")	16351	10	
			31mm (1 1/4")	16347	10	
			35mm (1 3/8")	16353	10	



SenSura® Mio Closed Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 18cm (7") 340 mL	Cut-to-Fit		Neutral grey with inspection window	10-40mm ($\frac{3}{8}$ - 1 $\frac{9}{16}$ "	18151	10
				10-50mm ($\frac{3}{8}$ - 2"	18160	10
	Pre-Cut			25mm (1"	18161	10
				30mm (1 $\frac{1}{8}$ "	18162	10
				35mm (1 $\frac{3}{8}$ "	18163	10
Maxi 22cm (8 $\frac{1}{2}$ ") 550 mL	Cut-to-Fit		Neutral grey with inspection window	10-40mm ($\frac{3}{8}$ - 1 $\frac{9}{16}$ "	18152	10
				10-50mm ($\frac{3}{8}$ - 2"	18170	10
			Transparent	10-60mm ($\frac{3}{8}$ - 2 $\frac{3}{8}$ "	18180	10
	Pre-Cut			10-60mm ($\frac{3}{8}$ - 2 $\frac{3}{8}$ "	18185	10
			25mm (1"	18171	10	
			30mm (1 $\frac{1}{8}$ "	18172	10	
			35mm (1 $\frac{3}{8}$ "	18173	10	
			40mm (1 $\frac{9}{16}$ "	18181	10	
			45mm (1 $\frac{3}{4}$ "	18182	10	



SenSura® Mio One Piece Urostomy

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 26cm ($10\frac{1}{4}$ ") 400 mL	Cut-to-Fit	Non-Convex	Neutral grey	10-45mm ($\frac{3}{8}$ - $1\frac{3}{4}$ "	10582	10
			Transparent	10-45mm ($\frac{3}{8}$ - $1\frac{3}{4}$ "	10586	10
				20mm ($\frac{3}{4}$ "	10583	10
	Pre-Cut	Non-Convex	Neutral grey	25mm (1")	10584	10
				30mm ($1\frac{1}{8}$ "	10585	10
Midi 24cm ($9\frac{1}{2}$ ") 330 mL	Cut-to-Fit	Convex Soft 	Neutral grey	10-33mm ($\frac{3}{8}$ - $1\frac{5}{16}$ "	16801	10
Maxi 26cm ($10\frac{1}{4}$ ") Cut-to-Fit 400 mL	Cut-to-Fit	Convex Soft 	Neutral grey	15-33mm ($\frac{5}{8}$ - $1\frac{1}{4}$ "	16805	10
				25 mm (1")	13680	10
			Transparent	10-50 mm ($\frac{3}{8}$ " - 2")	13678	10
				10-33 mm ($\frac{3}{8}$ " - $1\frac{5}{16}$ "	13679	10
				10-50mm ($\frac{3}{8}$ - 2")	16810	10
Maxi 26cm ($10\frac{1}{4}$ ") Cut-to-Fit 400 mL	Cut-to-Fit	Convex Light 	Neutral grey	10-23mm ($\frac{3}{8}$ - $\frac{7}{8}$ "	16825	10
				15-33mm ($\frac{5}{8}$ - $1\frac{1}{4}$ "	16826	10
			Transparent	10-23mm ($\frac{3}{8}$ - $\frac{7}{8}$ "	16835	10
				15-33mm ($\frac{5}{8}$ - $1\frac{1}{4}$ "	16836	10
				10-43mm ($\frac{3}{8}$ - $1\frac{3}{4}$ "	16837	10
			Neutral grey	31mm ($1\frac{1}{4}$ "	16827	10
	Pre-Cut	Convex Light 	Transparent	21mm ($\frac{7}{8}$ "	16840	10
				25mm (1")	16841	10
				28mm ($1\frac{1}{8}$ "	16842	10
Maxi 26cm ($10\frac{1}{4}$ ") 400 mL	Cut-to-Fit	Convex Deep 	Neutral grey	10-23mm ($\frac{3}{8}$ - $\frac{7}{8}$ "	16855	10
				15-33mm ($\frac{5}{8}$ - $1\frac{1}{4}$ "	16856	10
			Transparent	10-23mm ($\frac{3}{8}$ - $\frac{7}{8}$ "	16865	10
				15-33mm ($\frac{5}{8}$ - $1\frac{1}{4}$ "	16866	10
	Pre-Cut	Convex Deep 	Transparent	10-43mm ($\frac{3}{8}$ - $1\frac{3}{4}$ "	16867	10
				31mm ($1\frac{1}{4}$ "	16857	10
				21mm ($\frac{7}{8}$ "	16870	10
			Transparent	25mm (1")	16871	10
				28mm ($1\frac{1}{8}$ "	16872	10



SenSura® Mio One Piece Urostomy

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 26cm (10 1/4") 400 mL	Cut-to-Fit	Concave 	Neutral grey	10-40mm (3/8 - 1 9/16")	18401	10
			Transparent	10-40mm (3/8 - 1 9/16")	18405	10
	Pre-Cut	Concave 		20mm (3/4")	18402	10
			Neutral grey	25mm (1")	18403	10
				30mm (1 1/8")	18404	10



SenSura® Mio One Piece High Output Pouch

Cut/Outlet	Convexity	Colour	Volume	Stoma Size	Product Code	Units
Cut-to-fit/ Soft Outlet Universal Connector 28cm (11 3/16") 630 mL	Non-Convex	Neutral grey with inspection window	630 mL	10-70mm (3/8 - 2 3/4")	18660	10
			725 mL	10-100mm (3/8 - 4")	18661	10
		Transparent	630 mL	10-70mm (3/8 - 2 3/4")	18665	10
			725 mL	10-100mm (3/8 - 4")	18666	10
	Convex Soft 	Neutral grey with inspection window	630 mL	10-50mm (3/8 - 2")	18670	10
		Transparent	630 mL	10-50mm (3/8 - 2")	18671	10
	Convex Light 	Neutral grey with inspection window	630 mL	10-43mm (3/8 - 1 3/4")	18675	10
		Transparent	630 mL	10-43mm (3/8 - 1 3/4")	18676	10



Soft outlet can connect to Coloplast's fistula drain bag #14010, or any standard night bag

SenSura® Mio One Piece Post-Op Pouch

Size	Cut	Convexity	Volume	Stoma Size	Product Code	Units
28cm (10 15/18")	Cut-to-Fit	Transparent	680 mL	10-70mm (3/8 - 2 3/4")	18680	Sterile 5
29cm (11 1/2")	Cut-to-Fit	Transparent	725 mL	10-100mm (3/8 - 4")	18681	Sterile 5
28cm (10 15/18")	Cut-to-Fit	Transparent	680 mL	10-70mm (3/8 - 2 3/4")	18690	5
29cm (11 1/2")	Cut-to-Fit	Transparent	725 mL	10-100mm (3/8 - 4")	18691	5



Soft outlet can connect to Coloplast's fistula drain bag #14010, or any standard night bag

SenSura[®] Mio Baby

Ostomy Barrier and Pouch - Two Piece System

SenSura[®] Mio Baby Flex Skin Barrier

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	35mm (1 3/8")	0-35mm	18705	5



SenSura[®] Mio Baby Flex Drainable Pouch

Size	Colour	Coupling Size	Product Code	Units
20cm (7 3/4")	Transparent	35mm (1 3/8")	18700	10



Brava[®] Baby Adhesive Remover

	Product Code	Units
Single use	12014	30



SenSura[®] Mio Kids

Ostomy Barrier and Pouch - Two Piece System

SenSura[®] Mio Kids Flex Skin Barrier

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	27 mm	0-27 mm	18735	5
		35 mm	0-35 mm	18747	5



SenSura[®] Mio Kids Flex Drainable Pouch

Cut	Colour	Coupling Size	Product Code	Units
Mini 16cm (6 1/8")	Neutral grey	27 mm	18720	10
	Transparent	27 mm	18721	10
Midi 17cm (6 3/4")	Neutral grey	27 mm	18722	10
	Transparent	27 mm	18723	10
	Neutral grey	35 mm	18724	10
	Transparent	35 mm	18725	10



SenSura[®] Mio Kids Flex Urostomy Pouch

Cut	Colour	Coupling Size	Stoma Size	Product Code	Units
Midi 17cm (6 3/4")	Neutral grey	27 mm		18730	10
		35 mm		18731	10



Ostomy Barrier and Pouch - One Piece System

SenSura® Mio Kids Drainable Pouch 1 Piece

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 17cm (6 ³ / ₄ ")	Cut-to-fit	Non-Convex	Neutral grey	0-35mm	18710	10
			Transparent	0-35mm	18711	10



SenSura® Mio Kids Pouch 1 Piece Urostomy

Cut	Convexity	Colour	Stoma Size	Product Code	Units
Cut-to-fit 19cm (7 ¹ / ₄ ")	Non-Convex	Transparent	10-35mm	18712	10
		Neutral grey	10-35mm	18713	10



SenSura® Click

Ostomy Barrier and Pouch- Two Piece System

Click mechanical coupling with a lock ring for security. Provides a tight, secure fit to individual body shapes. All 2pc SenSura Click Flanges have 2 belt attachment points.

SenSura® Click Skin Barrier

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	40mm (1 ⁹ / ₁₆ ")	10-35mm (³ / ₈ - 1 ³ / ₈ ")	10011	5
		50mm (2")	10-45mm (³ / ₈ - 1 ³ / ₄ ")	10021	5
		60mm (2 ³ / ₈ ")	10-55mm (³ / ₈ - 2 ¹ / ₈ ")	10031	5
		70mm (2 ³ / ₄ ")	10-65mm (³ / ₈ - 2 ¹ / ₂ ")	10041	5
	Convex Light	40mm (1 ⁹ / ₁₆ ")	15-23mm (⁵ / ₈ - ⁷ / ₈ ")	11011	5
		50mm (2")	15-33mm (⁵ / ₈ - 1 ¹ / ₄ ")	11021	5
		60mm (2 ³ / ₈ ")	15-43mm (⁵ / ₈ - 1 ³ / ₄ ")	11031	5
		70mm (2 ³ / ₄ ")	15-53mm (⁵ / ₈ - 2 ¹ / ₁₆ ")	11041	5
Pre-Cut	Non-Convex	40mm (1 ⁹ / ₁₆ ")	20mm (³ / ₄ ")	10012	5
			25mm (1")	10013	5
		50mm (2")	30mm (1 ¹ / ₈ ")	10022	5
			35mm (1 ³ / ₈ ")	10023	5
		60mm (2 ³ / ₈ ")	40mm (1 ⁹ / ₁₆ ")	10032	5
			45mm (1 ¹ / ₂ ")	10033	5
	Convex Light	40mm (1 ⁹ / ₁₆ ")	21mm (⁷ / ₈ ")	11012	5
			25mm (1")	11022	5
		50mm (2")	28mm (1 ¹ / ₈ ")	11023	5
			31mm (1 ¹ / ₄ ")	11024	5
		60mm (2 ³ / ₈ ")	35mm (1 ³ / ₈ ")	11032	5
			40mm (1 ¹ / ₂ ")	11033	5



SenSura® Click Xpro Skin Barrier

(Engineered to offer extra protection for liquid ileostomies and urostomies)

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	40mm (1 9/16")	10-35mm (3/8 - 1 3/8")	10015	5
		50mm (2")	10-45mm (3/8 - 1 3/4")	10025	5
		60mm (2 3/8")	10-55mm (3/8 - 2 1/8")	10035	5
		70mm (2 3/4")	10-65mm (3/8 - 2 1/2")	10045	5
	Convex Light	40mm (1 9/16")	15-23mm (5/8 - 7/8")	11015	5
		50mm (2")	15-33mm (5/8 - 1 1/4")	11025	5
		60mm (2 3/8")	15-43mm (5/8 - 1 3/4")	11035	5
Pre-Cut	Non-Convex	70mm (2 3/4")	15-53mm (5/8 - 2 1/16")	11045	5
		40mm (1 9/16")	20mm (3/4")	10016	5
			25mm (1")	10017	5
		50mm (2")	30mm (1 1/8")	10026	5
			35mm (1 3/8")	10027	5
	Convex Light	60mm (2 3/8")	40mm (1 9/16")	10036	5
		40mm (1 9/16")	21mm (7/8")	11016	5
			25mm (1")	11026	5
		50mm (2")	28mm (1 1/8")	11027	5
			31mm (1 1/4")	11028	5
		60mm (2 3/8")	35mm (1 3/8")	11036	5



SenSura® Click Drainable Pouch with Easyclose™ WIDE-Outlet

Pouch Size	Pouch Colour	Coupling Size		Product Code	Units
Midi 27cm (10 1/2")	Opaque	40mm (1 9/16")		11114	10
		50mm (2")		11115	10
Maxi 30cm (12")	Opaque	40mm (1 9/16")		11124	10
		50mm (2")		11125	10
		50mm (2")	 No filter	11135	10
		60mm (2 3/8")		11126	10
		60mm (2 3/8")	 No filter	11136	10
		70mm -2 3/4"		11127	10
	Transparent	40mm (1 9/16")		11194	10
		50mm (2")		11195	10
		50mm (2")	 No filter	11185	10
		60mm (2 3/8")		11196	10
60mm (2 3/8")		 No filter	11186	10	
70mm -2 3/4"		11197	10		



SenSura® Click Closed Pouch

Pouch Size	Pouch Colour	Coupling Size	Product Code	Units
Midi 18cm (7")	Opaque	40mm (1 9/16")	10154	30
		50mm (2")	10155	30
		60mm (2 3/8")	10156	30
Maxi 22cm (8 1/2")	Opaque	40mm (1 9/16")	10164	30
		50mm (2")	10165	30
		60mm (2 3/8")	10166	30
		70mm (2 3/4")	10167	30



SenSura® Click Urostomy Pouch

Pouch Size	Pouch Colour	Coupling Size	Product Code	Units
Midi 24cm (9 1/4")	Opaque	40mm (1 9/16") ●	11841	10
		50mm (2") ●	11842	10
	Transparent	40mm (1 9/16") ●	11851	10
		50mm (2") ●	11852	10
Maxi 26cm (10 1/4")	Opaque	40mm (1 9/16") ●	11844	10
		50mm (2") ●	11845	10
		60mm (2 3/8") ●	11846	10
	Transparent	40mm (1 9/16") ●	11854	10
		50mm (2") ●	11855	10
		60mm (2 3/8") ●	11856	10



SenSura® High Output Pouch

Colour	Coupling Size	Volume	Product Code	Units
Opaque	40mm (1 9/16") ●	1160 mL	19040	10
	50mm (2") ●	1120 mL	19041	10
	60mm (2 3/8") ●	1080 mL	19042	10
Transparent	40mm (1 9/16") ●	1160 mL	19050	10
	50mm (2") ●	1120 mL	19051	10
	60mm (2 3/8") ●	1080 mL	19052	10



Will connect to
Coloplast Drainage
Bag 14010

SenSura® Flex

Ostomy Barrier and Pouch - Two Piece System

Low profile adhesive coupling system that provide flexibility and freedom of movement. Double layer adhesive for secure, safe, long lasting wear. Drainable pouch with velcro closure system. Only SenSura Flex Convex Flanges have 2 belt attachment points.

SenSura® Flex Skin Barrier

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	35mm (1 3/8") ●	10-33mm (3/8 - 1 1/4")	10101	10
		50mm (2") ●	10-48mm (3/8 - 2")	10102	10
		70mm (2 3/4") ●	10-68mm (3/8 - 2 3/4")	10103	10
		90mm - 3 1/2" ●	10-88mm (3/8 - 3 1/2")	10108	5
Pre-Cut	Convex light	35mm (1 3/8") ●	15-23mm (5/8 - 7/8")	11301	5
		50mm (2") ●	15-33mm (5/8 - 1 1/4")	11302	5
			15-43mm (5/8 - 1 3/4")	11303	5
		70mm (2 3/4") ●	15-56mm (5/8 - 2 1/8")	11307	5
	Non-Convex	35mm (1 3/8") ●	21mm (7/8")	10111	10
			28mm (1 1/8")	10113	10
		50mm (2") ●	25mm (1")	10114	10
			31mm (1 1/4")	10116	10
			35mm (1 3/8")	10117	10
			21mm (7/8")	11312	5
			25mm (1")	11313	5
			28mm (1 1/8")	11314	5
			31mm (1 1/4")	11315	5
			35mm (1 3/8")	11316	5



SenSura® Flex Xpro Skin Barrier

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	35mm (1 3/8")	10-33mm (3/8 - 1 1/4")	10104	10
		50mm (2")	10-48mm (3/8 - 2")	10105	10
		70mm (2 3/4")	10-68mm (3/8 - 2 3/4")	10106	10
	Convex light	35mm (1 3/8")	15-23mm (5/8 - 7/8")	11304	5
		50mm (2")	15-33mm (5/8 - 1 1/4")	11305	5
		70mm (2 3/4")	15-43mm (5/8 - 1 3/4")	11306	5
			15-56mm (5/8 - 2 1/8")	11308	5



SenSura® Flex Drainable Pouch with Easyclose™ WIDE-Outlet

Size	Colour	Coupling Size			Product Code	Units
Midi 27cm (10 1/2")	Opaque	50mm (2")			11505	20
Maxi 30cm (12")	Opaque	35mm (1 3/8")			11514	20
		50mm (2")			11515	20
		50mm (2")		No filter	11615	20
		70mm (2 3/4")			11518	20
	Transparent	50mm (2")		No filter	11612	20
		90mm (3 1/2")		No filter	11690	20



SenSura® Flex Closed Pouch

Size	Colour	Coupling Size	Product Code	Units
Midi 18cm (7")	Opaque	35mm (1 3/8")	10904	30
		50mm (2")	10905	30
Maxi 22cm (8 1/2")	Opaque	35mm (1 3/8")	10914	30
		50mm (2")	10915	30
		70mm (2 3/4")	10918	30



SenSura® Flex Urostomy Pouch

Size	Colour	Coupling Size	Product Code	Units
Midi 24cm (9 1/2")	Opaque	35mm (1 3/8")	11821	20
		50mm (2")	11822	20
Maxi 26cm (10 1/4")	Opaque	35mm (1 3/8")	11824	20
		50mm (2")	11825	20
	Transparent	35mm (1 3/8")	11827	20
		50mm (2")	11828	20



Ostomy Barrier and Pouch - One Piece System

Easiclose™ system with Velcro® to ensure security. Double layer adhesive provides secure adhesion, flexibility and ease of removal. Only SenSura 1pc Convex Systems have 2 belt attachment points.

SenSura® Drainable Pouch with Easiclose™ WIDE-Outlet						
Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 30cm (12")	Cut-to-Fit	Non-Convex	Opaque	10-76mm ($\frac{3}{8}$ -3")	15531	10
			Transparent	10-76mm ($\frac{3}{8}$ -3")	15521	10
				10-76mm ($\frac{3}{8}$ -3")	15821 <i>No filter</i>	10
		Convex Light	Opaque	15-33mm ($\frac{5}{8}$ -1 $\frac{1}{4}$ ")	15607	10
				15-43mm ($\frac{5}{8}$ -1 $\frac{3}{4}$ ")	15608	10
			Transparent	15-33mm ($\frac{5}{8}$ -1 $\frac{1}{4}$ ")	15605	10
				15-43mm ($\frac{5}{8}$ -1 $\frac{3}{4}$ ")	15606	10
	Pre-Cut	Non-Convex	Opaque	25mm (1")	15532	10
				30mm (1 $\frac{1}{8}$ ")	15533	10
				35mm (1 $\frac{3}{8}$ ")	15534	10
		Convex Light	Opaque	25mm (1")	15633	10
				28mm (1 $\frac{1}{8}$ ")	15634	10
				31mm (1 $\frac{1}{4}$ ")	15635	10
				25mm (1")	15522	10
			Transparent	30mm (1 $\frac{1}{8}$ ")	15523	10
				35mm (1 $\frac{3}{8}$ ")	15524	10



SenSura® Xpro Drainable Pouch with Easiclose™ WIDE-Outlet						
Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 30cm (12")	Cut-to-Fit	Non-Convex	Transparent	10-76mm ($\frac{3}{8}$ - 3")	15981	20
		Convex Light	Transparent	15-33mm ($\frac{5}{8}$ - 1 $\frac{1}{4}$ ")	15695	10
				15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	15696	10

SenSura® Closed Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 18cm (7")	Cut-to-Fit	Non-Convex	Opaque	10-66mm ($\frac{3}{8}$ - 2 $\frac{5}{8}$ ")	15450	30
			Transparent	10-66mm ($\frac{3}{8}$ - 2 $\frac{5}{8}$ ")	15441	30
		Convex Light	Opaque	15-33mm ($\frac{5}{8}$ - 1 $\frac{1}{4}$ ")	15403	10
				15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	15404	10
	Pre-Cut	Convex Light	Opaque	25mm (1")	15413	10
				28mm (1 $\frac{1}{8}$ ")	15414	10
				31mm (1 $\frac{1}{4}$ ")	15415	10
	Cut-to-Fit	Non-Convex	Opaque	10-76mm ($\frac{3}{8}$ - 3")	15480	30
				10-76mm ($\frac{3}{8}$ - 3")	15470	30
		Convex Light	Opaque	15-33mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	15407	10
				15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	15408	10
			Transparent	15-33mm ($\frac{5}{8}$ - 1 $\frac{1}{4}$ ")	15405	10
				15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	15406	10
Maxi 22cm (8 $\frac{1}{2}$ ")	Pre-Cut	Non-Convex	Opaque	25mm (1")	15482	30
				30mm (1 $\frac{1}{8}$ ")	15483	30
				35mm (1 $\frac{3}{8}$ ")	15484	30
				25mm (1")	15472	30
			Transparent	30mm (1 $\frac{1}{8}$ ")	15473	30
				35mm (1 $\frac{3}{8}$ ")	15474	30



SenSura® Urostomy Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 24cm (9 $\frac{1}{2}$ ")	Cut-to-Fit	Non-Convex	Opaque	10-66mm ($\frac{3}{8}$ - 2 $\frac{5}{8}$ ")	11801	10
			Transparent	10-66mm ($\frac{3}{8}$ - 2 $\frac{5}{8}$ ")	11803	10
		Non-Convex	Opaque	10-76mm ($\frac{3}{8}$ - 3")	11802	10
			Transparent	10-76mm ($\frac{3}{8}$ - 3")	11804	10
Maxi 26cm (10 $\frac{1}{4}$ ")	Cut-to-Fit	Convex Light	Opaque	15-33mm ($\frac{5}{8}$ - 1 $\frac{1}{4}$ ")	11811	10
				15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	11812	10
			Transparent	15-33mm ($\frac{5}{8}$ - 1 $\frac{1}{4}$ ")	11814	10
				15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	11815	10



SenSura® One-Piece High Output Pouch

Cut	Sterile	Stoma Size	Volume	Product Code	Units
Cut-to-Fit	Sterile	10-76mm ($\frac{3}{8}$ -3")	630 mL	19001	5
	Non-Sterile	10-76mm ($\frac{3}{8}$ -3")	630 mL	19003	6



SenSura® Post-Op Pouch

Cut	Sterile	Stoma Size	Window Size	Volume	Product Code	Units
Cut-to-Fit	Sterile	10-76mm ($\frac{3}{8}$ -3")	70mm	680 mL	19010	5
		10-115mm ($\frac{3}{8}$ - 4 $\frac{1}{2}$ ")	100mm	725 mL	19011	5
	Non-Sterile	10-76mm ($\frac{3}{8}$ -3")	70mm	680 mL	19020	6
		10-115mm ($\frac{3}{8}$ - 4 $\frac{1}{2}$ ")	100mm	725 mL	19021	6
		10-76mm ($\frac{3}{8}$ -3")	No window	680 mL	19030	6



Ostomy Barrier and Pouch - Two Piece System


Click mechanical coupling with a lock ring for security. Skin barrier with spiral technology for secure, long lasting wear. All 2pc Assura Click Flanges have 2 belt attachment points.

Assura® Skin Barrier						
Cut	Convexity	Coupling Size		Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	40mm (1 9/16")	●	10-35mm (3/8 - 1 3/8")	2881	5
		50mm (2")	●	10-45mm (3/8 - 1 3/4")	2882	5
		60mm (2 3/8")	●	10-55mm (3/8 - 2 1/8")	2883	5
	Convex Light	40mm (1 9/16")	●	15-23mm (5/8 - 7/8")	14261	5
		50mm (2")	●	15-33mm (5/8 - 1 1/4")	14262	5
		60mm (2 3/8")	●	15-43mm (5/8 - 1 3/4")	14263	5
	Convex Deep	50mm (2")	●	15-33mm (5/8 - 1 1/4")	12716	5
		60mm (2 3/8")	●	15-43mm (5/8 - 1 3/4")	12719	5
Pre-Cut	Non-Convex	40mm (1 9/16")	●	20mm (3/4")	12841	5
				25mm (1")	12842	5
		50mm (2")	●	30mm (1 1/8")	12844	5
				35mm (1 3/8")	12845	5
		60mm (2 3/8")	●	40mm (1 5/8")	12847	5
		60mm (2 3/8")	●	45mm (1 3/4")	12848	5
	Convex Light	40mm (1 9/16")	●	18mm (3/4")	14271	5
				21mm (7/8")	14272	5
		50mm (2")	●	25mm (1")	14273	5
				28mm (1 1/8")	14274	5
				31mm (1 1/4")	14275	5
		60mm (2 3/8")	●	35mm (1 3/8")	14276	5
	38mm (1 1/2")	14277		5		
	Convex Deep	40mm (1 9/16")	●	41mm (1 5/8")	14278	5
				15mm (5/8")	12701	5
				18mm (3/4")	12702	5
		50mm (2")	●	21mm (7/8")	12703	5
25mm (1")				12704	5	
28mm (1 1/8")				12705	5	
60mm (2 3/8")		●	31mm (1 1/4")	12706	5	
			35mm (1 3/8")	12707	5	
			38mm (1 1/2")	12708	5	
			41mm (1 5/8")	12709	5	



Assura® Skin Barrier - Extra Extended Wear

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	40mm (1 9/16")	10-35mm (3/8 - 1 3/8")	2831	5
		50mm (2")	10-45mm (3/8 - 1 3/4")	2832	5
		60mm (2 3/8")	10-55mm (3/8 - 2 1/8")	2833	5
	Convex Light	40mm (1 9/16")	15-23mm (5/8 - 7/8")	14281	5
		50mm (2")	15-33mm (5/8 - 1 1/4")	14282	5
		60mm (2 3/8")	15-43mm (5/8 - 1 3/4")	14283	5
	Convex Deep	40mm - (1 9/16")	15-23mm (5/8 - 7/8")	14243	5
		50mm (2")	15-33mm (5/8 - 1 1/4")	14246	5
		60mm (2 3/8")	15-43mm (5/8 - 1 3/4")	14249	5
Pre-Cut	Non-Convex	40mm (1 9/16")	20mm (3/4")	14251	5
			25mm (1")	14252	5
		50mm (2")	30mm (1 1/8")	14254	5
			35mm (1 3/8")	14255	5
		60mm (2 3/8")	45mm - 1 3/4"	14258	5
	Convex Light	40mm (1 9/16")	18mm (3/4")	14291	5
			21mm (7/8")	14292	5
		50mm (2")	25mm (1")	14293	5
			28mm (1 1/8")	14294	5
		60mm (2 3/8")	31mm (1 1/4")	14295	5
			35mm (1 3/8")	14296	5
	Convex Deep	40mm (1 9/16")	18mm (3/4")	14232	5
			21mm (7/8")	14233	5
		50mm (2")	25mm (1")	14234	5
			28mm (1 1/8")	14235	5
		60mm (2 3/8")	31mm (1 1/4")	14236	5
			35mm (1 3/8")	14237	5
		60mm (2 3/8")	38mm (1 1/2")	14238	5
			41mm (1 5/8")	14239	5



Assura® 2 pc. High Output Drainable Pouch

Size	Colour	Coupling Size	Product Code	Units
Maxi 32cm (12 1/2")	Opaque	50mm (2")	2846	10
		60mm (2 3/8")	2847	10



Assura® 2 pc. High Output Drainable System, Ileo Night Pouch 2L

Size	Colour	Coupling Size	Product Code	Units
2 L	Opaque	50mm (2")	2836	5
		60mm (2 3/8")	2837	5



Assura® Drainable Pouch Clamp Closure

Size	Colour	Coupling Size	Product Code	Units
Maxi 30cm (12 1/2")	Transparent	50mm (2") ●	2792	10
		60mm (2 3/8") ●	2793	10



Assura® Drainable Pouch with EasiClose™

Size	Colour	Coupling Size	Product Code	Units
Mini 23cm (9")	Opaque	40mm (1 9/16") ●	13924	10
		50mm (2") ●	13925	10
Midi 26cm (10 1/4")	Opaque with Graphic	40mm (1 9/16") ●	13964	10
		50mm (2") ●	13965	10
		60mm (2 3/8") ●	13966	10
Maxi 28cm (11")	Opaque	40mm (1 9/16") ●	13984	10
		40mm (1 9/16") ● No filter	14494	10
		50mm (2") ●	13985	10
		60mm (2 3/8") ●	13986	10
	Transparent	60mm (2 3/8") ● No filter	14498	10
		40mm (1 9/16") ●	13974	10
		50mm (2") ●	13975	10
		60mm (2 3/8") ●	13976	10



Assura® Closed Pouch

Size	Colour	Coupling Size	Product Code	Units
Midi 18cm (7")	Opaque	40mm (1 9/16") ●	12354	30
		50mm (2") ●	12355	30
		60mm (2 3/8") ●	12356	30
	Transparent	40mm (1 9/16") ●	12344	30
		60mm (2 3/8") ●	12346	30
Maxi 22cm (8 1/2")	Opaque	40mm (1 9/16") ●	12384	30
		50mm (2") ●	12385	30
		60mm (2 3/8") ●	12386	30
	Transparent	40mm (1 9/16") ●	12374	30
		50mm (2") ●	12375	30
		60mm (2 3/8") ●	12376	30



Assura® Urostomy Pouch

Size	Colour	Coupling Size	Product Code	Units
Maxi 25cm (10")	Opaque	40mm (1 9/16") ●	14224	10
		50mm (2") ●	14225	10
		60mm (2 3/8") ●	14226	10
	Transparent	40mm (1 9/16") ●	14227	10
		50mm (2") ●	14228	10
		60mm (2 3/8") ●	14229	10



Ostomy Barrier and Pouch - One Piece System

Skin Barrier with spiral technology for secure, long lasting wear. Different convexity and pouch shape. All Assura 1pc Convex Systems have 2 belt attachment points.

Assura® Drainable Pouch Clamp Closure							
Size	Cut	Convexity	Colour	Stoma Size	Product Code		Units
Maxi 30cm (12")	Cut-to-Fit	Non-Convex	Opaque	10-70mm ($\frac{3}{8}$ - 2 $\frac{3}{4}$ ")	2520	No filter	10
			Transparent	10-70mm ($\frac{3}{8}$ - 2 $\frac{3}{4}$ ")	2490	No filter	10
		Convex Deep	Opaque	15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	12533	No filter	10
			Transparent	15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	13706	No filter	10



Assura® Drainable Pouch with EasiClose™							
Size	Cut	Convexity	Colour	Stoma Size	Product Code		Units
Midi 26cm (10 1/4")	Cut-to-Fit	Non-Convex	Opaque	10-55mm ($\frac{3}{8}$ - 2 $\frac{1}{8}$ ")	13840		10
		Convex Light	Opaque	15-33mm ($\frac{5}{8}$ - 1 $\frac{1}{4}$ ")	14404		10
				15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	14405		10
		Convex Deep	Opaque	15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	14196		10
	Pre-Cut	Non-Convex	Opaque	25mm (1")	13844		10
				30mm (1 $\frac{1}{8}$ ")	13845		10
				35mm (1 $\frac{3}{8}$ ")	13846		10
Maxi 28cm (11")	Cut-to-Fit	Non-Convex	Opaque	10-70mm ($\frac{3}{8}$ - 2 $\frac{3}{4}$ ")	13870		10
			Transparent	10-70mm ($\frac{3}{8}$ - 2 $\frac{3}{4}$ ")	13860		10
		Convex Light	Opaque	15-33mm ($\frac{5}{8}$ - 1 $\frac{1}{4}$ ")	14414		10
				15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	14415		10
			Transparent	15-33mm ($\frac{5}{8}$ - 1 $\frac{1}{4}$ ")	14411		10
				15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	14412		10



Assura® Drainable Pouch with EasiClose™								
Size	Cut	Convexity	Colour	Stoma Size	Product Code		Units	
Maxi 28cm (11")	Cut-to-Fit	Convex Deep	Opaque	15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	14106		10	
			Transparent	15-43mm ($\frac{5}{8}$ - 1 $\frac{3}{4}$ ")	14103		10	
	Pre-Cut	Non-Convex	Opaque	25mm (1")	13874		10	
				30mm (1 $\frac{1}{8}$ ")	13875		10	
				35mm (1 $\frac{3}{8}$ ")	13876		10	
		Convex Light	Opaque	25mm (1")	14422		10	
				28mm (1 $\frac{1}{8}$ ")	14423		10	
				31mm (1 $\frac{1}{4}$ ")	14424		10	
	Pre-Cut	Convex Deep	Opaque	21mm ($\frac{7}{8}$ ")	14173		10	
				25mm (1")	14174		10	
				28mm (1 $\frac{1}{8}$ ")	14175		10	
			Transparent	31mm (1 $\frac{1}{4}$ ")	14176		10	
				21mm ($\frac{7}{8}$ ")	14163		10	
				25mm (1")	14164		10	
		Pre-Cut	Convex Deep	Transparent	28mm (1 $\frac{1}{8}$ ")	14165		10
					31mm (1 $\frac{1}{4}$ ")	14166		10



Assura® Closed Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Mini 12cm (4 3/4")	Cut-to-Fit	Non-Convex	Opaque	20-55mm (3/8 - 2 1/8")	12110	30
			Opaque	20-55mm (3/8 - 2 1/8")	12140	30
	Cut-to-Fit	Non-Convex	Transparent	20-55mm (3/8 - 2 1/8")	12130	30
				25mm (1")	12144	30
Midi 18cm (7")	Pre-Cut	Non-Convex	Opaque	30mm (1 1/8")	12145	30
				35mm (1 3/8")	12146	30
				40mm (1 9/16")	12147	30
				50mm (2")	12149	30



Assura® 1-Piece Post-Op Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Mini 31cm (12 1/4")	Cut-to-Fit	Non-Convex	Transparent	13-70mm (1/2 - 2 3/4")	12805	5
				13-70mm (1/2 - 2 3/4")	12800	5
Midi 33cm (13")	Cut-to-Fit	Non-Convex	Transparent	13-100mm (1/2 - 4")	12802	5
				13-70mm (1/2 - 2 3/4")	12808 Sterile	5
				13-100mm (1/2 - 4")	12810 Sterile	5



Assura® Closed Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 18cm (7")	Pre-Cut	Convex Light	Opaque	25mm (1")	14452	10
				28mm (1 1/8")	14453	10
				31mm (1 1/4")	14454	10
		Non-Convex	Opaque	20-70mm (3/4 - 2 3/4")	12170	30
Maxi 22cm (8 3/4")	Cut-to-Fit	Convex Light	Opaque	15-43mm (5/8 - 1 3/4")	14445	10
			Transparent	15-33mm (5/8 - 1 1/4")	14441	10
	Pre-Cut	Non-Convex	Opaque	25mm (1")	12174	30
				30mm (1 1/8")	12175	30
				35mm (1 3/8")	12176	30
				40mm (1 9/16")	12177	30



Assura® Urostomy Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 27cm (10 1/2")	Cut-to-Fit	Non-Convex	Opaque	10-55mm (3/8 - 2 1/8")	5580	10
			Transparent	10-55mm (3/8 - 2 1/8")	5585	10
		Convex Deep	Transparent	15-33mm (5/8 - 1 1/4")	12595	10
				15-43mm (5/8 - 1 3/4")	12596	10
	Pre-Cut	Convex Deep	Transparent	15mm (5/8")	12991	10
				18mm (3/4")	12992	10
				21mm (7/8")	12993	10
				25mm (1")	12994	10
				28mm (1 1/8")	12995	10
				31mm (1 1/4")	12996	10
				35mm (1 3/8")	12997	10



Assura® Urostomy Pouch - Extra Extended Wear

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 25cm (10")	Cut-to-Fit	Non-Convex	Opaque	10-55mm (3/8 - 2 1/8")	14201	10
			Transparent	10-55mm (3/8 - 2 1/8")	14202	10
		Convex Light	Opaque	15-43mm (5/8 - 1 3/4")	14701	10
			Transparent	15-43mm (5/8 - 1 3/4")	14702	10
Maxi 27cm (10 1/2")	Cut-to-Fit	Non-Convex	Opaque	10-55mm (3/8 - 2 1/8")	12471	10
			Transparent	10-55mm (3/8 - 2 1/8")	12472	10
		Convex Deep	Transparent	15-33mm (5/8 - 1 1/4")	12432	10
				15-43mm (5/8 - 1 3/4")	12433	10



Coloplast Urostomy Night Bag

Description	Product Code	Units
Includes a fixed connector that fits on to the urostomy pouch. The bag holds up to 2 litres and has a 120cm long flexible tube with anti-reflux valve. Not Compatible with External Catheters or Leg Bags	21365	10



Assura®

Ostomy Barrier and Pouch - Two Piece System - Pediatric

Click mechanical coupling.

Assura® Pediatric Skin Barrier

Cut	Convexity	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	10-35mm (3/8 - 1 3/8")	2181	5



Assura® Pediatric Drainable Pouch

Size	Colour	Product Code	Units
13cm (5")	Opaque	2151	10



Assura® Pediatric Closed Pouch

Size	Colour	Product Code	Units
13cm (5")	Transparent	2166	30

EasiFlex[®]

Ostomy Barrier and Pouch - Two Piece System - Pediatric

Flex adhesive coupling.

Easiflex [®] Pediatric Skin Barrier					
Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	17mm	0-15mm (0-5/8")	14307	10
		27mm	0-25mm (0-1")	14308	10



Easiflex [®] Pediatric Drainable Pouch					
Size	Colour	Coupling Size	Product Code	Units	
105 mL	Teddy Print	17mm	14691	10	
	Transparent	17mm	14681	10	
150 mL	Teddy Print	27mm	14692	10	
	Transparent	27mm	14682	10	

Assura[®]

Ostomy Barrier and Pouch - One Piece System - Pediatric

Assura [®] Pediatric Drainable Pouch						
Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
15cm (6")	Cut-to-Fit	Non-Convex	Transparent	10-35mm (3/8 - 1 3/8")	2115	10



Assura [®] Pediatric Closed Pouch						
Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
13cm (5")	Cut-to-Fit	Non-Convex	Transparent	10-35mm (3/8 - 1 3/8")	2125	30



Assura [®] Pediatric Urostomy Pouch						
Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
15cm (6")	Cut-to-Fit	Non-Convex	Transparent	10-35mm (3/8 - 1 3/8")	8009	10

The Brava® range covers three distinct needs

Fit to body

- *Elastic Barrier Strips* provide extra security by helping prevent baseplate edges from lifting or rolling up.
- *Protective Ring, Protective Ring Convex, Paste & Strip Paste* fill in cavities/folds in the skin around the stoma, providing a better seal to the ostomy baseplate.
- *Belt* holds appliances tight to the skin and ensures that the baseplate stays in place.



Elastic Barrier Strips

Protect the skin

- *Skin Barrier Wipe & Spray* create a thin breathable film on the skin and reduce skin problems associated with body waste and adhesive
- *Adhesive Remover Spray & Wipe* provide a skin friendly way to remove appliances more easily.
- *Protective Ring, Protective Ring Convex, Powder & Protective Sheet* absorb moisture, helping reduce moisture-associated skin damage and irritation
- *Barrier Cream* creates a long lasting barrier effect that protects the skin against body fluids.



Skin Barrier

Adhesive Remover

Facilitate routine

- *Adhesive Remover Spray & Wipe* help facilitate easier removal of appliances.
- *Lubricating Deodorant's* multi-lubricating effect allows easier emptying of the bag and helps reduce odour when changing or emptying an appliance.



Adhesive Remover



Protective Ring



Protective Ring
Convex



Paste



Strip Paste



Belt



Protective Ring



Protective Ring
Convex



Protective Sheet



Powder



Barrier Cream



Lubricating Deodorant

Ostomy Accessories

Brava® Protective Sheets

Size (cms)	Size (inches)	Product Code	Units
10 x 10 cm	4 x 4"	3210	10
15 x 15 cm	6 x 6"	3215	5
20 x 20 cm	8 x 8"	3220	5



Brava® Protective Rings

Size (mm) Internal diameter	Size (mm) External diameter	Size (mm) thickness	Product Code	Units
18 mm	48 mm	2.5 mm	12035	10
18 mm	48 mm	4.2 mm	12045	10
27 mm	57 mm	2.5 mm	12037	10
27 mm	57 mm	4.2 mm	12047	10
34 mm	64 mm	2.5 mm	12039	10
34 mm	64 mm	4.2 mm	12049	10



Brava® Protective Rings Wide

Size (mm) Internal diameter	Size (mm) External diameter	Size (mm) thickness	Product Code	Units
18 mm	57 mm	2.5 mm	12036	10
18 mm	64 mm	2.5 mm	12038	10
18 mm	57 mm	4.2 mm	12046	10
18 mm	64 mm	4.2 mm	12048	10



Brava® Protective Rings Extra Wide

Size (mm) Internal diameter	Size (mm) External diameter	Size (mm) thickness	Product Code	Units
18 mm	76 mm	2.5 mm	12032	10
27 mm	76 mm	2.5 mm	12033	10



Brava® Protective Ring Convex

Size (mm) Internal diameter	Size (mm) External diameter	Size (mm) thickness	Product Code	Units
20 mm	60 mm	8.0 mm	12090	10
25 mm	65 mm	8.0 mm	12091	10
30 mm	70 mm	8.0 mm	12092	10
35 mm	75 mm	8.0 mm	12093	10
40 mm	80 mm	8.0 mm	12094	10
45 mm	85 mm	8.0 mm	12095	10



Brava® Mouldable Rings

Size (mm) Internal diameter	Size (mm) thickness	Product Code	Units
18 mm	2.0 mm	12030	10
18 mm	4.2 mm	12042	10



Brava® Skin Barrier Spray and Wipe

Size	Product Code	Units
50 mL Spray	12020	1
Single unit wipe	12021	30



Brava® Adhesive Remover Spray and Wipe

Size	Product Code	Units
50 mL Spray	12010	1
Single unit wipe	12011	30

**Brava® Barrier Cream**

Size	Product Code	Units
60g Tube	12000	1

**Brava® Ostomy Strip Paste**

Size	Product Code	Units
57g / 2oz	2655	10

**Brava® Elastic Barrier Strips**

Size	Product Code	Units
Brava Elastic Barrier Strip Original (a.)	12070	20
Brava Elastic Barrier Strip Straight (b.)	12074	20
Brava Elastic Barrier Strip Belt (c.)	12072	20
Brava Elastic Barrier Strip XL (d.)	12076	20

**Brava® Ostomy Powder**

Size	Product Code	Units
25g / 1oz bottle	1907	16

**Brava® Mio Belt for All SenSura Mio Products**

New design for easier attachment and improved security

Size	Product Code	Units
Standard - 100 cm - 40"	0423	15
XXL - 155cm - 61"	0424	15

**Brava® Belt for SenSura Original + Assura Products**

Size	Product Code	Units
Standard - 100 cm - 40"	0421	15
XL - 135cm - 54"	0422	15

**Brava® Lubricating Deodorant**

Size	Product Code	Units
7.5 mL sachet	12060	20
240 mL bottle	12061	1

**Brava® Alcohol Free Paste**

Size	Product Code	Units
60 g tube	12050	1

**Ostomy Paste (contains alcohol)**

Size	Product Code	Units
60 g tube	2650	12



Skin Barrier Rings		
Pre-cut Size	Product Code	Units
20mm ($\frac{3}{4}$ "	2320	30
25mm (1"	2325	30
30mm ($1\frac{1}{8}$ "	2330	30
40mm ($1\frac{9}{16}$ "	2340	30
50mm (2"	2350	30



Coloplast Prep™ Protective Liquid Skin Barrier (contains alcohol)		
Size	Product Code	Units
Single Application Wipe	2041	12 boxes/54
59 mL - 2oz. dab bottle	0925	12



Coloplast Filtrodor® Pouch Filter		
	Product Code	Units
Charcoal-activated deodorizing filter with self-adhering foam center, adheres to front of Assura & SenSura Pouches	0509	5 pouches (10 units per pouch)



Coloplast Pouch Clamp		
	Product Code	Units
Simple closure clamp for Assura Drainable Pouches	9500	20



Ostomy Scissors		
	Product Code	Units
Stainless steel medical grade ostomy scissors with curved blades and rounded ends - right handed	9505	30



Hex-On®

Odour antagonist

Eliminates odour with a light, mild scent.

Hex-On®				
Size (ounces)	Size (mL)	Product Code	Units	
2 oz.	59 mL	7583	12	




Ostomy Cap

Designed for discreet situations. 2 piece compatible with SenSura Click and Assura 2 piece skin barriers.

**Note Uro Minicap and Coloplast Stoma cap do not fit SenSura Mio products*

Coloplast Uro Minicap (100 ml)						
Size	Colour	Coupling Size		Product Code	Units	
Mini	Opaque	40mm (1 9/16")	●	2807	30	
		50mm (2")	●	2808	30	

Coloplast Stoma Cap						
Size	Colour	Coupling Size		Product Code	Units	
Mini	Opaque	40mm (1 9/16")	●	2801	30	
		50mm (2")	●	2802	30	
		60mm (2 3/8")	●	2803	30	

Coloplast Stoma Cap - One Piece						
Size	Colour	Cut to Fit - Stoma size		Product Code	Units	
Mini	Opaque	20-55mm (3/4 - 2 1/8")		2501	30	

Mio Irrigation Products

Mio Irrigation One Piece Sleeve							
Cut	Convexity	Color	Description	Stoma Size	Product Code	Units	
Cut-to-Fit	Non-Convex	Transparent	Can be used with the Assura Irrigation Kit, but is non compatible with included 2-piece sleeve/ flanges. Mio Irrigation One Piece Sleeve is meant to be used on its own to support through the irrigation process.	15mm-55mm	11750	10	

Assura® Irrigation Products

Assura® Irrigation Set			
Description	Product Code	Units	
Includes: - 1 Water Bag/Regulator - 1 Cone	1500	1	
Assura® Irrigation Accessories			
Size	Product Code	Units	
Stoma Cone - Soft, flexible stoma cone with tubing	1110	4	
Irrigation Faceplate. For use with ostomy belt 0421 or 0422. Use with #12836 Irrigation Sleeve	12820	1	

The irrigation sleeves below must be used with their corresponding colour coded flange or faceplate. Consult this chart to order the correct transparent sleeves.

Description	Coupling Size	Use with Assura and SenSura Barrier Product Codes (will not work with SenSura Mio)	Product Code	Units	
Irrigation sleeves	40mm (1 ⁹ / ₁₆ ") ●	Will couple to any green 40 mm barriers within the Assura and SenSura lines	12834	5	
	50mm (2") ●	Will couple to any red 50 mm barriers within the Assura and SenSura lines	12835	5	
	60mm (2 ³ / ₈ ") ●	Will couple to any blue 60 mm barriers within the Assura and SenSura lines	12836	5	

Coloplast Fistula & Wound Management System

Pouch				
Size	Cutting Surface	Product Code	Units	
Mini	104-159mm (4 1/8 - 6 1/4")	14050	6	
Midi	156-228mm (6 3/16 - 9")	14060	3	
Maxi	208-297mm (8 1/4 - 11 3/4")	14070		
Flexible Lids				
Size	For	Product Code	Units	
Mini	Mini Fistula & Wound Management System	14021	10	
Midi	Midi Fistula & Wound Management System	14031		
Maxi	Maxi Fistula & Wound Management System	14041		
Accessories				
Size	For	Product Code	Units	
Bed Drainage Bag	All Fistula & Wound Management systems (Capacity 2000 mL)	14010	10	
Drain Port	All Fistula & Wound Management systems	14015	12	

Continence

The truth is:


Not all coatings are the same

To protect the mucosa from risk of damage, complications,* and UTIs, try a coating that stays bonded, stays smooth and stays hydrated.

SpeediCath®




The coating plays an important role in protecting the urethra
But unfortunately, not all hydrophilic catheter coatings are the same.
Some of the common pitfalls can be:



Dry spots

Some have dry spots where the coating doesn't adhere to the catheter during insertion or withdrawal.



Short dry-out time

Insufficient hydration may cause the coating to dry out over time and become rough or sticky.



The truth is:
**Not all coatings
are the same.**

Anthony | User, Continence Care

All **SpeediCath®** catheters come with a unique hydrophilic coating, featuring our **Triple Action Coating Technology™**. It consists of three strongly bonded elements, designed to protect the urethra in three ways:



Stays bonded for an even coating without dry spots.



Stays smooth, reducing friction on the urethra¹.



Stays hydrated so it's instantly ready to use.



Scan this **QR code** to learn more about hydrophilic catheters:

^{*}Limitations apply

1. Compared to uncoated catheters.

Anthony has received compensation from Coloplast to share his information. Each person's situation is unique so your experience may not be the same.

Talk to your healthcare provider about whether this product is right for you.

Prior to use, please refer to product 'Instructions for Use' for intended use and relevant safety information. PM-24135



SpeediCath®

SpeediCath® Flex

Hydrophilic Intermittent Catheter

A soft catheter for men with a dry-sleeve and a flexible tip with all the benefits of the SpeediCath® catheter. Offered with a loop system that fits into a discreet packaging which is 100% aluminum free. Latex, phthalate and PVC free. Sterile.

SpeediCath® Flex					
Effective Length (cm)	Effective Length (inches)	Size	Product Code	Units	
33 cm	13"	10 FR ●	28920	30	
		12 FR ○	28922	30	
		14 FR ●	28924	30	
		16 FR ●	28916	30	

SpeediCath® Standard

Hydrophilic Intermittent Catheter

Instantly ready to use catheter. Featuring a smooth and uniform hydrophilic coating and polished eyelets to ensure a convenient and simple catheterization. Latex, phthalate and PVC free. Sterile.

SpeediCath® Standard Male Straight					
Effective Length (cm)	Effective Length (inches)	Size	Product Code	Units	
35 cm	13.9"	8 FR ●	27408	30	
		10 FR ●	27410	30	
		12 FR ○	27412	30	
		14 FR ●	27414	30	
		16 FR ●	27416	30	
		18 FR ●	27418	30	

SpeediCath® Standard Male Coude					
Effective Length (cm)	Effective Length (inches)	Size	Product Code	Units	
35 cm	13.9"	10 FR ●	27490	30	
		12 FR ○	27492	30	
		14 FR ●	27494	30	
		16 FR ●	27496	30	

SpeediCath® Standard Female Straight					
Effective Length (cm)	Effective Length (inches)	Size	Product Code	Units	
16 cm	6.3"	8 FR ●	27508	30	
		10 FR ●	27510	30	
		12 FR ○	27512	30	
		14 FR ●	27514	30	
		16 FR ●	27516	30	

SpeediCath® Standard Pediatric Straight					
Effective Length (cm)	Effective Length (inches)	Size	Product Code	Units	
16 cm	6.3"	6 FR ●	27706	30	
		8 FR ●	27708	30	
		10 FR ●	27710	30	

SpeediCath® Compact Male

Hydrophilic Intermittent Catheter

Developed to offer discretion and practicality with all the benefits of the SpeediCath® catheter and more. Telescopic design makes the packaging half the size of a standard male catheter. Latex, phthalate and PVC free. Sterile.

SpeediCath® Compact Male					
Effective Length (cm)	Effective Length (inches)	Size	Product Code	Units	
30 cm	11.8"	12/18 FR (one size)	28692	30	



SpeediCath® Compact Eve

Hydrophilic Intermittent Catheter

Stylish and discreet compact catheter for women. Designed for modern life, delivering a one-step easy opening and secure reclosing. Triangular shape ensures easy handling and grip. Latex, phthalate and PVC free. Sterile.

SpeediCath® Compact Eve					
Effective Length (cm)	Effective Length (inches)	Size	Product Code	Units	
8 cm	3.1"	10 FR ●	28110	30	
		12 FR ○	28112	30	
		14 FR ●	28114	30	



SpeediCath® Compact Female

Hydrophilic Intermittent Catheter

Compact and discreet solution that is instantly ready to use. The catheter can easily be kept in a pocket or handbag. Latex, phthalate and PVC free. Sterile.

SpeediCath® Compact Female					
Effective Length (cm)	Effective Length (inches)	Size	Product Code	Units	
7 cm	2.7"	8 FR ●	28578	30	
		10 FR ●	28580	30	
		12 FR ○	28582	30	
		14 FR ●	28584	30	



SpeediCath® Compact Female Plus					
Effective Length (cm)	Effective Length (inches)	Size	Product Code	Units	
9 cm	3.5"	10 FR ●	28810	30	
		12 FR ○	28812	30	
		14 FR ●	28814	30	



SpeediCath® Compact Set

Hydrophilic Intermittent Catheter

Compact and discreet catheter and bag solution that is instantly ready to use. Can easily be kept in a pocket or handbag. Latex, phthalate and PVC free. Sterile.

SpeediCath® Compact Set Male

Effective Length (cm)	Effective Length (inches)	Size		Product Code	Units
30 cm	11.8"	12/18 FR	(one size)	28422	20



SpeediCath® Compact Set Female

Effective Length (cm)	Effective Length (inches)	Size		Product Code	Units
9 cm	3.5"	10 FR	●	28520	20
		12 FR	○	28522	20
		14 FR	●	28524	20



Self-Cath®

Intermittent Catheter

Siliconized surface for smooth insertion. Made of flexible medical grade PVC with smooth fire polished eyelets for patient comfort and ease of use. Latex free. Sterile.

Self-Cath® Female Straight

Effective Length (cm)	Effective Length (inches)	Size	Reference Code	Product Code	Units
12 cm	4.8"	8 FR	●	504600	208
		10 FR	●	504610	210
		12 FR	○	504620	212
		14 FR	●	504630	214
		14 FR	●	501960 without connector	240



Self-Cath® Pediatric

Effective Length (cm)	Effective Length (inches)	Size	Reference Code	Product Code	Units
20 cm	8"	5 FR	○	504400	305
		6 FR	●	504410	306
		8 FR	●	504420 without connector	308
		10 FR	●	504430 without connector	310



Self-Cath® Extension Tube

Total Length (cm)	Total Length (inches)	Size	Reference Code	Product Code	Units
61 cm	24"	n/a	507850	475	12



Self-Cath® Male Straight

Effective Length (cm)	Effective Length (inches)	Size	Reference Code	Product Code	Units
38 cm	14.9"	8 FR	● 504500	408	30
		10 FR	● 504510	410	30
		12 FR	○ 504520	412	30
		14 FR	● 504530	414	30
		16 FR	● 504540	416	30
		18 FR	● 504550	418	30



Self-Cath® Male Soft Straight

Effective Length (cm)	Effective Length (inches)	Size	Reference Code	Product Code	Units
38 cm	14.9"	12 FR	○ 501810	112	30
		14 FR	● 501820	114	30
		16 FR	● 501830	116	30

Self-Cath® Male Straight Curved Packaging

Effective Length (cm)	Effective Length (inches)	Size	Reference Code	Product Code	Units
38 cm	14.9"	12 FR	○ 502200	460	50
		14 FR	● 502210	450	50

Self-Cath® Male Coude Tapered with Guide Stripe

Effective Length (cm)	Effective Length (inches)	Size	Reference Code	Product Code	Units
38 cm	14.9"	8 FR	● 504650	608	30
		10 FR	● 504660	610	30
		12 FR	○ 504670	612	30
		14 FR	● 504680	614	30



Self-Cath® Male Coude Olive with Guide Stripe

Effective Length (cm)	Effective Length (inches)	Size	Reference Code	Product Code	Units
38 cm	14.9"	8 FR	● 504710	808	30
		10 FR	● 504720	810	30
		12 FR	○ 504730	812	30
		14 FR	● 504740	814	30
		16 FR	● 504750	816	30
		18 FR	● 504760	818	30



Self-Cath® Male Closed System Straight

Effective Length (cm)	Effective Length (inches)	Size	Reference Code	Product Code	Units
31 cm	12.2"	14 FR	● 502730	1114	50
		16 FR	● 502740	1116	50



Self-Cath®

Designed for patients
comfort and ease of use.

Brianna
Winnipeg, MB


Features & Benefits

Material:

Self-Cath® catheters are made with a flexible
polyvinyl chloride (PVC).


Smooth and comfortable insertion:

Self-Cath® offers *fire polished eyelets* that allow the urethral mucosa to slide over the eyelets reducing the risk of being drawn into the lumen during insertion.



Easy guidance:

Self-Cath® Coudé catheters have a patented *Guide Stripe* which have a blue stripe indicating tip orientation. Therefore, patients can be sure of where the tip faces for the entire cathing experience.



Easy navigation:

Self-Cath® is available in *Straight tip*, *Tapered Tip* and *Olive Tip* option. Whichever your patients use, they may experience easier navigation through the curves and bends of their anatomy.



Self-Cath® catheters offer a range of options for many patient needs:

- Male straight and coude tip
- Female straight
- Pediatric (including 5 FR)
- Tapered coude and olive coude tips
- Available in a soft format and closed system



Please scan the
QR code to view
the product guide:





Brad
User | Ostomy Care

Mental Health Support through Coloplast® Care for your patients

At Coloplast, we recognize your importance as a healthcare professional in the lives of your patients. We also recognize that all your patients are different.



We know your time is in short supply and personalized care tailored to patients is vital. That's why we have teamed up with **Teladoc Health** to add mental health support to our

care program. This is available to your patients who are referred to the Coloplast® Care program by you. This support provides three virtual counselling sessions and full access to the **myStrength** app for guidance and support.

We are here for you, supporting your care for every patient.

To learn more about Coloplast® Care contact your Territory Manager or contact us at:
1-866-293-6349 | CProfessional_CA@Coloplast.com

To sign your patients up for mental health support through Coloplast® Care scan the **QR Code** or visit:
visit.coloplast.ca/MentalHealthSupport



Emotional support

Patients have access to 3 sessions with a mental health specialist and unlimited use of the 'myStrength' app through **Teladoc Health**.



Reimbursement and product guidance

Patients can order product samples, receive reimbursement support, and manage their daily needs with a Care Advisor.



Resources and tools

Your patients can browse and download our digital educational resources, eBooks, and tools.

Conveen® Optima

Male External Catheter - Silicone

Discreet, easy to use, reliable all day self-adhering external catheter. Double-grip strip ensures easy and smooth roll-out for a secure application. Push ring with sure grip ribs provides secure connection to the collecting bag. Anti-kink bellow stop kinks ensuring a free flow. Compact easy-to-open packaging. Latex free.

Conveen® Optima Standard

Length (cm)	Length (inches)	Size	Product Code	Units
7.8 cm	3"	25 mm	22025	30
		28 mm	22028	30
		30 mm	22030	30
		35 mm	22035	30
		40 mm	22040	30

Conveen® Optima Sport

Length (cm)	Length (inches)	Size	Product Code	Units
4.6 cm	2"	21 mm	22121	30
		25 mm	22125	30
		30 mm	22130	30
		35 mm	22135	30



Conveen® Security+

Male External Catheter - SEBS*

Easy to use, reliable all day self-sealing external catheter. Applicator grip and single-grip strip ensures easy and smooth roll-out for a secure application. Push ring with sure grip ribs provides secure connection to the collecting bag. Anti-kink bellow stop kinks ensuring a free flow. Latex free.

Conveen® Security+ Self Sealing

Length (cm)	Length (inches)	Size	Product Code	Units
8.2 cm	3"	21 mm	5221	35
		25 mm	5225	35
		25 mm	5225H	100
		30 mm	5230	35
		30 mm	5230H	100
		35 mm	5235	35
		35 mm	5235H	100



Conveen® Security+ Adhesive Liner

Length (cm)	Length (inches)	Size	Product Code	Units
8.2 cm	3"	30 mm	5030	35
		30 mm	5030H	100
		35 mm	5035	35
		35 mm	5035H	100



*Styrene Ethylene Butylene Styrene

Conveen®

Male External Catheter - Latex

Conveen® Self-Sealing					
Length (cm)	Length (inches)	Size		Product Code	Units
10.5 cm	4"	20 mm	●	5212	35
		25 mm	●	5200	35
		30 mm	●	5205	35
		35 mm	●	5210	35
		40 mm	●	5215	35



Conveen® Adhesive Liner					
Length (cm)	Length (inches)	Size		Product Code	Units
10.5 cm	4"	30 mm	●	5130	35
		35 mm	●	5135	35



Active Cath®

Male External Catheter - Latex

Self-sealing male external catheter.

Active Cath®					
Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
7.9 cm	3"	23 mm	●	506220	100
		28 mm	●	506200	100
		31 mm	●	506240	100
		35 mm	●	506260	100



Freedom Cath®

Male External Catheter - Latex

Self-sealing male external catheter.

Freedom Clear®					
Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
7.9 cm	3"	23 mm	●	506110	100
		28 mm	●	506130	100
		31 mm	●	506140	100
		35 mm	●	506150	100



Conveen® Active

Urine Leg Bag

Developed for an active lifestyle. Integrated in a comfortable elastic band that wraps around the thigh for ease of use, mobility and discretion. Hide-away outlet tap concealed in elastic band when not in use. Soft fabric band for increased wearing comfort and a no-slip silicone feature inside straps for security. Latex and PVC free.

Conveen® Active Leg Bag, Non Sterile

Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
9 oz	250 mL	3"	25501	10



Conveen® Security+

Urine Leg bag

Secure and reliable urine leg bag with anti-kink tubing to reduce the risk of accidental leakage. Soft backing for increased wearing comfort (except codes 21034 and 21054). Straps Included. Latex and PVC free.

Conveen® Security+ Leg Bag, Clamp Outlet, Non Sterile

Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
17 oz	500 mL	20"	5161	10
26 oz	750 mL	20"	5167	10

Conveen® Security+ Leg Bag, Clamp Outlet, Sample Port, Sterile

Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
17 oz	500 mL	20"	5163	10

Conveen® Security+ Leg Bag, Lever Outlet, Sterile

Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
17 oz	500 mL	20"	21034	10
34 oz	1000 mL	20"	21054	10



Conveen® Standard

Day / Night Urine Bag

Secure and reliable urine bag for bedside drainage or large capacity leg bag. Clamp outlet makes it easy to verify proper closure and reduce the risk of accidental leakage. Standard tubing. Leg bag straps sold separately. Soft backing for increased wearing comfort. Latex free.

Conveen® Standard Day / Night Bag, Clamp Outlet, Non Sterile

Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
51 oz	1500 mL	36"	5062	10

Conveen® Standard Day / Night Bag, Clamp Outlet, Sample Port, Sterile

Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
51 oz	1500 mL	36"	5063	10



Conveen® Contour

Urine Leg Bag

Secure and reliable urine leg bag. Enhanced chambered design for added discretion. Soft backing for increased wearing comfort. Straps Included. Latex free.

Conveen® Contour Leg Bag, Clamp Outlet, Non Sterile				
Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
21 oz	600 mL	18"	5170	10
28 oz	800 mL	18"	5174	10



Conveen® Night Bag

Urine Night Bag

Secure and reliable urine night bag with anti-kink tubing to prevent accidental leakage. Large capacity allows the user to feel secure all night for a full night's sleep. Latex and PVC free. Not Compatible with Urostomy Pouches.

Conveen® Bedside Night Bag, Non-Sterile			
Capacity (oz)	Capacity (L)	Product Code	Units
68 oz	2L	21346	30

Conveen® Bedside Night Bag, Sterile			
Capacity (oz)	Capacity (L)	Product Code	Units
68 oz	2L	21356	10



Accessories

Conveen® Bag Hanger			
Length (cm)	Length (inches)	Product Code	Units
18 cm	7"	5070	1

Conveen® Leg Bag Strap			
Length (cm)	Length (inches)	Product Code	Units
63 cm	25"	Sold in pairs 5050	10

Conveen® Catheter retaining Strap			
Length (cm)	Length (inches)	Product Code	Units
90 cm	35"	383004	2



Peristeen Plus® Anal Irrigation

Transanal irrigation system for fecal incontinence and constipation. Minimizes the likelihood of involuntary bowel leakage and/or constipation.

Peristeen Plus® Anal Irrigation System

Items	Description	Product Code
- 1 control unit		
- 2 sterile rectal catheters		
- 1 bag (with lid)	** 90 irrigations = 6 months	29140
- 2 straps		
- 1 toiletries bag		



Peristeen Plus® Anal Irrigation System - Small

Items	Description	Product Code
- 1 control unit		
- 2 sterile rectal catheters - small		
- 1 bag (with lid)	** 90 irrigations = 6 months	29147
- 2 straps		
- 1 toiletries bag		

Peristeen Plus® Anal Irrigation Accessory unit

Items	Description	Product Code
- 15 sterile rectal catheters	**15 irrigations = 1 month	29142
- 1 bag (no lid)		



Peristeen Plus® Anal Irrigation Accessory unit - Small

Items	Description	Product Code
- 15 sterile rectal catheters - small	**15 irrigations = 1 month	29149
- 1 bag (no lid)		

Peristeen Plus® Anal Irrigation Rectal Catheter

Items	Product Code
- 10 sterile rectal catheters	29143



Peristeen Plus® Anal Irrigation Leg Strap

Items	Product Code
- 1 pair	29145



Peristeen Plus® Anal Irrigation Tube

Items	Product Code
- 2 tubes with blue connectors	29146



** if used every other day

Advanced Wound Dressings

Biatain® Silicone

Foam Dressing with Silicone Adhesive

Soft, flexible and absorbent foam dressing with gentle silicone adhesive for use on a broad range of exuding wounds, making it the obvious choice for moist wound healing of commonly occurring wounds. For more information please visit www.biatainsiliconeca.com.



reddot award 2017
winner



Biatain® Silicone			
Size (cm)	Size (inches)	Product Code	Units
7.5 x 7.5 cm	3 x 3"	33434	10
10 x 10 cm	4 x 4"	33435	10
12.5 x 12.5 cm	5 x 5"	33436	10
15 x 15 cm	6 x 6"	33437	5
17.5 x 17.5 cm	7 x 7"	33438	5
10 x 20 cm	3.9 x 7.9"	33400	5
10 x 30 cm	3.9 x 11.8"	33401	5
18 x 18 cm	7.1 x 7.1" Heel	33406	5
25 x 25 cm	9.8 x 9.8" Large Sacral	39000*	5
15 x 19 cm	5.9 x 7.5" Small Sacral	39001*	5
14 x 19.5 cm	5.5 x 7.7" Multi-shape	39010*	5



* This product has the Pressure Injury Prevention indication.

Biatain® Silicone Non-Border

Designed for high absorption and retention under compression to minimise the risk of leakage and maceration of fragile periwound skin. Can also be cut to fit areas that are difficult to dress, e.g. on the heel or in between toes and fingers. For high absorption and retention under compression. 3DFit Technology minimizes risk of leakage and maceration.

Biatain® Silicone			
Size (cm)	Size (inches)	Product Code	Units
5 x 7.5 cm	2 x 3"	39020	10
7.5 x 7.5 cm	3 x 3"	39021	10
10 x 10 cm	4 x 4"	39022	10
12.5 x 12.5 cm	5 x 5"	39023	10
10 x 20 cm	4 x 8"	39024	5
15 x 15 cm	6 x 6"	39025	5
20 x 20 cm	8 x 8"	39026	5



Biatain® Silicone Lite

Thin Foam Dressing with Silicone Adhesive

Thin and ultra flexible foam dressing with gentle silicone adhesive for use on a wide range of low to moderate exuding wounds.



Biatain® Silicone Lite			
Size (cm)	Size (inches)	Product Code	Units
5 x 5 cm	2 x 2"	33452	5
5 x 12.5 cm	2 x 5"	33453	5
7.5 x 7.5 cm	3 x 3"	33444	10
10 x 10 cm	4 x 4"	33445	10
12.5 x 12.5 cm	5 x 5"	33446	10



Biatain® Silicone Ag

Foam Dressing with Silicone Adhesive and Silver

Flexible and absorbent foam dressing with gentle silicone adhesive and a patented silver complex for use on a broad range of exuding infected wounds or wounds at risk of infection.



Biatain® Silicone Ag				
Size (cm)	Size (inches)	Product Code	Units	
7.5 x 7.5 cm	3 x 3"	39636	5	     
10 x 10 cm	4 x 4"	39637	5	
12.5 x 12.5 cm	5 x 5"	39638	5	
15 x 15 cm	6 x 6"	39639	5	
17.5 x 17.5 cm	7 x 7"	39640	5	
10 x 20 cm	3.9 x 7.9"	39644	5	
10 x 30 cm	3.9 x 11.8"	39645	5	
15 x 19 cm	5.9 x 3.5" Small Sacral	39650	5	
25 x 25 cm	9.8 x 9.8" Large Sacral	39651	5	
18 x 18 cm	7.1 x 7.1" Heel	39652	5	

Biatain®

Foam Dressing

Soft and conformable polyurethane foam dressing that effectively absorbs and retains wound exudate. This ensures a balance that is optimal for management of exuding wounds. Available in adhesive and non-adhesive. The adhesive has a hydrocolloid border – for a secure adhesion. This product is also used in high friction and shear areas such as the coccyx.



Biatain® Adhesive				
Size (cm)	Size (inches)		Product Code	Units
7.5 x 7.5 cm	3 x 3"		3462	10
10 x 10 cm	4 x 4"		3430	10
12.5 x 12.5 cm	5 x 5"		3420	10
17 x 17 cm	6.7 x 6.7"	Small Sacral	3483	5
18 x 18 cm	7 x 7"		3423	5
19 x 20 cm	7.5 x 8"	Heel	3488	5
23 x 23 cm	9 x 9"	Large Sacral	3485	5



Biatain® Non-Adhesive				
Size (cm)	Size (inches)		Product Code	Units
5 x 7 cm	2 x 2.8"		6105	10
10 x 10 cm	4 x 4"		3410	10
10 x 20 cm	4 x 8"		3412	5
15 x 15 cm	6 x 6"		3413	5
20 x 20 cm	8 x 8"		3416	5



Biatain® Ag

Foam Dressing with Silver

Soft and conformable polyurethane foam dressing with a patented silver complex for use on infected wounds or wounds at risk of infection. Available in non-adhesive. This product is also used in high friction and shear areas such as the coccyx.



Biatain® Ag Non-Adhesive

Size (cm)	Size (inches)	Product Code	Units
5 x 7 cm	2 x 2.8"	5105	5
10 x 10 cm	4 x 4"	9622	5
10 x 20 cm	4 x 8"	9623	5
15 x 15 cm	6 x 6"	9625	5
20 x 20 cm	8 x 8"	9626	5



Biatain® Fiber

Suitable for highly exuding wounds and locks in exudate and bacteria. Effectively absorbs and retains exudate and bacteria, thereby reducing exudate pooling. HexaLock Technology is the new gelling fiber for all exuding wounds.



reddot award 2017
winner

Biatain® Fiber

Size (cm)	Size (inches)	Product Code	Units
5 x 5 cm	2 x 2"	33501	10
10 x 10 cm	4 x 4"	33503	10
15 x 15 cm	6 x 6"	33506	10
19 x 25 cm	7.5 x 9.8"	33507	10
2.5 x 46 cm	1.5 x 18"	33508	10



Biatain® Alginate

Calcium Alginate Dressing

Highly absorbent alginate dressing for use on moderate to heavily exuding wounds of any size, shape and depth.

Biatain® Alginate

Size (cm)	Size (inches)	Product Code	Units
5 x 5 cm	2 x 2"	3705	30
10 x 10 cm	4 x 4"	3710	10
15 x 15 cm	6 x 6"	3715	10
3 x 44 cm	1 x 17.5" Ribbon	3740	6



Biatain® Alginate Ag

Calcium Alginate Dressing with Silver

Highly absorbent alginate dressing with silver for use on moderate to heavily exuding infected wounds or wounds at risk of infection of any size, shape and depth.

Biatain® Alginate Ag			
Size (cm)	Size (inches)	Product Code	Units
5 x 5 cm	2 x 2"	3755	30
10 x 10 cm	4 x 4"	3760	10
15 x 15 cm	6 x 6"	3765	10
3 x 44 cm	1 x 17.5" Ribbon	3780	10



Biatain® Ibu

Foam Dressing with Ibuprofen

Soft and flexible absorbent foam dressing containing ibuprofen for use on painful wounds. Available in non-adhesive and soft hold.



Biatain® Ibu Non-Adhesive			
Size (cm)	Size (inches)	Product Code	Units
10 x 10 cm	4 x 4"	4110	5
10 x 20 cm	4 x 8"	4112	5
15 x 15 cm	6 x 6"	4115	5
20 x 20 cm	8 x 8"	4120	5

Biatain® Ibu Soft-Hold			
Size (cm)	Size (inches)	Product Code	Units
10 x 10 cm	4 x 4"	4140	5
10 x 20 cm	4 x 8"	4142	5



Biatain® Super

Hydrocapillary Pad Dressing

Highly absorbent hydrocapillary pad dressing for use on heavily exuding wounds. Available in non-adhesive and adhesive.

Biatain® Super Adhesive			
Size (cm)	Size (inches)	Product Code	Units
10 x 10 cm	4 x 4"	4610	10
12.5 x 12.5 cm	5 x 5"	4612	10
12 x 20 cm	4.5 x 8"	4625	10
15 x 15 cm	6 x 6"	4615	10
20 x 20 cm	8 x 8"	4620	10

Biatain® Super Non-Adhesive			
Size (cm)	Size (inches)	Product Code	Units
10 x 10 cm	4 x 4"	4630	10
12.5 x 12.5 cm	5 x 5"	4632	10
12 x 20 cm	4.5 x 8"	4645	10
15 x 15 cm	6 x 6"	4635	10
20 x 20 cm	8 x 8"	4639	10



Comfeel® Plus

Hydrocolloid Dressing

Hydrocolloid dressing for use on low to moderate exuding wounds. Provides faster wound healing by sealing and protecting the wounds.

Comfeel® Plus				
Size (cm)	Size (inches)		Product Code	Units
4 x 6 cm	1.5 x 2.5"		33146	30
10 x 10 cm	4 x 4"		33110	10
15 x 15 cm	6 x 6"		33115	5
20 x 20 cm	8 x 8"		33120	5
18 x 20 cm	7 x 8"	Sacral	33285	5
9 x 11 cm	3 x 4"	Contour	33283	5



Comfeel® Plus Transparent

Hydrocolloid Dressing

Thin and flexible hydrocolloid dressing for use on dry or low exuding wounds. Provides faster wound healing by sealing and protecting the wound. Transparency allows for visual monitoring of the wound.

Comfeel® Plus Transparent				
Size (cm)	Size (inches)		Product Code	Units
5 x 7 cm	2 x 2.75"		33530	10
10 x 10 cm	4 x 4"		33533	10
9 x 14 cm	3.5 x 5.5"		33536	10
15 x 15 cm	6 x 6"		33539	5
9 x 25 cm	3.5 x 10"		33537	5
5 x 15 cm	2 x 6"		33547	10
5 x 25 cm	2 x 10"		33548	10



Triad®

Hydrophilic Wound Dressing

Hydrophilic wound dressing. Used for autolytic debridement. The ideal alternative for difficult-to-dress areas. Sterile on first application. For more information please visit triad.coloplast.ca.

Triad®				
Size (ounces)	Size (grams)		Product Code	Units
2.5 oz.	71 g tube		1964	12
6 oz.	170 g tube		1967	12



Purilon® Gel

Hydrogel with natural ingredients

Effective and gentle debridement.

Purilon® Gel			
Size (ounces)	Size (grams)	Product Code	Units
28 oz.	8 g tube	3906	10
5 oz.	15 g tube	3900	10



InterDry®

Moisture-wicking textile with antimicrobial silver complex

Effectively manage complications associated with skin folds. Appropriate for a variety of skin-to-skin and skin-to-device application sites. For more information please visit www.interdry.ca.

InterDry®			
Size (cm)	Size (Inches)	Product Code	Units
25 x 46 cm	10 x 18"	7915	40
25 x 92 cm	10 x 36"	7912	10
25 x 366 cm	10 x 144" textile on a roll	7910	10



Brief Statements

Biatain® Family of Dressings

Biatain® dressings are single-use only. See the device manual for detailed information regarding complete instructions, contraindications, warnings, precautions and potential adverse events. For further questions, call Coloplast Canada Corporation at 1-888-880-8605 and/or consult the company website at www.coloplast.ca.

InterDry®

InterDry® is a skin protectant for the management of skin folds and other skin-to-device contact areas. It reduces microbial colonization in the fabric. This single use fabric provides moisture transportation to keep skin dry while the antimicrobial in the fabric reduces odor. The safe use of InterDry fabric during pregnancy, lactation and on children has not been demonstrated. Although many patients benefit from the use of this device, results may vary. See the device labeling for detailed information including the indication, contraindications, warnings, precautions, and potential complications/adverse events. For further questions, call Coloplast Canada Corporation at 1-888-880-8605 and/or consult the company website at www.coloplast.ca.

Peristeen Plus®

The Peristeen Plus® Anal Irrigation System is intended to instill water into the colon through a rectal catheter, which incorporates an inflatable balloon, inserted into the rectum to promote evacuation of the contents of the lower colon. The Peristeen Plus® Anal Irrigation System is indicated for use by children (3 years - <12 years old), adolescent (12 years - < 18 years old), transitional adolescent (18 - <21 years old) and adult patients with neurogenic bowel dysfunction who suffer from fecal incontinence, chronic constipation, and/or time-consuming bowel management procedures.

Contraindications: Peristeen Plus® Anal Irrigation must not be used in the following situations: known obstruction of the large bowel due to strictures or tumors, acute inflammatory bowel disease, diverticulitis, complex diverticular disease, abdominal or anal surgery within the last 3 months, spinal shock phase, and in patients who are pregnant and have not used the system before (if the patient is pregnant and has never used anal irrigation before, they should not start the irrigation procedure during pregnancy).

Warning: Anal irrigation procedure should always be carried out with care. Bowel perforation is an extremely rare, but serious and potentially lethal complication to anal irrigation and will require immediate admission to a hospital, often requiring surgery.

See the device manual for detailed information regarding the procedure, contraindications, warnings, precautions, and potential complications/adverse events. For further information, call Coloplast Canada Corporation at 1-866-293-6349 and/or consult the company website at www.coloplast.ca.

SpeediCath®

SpeediCath® catheters are indicated for use by patients with chronic urine retention and patients with a post void residual volume (PVR) due to neurogenic and non-neurogenic voiding dysfunction. The catheter is inserted into the urethra to reach the bladder allowing urine to drain. There is a separate SpeediCath device intended for either males or females. This is a single-use only product. Patients are to contact their healthcare professional before performing self-catheterization for the first time and follow the instructions carefully. See the device manual for further information or, call Coloplast Canada Corporation at 1-866-293-6349 and/or consult the company website at www.coloplast.ca.

Self-Cath®

Self-Cath® Catheters are intended for use in male, female, and pediatric patients (neonates, infants, children, adolescents, and transitional adolescents) requiring bladder drainage as determined by their physician. The catheter is indicated for those individuals unable to promote a natural urine flow or for those individuals who have a significant volume of residual urine following a natural bladder-voiding episode. Single-use only. See the device manual for detailed information regarding the procedure, contraindications, warnings, precautions, and potential complications/adverse events. For further information, call Coloplast Canada Corporation at 1-866-293-6349 and/or consult the company website at www.coloplast.ca.

Triad®

Triad® hydrophilic wound dressing is indicated for the local management of pressure and venous stasis ulcers, dermal lesions/injuries, superficial wounds, scrapes, first-and second-degree burns, partial-and full-thickness wounds. Suitable for use on neonates.

Contraindications: Triad® is contraindicated for third degree burns or infected wounds. See device label for complete information including warnings and precautions. For further questions, call Coloplast Canada Corporation at 1-888-880-8605 and/or consult the company website at www.coloplast.ca.

Product Code	Page
112	50
114	50
116	50
208	49
210	49
212	49
214	49
240	49
305	49
306	49
308	49
310	49
408	50
410	50
412	50
414	50
416	50
418	50
0421	39
0422	39
0423	39
0424	39
450	50
460	50
475	49
0509	40
608	50
610	50
612	50
614	50
808	50
810	50
812	50
814	50
816	50
818	50
0925	40
1110	41
1114	50
1116	50
1500	41
1907	39
1964	65
1967	65
2041	40
2115	35
2125	35
2151	34
2166	34
2181	34
2320	40
2325	40
2330	40
2340	40
2350	40
2490	32

Product Code	Page
2501	41
2520	32
2650	39
2655	39
2792	31
2793	31
2801	41
2802	41
2803	41
2807	41
2808	41
2831	30
2832	30
2833	30
2836	30
2837	30
2846	30
2847	30
2881	29
2882	29
2883	29
3210	38
3215	38
3220	38
3410	62
3412	62
3413	62
3416	62
3420	62
3423	62
3430	62
3462	62
3483	62
3485	62
3488	62
3705	63
3710	63
3715	63
3740	63
3755	64
3760	64
3765	64
3780	64
3900	66
3906	66
4110	64
4112	64
4115	64
4120	64
4140	64
4142	64
4610	64
4612	64
4615	64
4620	64
4625	64

Product Code	Page
4630	64
4632	64
4635	64
4639	64
4645	64
5030	53
5030H	53
5035	53
5035H	53
5050	56
5062	55
5063	55
5070	56
5105	63
5130	54
5135	54
5161	55
5163	55
5167	55
5170	56
5174	56
5200	54
5205	54
5210	54
5212	54
5215	54
5221	53
5225	53
5225H	53
5230	53
5230H	53
5235	53
5235H	53
5580	34
5585	34
6105	62
7583	40
7910	66
7912	66
7915	66
8000	54
8009	35
8100	54
8200	54
8205	54
8300	54
8305	54
8400	54
8500	54
9500	40
9505	40
9622	63
9623	63
9625	63
9626	63
10011	23

Product Code	Page
10012	23
10013	23
10015	24
10016	24
10017	24
10021	23
10022	23
10023	23
10025	24
10026	24
10027	24
10031	23
10032	23
10035	24
10036	24
10041	23
10045	24
10101	25
10102	25
10103	25
10104	26
10105	26
10106	26
10108	25
10111	25
10113	25
10114	25
10116	25
10117	25
10154	24
10155	24
10156	24
10164	24
10165	24
10166	24
10167	24
10451	16
10461	16
10471	16
10472	16
10473	16
10474	16
10475	16
10476	16
10477	16
10478	16
10479	16
10480	16
10481	16
10482	17
10483	17
10484	17
10487	17
10489	16
10502	9
10503	9

Product Code	Page
10504	9
10505	9
10506	9
10512	9
10513	9
10514	9
10515	9
10516	9
10522	9
10523	9
10524	9
10525	9
10531	9
10551	14
10561	14
10562	14
10563	14
10564	14
10565	14
10566	14
10567	14
10568	14
10571	14
10572	14
10573	14
10582	20
10583	20
10584	20
10585	20
10586	20
10870	19
10871	19
10872	19
10873	19
10874	19
10881	19
10883	19
10884	19
10885	19
10904	26
10905	26
10914	26
10915	26
10918	26
11011	23
11012	23
11015	24
11016	24
11021	23
11022	23
11023	23
11024	23
11025	24
11026	24
11027	24
11028	24

Product Code	Page
11031	23
11032	23
11035	24
11036	24
11041	23
11045	24
11114	24
11115	24
11124	24
11125	24
11126	24
11127	24
11135	24
11136	24
11185	24
11186	24
11194	24
11195	24
11196	24
11197	24
11301	25
11302	25
11303	25
11304	26
11305	26
11306	26
11307	25
11308	26
11312	25
11313	25
11314	25
11315	25
11316	25
11401	10
11402	10
11412	10
11413	10
11421	10
11422	10
11425	10
11432	10
11443	10
11461	10
11462	10
11463	10
11471	10
11472	10
11473	10
11474	10
11482	10
11483	10
11484	10
11486	10
11487	10
11491	11
11492	11

Product Code	Page
11496	11
11497	11
11498	11
11499	11
11500	11
11505	26
11514	26
11515	26
11518	26
11612	26
11615	26
11690	26
11750	41
11801	28
11802	28
11803	28
11804	28
11811	28
11812	28
11814	28
11815	28
11821	26
11822	26
11824	26
11825	26
11827	26
11828	26
11841	25
11842	25
11844	25
11845	25
11846	25
11851	25
11852	25
11854	25
11855	25
11856	25
12000	39
12010	39
12011	39
12014	22
12020	38
12021	38
12030	38
12032	38
12033	38
12035	38
12036	38
12037	38
12038	38
12039	38
12042	38
12045	38
12046	38
12047	38
12048	38

Product Code	Page
12049	38
12050	39
12060	39
12061	39
12070	39
12072	39
12074	39
12076	39
12090	38
12091	38
12092	38
12093	38
12094	38
12095	38
12110	33
12130	33
12140	33
12144	33
12145	33
12146	33
12147	33
12149	33
12170	33
12174	33
12175	33
12176	33
12177	33
12201	15
12202	15
12212	15
12213	15
12222	15
12232	15
12243	15
12262	15
12271	15
12272	15
12274	15
12282	15
12284	15
12291	15
12292	15
12297	15
12298	15
12344	31
12346	31
12354	31
12355	31
12356	31
12374	31
12375	31
12376	31
12384	31
12385	31
12386	31
12432	34

Product Code	Page
12433	34
12471	34
12472	34
12533	32
12595	34
12596	34
12701	29
12702	29
12703	29
12704	29
12705	29
12706	29
12707	29
12708	29
12709	29
12716	29
12719	29
12800	33
12802	33
12805	33
12808	33
12810	33
12820	41
12834	42
12835	42
12836	42
12841	29
12842	29
12844	29
12845	29
12847	29
12848	29
12991	34
12992	34
12993	34
12994	34
12995	34
12996	34
12997	34
13671	17
13672	17
13673	17
13674	17
13675	17
13676	17
13678	20
13679	20
13680	20
13683	19
13684	19
13685	19
13706	32
13840	32
13844	32
13845	32
13846	32

Product Code	Page
13860	32
13870	32
13874	32
13875	32
13876	32
13924	31
13925	31
13964	31
13965	31
13966	31
13974	31
13975	31
13976	31
13984	31
13985	31
13986	31
14010	42
14015	42
14021	42
14031	42
14041	42
14050	42
14060	42
14070	42
14103	32
14106	32
14163	32
14164	32
14165	32
14166	32
14173	32
14174	32
14175	32
14176	32
14196	32
14201	34
14202	34
14224	31
14225	31
14226	31
14227	31
14228	31
14229	31
14232	30
14233	30
14234	30
14235	30
14236	30
14237	30
14238	30
14239	30
14243	30
14246	30
14249	30
14251	30
14252	30

Product Code	Page
14254	30
14255	30
14258	30
14261	29
14262	29
14263	29
14271	29
14272	29
14273	29
14274	29
14275	29
14276	29
14277	29
14278	29
14281	30
14282	30
14283	30
14291	30
14292	30
14293	30
14294	30
14295	30
14296	30
14307	35
14308	35
14404	32
14405	32
14411	32
14412	32
14414	32
14415	32
14422	32
14423	32
14424	32
14441	33
14445	33
14452	33
14453	33
14454	33
14494	31
14498	31
14681	35
14682	35
14691	35
14692	35
14701	34
14702	34
15403	28
15404	28
15405	28
15406	28
15407	28
15408	28
15413	28
15414	28
15415	28

Product Code	Page
15441	28
15450	28
15470	28
15472	28
15473	28
15474	28
15480	28
15482	28
15483	28
15484	28
15521	27
15522	27
15523	27
15524	27
15531	27
15532	27
15533	27
15534	27
15605	27
15606	27
15607	27
15608	27
15633	27
15634	27
15635	27
15695	27
15696	27
15821	27
15981	27
16301	19
16305	19
16306	19
16310	19
16311	19
16312	19
16313	19
16315	19
16322	19
16323	19
16325	19
16326	19
16327	19
16330	19
16331	19
16332	19
16333	19
16335	19
16345	19
16346	19
16347	19
16350	19
16351	19
16352	19
16353	19
16471	14
16472	14

Product Code	Page
16481	14
16483	14
16484	14
16485	14
16486	14
16487	14
16488	14
16491	14
16493	14
16494	14
16705	17
16706	17
16715	17
16716	17
16725	17
16726	17
16727	17
16730	17
16731	17
16732	17
16733	17
16734	17
16735	17
16737	17
16738	17
16739	17
16748	17
16755	17
16756	17
16757	17
16760	17
16761	17
16762	17
16763	17
16764	17
16767	17
16768	17
16770	17
16771	17
16772	17
16773	17
16778	17
16801	20
16805	20
16810	20
16825	20
16826	20
16827	20
16835	20
16836	20
16837	20
16840	20
16841	20
16842	20
16855	20
16856	20

Product Code	Page
16857	20
16865	20
16866	20
16867	20
16870	20
16871	20
16872	20
16880	14
16881	14
16882	14
16883	14
16884	14
16885	14
16886	14
16887	14
16888	14
16889	14
16890	14
16891	14
16901	9
16903	9
16911	9
16913	9
16914	9
16915	9
16921	9
16923	9
16924	9
16931	9
16933	9
16934	9
16941	9
16942	9
16951	9
16953	9
16954	9
16955	9
16961	9
16963	9
16964	9
16971	9
16973	9
16974	9
18151	20
18152	20
18160	20
18161	20
18162	20
18163	20
18170	20
18171	20
18172	20
18173	20
18180	20
18181	20
18182	20

Product Code	Page
18185	20
18300	18
18301	18
18302	18
18305	18
18310	18
18312	18
18313	18
18314	18
18315	18
18320	18
18321	18
18322	18
18323	18
18324	18
18330	18
18331	18
18332	18
18333	18
18340	18
18341	18
18342	18
18343	18
18344	18
18371	16
18373	16
18382	16
18383	16
18386	15
18389	15
18390	14
18401	21
18402	21
18403	21
18404	21
18405	21
18500	10
18501	10
18510	10
18511	10
18512	10
18513	10
18520	10
18521	10
18522	10
18523	10
18524	10
18530	10
18531	10
18532	10
18600	11
18601	11
18602	11
18603	11
18610	11
18611	11

Product Code	Page
18612	11
18613	11
18620	11
18621	11
18622	11
18630	11
18631	11
18632	11
18640	11
18641	11
18642	11
18650	15
18651	15
18652	15
18655	15
18656	15
18657	15
18660	21
18661	21
18665	21
18666	21
18670	21
18671	21
18675	21
18676	21
18680	21
18681	21
18690	21
18691	21
18700	22
18705	22
18710	23
18711	23
18712	23
18713	23
18720	22
18721	22
18722	22
18723	22
18724	22
18725	22
18730	22
18731	22
18735	22
18747	22
19001	28
19003	28
19010	28
19011	28
19020	28
19021	28
19030	28
19040	25


Product Code	Page
19041	25
19042	25
19050	25
19051	25
19052	25
21034	55
21054	55
21346	56
21356	56
21365	34
22025	53
22028	53
22030	53
22035	53
22040	53
22121	53
22125	53
22130	53
22135	53
25501	55
27408	47
27410	47
27412	47
27414	47
27416	47
27418	47
27490	47
27492	47
27494	47
27496	47
27508	47
27510	47
27512	47
27514	47
27516	47
27706	47
27708	47
27710	47
28110	48
28112	48
28114	48
28422	49
28520	49
28522	49
28524	49
28578	48
28580	48
28582	48
28584	48
28692	48
28810	48
28812	48
28814	48

Product Code	Page
28916	47
28920	47
28922	47
28924	47
29140	57
29142	57
29143	57
29145	57
29146	57
29147	57
29149	57
33110	65
33115	65
33120	65
33146	65
33283	65
33285	65
33400	61
33401	61
33406	61
33434	61
33435	61
33436	61
33437	61
33438	61
33444	61
33445	61
33446	61
33452	61
33453	61
33501	63
33503	63
33506	63
33507	63
33508	63
33530	65
33533	65
33536	65
33537	65
33539	65
33547	65
33548	65
39000*	61
39001*	61
39010*	61
39020	61
39021	61
39022	61
39023	61
39024	61
39025	61
39026	61
39636	62

Product Code	Page
39637	62
39638	62
39639	62
39640	62
39644	62
39645	62
39650	62
39651	62
39652	62
383004	56
506110	54
506130	54
506140	54
506150	54
506200	54
506220	54
506240	54
506260	54

Coloplast[®] Care

Support for your ostomy or continence care needs



Website & Emails

A dedicated website and regular emails with articles and reliable advice.




Self-assessment Tools

Helping users keep track of their ostomy or continence health through regular self-assessment.



Phone Support

Advisors ready to answer questions and provide support by phone.



Product Guidance

Access to product samples on request.

Enrolling in Care is easy!

Online at: www.coloplast.ca/careform

By phone at: [1-866-293-6349](tel:1-866-293-6349)

Contact your [Coloplast Territory Manager](#) for more information and Coloplast[®] Care enrollment forms.

Velcro[®] is a registered trademark of Velcro Industries B.V.

Coloplast Canada, 2030 Bristol Circle, Unit 117 Oakville, ON L6H 0H2 www.coloplast.ca
The Coloplast logo is a registered trademark of Coloplast A/S. © 2024-06.
All rights reserved Coloplast A/S, 3050 Humlebaek, Denmark.

