Coloplast Product Guide



Ostomy, Continence, Wound and Skin Care

Effective November, 2022

Coloplast





”I was so excited to try this technique because it would give me the independence I desperately needed. It was so easy I could hardly believe it. As long as I have my standard SpeediCath® or my SpeediCath® Compact Set (which is a catheter and bag all in one) I am good to go. Both of these products are small, require no lubrication, and are easy to use!”

Brittney
Alberta, Canada
Speedicath® User

OSTOMY

Appliances

SenSura® Mio Click

Ostomy Barrier and Pouch - Two Piece System

SenSura® Mio Click Skin Barrier 9
SenSura® Mio Click Drainable Pouch 10
SenSura® Mio Click Drainable Pouch 10
SenSura® Mio Click Closed Pouch 10
SenSura® Mio Click Urostomy Pouch 11
SenSura® Mio Click High Output Pouch 11

SenSura® Mio Flex

Ostomy Barrier and Pouch - Two Piece System

SenSura® Mio Flex Skin Barrier 12
SenSura® Mio Flex Drainable Pouch 12
SenSura® Mio Flex Closed Pouch 13
SenSura® Mio Flex Urostomy Pouch 13
SenSura® Mio Flex High Output Pouch 13

SenSura® Mio

Ostomy Barrier and Pouch - One Piece System

SenSura® Mio Drainable Pouch 13
SenSura® Mio Drainable Pouch 14
SenSura® Mio Drainable Pouch 15
SenSura® Mio Closed Pouch 16
SenSura® Mio Closed Pouch 17
SenSura® Mio One Piece Urostomy 18
SenSura® Mio One Piece High Output Pouch 18
SenSura® Mio One Piece Post-Op Pouch 18

SenSura® Mio Baby

Ostomy Barrier and Pouch - Two Piece System

SenSura® Mio Baby Flex Skin Barrier 19
SenSura® Mio Baby Flex Drainable Pouch 19
Brava® Baby Adhesive Remover Wipe 19

SenSura® Mio Kids

Ostomy Barrier and Pouch - Two Piece System

SenSura® Mio Kids **Flex** Skin Barrier 19
SenSura® Mio Kids **Flex** Drainable Pouch 19
SenSura® Mio Kids Flex Urostomy Pouch 19

Ostomy Barrier and Pouch - One Piece System

SenSura® Mio Kids Drainable Pouch 1 Piece 20
SenSura® Mio Kids Pouch 1 Piece Urostomy 20

SenSura® Click

Ostomy Barrier and Pouch- Two Piece System

SenSura® Click Skin Barrier 20
SenSura® Click Xpro Skin Barrier 21
SenSura® Click Drainable Pouch 21
SenSura® Click Closed Pouch 21
SenSura® Click Urostomy Pouch 22
SenSura® High Output Pouch 22

SenSura® Flex

Ostomy Barrier and Pouch - Two Piece System

SenSura® Flex Skin Barrier 22
SenSura® Flex Xpro Skin Barrier 23
SenSura® Flex Drainable Pouch 23
SenSura® Flex Closed Pouch 23
SenSura® Flex Urostomy Pouch 23

SenSura® - One Piece

Ostomy Barrier and Pouch - One Piece System

SenSura® Drainable Pouch with Easiclose™ 24
SenSura® Xpro Drainable Pouch 24
SenSura® Closed Pouch 25
SenSura® Urostomy Pouch 25
SenSura® One-Piece High Output Pouch 25
SenSura® Post-Op Pouch 25

Assura®

Ostomy Barrier and Pouch - Two Piece System

Assura® Skin Barrier 26
Assura® Skin Barrier - Extra Extended Wear 27
Assura® 2 pc. High Output Drainable Pouch 27
Assura® 2 pc. High Output Drainable System 27
Assura® Drainable Pouch Clamp Closure 28
Assura® Drainable Pouch with EasiClose™ 28
Assura® Closed Pouch 28
Assura® Urostomy Pouch 28

Assura® - One Piece

Ostomy Barrier and Pouch - One Piece System

Assura® Drainable Pouch Clamp Closure 29
Assura® Drainable Pouch with EasiClose™ 29
Assura® Drainable Pouch with EasiClose™ 29
Assura® Closed Pouch 30
Assura® 1-Piece Post-Op Pouch 30
Assura® Closed Pouch 30
Assura® Urostomy Pouch 31
Assura® Urostomy Pouch - Extra Extended Wear 31

Assura® - Pediatric - Two Piece

Ostomy Barrier and Pouch - Two Piece System - Pediatric

Assura® Pediatric Skin Barrier 31
Assura® Pediatric Drainable Pouch 31
Assura® Pediatric Closed Pouch 31

EasiFlex® - Pediatric

Ostomy Barrier and Pouch - Two Piece System - Pediatric

Easiflex® Pediatric Skin Barrier 32
Easiflex® Pediatric Drainable Pouch 32

Assura® - Pediatric - One Piece

Ostomy Barrier and Pouch - One Piece System - Pediatric
Assura® Pediatric Drainable Pouch 32
Assura® Pediatric Closed Pouch 32
Assura® Pediatric Urostomy Pouch 32

Accessories

Ostomy Accessories

Brava® Protective Sheets 32
Brava® Protective Rings 32
Brava® Protective Rings Wide 33
Brava® Protective Rings Extra Wide 33
Brava® Protective Ring Convex 33
Brava® Mouldable Rings 33
Brava® Skin Barrier Spray and Wipe 33
Brava® Adhesive Remover Spray and Wipe 33
Brava® Barrier Cream 33
Brava® Ostomy Strip Paste 33
Brava® Elastic Barrier Strips 33
Brava® Ostomy Powder 34
Brava® Mio Belt for All SenSura Mio Products 34
Brava® Belt for SenSura Original + Assura Products 34
Brava® Lubricating Deodorant 34
Brava® Alcohol Free Paste 34
Ostomy Paste (contains alcohol) 34
Skin Barrier Rings 34
Coloplast Prep™ Protective Liquid Skin Barrier 34
Coloplast Filtrodor® Pouch Filter 35
Coloplast Pouch Clamp 35
Ostomy Scissors 35

Hex-On®

Odour antagonist
Hex-On® 35

Assura® - Ostomy Cap

Coloplast Uro Minicap (100 ml) 35
Coloplast Stoma Cap 35
Coloplast Stoma Cap - One Piece 35

Assura® Irrigation Products

Assura® Irrigation Set 36
Assura® Irrigation Accessories 36
Assura® Irrigation Sleeves 36
Coloplast Urostomy Night Bag 37

Coloplast Fistula & Wound Management System

Pouch 37
Flexible Lids 37
Accessories 37

CONTINENCE

Hydro Intermittent Catheters

SpeediCath® Flex

Hydrophilic Intermittent Catheter
SpeediCath® Flex 41

SpeediCath® Standard

Hydrophilic Intermittent Catheter
SpeediCath® Standard Male Straight 41
SpeediCath® Standard Male Coude 41
SpeediCath® Standard Female Straight 41
SpeediCath® Standard Pediatric Straight 41

SpeediCath® Compact Male

Hydrophilic Intermittent Catheter
SpeediCath® Compact Male 42

SpeediCath® Compact Eve

Hydrophilic Intermittent Catheter
SpeediCath® Compact Eve 42

SpeediCath® Compact Female

Hydrophilic Intermittent Catheter
SpeediCath® Compact Female 42
SpeediCath® Compact Female Plus 42

SpeediCath® Compact Set

Hydrophilic Intermittent Catheter
SpeediCath® Compact Set Male 43
SpeediCath® Compact Set Female 43

Self-Cath®

Intermittent Catheter
Self-Cath® Female Straight 43
Self-Cath® Pediatric 43
Self-Cath® Extension Tube 43
Self-Cath® Male Straight 44
Self-Cath® Male Soft Straight 44
Self-Cath® Male Straight Curved Packaging 44
Self-Cath® Male Coude Tapered with Guide Stripe 44
Self-Cath® Male Coude Olive with Guide Stripe 44
Self-Cath® Male Closed System Straight 44

External

Conveen® Optima

Male External Catheter - Silicone
Conveen® Optima Standard 45
Conveen® Optima Short 45

Conveen® Security+

Male External Catheter - SEBS*
Conveen® Security+ Self Sealing 45
Conveen® Security+ Adhesive Liner 45

Clear Advantage® / Freedom Clear™

Male External Catheter - Silicone

Clear Advantage® with Aloe Vera 46
Freedom Clear™ 46
Freedom Clear™ "LS" Long Seal 46
Freedom Clear™ "SS" Sport Sheath 46

Conveen®

Male External Catheter - Latex

Conveen® Self-Sealing 46
Conveen® Adhesive Liner 46

Active Cath®

Male External Catheter - Latex

Active Cath® 47

Freedom Cath®

Male External Catheter - Latex

Freedom Clear® 47

Conveen® Active

Urine Leg Bag

Conveen® Active Leg Bag, Non Sterile 47

Conveen® Security+

Urine Leg bag

Conveen® Security+ Leg Bag, Clamp Outlet, Straps, Non Sterile 47
Conveen® Security+ Leg Bag, Clamp Outlet, Sample Port, Straps, Sterile 47
Conveen® Security+ Leg Bag, Lever Outlet, Straps, Sterile 47

Conveen® Standard

Day / Night Urine Bag

Conveen® Standard Day / Night Bag, Clamp Outlet, Non Sterile 48
Conveen® Standard Day / Night Bag, Clamp Outlet, Sample Port, Sterile 48

Conveen® Contour

Urine Leg Bag

Conveen® Contour Leg Bag, Clamp Outlet, Straps, Non Sterile 48

Conveen® Night Bag

Urine Night Bag

Conveen® Bedside Night Bag, Non-Sterile 48
Conveen® Bedside Night Bag, Sterile 48

Accessories

Conveen® Bag Hanger 48
Conveen® Leg Bag Strap 48
Conveen® Catheter retaining Strap 48

Peristeen Plus® Anal Irrigation

Peristeen Plus® Anal Irrigation System 49
Peristeen Plus® Anal Irrigation System - Small 49
Peristeen Plus® Anal Irrigation Accessory unit 49
Peristeen Plus® Anal Irrigation Accessory unit - Small 49
Peristeen Plus® Anal Irrigation Rectal Catheter 49
Peristeen Plus® Anal Irrigation Leg Strap 49
Peristeen Plus® Anal Irrigation Tube 49

WOUND

Biatain® Silicone

Foam Dressing with Silicone Adhesive

Biatain® Silicone 53

Biatain® Silicone Non-Border

Biatain® Silicone Non-Border 53

Biatain® Silicone Lite

Thin Foam Dressing with Silicone Adhesive

Biatain® Silicone Lite 53

Biatain® Silicone Ag

Foam Dressing with Silicone Adhesive and Silver

Biatain® Silicone Ag 54

Biatain®

Foam Dressing

Biatain® Adhesive 54
Biatain® Non-Adhesive 54

Biatain® Soft Hold

Foam Dressing with Gentle Adhesion

Biatain® Soft Hold 55

Biatain® Ag

Foam Dressing with Silver

Biatain® Ag Adhesive 55
Biatain® Ag Non-Adhesive 55

Biatain® Fiber

Biatain® Fiber 55

Biatain® Alginate

Calcium Alginate Dressing

Biatain® Alginate 56

Biatain® Alginate Ag

Calcium Alginate Dressing with Silver

Biatain® Alginate Ag 56

Biatain® Ibu

Foam Dressing with Ibuprofen

Biatain® Ibu Non-Adhesive 56
Biatain® Ibu Soft-Hold 56

Biatain® Super

Hydrocapillary Pad Dressing
Biatain® Super Adhesive 57
Biatain® Super Non-Adhesive 57

Comfeel® Plus

Hydrocolloid Dressing
Comfeel® Plus 57

Comfeel® Plus Transparent

Hydrocolloid Dressing
Comfeel® Plus Transparent 57

Triad™

Hydrophilic Wound Dressing
Triad™ 58

Purilon® Gel

Hydrogel with natural ingredients
Purilon® Gel 58

Woun'Dres®

Collagen Hydrogel
Woun'Dres® Collagen Hydrogel 58

Sea-Clens®

Wound Cleanser
Sea-Clens® Wound Cleanser 58

SKIN

Sproam®

No-Rinse Spray or Foam Cleanser
Sproam® 61

Bedside-Care® Sensitive Skin

No-Rinse Foam Cleanser
Bedside-Care® Sensitive Skin 61

Bedside-Care® Foam

No-Rinse Foam Cleanser
Bedside-Care® 61

Bedside-Care®

No-Rinse Spray Cleanser
Bedside-Care® 61

EasiCleanse® Bath

No-Rinse Wash Cloth
EasiCleanse® Bath 62

Gentle Rain® Extra Mild

Rinse Cleanser
Gentle Rain® Extra Mild 62

Gentle Rain® Extra Mild Fragrance Free

Rinse Cleanser
Gentle Rain® Extra Mild Fragrance Free 62

Atrac-Tain®

Moisturizer
Atrac-Tain® 62

Sween® 24

Moisturizer
Sween® 24 63

Sween Cream®

Moisturizer
Sween Cream® 63

Sween® Lotion Unscented

Moisturizer
Sween® Lotion Unscented 63

Critic-Aid® Clear

Skin Protectant
Critic-Aid® Clear 64

Baza® Protect II

Skin Protectant
Baza® Protect II 64

Baza Cleanse and Protect®

Skin Protectant
Baza Cleanse and Protect® 64

InterDry®

Moisture-wicking textile with antimicrobial silver complex
InterDry® 65

Isagel®

Hand Sanitizer Gel
Isagel® 65

Sween® Body Powder

All-body powder
Sween® Body Powder 65

Accessories

Coloplast® Wall Dispenser for Bag 66
Coloplast® Wall Dispenser Drip Tray 66
Embrace™ Dispenser Bracket 66

Ostomy

SenSura® Mio Click

Ostomy Barrier and Pouch - Two Piece System

Click mechanical coupling with a lock ring for security. Provides a tight, secure fit to individual body shapes with optimal discretion. Elastic adhesive, neutral grey textile and full-circle filter. All 2pc Mio Click Flanges have 4 belt attachment points.

SenSura® Mio Click Skin Barrier


Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	40mm (1 9/16")	10-35mm (3/8 - 1 3/8")	10502	5
		50mm (2")	10-45mm (3/8 - 1 3/4")	10512	5
		60mm (2 3/8")	10-55mm (3/8 - 2 1/8")	10522	5
		70mm (2 3/4")	10-65mm (3/8 - 2 1/2")	10531	5
Pre-Cut	Non-Convex	40mm (1 9/16")	20mm (3/4")	10503	5
			22mm (7/8")	10506 *New!*	5
			25mm (1")	10504	5
			28mm (1 1/8")	10505 *New!*	5
		50mm (2")	30mm (1 1/8")	10513	5
			32mm (1 1/4")	10514 *New!*	5
			35mm (1 3/8")	10515	5
			38mm (1 1/2")	10516 *New!*	5
		60mm (2 3/8")	40mm (1 9/16")	10523 *New!*	5
			45mm (1 3/4")	10524 *New!*	5
			50mm (2")	10525 *New!*	5
Cut-to-Fit	Convex Light	40mm (1 9/16")	10-20mm (3/8 - 3/4")	16901	5
		50mm (2")	15-30mm (5/8 - 1 1/8")	16911	5
		60mm (2 3/8")	15-40mm (5/8 - 1 9/16")	16921	5
		70mm (2 3/4")	15-53mm (5/8"- 2 1/16")	16931	5
Pre-Cut	Convex Light	40mm (1 9/16")	18mm (3/4")	16903	5
		50mm (2")	21mm (13/16")	16913 *New!*	5
			25mm (1")	16914	5
			28mm (1 1/8")	16915	5
		60mm (2 3/8)	31mm (1 1/4")	16923	5
			35mm (1 3/8")	16924	5
		70mm (2 3/4")	38mm (1 1/2")	16933 *New!*	5
			41mm (1 5/8")	16934 *New!*	5
Cut-to-Fit	Convex Deep	40mm (1 9/16")	10-20mm (3/8 - 3/4")	16941	5
		50mm (2")	15-30mm (5/8 - 1 1/8")	16951	5
		60mm (2 3/8")	15-40mm (5/8 - 1 9/16")	16961	5
		70mm (2 3/4")	15-53mm (5/8 - 2 1/16")	16971	5
Pre-Cut	Convex Deep	40mm (1 9/16")	18mm (3/4")	16942	5
		50mm (2")	21mm (13/16")	16953 *New!*	5
			25mm (1")	16954	5
			28mm (1 1/8")	16955	5
		60mm (2 3/8")	31mm (1 1/4")	16963 *New!*	5
			35mm (1 3/8")	16964 *New!*	5
		70mm (2 3/4")	38mm (1 1/2")	16973 *New!*	5
			41mm (1 5/8")	16974 *New!*	5





SenSura® Mio Click Skin Barrier

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Concave	40mm (1 9/16")	10-35mm (3/8 - 1 3/8")	18500	5
		50mm (2")	10-45mm (3/8 - 1 3/4")	18510	5
		60mm (2 3/8")	10-55mm (3/8 - 2 1/8")	18520	5
		70mm (2 3/4")	10-65mm (3/8 - 2 1/2")	18530	5
Pre-Cut	Concave	40mm (1 9/16")	20mm (3/4")	18501 *New!*	5
		50mm (2")	25mm (1")	18511	5
			30mm (1 1/8")	18512	5
			35mm (1 3/8")	18513	5
		60mm (2 3/8")	30mm (1 1/8")	18521	5
			35mm (1 3/8")	18522	5
			40mm (1 9/16")	18523	5
			45mm (1 3/4")	18524	5
		70mm (2 3/4")	45mm (1 3/4")	18531	5
			50mm (2")	18532	5



SenSura® Mio Click Drainable Pouch with Easiclose™ WIDE-Outlet

Size	Colour	Coupling Size		Product Code	Units
Midi 26cm (10 1/4")	Neutral grey with inspection window	40mm (1 9/16")		11461	10
		50mm (2")		11471	10
Maxi 28cm (11")	Neutral grey with inspection window	40mm (1 9/16")		11462	10
		50mm (2")		11472	10
		60mm (2 3/8")		11482	10
		70mm (2 3/4")		11486	10
		40mm (1 9/16")	Narrow outlet with Ileo filter	11432	10
		50mm (2")	Narrow outlet with Ileo filter	11443	10
Maxi 28cm (11")	Transparent	40mm (1 9/16")		11463	10
		50mm (2")		11473	10
		50mm (2")	No filter	11474	10
		60mm (2 3/8")		11483	10
		60mm (2 3/8")	No filter	11484	10
		70mm (2 3/4")		11487	10





SenSura® Mio Click Closed Pouch

Size	Colour	Coupling Size	Product Code	Units
Midi 18cm (7")	Neutral grey with inspection window	40mm (1 9/16")	11401	30
		50mm (2")	11412	30
		60mm (2 3/8")	11421	30
Maxi 22cm (8 1/2")	Neutral grey with inspection window	40mm (1 9/16")	11402	30
		50mm (2")	11413	30
		60mm (2 3/8")	11422	30
		70mm (2 3/4")	11425	30



SenSura® Mio Click Urostomy Pouch

Size	Colour	Coupling Size	Product Code	Units
Midi 24cm (9 1/2")	Neutral grey	50mm (2")	11496	10
Maxi 26cm (10 1/4")	Neutral grey	40mm (1 9/16")	11491	10
		50mm (2")	11497	10
		60mm (2 3/8")	11499	10
	Transparent	40mm (1 9/16")	11492	10
		50mm (2")	11498	10
		60mm (2 3/8")	11500	10



SenSura® Mio Click High Output Pouch

Outlet	Colour	Coupling Size	Volume	Product Code	Units
Wide-lumen Hard Outlet 32cm (12 11/16")	Neutral grey with inspection window	40mm (1 9/16")	1160 mL	18600	10
		50mm (2")	1120 mL	18601	10
		60mm (2 3/8")	1080 mL	18602	10
		70mm (2 3/4")	1040 mL	18603	10
	Transparent	40mm (1 9/16")	1160 mL	18610	10
		50mm (2")	1120 mL	18611	10
		60mm (2 3/8")	1080 mL	18612	10
		70mm (2 3/4")	1040 mL	18613	10
Mini Pouch Wide-lumen Hard Outlet 24cm (9 1/2")	Neutral grey with inspection window	40mm (1 9/16")	650 mL	18620	10
		50mm (2")	581 mL	18621	10
		60mm (2 3/8")	546 mL	18622	10
Soft Outlet Universal Connector 32cm (12 11/16")	Neutral grey with inspection window	50mm (2")	1120 mL	18630	10
		60mm (2 3/8")	1080 mL	18631	10
		70mm (2 3/4")	1040 mL	18632	10
	Transparent	50mm (2")	1120 mL	18640	10
		60mm (2 3/8")	1080 mL	18641	10
		70mm (2 3/4")	1040 mL	18642	10



Pouches with Wide Lumen Hard Outlet will connect to Coloplast drainage bag 14010

SenSura® Mio **Flex**

Ostomy Barrier and Pouch - Two Piece System

Flex adhesive coupling allows flexibility and freedom of movement. Provides a tight, secure fit to individual body shapes with optimal discretion. Elastic adhesive, neutral grey textile and full-circle filter. Only Mio Flex Convex Flanges have 4 belt attachment points.

SenSura® Mio Flex Skin Barrier

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	35mm (1 3/8")	10-33mm (3/8 - 1 1/4")	10551	10
		50mm (2")	10-48mm (3/8 - 2")	10561	10
		70mm (2 3/4")	10-68mm (3/8 - 2 3/4")	10571	10
		90mm (3 1/2")	10-88mm (3/8" - 3 1/2")	18390 *New!*	10
Pre-Cut	Non-Convex	50mm (2")	25mm (1")	10562	10
			28mm (1 1/8")	10566 *New!*	10
			30mm (1 1/4")	10563	10
			32mm (1 1/4")	10567 *New!*	10
			35mm (1 3/8")	10564	10
			38mm (1 1/2")	10568 *New!*	10
			40mm (1 9/16")	10565 *New!*	10
		70mm (2 3/4")	45mm (1 3/4")	10572 *New!*	5
			50mm (2")	10573 *New!*	5
Cut-to-Fit	Convex Light	35mm (1 3/8")	10-20mm (3/8 - 3/4")	16471	5
		50mm (2")	15-30mm (5/8 - 1 1/8")	16481	5
			15-40mm (5/8 - 1 9/16")	16483	5
		70mm (2 3/4")	15-53mm (5/8 - 2 1/16")	16491	5
Pre-Cut	Convex Light	35mm (1 3/8")	18mm (3/4")	16472	5
		50mm (2")	21mm (13/16")	16484 *New!*	5
			25mm (1")	16485	5
			28mm (1 1/8")	16486	5
			31mm (1 1/4")	16487	5
			35mm (1 3/8")	16488	5
		70mm (2 3/4")	38mm (1 1/2")	16493 *New!*	5
			41mm (1 5/8")	16494 *New!*	5



SenSura® Mio Flex Drainable Pouch with Easiclose™ WIDE-Outlet

Size	Colour	Coupling Size		Product Code	Units
Midi 26cm (10 1/4")	Neutral grey with inspection window	50mm (2")		12271	10
Maxi 28cm (11")	Neutral grey with inspection window	35mm (1 3/8")	*Narrow outlet with Ileo filter*	12232	10
		35mm (1 3/8")		12262	10
		50mm (2")	*Narrow outlet with Ileo filter*	12243	10
		50mm (2")		12272	10
		70mm (2 3/4")		12282	10
	Transparent	50mm (2")	*No filter*	12274	10
		70mm (2 3/4")	*No filter*	12284	10
		90mm (3 1/2")		18389 *New!*	10



SenSura® Mio Flex Closed Pouch

Size	Colour	Coupling Size	Product Code	Units
Midi 18cm (7")	Neutral grey with inspection window	35mm (1 3/8")	12201	30
		50mm (2")	12212	30
Maxi 22cm (8 1/2")	Neutral grey with inspection window	35mm (1 3/8")	12202	30
		50mm (2")	12213	30
		70mm (2 3/4")	12222	30
		90mm (3 1/2")	18386 *New!*	10



SenSura® Mio Flex Urostomy Pouch

Size	Colour	Coupling Size	Product Code	Units
Maxi 26cm (10 1/4")	Neutral grey	35mm (1 3/8")	12291	10
		50mm (2")	12297	10
	Transparent	35mm (1 3/8")	12292	10
		50mm (2")	12298	10



SenSura® Mio Flex High Output Pouch

Outlet	Colour	Coupling Size	Volume	Product Code	Units
Wide-lumen Hard Outlet 32cm (12 11/16")	Neutral grey with inspection window	50mm (2")	980 mL	18650	10
		70mm (2 3/4")	906 mL	18651	10
	Transparent	50mm (2")	980 mL	18655	10
		70mm (2 3/4")	906 mL	18656	10
Soft outlet 32 cm (12 11/16")	Neutral grey with inspection window	50mm (2")	1080 mL	18652 *New!*	10
	Transparent	50mm (2")	1080 mL	18657 *New!*	10



Pouches with Wide Lumen Hard Outlet will connect to Coloplast drainage bag 14010

SenSura® Mio

Ostomy Barrier and Pouch - One Piece System

Elastic adhesive, neutral grey textile and full-circle filter. Provides a tight, secure fit to individual body shapes with optimal discretion. Only 1pc Convex systems have 4 belt attachment points.

SenSura® Mio Drainable Pouch with Easiclose™ WIDE-Outlet

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 26cm (10 1/4")	Cut-to-Fit	Non-Convex	Neutral grey with inspection window	10-45mm (3/8 - 1 3/4")	10461	10
Maxi 28cm (11")	Cut-to-Fit	Non-Convex	Neutral grey with inspection window	10-55mm (3/8 - 2 1/8") *Narrow outlet with Ileo filter*	10451	10
				10-55mm (3/8 - 2 1/8")	10471	10
			Transparent	10-55mm (3/8 - 2 1/8")	10481	10
				10-55mm (3/8 - 2 1/8")	10489 *No filter*	10
	Pre-Cut	Non-Convex	Neutral grey with inspection window	25mm (1")	10472	10
				28mm (1 1/8")	10473 *New!*	10
				30mm (1 1/8")	10474	10
				32mm (1 1/4")	10475 *New!*	10
				35mm (1 3/8")	10476	10
				38mm (1 1/2")	10480 *New!*	10
				40mm (1 9/16")	10479	10
				45mm (1 3/4")	10477	10
				50mm (2")	10478	10



SenSura® Mio Drainable Pouch with Easiclose™ WIDE-Outlet

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 28cm (11")	Pre-Cut	Non-Convex	Transparent	25mm (1")	10482	10
				30mm (1 1/8")	10483	10
				35mm (1 3/8")	10484	10
				40mm (1 9/16")	10487	10
Maxi 28cm (11")	Cut-to-Fit	Convex Soft	Neutral grey with inspection window	15-33mm (5/8 - 1 1/4")	16705	10
				10-50mm (3/8 - 2")	16706	10
			Transparent	15-33mm (5/8 - 1 1/4")	16715	10
				10-50mm (3/8 - 2")	16716	10
	Pre-Cut	Convex Soft	Neutral grey with inspection window	21mm (7/8")	13671 *New!*	10
				25 mm (1")	13672 *New!*	10
				28 mm (1 1/8")	13673 *New!*	10
				31mm (1 1/4")	13674 *New!*	10
				35 mm (1 3/8")	13675 *New!*	10
				38 mm (1 1/2")	13676 *New!*	10
	Cut-to-Fit	Convex Light	Neutral grey with inspection window	10-23mm (3/8 - 7/8")	16725	10
				15-33mm (5/8 - 1 1/4")	16726	10
				10-43mm (3/8 - 1 3/4")	16733	10
			Transparent	10-23mm (3/8 - 7/8")	16737	10
				15-33mm (5/8 - 1 1/4")	16738	10
				10-43mm (3/8 - 1 3/4")	16739	10
				10-43mm (3/8 - 1 3/4")	16748 *No filter*	10
	Pre-Cut	Convex Light	Neutral grey with inspection window	21mm (13/16")	16730 *New!*	10
				25mm (1")	16731	10
				28mm (1 1/8")	16732	10
				31mm (1 1/4")	16727	10
				35mm (1 3/8")	16734	10
				38mm (1 1/2")	16735 *New!*	10
	Cut-to-Fit	Convex Deep	Neutral grey with inspection window	10-23mm (3/8 - 7/8")	16755	10
				15-33mm (5/8 - 1 1/4")	16756	10
				10-43mm (3/8 - 1 3/4")	16768	10
			Transparent	10-43mm (3/8 - 1 3/4")	16767	10
				10-43mm (3/8 - 1 3/4")	16778 *No filter*	10
	Pre-Cut	Convex Deep	Neutral grey with inspection window	21mm (7/8")	16760	10
				25mm (1")	16761	10
				28mm (1 1/8")	16762	10
				31mm (1 1/4")	16757	10
				35mm (1 3/8")	16763	10
				38mm (1 1/2")	16764 *New!*	10
			Transparent	21mm (13/16")	16770 *New!*	10
				25mm (1")	16771	10
				28mm (1 1/8")	16772	10
				31mm (1 1/4")	16773	10



SenSura® Mio Drainable Pouch with Easiclose™ WIDE-Outlet

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 26cm (10 1/4")	Cut-to-fit	Concave	Neutral grey with inspection window	10-40mm (3/8 - 1 9/16")	18300	10
				10-50mm (3/8 - 2")	18310	10
			Neutral grey with inspection window	10-50mm (3/8 - 2")	18312	10
				10-60mm (3/8 - 2 3/8")	18330	10
			Transparent	10-40mm (3/8 - 1 9/16")	18305	10
	Cut-to-fit	Concave		10-50mm (3/8 - 2")	18320 *No filter*	10
				10-50mm (3/8 - 2")	18321	10
				10-60mm (3/8 - 2 3/8")	18340 *No filter*	10
				10-60mm (3/8 - 2 3/8")	18341	10
Maxi 28cm (11")	Pre-Cut	Concave	Neutral grey with inspection window	25-40mm (1 - 1 9/16")	18301	10
				30-40mm (1 1/8 - 1 9/16")	18302	10
				25-50mm (1 - 2")	18313	10
				30-50mm (1 1/8 - 2")	18314	10
				35-50mm (1 3/8 - 2")	18315	10
				40-60mm (1 9/16 - 2 3/8")	18331	10
				45-60mm (1 3/4 - 2 3/8")	18332	10
				50-60mm (2" - 2 3/8")	18333	10
			Transparent	25-50mm (1 - 2")	18322	10
				30-50mm (1 1/8 - 2")	18323	10
				35-50mm (1 3/8 - 2")	18324	10
				40-60mm (1 9/16 - 2 3/8")	18342	10
				45-60mm (1 3/4 - 2 3/8")	18343	10
				50-60mm (2" - 2 3/8")	18344	10



SenSura® Mio Closed Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Mini 14cm (5 1/2")	Cut-to-fit	Non-Convex	Neutral grey with inspection window	15-45mm (5/8 - 1 3/4")	10870	30
Midi 18cm (7")	Cut-to-fit	Non-Convex	Neutral grey with inspection window	15-45mm (5/8 - 1 3/4")	10871	30
	Pre-Cut	Non-Convex	Neutral grey with inspection window	25mm (1")	10872	30
				30mm (1 1/8")	10873	30
				35mm (1 3/8")	10874	30
	Cut-to-Fit	Convex Soft	Neutral grey with inspection window	15-33mm (5/8 - 1 1/4")	16301	10
		Convex Light	Neutral grey with inspection window	15-33mm (5/8 - 1 1/4")	16322	10
	Pre-Cut	Convex Light	Neutral grey with inspection window	31mm (1 1/4")	16323	10
Maxi 22cm (8 1/2")	Cut-to-Fit	Non-Convex	Neutral grey with inspection window	15-45mm (5/8 - 1 3/4")	10881	30
				40-55mm (1 9/16 - 2 1/8")	10885	30
	Pre-Cut	Non-Convex	Neutral grey with inspection window	30mm (1 1/8")	10883	30
				35mm (1 3/8")	10884	30
	Cut-to-Fit	Convex Soft	Neutral grey with inspection window	15-33mm (5/8 - 1 1/4")	16305	10
				10-50mm (3/8 - 2")	16306	10
			Transparent	10-50mm (3/8 - 2")	16315	10
	Pre-Cut	Convex Soft	Neutral grey with inspection window	25mm (1")	16310	10
				28mm (1 1/8")	16311	10
				31mm (1 1/4")	16312	10
				35mm (1 3/8")	16313	10
				38 mm (1 1/2")	13684 *New!*	10
				41 mm (1 5/8")	13685 *New!*	10
	Cut-to-Fit	Convex Light	Neutral grey with inspection window	10-23mm (3/8 - 7/8")	16325	10
				15-33mm (5/8 - 1 1/4")	16326	10
				10-43mm (3/8 - 1 3/4")	16332	10
			Transparent	10-43mm (3/8 - 1 3/4")	16335	10
	Pre-Cut	Convex Light	Neutral grey with inspection window	25mm (1")	16330	10
				28mm (1 1/8")	16331	10
				31mm (1 1/4")	16327	10
				35mm (1 3/8")	16333	10
				41 mm (1 5/8")	13683 *New!*	10
	Cut-to-Fit	Convex Deep	Neutral grey with inspection window	10-23mm (3/8 - 7/8")	16345	10
				15-33mm (5/8 - 1 1/4")	16346	10
				10-43mm (3/8 - 1 3/4")	16352	10
	Pre-Cut	Convex Deep	Neutral grey with inspection window	25mm (1")	16350	10
				28mm (1 1/8")	16351	10
				31mm (1 1/4")	16347	10
				35mm (1 3/8")	16353	10



SenSura® Mio Closed Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 18cm (7")	Cut-to-Fit	Concave	Neutral grey with inspection window	10-40mm (3/8 - 1 9/16")	18151	10
				10-50mm (3/8 - 2")	18160	10
	Pre-Cut	Concave		25mm (1")	18161	10
				30mm (1 1/8")	18162	10
				35mm (1 3/8")	18163	10
Maxi 22cm (8 1/2")	Cut-to-Fit	Concave	Neutral grey with inspection window	10-40mm (3/8 - 1 9/16")	18152	10
				10-50mm (3/8 - 2")	18170	10
				10-60mm (3/8 - 2 3/8")	18180	10
			Transparent	10-60mm (3/8 - 2 3/8")	18185	10
	Pre-Cut	Concave	Neutral grey with inspection window	25mm (1")	18171	10
				30mm (1 1/8")	18172	10
				35mm (1 3/8")	18173	10
				40mm (1 9/16")	18181	10
				45mm (1 3/4")	18182	10



SenSura® Mio One Piece Urostomy

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 26cm (10 1/4")	Cut-to-Fit	Non-Convex	Neutral grey	10-45mm (3/8 - 1 3/4")	10582	10
			Transparent	10-45mm (3/8 - 1 3/4")	10586	10
	Pre-Cut	Non-Convex	Neutral grey	20mm (3/4")	10583	10
				25mm (1")	10584	10
				30mm (1 1/8")	10585	10
Midi 24cm (9 1/2")	Cut-to-Fit	Convex Soft	Neutral grey	10-33mm (3/8 - 1 5/16")	16801 *New!*	10
Maxi 26cm (10 1/4")	Cut-to-Fit	Convex Soft	Neutral grey	15-33mm (5/8 - 1 1/4")	16805	10
				25 mm (1")	13680 *New!*	10
				10-50 mm (3/8" - 2")	13678 *New!*	10
			Transparent	10-33 mm (3/8" - 1 5/16")	13679 *New!*	10
				10-50mm (3/8 - 2")	16810	10
Maxi 26cm (10 1/4")	Cut-to-Fit	Convex Light	Neutral grey	10-23mm (3/8 - 7/8")	16825	10
				15-33mm (5/8 - 1 1/4")	16826	10
			Transparent	10-23mm (3/8 - 7/8")	16835	10
				15-33mm (5/8 - 1 1/4")	16836	10
				10-43mm (3/8 - 1 3/4")	16837	10
Maxi 26cm (10 1/4")	Pre-Cut	Convex Light	Neutral grey	31mm (1 1/4")	16827	10
			Transparent	21mm (7/8")	16840	10
				25mm (1")	16841	10
				28mm (1 1/8")	16842	10
	Cut-to-Fit	Convex Deep	Neutral grey	10-23mm (3/8 - 7/8")	16855	10
				15-33mm (5/8 - 1 1/4")	16856	10
			Transparent	10-23mm (3/8 - 7/8")	16865	10
				15-33mm (5/8 - 1 1/4")	16866	10
				10-43mm (3/8 - 1 3/4")	16867	10
	Pre-Cut	Convex Deep	Neutral grey	31mm (1 1/4")	16857	10
			Transparent	21mm (7/8")	16870	10
				25mm (1")	16871	10
				28mm (1 1/8")	16872	10



SenSura® Mio One Piece Urostomy

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 26cm (10 1/4")	Cut-to-Fit	Concave	Neutral grey	10-40mm (3/8- 1 9/16")	18401	10
			Transparent	10-40mm (3/8- 1 9/16")	18405	10
	Pre-Cut	Concave	Neutral grey	20mm (3/4")	18402	10
				25mm (1")	18403	10
				30mm (1 1/8")	18404	10



SenSura® Mio One Piece High Output Pouch

Cut/Outlet	Convexity	Colour	Volume	Stoma Size	Product Code	Units
Cut-to-fit/ Soft Outlet Universal Connector 28cm (11 3/16")	Non-Convex	Neutral grey with inspection window	630 mL	10-70mm (3/8- 2 3/4")	18660	10
			725 mL	10-100mm (3/8- 4")	18661	10
		Transparent	630 mL	10-70mm (3/8- 2 3/4")	18665	10
			725 mL	10-100mm (3/8- 4")	18666	10
	Convex Soft	Neutral grey with inspection window	630 mL	10-50mm (3/8- 2")	18670	10
		Transparent	630 mL	10-50mm (3/8- 2")	18671	10
	Convex Light	Neutral grey with inspection window	630 mL	10-43mm (3/8- 1 3/4")	18675	10
		Transparent	630 mL	10-43mm (3/8- 1 3/4")	18676	10



Does not connect to Coloplast drainage bag 14010

SenSura® Mio One Piece Post-Op Pouch

Size	Cut	Convexity	Volume	Stoma Size	Product Code		Units
28cm (10 15/18")	Cut-to-Fit	Transparent	680 mL	10-70mm (3/8- 2 3/4")	18680	*Sterile*	5
29cm (11 1/2")	Cut-to-Fit	Transparent	725 mL	10-100mm (3/8- 4")	18681	*Sterile*	5
28cm (10 15/18")	Cut-to-Fit	Transparent	680 mL	10-70mm (3/8- 2 3/4")	18690		5
29cm (11 1/2")	Cut-to-Fit	Transparent	725 mL	10-100mm (3/8- 4")	18691		5



SenSura® Mio Baby

Ostomy Barrier and Pouch - Two Piece System

SenSura® Mio Baby Flex Skin Barrier

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	35mm (1 3/8")	0-35mm	18705	5



SenSura® Mio Baby Flex Drainable Pouch

Size	Colour	Coupling Size	Product Code	Units
20cm (7 3/4")	Transparent	35mm (1 3/8")	18700	10



Brava® Baby Adhesive Remover

	Product Code	Units
Single use	12014	30



SenSura® Mio Kids

Ostomy Barrier and Pouch - Two Piece System

SenSura® Mio Kids Flex Skin Barrier

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	27 mm	0-27 mm	18735	5
		35 mm	0-35 mm	18747	5



SenSura® Mio Kids Flex Drainable Pouch

Cut	Colour	Coupling Size	Product Code	Units
Mini 16cm (6 1/8")	Neutral grey	27 mm	18720	10
	Transparent	27 mm	18721	10
Midi 17cm (6 3/4")	Neutral grey	27 mm	18722	10
	Transparent	27 mm	18723	10
	Neutral grey	35 mm	18724	10
	Transparent	35 mm	18725	10





SenSura® Mio Kids Flex Urostomy Pouch

Cut	Colour	Coupling Size	Stoma Size	Product Code	Units
Midi 17cm (6 3/4")	Neutral grey	27 mm		18730	10
		35 mm		18731	10



Ostomy Barrier and Pouch - One Piece System

SenSura® Mio Kids Drainable Pouch 1 Piece

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 17cm (6 3/4")	Cut-to-fit	Non-Convex	Neutral grey	0-35mm	18710	10
			Transparent	0-35mm	18711	10



SenSura® Mio Kids Pouch 1 Piece Urostomy

Cut	Convexity	Colour	Stoma Size	Product Code	Units
Cut-to-fit 19cm (7 1/4")	Non-Convex	Transparent	10-35mm	18712	10
		Neutral grey	10-35mm	18713	10




SenSura® Click

Ostomy Barrier and Pouch- Two Piece System

Click mechanical coupling with a lock ring for security. Provides a tight, secure fit to individual body shapes. All 2pc SenSura Click Flanges have 2 belt attachment points.

SenSura® Click Skin Barrier

Cut	Convexity	Coupling Size		Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	40mm (1 9/16")	●	10-35mm (3/8 - 1 3/8")	10011	5
		50mm (2")	●	10-45mm (3/8 - 1 3/4")	10021	5
		60mm (2 3/8")	●	10-55mm (3/8 - 2 1/8")	10031	5
		70mm (2 3/4")	●	10-65mm (3/8 - 2 1/2")	10041	5
	Convex Light	40mm (1 9/16")	●	15-23mm (5/8 - 7/8")	11011	5
		50mm (2")	●	15-33mm (5/8 - 1 1/4")	11021	5
		60mm (2 3/8")	●	15-43mm (5/8 - 1 3/4")	11031	5
		70mm (2 3/4")	●	15-53mm (5/8 - 2 1/16")	11041	5
Pre-Cut	Non-Convex	40mm (1 9/16")	●	20mm (3/4")	10012	5
				25mm (1")	10013	5
		50mm (2")	●	30mm (1 1/8")	10022	5
				35mm (1 3/8")	10023	5
		60mm (2 3/8")	●	40mm (1 9/16")	10032	5
	Convex Light	40mm (1 9/16")	●	21mm (7/8")	11012	5
		50mm (2")	●	25mm (1")	11022	5
				28mm (1 1/8")	11023	5
				31mm (1 1/4")	11024	5
		60mm (2 3/8")	●	35mm (1 3/8")	11032	5



SenSura® Click Xpro Skin Barrier
(Engineered to offer extra protection for liquid ileostomies and urostomies)

Cut	Convexity	Coupling Size		Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	40mm (1 9/16")		10-35mm (3/8 - 1 3/8")	10015	5
		50mm (2")		10-45mm (3/8 - 1 3/4")	10025	5
		60mm (2 3/8")		10-55mm (3/8 - 2 1/8")	10035	5
		70mm (2 3/4")		10-65mm (3/8 - 2 1/2")	10045	5
	Convex Light	40mm (1 9/16")		15-23mm (5/8 - 7/8")	11015	5
		50mm (2")		15-33mm (5/8 - 1 1/4")	11025	5
		60mm (2 3/8")		15-43mm (5/8 - 1 3/4")	11035	5
		70mm (2 3/4")		15-53mm (5/8 - 2 1/16")	11045	5
Pre-Cut	Non-Convex	40mm (1 9/16")		20mm (3/4")	10016	5
				25mm (1")	10017	5
		50mm (2")		30mm (1 1/8")	10026	5
				35mm (1 3/8")	10027	5
		60mm (2 3/8")		40mm (1 9/16")	10036	5
	Convex Light	40mm (1 9/16")		21mm (7/8")	11016	5
		50mm (2")		25mm (1")	11026	5
				28mm (1 1/8")	11027	5
				31mm (1 1/4")	11028	5
		60mm (2 3/8")		35mm (1 3/8")	11036	5



SenSura® Click Drainable Pouch with Easiclose™ WIDE-Outlet

Pouch Size	Pouch Colour	Coupling Size			Product Code	Units
Midi 27cm (10 1/2")	Opaque	40mm (1 9/16")			11114	10
		50mm (2")			11115	10
Maxi 30cm (12")	Opaque	40mm (1 9/16")			11124	10
		50mm (2")			11125	10
		50mm (2")		No filter	11135	10
		60mm (2 3/8")			11126	10
		60mm (2 3/8")		No filter	11136	10
		70mm -2 3/4"			11127	10
	Transparent	40mm (1 9/16")			11194	10
		50mm (2")			11195	10
		50mm (2")		No filter	11185	10
		60mm (2 3/8")			11196	10
		60mm (2 3/8")		No filter	11186	10
		70mm -2 3/4"			11197	10



SenSura® Click Closed Pouch

Pouch Size	Pouch Colour	Coupling Size			Product Code	Units
Midi 18cm (7")	Opaque	40mm (1 9/16")			10154	30
		50mm (2")			10155	30
		60mm (2 3/8")			10156	30
Maxi 22cm (8 1/2")	Opaque	40mm (1 9/16")			10164	30
		50mm (2")			10165	30
		60mm (2 3/8")			10166	30
		70mm (2 3/4")			10167	30



SenSura® Click Urostomy Pouch

Pouch Size	Pouch Colour	Coupling Size		Product Code	Units
Midi 24cm (9 1/4")	Opaque	40mm (1 9/16")	●	11841	10
		50mm (2")	●	11842	10
	Transparent	40mm (1 9/16")	●	11851	10
		50mm (2")	●	11852	10
Maxi 26cm (10 1/4")	Opaque	40mm (1 9/16")	●	11844	10
		50mm (2")	●	11845	10
		60mm (2 3/8")	●	11846	10
	Transparent	40mm (1 9/16")	●	11854	10
		50mm (2")	●	11855	10
		60mm (2 3/8")	●	11856	10



SenSura® High Output Pouch

Colour	Coupling Size		Volume	Product Code	Units
Opaque	40mm (1 9/16")	●	1160 mL	19040	10
	50mm (2")	●	1120 mL	19041	10
	60mm (2 3/8")	●	1080 mL	19042	10
Transparent	40mm (1 9/16")	●	1160 mL	19050	10
	50mm (2")	●	1120 mL	19051	10
	60mm (2 3/8")	●	1080 mL	19052	10



Will connect to Coloplast Drainage Bag 14010

SenSura® Flex

Ostomy Barrier and Pouch - Two Piece System

Low profile adhesive coupling system that provide flexibility and freedom of movement. Double layer adhesive for secure, safe, long lasting wear. Drainable pouch with velcro closure system. Only SenSura Flex Convex Flanges have 2 belt attachment points.

SenSura® Flex Skin Barrier

Cut	Convexity	Coupling Size		Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	35mm (1 3/8")	●	10-33mm (3/8 - 1 1/4")	10101	10
		50mm (2")	●	10-48mm (3/8 - 2")	10102	10
		70mm (2 3/4")	●	10-68mm (3/8 - 2 3/4")	10103	10
		90mm - 3 1/2"	●	10-88mm (3/8 - 3 1/2")	10108	5
	Convex light	35mm (1 3/8")	●	15-23mm (5/8 - 7/8")	11301	5
		50mm (2")	●	15-33mm (5/8 - 1 1/4")	11302	5
				15-43mm (5/8 - 1 3/4")	11303	5
		70mm (2 3/4")	●	15-56mm (5/8 - 2 1/8")	11307	5
Pre-Cut	Non-Convex	35mm (1 3/8")	●	21mm (7/8")	10111	10
				28mm (1 1/8")	10113	10
		50mm (2")	●	25mm (1")	10114	10
				31mm (1 1/4")	10116	10
				35mm (1 3/8")	10117	10
	Convex light	35mm (1 3/8")	●	21mm (7/8")	11312	5
		50mm (2")	●	25mm (1")	11313	5
				28mm (1 1/8")	11314	5
				31mm (1 1/4")	11315	5
				35mm (1 3/8")	11316	5



SenSura® Flex Xpro Skin Barrier

Cut	Convexity	Coupling Size		Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	35mm (1 3/8")		10-33mm (3/8 - 1 1/4")	10104	10
		50mm (2")		10-48mm (3/8 - 2")	10105	10
		70mm (2 3/4")		10-68mm (3/8 - 2 3/4")	10106	10
	Convex light	35mm (1 3/8")		15-23mm (5/8 - 7/8")	11304	5
		50mm (2")		15-33mm (5/8 - 1 1/4")	11305	5
				15-43mm (5/8 - 1 3/4")	11306	5
		70mm (2 3/4")		15-56mm (5/8 - 2 1/8")	11308	5



SenSura® Flex Drainable Pouch with Easiclose™ WIDE-Outlet

Size	Colour	Coupling Size			Product Code	Units
Midi 27cm (10 1/2")	Opaque	50mm (2")			11505	20
Maxi 30cm (12")	Opaque	35mm (1 3/8")			11514	20
		50mm (2")			11515	20
		50mm (2")		No filter	11615	20
		70mm (2 3/4")			11518	20
	Transparent	50mm (2")		No filter	11612	20
		90mm (3 1/2")		No filter	11690	20



SenSura® Flex Closed Pouch

Size	Colour	Coupling Size			Product Code	Units
Midi 18cm (7")	Opaque	35mm (1 3/8")			10904	30
		50mm (2")			10905	30
Maxi 22cm (8 1/2")	Opaque	35mm (1 3/8")			10914	30
		50mm (2")			10915	30
		70mm (2 3/4")			10918	30



SenSura® Flex Urostomy Pouch

Size	Colour	Coupling Size			Product Code	Units
Midi 24mm (9 1/2")	Opaque	35mm (1 3/8")			11821	20
		50mm (2")			11822	20
Maxi 26cm (10 1/4")	Opaque	35mm (1 3/8")			11824	20
		50mm (2")			11825	20
	Transparent	35mm (1 3/8")			11827	20
		50mm (2")			11828	20



SenSura®

Ostomy Barrier and Pouch - One Piece System

Easiclose™ system with Velcro® to ensure security. Double layer adhesive provides secure adhesion, flexibility and ease of removal. Only SenSura 1pc Convex Systems have 2 belt attachment points.

SenSura® Drainable Pouch with Easiclose™ WIDE-Outlet

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 30cm (12")	Cut-to-Fit	Non-Convex	Opaque	10-76mm (3/8-3")	15531	10
			Transparent	10-76mm (3/8-3")	15521	10
				10-76mm (3/8-3")	15821 *No filter*	10
		Convex Light	Opaque	15-33mm (5/8 - 1 1/4")	15607	10
				15-43mm (5/8 - 1 3/4")	15608	10
			Transparent	15-33mm (5/8 - 1 1/4")	15605	10
				15-43mm (5/8 - 1 3/4")	15606	10
	Pre-Cut	Non-Convex	Opaque	25mm (1")	15532	10
				30mm (1 1/8")	15533	10
				35mm (1 3/8")	15534	10
		Convex Light	Opaque	25mm (1")	15633	10
				28mm (1 1/8")	15634	10
				31mm (1 1/4")	15635	10
		Non-Convex	Transparent	25mm (1")	15522	10
				30mm (1 1/8")	15523	10
				35mm (1 3/8")	15524	10





SenSura® Xpro Drainable Pouch with Easiclose™ WIDE-Outlet

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 30cm (12")	Cut-to-Fit	Non-Convex	Transparent	10-76mm (3/8- 3")	15981	20
		Convex Light	Transparent	15-33mm (5/8 - 1 1/4")	15695	10
				15-43mm (5/8 - 1 3/4")	15696	10

SenSura® Closed Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 18cm (7″)	Cut-to-Fit	Non-Convex	Opaque	10-66mm (3/8 - 2 5/8″)	15450	30
			Transparent	10-66mm (3/8 - 2 5/8″)	15441	30
		Convex Light	Opaque	15-33mm (5/8 - 1 1/4″)	15403	10
				15-43mm (5/8 - 1 3/4″)	15404	10
	Pre-Cut	Convex Light	Opaque	25mm (1″)	15413	10
				28mm (1 1/8″)	15414	10
				31mm (1 1/4″)	15415	10
Maxi 22cm (8 1/2″)	Cut-to-Fit	Non-Convex	Opaque	10-76mm (3/8 - 3″)	15480	30
			Transparent	10-76mm (3/8 - 3″)	15470	30
		Convex Light	Opaque	15-33mm (5/8 - 1 3/4″)	15407	10
				15-43mm (5/8 - 1 3/4″)	15408	10
			Transparent	15-33mm (5/8 - 1 1/4″)	15405	10
				15-43mm (5/8 - 1 3/4″)	15406	10
	Pre-Cut	Non-Convex	Opaque	25mm (1″)	15482	30
				30mm (1 1/8″)	15483	30
				35mm (1 3/8″)	15484	30
			Transparent	25mm (1″)	15472	30
				30mm (1 1/8″)	15473	30
				35mm (1 3/8″)	15474	30



SenSura® Urostomy Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 24cm (9 1/2″)	Cut-to-Fit	Non-Convex	Opaque	10-66mm (3/8 - 2 5/8″)	11801	10
			Transparent	10-66mm (3/8 - 2 5/8″)	11803	10
Maxi 26cm (10 1/4″)	Cut-to-Fit	Non-Convex	Opaque	10-76mm (3/8 - 3″)	11802	10
			Transparent	10-76mm (3/8 - 3″)	11804	10
		Convex Light	Opaque	15-33mm (5/8 - 1 1/4″)	11811	10
				15-43mm (5/8 - 1 3/4″)	11812	10
			Transparent	15-33mm (5/8 - 1 1/4″)	11814	10
				15-43mm (5/8 - 1 3/4″)	11815	10



SenSura® One-Piece High Output Pouch

Cut	Sterile	Stoma Size	Volume	Product Code	Units
Cut-to-Fit	Sterile	10-76mm (3/8-3″)	630 mL	19001	5
	Non-Sterile	10-76mm (3/8-3″)	630 mL	19003	6



SenSura® Post-Op Pouch

Cut	Sterile	Stoma Size	Window Size	Volume	Product Code	Units
Cut-to-Fit	Sterile	10-76mm (3/8-3″)	70mm	680 mL	19010	5
		10-115mm (3/8 - 4 1/2″)	100mm	725 mL	19011	5
	Non-Sterile	10-76mm (3/8-3″)	70mm	680 mL	19020	6
		10-115mm (3/8 - 4 1/2″)	100mm	725 mL	19021	6
		10-76mm (3/8-3″)	No window	680 mL	19030	6



Assura®

Ostomy Barrier and Pouch - Two Piece System

Click mechanical coupling with a lock ring for security. Skin barrier with spiral technology for secure, long lasting wear. All 2pc Assura Click Flanges have 2 belt attachment points.


Assura® Skin Barrier

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	40mm (1 9/16")	10-35mm (3/8 - 1 3/8")	2881	5
		50mm (2")	10-45mm (3/8 - 1 3/4")	2882	5
		60mm (2 3/8")	10-55mm (3/8 - 2 1/8")	2883	5
	Convex Light	40mm (1 9/16")	15-23mm (5/8 - 7/8")	14261	5
		50mm (2")	15-33mm (5/8 - 1 1/4")	14262	5
		60mm (2 3/8")	15-43mm (5/8 - 1 3/4")	14263	5
	Convex Deep	50mm (2")	15-33mm (5/8 - 1 1/4")	12716	5
		60mm (2 3/8")	15-43mm (5/8 - 1 3/4")	12719	5
Pre-Cut	Non-Convex	40mm (1 9/16")	20mm (3/4")	12841	5
			25mm (1")	12842	5
		50mm (2")	30mm (1 1/8")	12844	5
			35mm (1 3/8")	12845	5
		60mm (2 3/8")	40mm (1 5/8")	12847	5
		60mm (2 3/8")	45mm (1 3/4")	12848	5
	Convex Light	40mm (1 9/16")	18mm (3/4")	14271	5
			21mm (7/8")	14272	5
		50mm (2")	25mm (1")	14273	5
			28mm (1 1/8")	14274	5
			31mm (1 1/4")	14275	5
		60mm (2 3/8")	35mm (1 3/8")	14276	5
			38mm (1 1/2")	14277	5
			41mm (1 5/8")	14278	5
	Convex Deep	40mm (1 9/16")	15mm (5/8")	12701	5
			18mm (3/4")	12702	5
			21mm (7/8")	12703	5
		50mm (2")	25mm (1")	12704	5
			28mm (1 1/8")	12705	5
			31mm (1 1/4")	12706	5
		60mm (2 3/8")	35mm (1 3/8")	12707	5
			38mm (1 1/2")	12708	5
			41mm (1 5/8")	12709	5



Assura® Skin Barrier - Extra Extended Wear

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	40mm (1 9/16")	10-35mm (3/8 - 1 3/8")	2831	5
		50mm (2")	10-45mm (3/8 - 1 3/4")	2832	5
		60mm (2 3/8")	10-55mm (3/8 - 2 1/8")	2833	5
	Convex Light	40mm (1 9/16")	15-23mm (5/8 - 7/8")	14281	5
		50mm (2")	15-33mm (5/8 - 1 1/4")	14282	5
		60mm (2 3/8")	15-43mm (5/8 - 1 3/4")	14283	5
	Convex Deep	40mm - (1 9/16")	15-23mm (5/8 - 7/8")	14243	5
		50mm (2")	15-33mm (5/8 - 1 1/4")	14246	5
		60mm (2 3/8")	15-43mm (5/8 - 1 3/4")	14249	5
Pre-Cut	Non-Convex	40mm (1 9/16")	20mm (3/4")	14251	5
			25mm (1")	14252	5
		50mm (2")	30mm (1 1/8")	14254	5
			35mm (1 3/8")	14255	5
		60mm (2 3/8")	45mm - 1 3/4"	14258	5
	Convex Light	40mm (1 9/16")	18mm (3/4")	14291	5
			21mm (7/8")	14292	5
		50mm (2")	25mm (1")	14293	5
			28mm (1 1/8")	14294	5
			31mm (1 1/4")	14295	5
		60mm (2 3/8")	35mm (1 3/8")	14296	5
	Convex Deep	40mm (1 9/16")	18mm (3/4")	14232	5
			21mm (7/8")	14233	5
		50mm (2")	25mm (1")	14234	5
			28mm (1 1/8")	14235	5
			31mm (1 1/4")	14236	5
		60mm (2 3/8")	35mm (1 3/8")	14237	5
			38mm (1 1/2")	14238	5
			41mm (1 5/8")	14239	5



Assura® 2 pc. High Output Drainable Pouch

Size	Colour	Coupling Size	Product Code	Units
Maxi 32cm (12 1/2")	Opaque	50mm (2")	2846	10
		60mm (2 3/8")	2847	10



Assura® 2 pc. High Output Drainable System, Ileo Night Pouch 2L

Size	Colour	Coupling Size	Product Code	Units
2 L	Opaque	50mm (2")	2836	5
		60mm (2 3/8")	2837	5



Assura® Drainable Pouch Clamp Closure

Size	Colour	Coupling Size	Product Code	Units
Maxi 30cm (12 1/2")	Transparent	50mm (2")	2792	10
		60mm (2 3/8")	2793	10



Assura® Drainable Pouch with EasiClose™

Size	Colour	Coupling Size	Product Code	Units
Mini 23cm (9")	Opaque	40mm (1 9/16")	13924	10
		50mm (2")	13925	10
Midi 26cm (10 1/4")	Opaque with Graphic	40mm (1 9/16")	13964	10
		50mm (2")	13965	10
		60mm (2 3/8")	13966	10
Maxi 28cm (11")	Opaque	40mm (1 9/16")	13984	10
		40mm (1 9/16") No filter	14494	10
		50mm (2")	13985	10
		60mm (2 3/8")	13986	10
		60mm (2 3/8") No filter	14498	10
	Transparent	40mm (1 9/16")	13974	10
		50mm (2")	13975	10
		60mm (2 3/8")	13976	10



Assura® Closed Pouch

Size	Colour	Coupling Size	Product Code	Units
Midi 18cm (7")	Opaque	40mm (1 9/16")	12354	30
		50mm (2")	12355	30
		60mm (2 3/8")	12356	30
	Transparent	40mm (1 9/16")	12344	30
		60mm (2 3/8")	12346	30
Maxi 22cm (8 1/2")	Opaque	40mm (1 9/16")	12384	30
		50mm (2")	12385	30
		60mm (2 3/8")	12386	30
	Transparent	40mm (1 9/16")	12374	30
		50mm (2")	12375	30
		60mm (2 3/8")	12376	30



Assura® Urostomy Pouch

Size	Colour	Coupling Size	Product Code	Units
Maxi 25cm (10")	Opaque	40mm (1 9/16")	14224	10
		50mm (2")	14225	10
		60mm (2 3/8")	14226	10
	Transparent	40mm (1 9/16")	14227	10
		50mm (2")	14228	10
		60mm (2 3/8")	14229	10



Assura®

Ostomy Barrier and Pouch - One Piece System

Skin Barrier with spiral technology for secure, long lasting wear. Different convexity and pouch shape. All Assura 1pc Convex Systems have 2 belt attachment points.

Assura® Drainable Pouch Clamp Closure

Size	Cut	Convexity	Colour	Stoma Size	Product Code		Units
Maxi 30cm (12")	Cut-to-Fit	Non-Convex	Opaque	10-70mm (3/8 - 2 3/4")	2520	No filter	10
			Transparent	10-70mm (3/8 - 2 3/4")	2490	No filter	10
		Convex Deep	Opaque	15-43mm (5/8 - 1 3/4")	12533	No filter	10
			Transparent	15-43mm (5/8 - 1 3/4")	13706	No filter	10



Assura® Drainable Pouch with EasiClose™

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 26cm (10 1/4")	Cut-to-Fit	Non-Convex	Opaque	10-55mm (3/8 - 2 1/8")	13840	10
		Convex Light	Opaque	15-33mm (5/8 - 1 1/4")	14404	10
				15-43mm (5/8 - 1 3/4")	14405	10
		Convex Deep	Opaque	15-43mm (5/8 - 1 3/4")	14196	10
	Pre-Cut	Non-Convex	Opaque	25mm (1")	13844	10
				30mm (1 1/8")	13845	10
				35mm (1 3/8")	13846	10
Maxi 28cm (11")	Cut-to-Fit	Non-Convex	Opaque	10-70mm (3/8 - 2 3/4")	13870	10
			Transparent	10-70mm (3/8 - 2 3/4")	13860	10
		Convex Light	Opaque	15-33mm (5/8 - 1 1/4")	14414	10
				15-43mm (5/8 - 1 3/4")	14415	10
			Transparent	15-33mm (5/8 - 1 1/4")	14411	10
				15-43mm (5/8 - 1 3/4")	14412	10



Assura® Drainable Pouch with EasiClose™

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 28cm (11")	Cut-to-Fit	Convex Deep	Opaque	15-43mm (5/8 - 1 3/4")	14106	10
			Transparent	15-43mm (5/8 - 1 3/4")	14103	10
	Pre-Cut	Non-Convex	Opaque	25mm (1")	13874	10
				30mm (1 1/8")	13875	10
				35mm (1 3/8")	13876	10
		Convex Light	Opaque	25mm (1")	14422	10
				28mm (1 1/8")	14423	10
				31mm (1 1/4")	14424	10
		Convex Deep	Opaque	21mm (7/8")	14173	10
				25mm (1")	14174	10
				28mm (1 1/8")	14175	10
				31mm (1 1/4")	14176	10
			Transparent	21mm (7/8")	14163	10
				25mm (1")	14164	10
				28mm (1 1/8")	14165	10
				31mm (1 1/4")	14166	10



Assura® Closed Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Mini 12cm (4 3/4")	Cut-to-Fit	Non-Convex	Opaque	20-55mm (3/8 - 2 1/8")	12110	30
Midi 18cm (7")	Cut-to-Fit	Non-Convex	Opaque	20-55mm (3/8 - 2 1/8")	12140	30
			Transparent	20-55mm (3/8 - 2 1/8")	12130	30
	Pre-Cut	Non-Convex	Opaque	25mm (1")	12144	30
				30mm (1 1/8")	12145	30
				35mm (1 3/8")	12146	30
				40mm (1 9/16")	12147	30
				50mm (2")	12149	30



Assura® 1-Piece Post-Op Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Mini 31cm (12 1/4")	Cut-to-Fit	Non-Convex	Transparent	13-70mm (1/2 - 2 3/4")	12805	5
Midi 33cm (13")	Cut-to-Fit	Non-Convex	Transparent	13-70mm (1/2 - 2 3/4")	12800	5
				13-100mm (1/2 - 4")	12802	5
				13-70mm (1/2 - 2 3/4")	12808 *Sterile*	5
				13-100mm (1/2 - 4")	12810 *Sterile*	5



Assura® Closed Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Midi 18cm (7")	Pre-Cut	Convex Light	Opaque	25mm (1")	14452	10
				28mm (1 1/8")	14453	10
				31mm (1 1/4")	14454	10
Maxi 22cm (8 3/4")	Cut-to-Fit	Non-Convex	Opaque	20-70mm (3/4-2 3/4")	12170	30
		Convex Light	Opaque	15-43mm (5/8 - 1 3/4")	14445	10
			Transparent	15-33mm (5/8 - 1 1/4")	14441	10
	Pre-Cut	Non-Convex	Opaque	25mm (1")	12174	30
				30mm (1 1/8")	12175	30
				35mm (1 3/8")	12176	30
				40mm (1 9/16")	12177	30



Assura® Urostomy Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 27cm (10 1/2")	Cut-to-Fit	Non-Convex	Opaque	10-55mm (3/8 - 2 1/8")	5580	10
			Transparent	10-55mm (3/8 - 2 1/8")	5585	10
		Convex Deep	Transparent	15-33mm (5/8 - 1 1/4")	12595	10
				15-43mm (5/8 - 1 3/4")	12596	10
	Pre-Cut	Convex Deep	Transparent	15mm (5/8")	12991	10
				18mm (3/4")	12992	10
				21mm (7/8")	12993	10
				25mm (1")	12994	10
				28mm (1 1/8")	12995	10
				31mm (1 1/4")	12996	10
				35mm (1 3/8")	12997	10



Assura® Urostomy Pouch - Extra Extended Wear

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
Maxi 25cm (10")	Cut-to-Fit	Non-Convex	Opaque	10-55mm (3/8 - 2 1/8")	14201	10
			Transparent	10-55mm (3/8 - 2 1/8")	14202	10
		Convex Light	Opaque	15-43mm (5/8 - 1 3/4")	14701	10
			Transparent	15-43mm (5/8 - 1 3/4")	14702	10
Maxi 27cm (10 1/2")	Cut-to-Fit	Non-Convex	Opaque	10-55mm (3/8 - 2 1/8")	12471	10
			Transparent	10-55mm (3/8 - 2 1/8")	12472	10
		Convex Deep	Transparent	15-33mm (5/8 - 1 1/4")	12432	10
				15-43mm (5/8 - 1 3/4")	12433	10



Assura®

Ostomy Barrier and Pouch - Two Piece System - Pediatric

Click mechanical coupling.

Assura® Pediatric Skin Barrier

Cut	Convexity	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	10-35mm (3/8 - 1 3/8")	2181	5



Assura® Pediatric Drainable Pouch

Size	Colour	Product Code	Units
13cm (5")	Opaque	2151	10

Assura® Pediatric Closed Pouch

Size	Colour	Product Code	Units
13cm (5")	Transparent	2166	30



EasiFlex®

Ostomy Barrier and Pouch - Two Piece System - Pediatric

Flex adhesive coupling.

Easiflex® Pediatric Skin Barrier

Cut	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	17mm	0-15mm (0-5/8")	14307	10
		27mm	0-25mm (0-1")	14308	10



Easiflex® Pediatric Drainable Pouch

Size	Colour	Coupling Size	Product Code	Units
105 mL	Teddy Print	17mm	14691	10
	Transparent	17mm	14681	10
150 mL	Teddy Print	27mm	14692	10
	Transparent	27mm	14682	10



Assura®

Ostomy Barrier and Pouch - One Piece System - Pediatric

Assura® Pediatric Drainable Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
15cm (6")	Cut-to-Fit	Non-Convex	Transparent	10-35mm (3/8 - 1 3/8")	2115	10

Assura® Pediatric Closed Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
13cm (5")	Cut-to-Fit	Non-Convex	Transparent	10-35mm (3/8 - 1 3/8")	2125	30

Assura® Pediatric Urostomy Pouch

Size	Cut	Convexity	Colour	Stoma Size	Product Code	Units
15cm (6")	Cut-to-Fit	Non-Convex	Transparent	10-35mm (3/8 - 1 3/8")	8009	10



Ostomy Accessories

Brava® Protective Sheets

Size (cms)	Size (inches)	Product Code	Units
10 x 10 cm	4 x 4"	3210	10
15 x 15 cm	6 x 6"	3215	5
20 x 20 cm	8 x 8"	3220	5



Brava® Protective Rings

Size (mm) Internal diameter	Size (mm) External diameter	Size (mm) thickness	Product Code	Units
18 mm	48 mm	2.5 mm	12035	10
18 mm	48 mm	4.2 mm	12045	10
27 mm	57 mm	2.5 mm	12037	10
27 mm	57 mm	4.2 mm	12047	10
34 mm	64 mm	2.5 mm	12039	10
34 mm	64 mm	4.2 mm	12049	10



Brava® Protective Rings Wide

Size (mm) Internal diameter	Size (mm) External diameter	Size (mm) thickness	Product Code	Units
18 mm	57 mm	2.5 mm	12036	10
18 mm	64 mm	2.5 mm	12038	10
18 mm	57 mm	4.2 mm	12046	10
18 mm	64 mm	4.2 mm	12048	10



Brava® Protective Rings Extra Wide

Size (mm) Internal diameter	Size (mm) External diameter	Size (mm) thickness	Product Code	Units
18 mm	76 mm	2.5 mm	12032	10
27 mm	76 mm	2.5 mm	12033	10



Brava® Protective Ring Convex

Size (mm) Internal diameter	Size (mm) External diameter	Size (mm) thickness	Product Code	Units
20 mm	60 mm	8.0 mm	12090	10
25 mm	65 mm	8.0 mm	12091	10
30 mm	70 mm	8.0 mm	12092	10
35 mm	75 mm	8.0 mm	12093	10
40 mm	80 mm	8.0 mm	12094	10
45 mm	85 mm	8.0 mm	12095	10



Brava® Mouldable Rings

Size (mm) Internal diameter	Size (mm) thickness	Product Code	Units
18 mm	2.0 mm	12030	10
18 mm	4.2 mm	12042	10



Brava® Skin Barrier Spray and Wipe

Size	Product Code	Units
50 mL Spray	12020	1
Single unit wipe	12021	30



Brava® Adhesive Remover Spray and Wipe

Size	Product Code	Units
50 mL Spray	12010	1
Single unit wipe	12011	30



Brava® Barrier Cream

Size	Product Code	Units
60g Tube	12000	1



Brava® Ostomy Strip Paste

Size	Product Code	Units
57g / 2oz	2655	10



Brava® Elastic Barrier Strips

Size	Product Code	Units
Brava Elastic Barrier Strip Original (a.)	12070	20
Brava Elastic Barrier Strip Straight (b.)	12074	20
Brava Elastic Barrier Strip Belt (c.)	12072	20
Brava Elastic Barrier Strip XL (d.)	12076	20



Brava® Ostomy Powder

Size	Product Code	Units
25g / 1oz bottle	1907	16



Brava® Mio Belt for All SenSura Mio Products
New design for easier attachment and improved security

Size	Product Code	Units
Standard - 100 cm - 40"	0423	15
XXL - 155cm - 61"	0424	15



Brava® Belt for SenSura Original + Assura Products

Size	Product Code	Units
Standard - 100 cm - 40"	0421	15
XL - 135cm - 54"	0422	15



Brava® Lubricating Deodorant

Size	Product Code	Units
7.5 mL sachet	12060	20
240 mL bottle	12061	1



Brava® Alcohol Free Paste

Size	Product Code	Units
60 g tube	12050	1



Ostomy Paste (contains alcohol)

Size	Product Code	Units
60 g tube	2650	12



Skin Barrier Rings

Pre-cut Size	Product Code	Units
20mm ($^{3}/_{4}$")	2320	30
25mm (1")	2325	30
30mm (1 $^{1}/_{8}$")	2330	30
40mm (1 $^{9}/_{16}$")	2340	30
50mm (2")	2350	30



Coloplast Prep™ Protective Liquid Skin Barrier
(contains alcohol)

Size	Product Code	Units
Single Application Wipe	2041	12 boxes/54
59 mL - 2oz. dab bottle	0925	12



Coloplast Filtrodor® Pouch Filter

	Product Code	Units
Charcoal-activated deodorizing filter with self-adhering foam center, adheres to front of Assura & SenSura Pouches	0509	5 pouches (10 units per pouch)



Coloplast Pouch Clamp

	Product Code	Units
Simple closure clamp for Assura Drainable Pouches	9500	20



Ostomy Scissors

	Product Code	Units
Stainless steel medical grade ostomy scissors with curved blades and rounded ends - right handed	9505	30



Hex-On®

Odour antagonist

Eliminates odour with a light, mild scent.

Hex-On®

Size (ounces)	Size (mL)	Product Code	Units
2 oz.	59 mL	7583	12



Ostomy Cap

Designed for discreet situations. 2 piece compatible with SenSura Click and Assura 2 piece skin barriers.

**Note Uro Minicap and Coloplast Stoma cap do not fit SenSura Mio products*

Coloplast Uro **Minicap** (100 ml)

Size	Colour	Coupling Size	Product Code	Units
Mini	Opaque	40mm (1 9/16")	2807	30
		50mm (2")	2808	30




Coloplast Stoma Cap

Size	Colour	Coupling Size	Product Code	Units
Mini	Opaque	40mm (1 9/16")	2801	30
		50mm (2")	2802	30
		60mm (2 3/8")	2803	30



Coloplast Stoma Cap - One Piece

Size	Colour	Cut to Fit - Stoma size	Product Code	Units
Mini	Opaque	20-55mm (3/4 - 2 1/8")	2501	30



Assura® Irrigation Products

Assura® Irrigation Set

Description	Product Code	Units
Includes: - 1 Water Bag/Regulator - 1 Cone	1500	1



Assura® Irrigation Accessories

Size	Product Code	Units
Stoma Cone - Soft, flexible stoma cone with tubing	1110	4



Irrigation Faceplate. For use with ostomy belt 0421 or 0422. Use with #12836 Irrigation Sleeve	12820	1



Assura® Irrigation Sleeves

Size	Product Code	Units
Transparent, disposable irrigation sleeve with adhesive back. - 1-pc sleeve	1003	10



The irrigation sleeves below must be used with their corresponding colour coded flange or faceplate. Consult this chart to order the correct transparent sleeves.

Description	Coupling Size	Use with Assura and SenSura Barrier Product Codes (will not work with SenSura Mio)	Product Code	Units
Irrigation sleeves	40mm (1 9/16")	Will couple to any green 40 mm barriers within the Assura and SenSura lines	12834	5
	50mm (2")	Will couple to any red 50 mm barriers within the Assura and SenSura lines	12835	5
	60mm (2 3/8")	Will couple to any blue 60 mm barriers within the Assura and SenSura lines	12836	5



Coloplast Urostomy Night Bag

Description	Product Code	Units
Includes a fixed connector that fits on to the urostomy pouch. The bag holds up to 2 litres and has a 120cm long flexible tube with anti-reflux valve. ***Not Compatible with External Catheters or Leg Bags***	21365	10



Coloplast Fistula & Wound Management System

Pouch

Size	Cutting Surface	Product Code	Units
Mini	104-159mm (4 1/8 - 6 1/4")	14050	6
Midi	156-228mm (6 3/16 - 9")	14060	3
Maxi	208-297mm (8 1/4 - 11 3/4")	14070	



Flexible Lids

Size	For	Product Code	Units
Mini	Mini Fistula & Wound Management System	14021	10
Midi	Midi Fistula & Wound Management System	14031	
Maxi	Maxi Fistula & Wound Management System	14041	



Accessories

Size	For	Product Code	Units
Bed Drainage Bag	All Fistula & Wound Management systems (Capacity 2000 mL)	14010	10
Drain Port	All Fistula & Wound Management systems	14015	12





Continence

SpeediCath® Flex

Hydrophilic Intermittent Catheter

A soft catheter for men with a dry-sleeve and a flexible tip with all the benefits of the SpeediCath® catheter. Offered with a loop system that fits into a discreet packaging which is 100% aluminum free. Latex, phthalate and PVC free. Sterile.

SpeediCath® Flex

Effective Length (cm)	Effective Length (inches)	Size	Product Code	Units
33 cm	13”	10 FR	28920	30
		12 FR	28922	30
		14 FR	28924	30
		16 FR	28916	30



SpeediCath® Standard

Hydrophilic Intermittent Catheter

Instantly ready to use catheter. Featuring a smooth and uniform hydrophilic coating and polished eyelets to ensure a convenient and simple catheterization. Latex, phthalate and PVC free. Sterile.

SpeediCath® Standard Male Straight

Effective Length (cm)	Effective Length (inches)	Size	Product Code	Units
35 cm	13.9”	8 FR	27408	30
		10 FR	27410	30
		12 FR	27412	30
		14 FR	27414	30
		16 FR	27416	30
		18 FR	27418	30



SpeediCath® Standard Male Coude

Effective Length (cm)	Effective Length (inches)	Size	Product Code	Units
35 cm	13.9”	10 FR	27490	30
		12 FR	27492	30
		14 FR	27494	30
		16 FR	27496	30



SpeediCath® Standard Female Straight

Effective Length (cm)	Effective Length (inches)	Size	Product Code	Units
16 cm	6.3”	8 FR	27508	30
		10 FR	27510	30
		12 FR	27512	30
		14 FR	27514	30
		16 FR	27516	30



SpeediCath® Standard Pediatric Straight

Effective Length (cm)	Effective Length (inches)	Size	Product Code	Units
16 cm	6.3”	6 FR	27706	30
		8 FR	27708	30
		10 FR	27710	30



SpeediCath® Compact Male

Hydrophilic Intermittent Catheter

Developed to offer discretion and practicality with all the benefits of the SpeediCath® catheter and more. Telescopic design makes the packaging half the size of a standard male catheter. Latex, phthalate and PVC free. Sterile.

SpeediCath® Compact Male

Effective Length (cm)	Effective Length (inches)	Size		Product Code	Units
30 cm	11.8"	12/18 FR	(one size)	28692	30



SpeediCath® Compact Eve

Hydrophilic Intermittent Catheter

Stylish and discreet compact catheter for women. Designed for modern life, delivering a one-step easy opening and secure reclosing. Triangular shape ensures easy handling and grip. Latex, phthalate and PVC free. Sterile.

SpeediCath® Compact Eve

Effective Length (cm)	Effective Length (inches)	Size	Product Code	Units
8 cm	3.1"	10 FR ●	28110	30
		12 FR ○	28112	30
		14 FR ●	28114	30



SpeediCath® Compact Female

Hydrophilic Intermittent Catheter

Compact and discreet solution that is instantly ready to use. The catheter can easily be kept in a pocket or handbag. Latex, phthalate and PVC free. Sterile.

SpeediCath® Compact Female

Effective Length (cm)	Effective Length (inches)	Size	Product Code	Units
7 cm	2.7"	8 FR ●	28578	30
		10 FR ●	28580	30
		12 FR ○	28582	30
		14 FR ●	28584	30

SpeediCath® Compact Female Plus

Effective Length (cm)	Effective Length (inches)	Size	Product Code	Units
9 cm	3.5"	10 FR ●	28810	30
		12 FR ○	28812	30
		14 FR ●	28814	30





SpeediCath® Compact Set

Hydrophilic Intermittent Catheter

Compact and discreet catheter and bag solution that is instantly ready to use. Can easily be kept in a pocket or handbag. Latex, phthalate and PVC free. Sterile.

SpeediCath® Compact Set Male

Effective Length (cm)	Effective Length (inches)	Size		Product Code	Units
30 cm	11.8"	12/18 FR	(one size)	28422	20



SpeediCath® Compact Set Female

Effective Length (cm)	Effective Length (inches)	Size		Product Code	Units
9 cm	3.5"	10 FR	●	28520	20
		12 FR	○	28522	20
		14 FR	●	28524	20



Self-Cath®

Intermittent Catheter

Siliconized surface for smooth insertion. Made of flexible medical grade PVC with smooth fire polished eyelets for patient comfort and ease of use. Latex free. Sterile.

Self-Cath® Female Straight

Effective Length (cm)	Effective Length (inches)	Size		Reference Code		Product Code	Units
12 cm	4.8"	8 FR	●	504600		208	30
		10 FR	●	504610		210	30
		12 FR	○	504620		212	30
		14 FR	●	504630		214	30
		14 FR	●	501960	without connector	240	30



Self-Cath® Pediatric

Effective Length (cm)	Effective Length (inches)	Size		Reference Code		Product Code	Units
20 cm	8"	5 FR	○	504400		305	30
		6 FR	●	504410		306	30
		8 FR	●	504420	without connector	308	30
		10 FR	●	504430	without connector	310	30



Self-Cath® Extension Tube

Total Length (cm)	Total Length (inches)	Size	Reference Code		Product Code	Units
61 cm	24"	n/a	507850		475	12



Self-Cath® Male Straight

Effective Length (cm)	Effective Length (inches)	Size	Reference Code	Product Code	Units
38 cm	14.9”	8 FR	504500	408	30
		10 FR	504510	410	30
		12 FR	504520	412	30
		14 FR	504530	414	30
		16 FR	504540	416	30
		18 FR	504550	418	30



Self-Cath® Male Soft Straight

Effective Length (cm)	Effective Length (inches)	Size	Reference Code	Product Code	Units
38 cm	14.9”	12 FR	501810	112	30
		14 FR	501820	114	30
		16 FR	501830	116	30

Self-Cath® Male Straight Curved Packaging

Effective Length (cm)	Effective Length (inches)	Size	Reference Code	Product Code	Units
38 cm	14.9”	12 FR	502200	460	50
		14 FR	502210	450	50

Self-Cath® Male Coude Tapered with Guide Stripe

Effective Length (cm)	Effective Length (inches)	Size	Reference Code	Product Code	Units
38 cm	14.9”	8 FR	504650	608	30
		10 FR	504660	610	30
		12 FR	504670	612	30
		14 FR	504680	614	30



Self-Cath® Male Coude Olive with Guide Stripe

Effective Length (cm)	Effective Length (inches)	Size	Reference Code	Product Code	Units
38 cm	14.9”	8 FR	504710	808	30
		10 FR	504720	810	30
		12 FR	504730	812	30
		14 FR	504740	814	30
		16 FR	504750	816	30
		18 FR	504760	818	30



Self-Cath® Male Closed System Straight

Effective Length (cm)	Effective Length (inches)	Size	Reference Code	Product Code	Units
31 cm	12.2”	14 FR	502730	1114	50
		16 FR	502740	1116	50



Conveen® Optima

Male External Catheter - Silicone

Discreet, easy to use, reliable all day self-adhering external catheter. Double-grip strip ensures easy and smooth roll-out for a secure application. Push ring with sure grip ribs provides secure connection to the collecting bag. Anti-kink bellow stop kinks ensuring a free flow. Compact easy-to-open packaging. Latex free.

Conveen® Optima Standard

Length (cm)	Length (inches)	Size	Product Code	Units
7.8 cm	3"	25 mm	22025	30
		28 mm	22028	30
		30 mm	22030	30
		35 mm	22035	30
		40 mm	22040	30

Conveen® Optima Short

Length (cm)	Length (inches)	Size	Product Code	Units
4.6 cm	2"	21 mm	22121	30
		25 mm	22125	30
		30 mm	22130	30
		35 mm	22135	30



Conveen® Security+

Male External Catheter - SEBS*

Easy to use, reliable all day self-sealing external catheter. Applicator grip and single-grip strip ensures easy and smooth roll-out for a secure application. Push ring with sure grip ribs provides secure connection to the collecting bag. Anti-kink bellow stop kinks ensuring a free flow. Latex free.

Conveen® Security+ Self Sealing

Length (cm)	Length (inches)	Size	Product Code	Units
8.2 cm	3"	21 mm	5221	35
		25 mm	5225	35
		25 mm	5225H	100
		30 mm	5230	35
		30 mm	5230H	100
		35 mm	5235	35
		35 mm	5235H	100



Conveen® Security+ Adhesive Liner

Length (cm)	Length (inches)	Size	Product Code	Units
8.2 cm	3"	30 mm	5030	35
		30 mm	5030H	100
		35 mm	5035	35
		35 mm	5035H	100



*Styrene Ethylene Butylene Styrene

Clear Advantage® / Freedom Clear™

Male External Catheter - Silicone

Self-adhering male external catheter. Latex-free.

Clear Advantage® with Aloe Vera

Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
7.9 cm	3"	23 mm	504311	C6100	100
		28 mm	504331	6200	100
		31 mm	504341	6300	100
		35 mm	504351	6400	100

Freedom Clear™

Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
7.9 cm	3"	23 mm	505711	5100	100
		28 mm	505731	C5200	100
		31 mm	505741	5300	100
		35 mm	505751	5400	100



Freedom Clear™ "LS" Long Seal

Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
8.9 cm	3.5"	28 mm	505421	5290	100
		31 mm	505431	5390	100
		35 mm	505341	5490	100
		40 mm	505451	5590	100

Freedom Clear™ "SS" Sport Sheath

Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
5.5 cm	2"	23 mm	505311	5110	100
		28 mm	505321	C5210	100
		31 mm	505331	C5310	100
		35 mm	505441	C5410	100

Conveen®

Male External Catheter - Latex

Conveen® Self-Sealing

Length (cm)	Length (inches)	Size	Product Code	Units
10.5 cm	4"	20 mm	5212	35
		25 mm	5200	35
		30 mm	5205	35
		35 mm	5210	35
		40 mm	5215	35



Conveen® Adhesive Liner

Length (cm)	Length (inches)	Size	Product Code	Units
10.5 cm	4"	30 mm	5130	35
		35 mm	5135	35



Active Cath®

Male External Catheter - Latex

Self-sealing male external catheter.

Active Cath®

Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
7.9 cm	3″	23 mm	506220	8100	100
		28 mm	506200	8300	100
		31 mm	506240	8305	100
		35 mm	506260	8500	100



Freedom Cath®

Male External Catheter - Latex

Self-sealing male external catheter.

Freedom Clear®

Length (cm)	Length (inches)	Size	Reference Code	Product Code	Units
7.9 cm	3″	23 mm	506110	8000	100
		28 mm	506130	8200	100
		31 mm	506140	8205	100
		35 mm	506150	8400	100



Conveen® Active

Urine Leg Bag

Developed for an active lifestyle. Integrated in a comfortable elastic band that wraps around the thigh for ease of use, mobility and discretion. Hide-away outlet tap concealed in elastic band when not in use. Soft fabric band for increased wearing comfort and a no-slip silicone feature inside straps for security. Latex and PVC free.

Conveen® Active Leg Bag, Non Sterile

Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
9 oz	250 mL	3″	25501	10



Conveen® Security+

Urine Leg bag

Secure and reliable urine leg bag with anti-kink tubing to reduce the risk of accidental leakage. Soft backing for increased wearing comfort (except codes 21034 and 21054). Latex and PVC free.

Conveen® Security+ Leg Bag, Clamp Outlet, Straps, Non Sterile

Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
17 oz	500 mL	20″	5161	10
26 oz	750 mL	20″	5167	10

Conveen® Security+ Leg Bag, Clamp Outlet, Sample Port, Straps, Sterile

Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
17 oz	500 mL	20″	5163	10

Conveen® Security+ Leg Bag, Lever Outlet, Straps, Sterile

Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
17 oz	500 mL	20″	21034	10
34 oz	1000 mL	20″	21054	10



Conveen® Standard

Day / Night Urine Bag

Secure and reliable urine bag for bedside drainage or large capacity leg bag. Clamp outlet makes it easy to verify proper closure and reduce the risk of accidental leakage. Standard tubing. Leg bag straps sold separately. Soft backing for increased wearing comfort. Latex free.

Conveen® Standard Day / Night Bag, Clamp Outlet, Non Sterile

Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
51 oz	1500 mL	36”	5062	10

Conveen® Standard Day / Night Bag, Clamp Outlet, Sample Port, Sterile

Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
51 oz	1500 mL	36”	5063	10



Conveen® Contour

Urine Leg Bag

Secure and reliable urine leg bag. Enhanced chambered design for added discretion. Soft backing for increased wearing comfort. Latex free.

Conveen® Contour Leg Bag, Clamp Outlet, Straps, Non Sterile

Capacity (oz)	Capacity (mL)	Tube Length (inches)	Product Code	Units
21 oz	600 mL	18”	5170	10
28 oz	800 mL	18”	5174	10



Conveen® Night Bag

Urine Night Bag

Secure and reliable urine night bag with anti-kink tubing to prevent accidental leakage. Large capacity allows the user to feel secure all night for a full night’s sleep. Latex and PVC free. Not Compatible with Urostomy Pouches.

Conveen® Bedside Night Bag, Non-Sterile

Capacity (oz)	Capacity (L)	Product Code	Units
68 oz	2L	21346	30

Conveen® Bedside Night Bag, Sterile

Capacity (oz)	Capacity (L)	Product Code	Units
68 oz	2L	21356	10



Accessories

Conveen® Bag Hanger

Length (cm)	Length (inches)	Product Code	Units
18 cm	7”	5070	1



Conveen® Leg Bag Strap

Length (cm)	Length (inches)		Product Code	Units
63 cm	25”	Sold in pairs	5050	10



Conveen® Catheter retaining Strap

Length (cm)	Length (inches)	Product Code	Units
90 cm	35”	383004	2



Peristeen Plus® Anal Irrigation

Transanal irrigation system for fecal incontinence and constipation. Minimizes the likelihood of involuntary bowel leakage and/or constipation.

Peristeen Plus® Anal Irrigation System

Items	Description	Product Code
- 1 control unit - 2 sterile rectal catheters - 1 bag (with lid) - 2 straps - 1 toiletries bag	** 90 irrigations = 6 months	29140



Peristeen Plus® Anal Irrigation System - Small

Items	Description	Product Code
- 1 control unit - 2 sterile rectal catheters - small - 1 bag (with lid) - 2 straps - 1 toiletries bag	** 90 irrigations = 6 months	29147

Peristeen Plus® Anal Irrigation Accessory unit

Items	Description	Product Code
- 15 sterile rectal catheters - 1 bag (no lid)	**15 irrigations = 1 month	29142



Peristeen Plus® Anal Irrigation Accessory unit - Small

Items	Description	Product Code
- 15 sterile rectal catheters - small - 1 bag (no lid)	**15 irrigations = 1 month	29149

Peristeen Plus® Anal Irrigation Rectal Catheter

Items	Product Code
- 10 sterile rectal catheters	29143



Peristeen Plus® Anal Irrigation Leg Strap

Items	Product Code
- 1 pair	29145



Peristeen Plus® Anal Irrigation Tube

Items	Product Code
- 2 tubes with blue connectors	29146



** if used every other day

Wound

WOUND

Biatain® Silicone

Foam Dressing with Silicone Adhesive

Soft, flexible and absorbent foam dressing with gentle silicone adhesive for use on a broad range of exuding wounds, making it the obvious choice for moist wound healing of commonly occurring wounds. For more information please visit www.biatainsiliconeca.com.

reddot award 2017 winner

3DFit™ Technology

Biatain® Silicone			
Size (cm)	Size (inches)	Product Code	Units
7.5 x 7.5 cm	3 x 3"	33434	10
10 x 10 cm	4 x 4"	33435	10
12.5 x 12.5 cm	5 x 5"	33436	10
15 x 15 cm	6 x 6"	33437	5
17.5 x 17.5 cm	7 x 7"	33438	5
10 x 20 cm	3.9 x 7.9"	33400	5
10 x 30 cm	3.9 x 11.8"	33401	5
18 x 18 cm	7.1 x 7.1" Heel	33406	5
25 x 25 cm	9.8 x 9.8" Large Sacral	39000*	5
15 x 19 cm	5.9 x 7.5" Small Sacral	39001*	5
14 x 19.5 cm	5.5 x 7.7" Multi-shape	39010*	5



* This product has the Pressure Injury Prevention indication.

Biatain® Silicone Non-Border

Designed for high absorption and retention under compression to minimise the risk of leakage and maceration of fragile periwound skin. Can also be cut to fit areas that are difficult to dress, e.g. on the heel or in between toes and fingers. For high absorption and retention under compression. 3DFit Technology minimizes risk of leakage and maceration.

Biatain® Silicone			
Size (cm)	Size (inches)	Product Code	Units
5 x 7.5 cm	2 x 3"	39020	10
7.5 x 7.5 cm	3 x 3"	39021	10
10 x 10 cm	4 x 4"	39022	10
12.5 x 12.5 cm	5 x 5"	39023	10
10 x 20 cm	4 x 8"	39024	10
15 x 15 cm	6 x 6"	39025	10
20 x 20 cm	8 x 8"	39026	10



Biatain® Silicone Lite

Thin Foam Dressing with Silicone Adhesive

Thin and ultra flexible foam dressing with gentle silicone adhesive for use on a wide range of low to moderate exuding wounds.

3DFit™ Technology

Biatain® Silicone Lite			
Size (cm)	Size (inches)	Product Code	Units
5 x 5 cm	2 x 2"	33452	5
5 x 12.5 cm	2 x 5"	33453	5
7.5 x 7.5 cm	3 x 3"	33444	10
10 x 10 cm	4 x 4"	33445	10
12.5 x 12.5 cm	5 x 5"	33446	10



Biatain® Silicone Ag

Foam Dressing with Silicone Adhesive and Silver

Flexible and absorbent foam dressing with gentle silicone adhesive and a patented silver complex for use on a broad range of exuding infected wounds or wounds at risk of infection.

3DFit™ Technology

Biatain® Silicone Ag

Size (cm)	Size (inches)		Product Code	Units
7.5 x 7.5 cm	3 x 3”		39636	5
10 x 10 cm	4 x 4”		39637	5
12.5 x 12.5 cm	5 x 5”		39638	5
15 x 15 cm	6 x 6”		39639	5
17.5 x 17.5 cm	7 x 7”		39640	5
10 x 20 cm	3.9 x 7.9”		39644	5
10 x 30 cm	3.9 x 11.8”		39645	5
15 x 19 cm	5.9 x 3.5”	Small Sacral	39650	5
25 x 25 cm	9.8 x 9.8”	Large Sacral	39651	5
18 x 18 cm	7.1 x 7.1”	Heel	39652	5



Biatain®

Foam Dressing

Soft and conformable polyurethane foam dressing that effectively absorbs and retains wound exudate. This ensures a balance that is optimal for management of exuding wounds. Available in adhesive and non-adhesive. The adhesive has a hydrocolloid border – for a secure adhesion. This product is also used in high friction and shear areas such as the coccyx.

3DFit™ Technology

Biatain® Adhesive

Size (cm)	Size (inches)		Product Code	Units
7.5 x 7.5 cm	3 x 3”		3462	10
10 x 10 cm	4 x 4”		3430	10
12.5 x 12.5 cm	5 x 5”		3420	10
17 x 17 cm	6.7 x 6.7”	Small Sacral	3483	5
18 x 18 cm	7 x 7”		3423	5
19 x 20 cm	7.5 x 8”	Heel	3488	5
23 x 23 cm	9 x 9”	Large Sacral	3485	5



Biatain® Non-Adhesive

Size (cm)	Size (inches)	Product Code	Units
5 x 7 cm	2 x 2.8”	6105	10
10 x 10 cm	4 x 4”	3410	10
10 x 20 cm	4 x 8”	3412	5
15 x 15 cm	6 x 6”	3413	5
20 x 20 cm	8 x 8”	3416	5



Biatain® Soft Hold

Foam Dressing with Gentle Adhesion

Soft and conformable polyurethane foam dressing with a gentle adherent layer covering less than 50% of the foam surface. This product is used under compression. This allows both hands to be free during dressing application and removal.

3DFit™ Technology

Biatain® Soft Hold

Size (cm)	Size (inches)	Product Code	Units
10 x 10 cm	4 x 4"	3470	5
10 x 20 cm	4 x 8"	3472	5
15 x 15 cm	6 x 6"	3475	5



Biatain® Ag

Foam Dressing with Silver

Soft and conformable polyurethane foam dressing with a patented silver complex for use on infected wounds or wounds at risk of infection. Available in adhesive and non-adhesive. The adhesive has a hydrocolloid border – for a secure adhesion. This product is also used in high friction and shear areas such as the coccyx.

3DFit™ Technology

Biatain® Ag Adhesive

Size (cm)	Size (inches)		Product Code	Units
7.5 x 7.5 cm	3 x 3"		9631	5
12.5 x12.5 cm	5 x 5"		9632	5
18 x 18 cm	7 x 7"		9635	5
23 x 23 cm	9 x 9"	Sacral	9641	5
19 x 20 cm	7.5 x 8"	Heel	9643	5



Biatain® Ag Non-Adhesive

Size (cm)	Size (inches)	Product Code	Units
5 x 7 cm	2 x 2.8"	5105	5
10 x 10 cm	4 x 4"	9622	5
10 x 20 cm	4 x 8"	9623	5
15 x 15 cm	6 x 6"	9625	5
20 x 20 cm	8 x 8"	9626	5



Biatain® Fiber

Suitable for highly exuding wounds and locks in exudate and bacteria. Effectively absorbs and retains exudate and bacteria, thereby reducing exudate pooling. HexaLock Technology is the new gelling fiber for all exuding wounds.

reddot award 2017 winner

Biatain® Fiber

Size (cm)	Size (inches)	Product Code	Units
5 x 5 cm	2 x 2"	33501	10
10 x 10 cm	4 x 4"	33503	10
15 x 15 cm	6 x 6"	33506	10
19 x 25 cm	7.5 x 9.8"	33507	10
2.5 x 46 cm	1.5 x 18"	33508	10



Biatain® Alginate

Calcium Alginate Dressing

Highly absorbent alginate dressing for use on moderate to heavily exuding wounds of any size, shape and depth.

Biatain® Alginate

Size (cm)	Size (inches)	Product Code	Units
5 x 5 cm	2 x 2"	3705	30
10 x 10 cm	4 x 4"	3710	10
15 x 15 cm	6 x 6"	3715	10
3 x 44 cm	1 x 17.5" Ribbon	3740	6



Biatain® Alginate Ag

Calcium Alginate Dressing with Silver

Highly absorbent alginate dressing with silver for use on moderate to heavily exuding infected wounds or wounds at risk of infection of any size, shape and depth.

Biatain® Alginate Ag

Size (cm)	Size (inches)	Product Code	Units
5 x 5 cm	2 x 2"	3755	30
10 x 10 cm	4 x 4"	3760	10
15 x 15 cm	6 x 6"	3765	10
3 x 44 cm	1 x 17.5" Ribbon	3780	10



Biatain® Ibu

Foam Dressing with Ibuprofen

Soft and flexible absorbent foam dressing containing Ibuprofen for use on painful wounds. Available in non-adhesive and soft hold.

3DFit™ Technology

Biatain® Ibu Non-Adhesive

Size (cm)	Size (inches)	Product Code	Units
10 x 10 cm	4 x 4"	4110	5
10 x 20 cm	4 x 8"	4112	5
15 x 15 cm	6 x 6"	4115	5
20 x 20 cm	8 x 8"	4120	5



Biatain® Ibu Soft-Hold

Size (cm)	Size (inches)	Product Code	Units
10 x 10 cm	4 x 4"	4140	5
10 x 20 cm	4 x 8"	4142	5



Biatain® Super

Hydrocapillary Pad Dressing

Highly absorbent hydrocapillary pad dressing for use on heavily exuding wounds. Available in non-adhesive and adhesive.

Biatain® Super Adhesive

Size (cm)	Size (inches)	Product Code	Units
10 x 10 cm	4 x 4"	4610	10
12.5 x 12.5 cm	5 x 5"	4612	10
12 x 20 cm	4.5 x 8"	4625	10
15 x 15 cm	6 x 6"	4615	10
20 x 20 cm	8 x 8"	4620	10



Biatain® Super Non-Adhesive

Size (cm)	Size (inches)	Product Code	Units
10 x 10 cm	4 x 4"	4630	10
12.5 x 12.5 cm	5 x 5"	4632	10
12 x 20 cm	4.5 x 8"	4645	10
15 x 15 cm	6 x 6"	4635	10
20 x 20 cm	8 x 8"	4639	10



Comfeel® Plus

Hydrocolloid Dressing

Hydrocolloid dressing for use on low to moderate exuding wounds. Provides faster wound healing by sealing and protecting the wounds.

Comfeel® Plus

Size (cm)	Size (inches)		Product Code	Units
4 x 6 cm	1.5 x 2.5"		33146	30
10 x 10 cm	4 x 4"		33110	10
15 x 15 cm	6 x 6"		33115	5
20 x 20 cm	8 x 8"		33120	5
18 x 20 cm	7 x 8"	Sacral	33285	5
9 x 11 cm	3 x 4"	Contour	33283	5



33285



Comfeel® Plus Transparent

Hydrocolloid Dressing

Thin and flexible hydrocolloid dressing for use on dry or low exuding wounds. Provides faster wound healing by sealing and protecting the wound. Transparency allows for visual monitoring of the wound.

Comfeel® Plus Transparent

Size (cm)	Size (inches)	Product Code	Units
5 x 7 cm	2 x 2.75"	33530	10
10 x 10 cm	4 x 4"	33533	10
9 x 14 cm	3.5 x 5.5"	33536	10
15 x 15 cm	6 x 6"	33539	5
9 x 25 cm	3.5 x 10"	33537	5
5 x 15 cm	2 x 6"	33547	10
5 x 25 cm	2 x 10"	33548	10



33537



Triad™

Hydrophilic Wound Dressing

Zinc-oxide based hydrophilic sterile coating. Used for autolytic debridement. The ideal alternative for difficult-to-dress areas. Sterile. For more information please visit https://triad.coloplsat.ca/.

Triad™			
Size (ounces)	Size (grams)	Product Code	Units
2.5 oz.	71 g tube	1964	12
6 oz.	170 g tube	1967	12



Purilon® Gel

Hydrogel with natural ingredients

Effective and gentle debridement.

Purilon® Gel			
Size (ounces)	Size (grams)	Product Code	Units
28 oz.	8 g tube	3906	10
5 oz.	15 g tube	3900	10



Woun'Dres®

Collagen Hydrogel

Rehydrates dry wounds and eschar.

Woun'Dres® Collagen Hydrogel			
Size (ounces)	Size (grams)	Product Code	Units
1 oz.	28 g tube	1166	36
3 oz.	84 g tube	7690	12



Sea-Clens®

Wound Cleanser

Saline based solution with a surfactant for cleansing and irrigating acute and chronic wounds.

Sea-Clens® Wound Cleanser			
Size (ounces)	Size (mL)	Product Code	Units
6 oz.	178 mL	1063	12
12 oz.	355 mL	1061	12



Skin

Sproam®

No-Rinse Spray or Foam Cleanser

Spray or foam body wash, shampoo and incontinence cleanser.

Sproam®		
Size	Product Code	Units
175 ml	0897	12
350 ml	0901	12



Bedside-Care® Sensitive Skin

No-Rinse Foam Cleanser

Unscented foaming body wash, shampoo and incontinence cleanser.

Bedside-Care® Sensitive Skin		
Size	Product Code	Units
125 mL	7300	36
240 mL	7301	12



Bedside-Care® Foam

No-Rinse Foam Cleanser

Scented foaming body wash, shampoo and incontinence cleanser.

Bedside-Care®		
Size	Product Code	Units
240 mL	67145	12



Bedside-Care®

No-Rinse Spray Cleanser

Spray body wash, shampoo and incontinence cleanser.

Bedside-Care®			
Size		Product Code	Units
120 mL	Scented	61768	36
240 mL	Scented	61762	12
240 mL	Unscented	61761	12



No-Rinse Cleanser · Rinse Cleanser · Moisturize · Protect · CHG Compatible · Neonate Safe

Easi-Cleanse® Bath

No-Rinse Wash Cloth.

Body wash, shampoo, and incontinence cleanser in a self-sudsing disposable wash cloth.

Easi-Cleanse		
Size	Product Code	Units
5 per package	7056	100 pkgs/case



Gentle Rain® Extra Mild

Rinse Cleanser

Body wash and shampoo.

Gentle Rain® Extra Mild		
Size	Product Code	Units
5 mL sachet	7231	300
50 mL	7230	36
600 mL (accessories available, see pg 62)	7233	12
1000 mL (accessories available, see pg 62)	7234	12



Gentle Rain® Extra Mild Fragrance Free

Rinse Cleanser

Body wash and shampoo.

Gentle Rain® Extra Mild Fragrance Free		
Size	Product Code	Units
5 mL sachet	7073	300
50 mL	7074	36
600 mL (accessories available, see pg 62)	7075	12



Atrac-Tain®

Moisturizer

Skin moisturizer for severely dry, scaly and fissured skin with 4% Alphahydroxyacid (AHA).

Atrac-Tain®		
Size	Product Code	Units
2 g sachet	1843	300
60 mL	1802	12
140 mL	1814	12



No-Rinse Cleanser · Rinse Cleanser · Moisturize · Protect · CHG Compatible · Neonate Safe

Sween® 24

Moisturizer

Once a day skin moisturizer for itchy, irritated, dry and flaky skin.

Sween® 24		
Size	Product Code	Units
4 g sachet	7090	300
90 g	7094	12



DIN 02250322

Sween Cream®

Moisturizer

Soothing body cream for irritated dry skin.

Sween Cream®		
Size	Product Code	Units
90 mL	7067	12
190 mL	7068	12



Sween® Lotion Unscented

Moisturizer

Everyday moisturizer. For all body use.

Sween® Lotion Unscented		
Size	Product Code	Units
60 mL	7222	36
120 mL	7223	36
620 mL	7224	12



No-Rinse Cleanser · Rinse Cleanser · Moisturize · Protect · CHG Compatible · Neonate Safe

Critic-Aid® Clear

Skin Protectant

Protects skin and treats irritation due to incontinence.

Critic-Aid® Clear		
Size	Product Code	Units
4 g sachet	7565	300
71 g	7573	12



DIN 02277778

Baza® Protect II

Skin Protectant

Protects skin and treats irritation due to incontinence. Helps seal out wetness & relieves sore, irritated skin.

Baza® Protect II		
Size	Product Code	Units
60 g	7086	12
140 g	7087	12



DIN 02315793

Baza Cleanse and Protect®

Skin Protectant

All-in-one cleanser, moisturizer and a breathable skin protectant for sore, irritated, chafed skin due to incontinence.

Baza Cleanse and Protect®			
Size		Product Code	Units
4 g sachet	(fragrance free)	7710	300
240 mL	(fragrance free)	7712	12
240 mL	(with fragrance)	7725	12



DIN 02243966

No-Rinse Cleanser | Rinse Cleanser | Moisturize | Protect | CHG Compatible | Neonate Safe

InterDry®

Moisture-wicking textile with antimicrobial silver complex

Effectively manage complications associated with skin folds. Appropriate for a variety of skin-to-skin and skin-to-device application sites. For more information please visit www.interdry.ca.

InterDry®			
Size (cm)	Size (Inches)	Product Code	Units
25 x 46 cm	10 x 18"	7915	40
25 x 92 cm	10 x 36"	7912	10
25 x 366 cm	10 x 144" textile on a roll	7910	10



Isagel®

compatible CHG

Hand Sanitizer Gel

Sanitizing hand gel with added moisturizer.

Isagel®		
Size	Product Code	Units
115 mL	1644	36
621 mL (*Accessories available. Please see page 62.*)	1645	12



Sween® Body Powder

All-body powder

Provides excellent absorption to keep patient dry.

Sween® Body Powder		
Size	Product Code	Units
90 g	0506	36
227 g	0505	36



No-Rinse Cleanser · Rinse Cleanser · Moisturize · Protect · CHG Compatible · Neonate Safe

Accessories

Coloplast® Wall Dispenser for Bag

Description	Product Code	Units
Compatible with 800 mL (Isagel®) and 1000 mL (Gentle Rain® Extra Mild) sized bags (7234 and 7041)*	7251	12



Coloplast® Wall Dispenser Drip Tray

Description	Product Code	Units
To be used with Wall Dispenser to prevent spillage (7251)*	7262	60



Embrace™ Dispenser Bracket

Description	Product Code	Units
Compatible with 600 mL (Gentle Rain® Extra Mild) and 621 mL (Isagel®) sized pump bottles (1645, 7224, 7233)*	0089	12



*Product not included.

Brief Statements

Biatain® Family of Dressings

Biatain® dressings are single-use only. See the device manual for detailed information regarding complete instructions, contraindications, warnings, precautions and potential adverse events. For further questions, call Coloplast Canada Corporation at 1-888-880-8605 and/or consult the company website at www.coloplast.ca.

InterDry®

InterDry® is a skin protectant for the management of skin folds and other skin-to-device contact areas. It reduces microbial colonization in the fabric. This single use fabric provides moisture transportation to keep skin dry while the antimicrobial in the fabric reduces odor. The safe use of InterDry fabric during pregnancy, lactation and on children has not been demonstrated. Although many patients benefit from the use of this device, results may vary. See the device labeling for detailed information including the indication, contraindications, warnings, precautions, and potential complications/adverse events. For further questions, call Coloplast Canada Corporation at 1-888-880-8605 and/or consult the company website at www.coloplast.ca.

Peristeen Anal Plug®

The Peristeen Anal Plug is intended as a barrier to faeces in the rectum. The anal plug can be used by people who are faecal incontinent or suffer from involuntary passing of air (flatus). It functions as an effective barrier to faeces in the rectum.

Contraindication: Peristeen Anal Plug must not be used, if you suffer from: hermorrhoids (3rd and 4th degree), and anal stenosis

Caution: This product should not be used without assessment by your health care professional. Contact your health care professional before you use Peristeen Anal Plug for the first time, as professional guidance in the use of the plug is vital. The anal plug should not be used if you are suffering from diarrhoea. If placed in the mouth, the Peristeen Anal Plug will expand and may cause choking. Remove the anal plug when you feel the urge to defecate. If you have no sensation in the rectal area and do not feel the urge to defecate, changing plug at regular intervals is recommended.

The Peristeen Anal Plug must be removed no later than 12 hours after insertion. If you experience difficulties removing the Peristeen Anal Plug, you must contact your health care professional. If you repeatedly experience discomfort when using the Peristeen Anal Plug, stop using it and contact your health care professional. Keep out of reach of children. See the device manual for detailed information regarding instructions, contraindications, warnings, precautions, and potential complications/adverse events or, call Coloplast Canada Corporation at 1-866-293-6349 and/or consult the company website at www.coloplast.ca.

Peristeen®

The Peristeen® Anal Irrigation System is intended to instill water into the colon through a rectal catheter, which incorporates an inflatable balloon, inserted into the rectum to promote evacuation of the contents of the lower colon. The Peristeen Anal Irrigation System is indicated for use by children (3 years - <12 years old), adolescent (12 years - < 18 years old), transitional adolescent (18 - <21 years old) and adult patients with neurogenic bowel dysfunction who suffer from fecal incontinence, chronic constipation, and/or time-consuming bowel management procedures.

Contraindications: Peristeen Anal Irrigation must not be used in the following situations: known obstruction of the large bowel due to strictures or tumors, acute inflammatory bowel disease, diverticulitis, complex diverticular disease, abdominal or anal surgery within the last 3 months, spinal shock phase, and in patients who are pregnant and have not used the system before (if the patient is pregnant and has never used anal irrigation before, they should not start the irrigation procedure during pregnancy).

Warning: Anal irrigation procedure should always be carried out with care. Bowel perforation is an extremely rare, but serious and potentially lethal complication to anal irrigation and will require immediate admission to a hospital, often requiring surgery.

See the device manual for detailed information regarding the procedure, contraindications, warnings, precautions, and potential complications/adverse events. For further information, call Coloplast Canada Corporation at 1-866-293-6349 and/or consult the company website at www.coloplast.ca.

SpeediCath®

SpeediCath® catheters are indicated for use by patients with chronic urine retention and patients with a post void residual volume (PVR) due to neurogenic and non-neurogenic voiding dysfunction. The catheter is inserted into the urethra to reach the bladder allowing urine to drain. There is a separate SpeediCath device intended for either males or females. This is a single-use only product. Patients are to contact their healthcare professional before performing self-catheterization for the first time and follow the instructions carefully. See the device manual for further information or, call Coloplast Canada Corporation at 1-866-293-6349 and/or consult the company website at www.coloplast.ca.

Self-Cath®

Self-Cath® Catheters are intended for use in male, female, and pediatric patients (neonates, infants, children, adolescents, and transitional adolescents) requiring bladder drainage as determined by their physician. The catheter is indicated for those individuals unable to promote a natural urine flow or for those individuals who have a significant volume of residual urine following a natural bladder-voiding episode. Single-use only. See the device manual for detailed information regarding the procedure, contraindications, warnings, precautions, and potential complications/adverse events. For further information, call Coloplast Canada Corporation at 1-866-293-6349 and/or consult the company website at www.coloplast.ca.

Triad™

Triad™ hydrophilic wound dressing is indicated for the local management of pressure and venous statis ulcers, dermal lesions/injuries, superficial wounds, scrapes, first-and second-degree burns, partial-and full-thickness wounds. Suitable for use on neonates.

Contraindications: Triad™ is contraindicated for third degree burns or infected wounds. See device label for complete information including warnings and precautions. For further questions, call Coloplast Canada Corporation at 1-888-880-8605 and/or consult the company website at www.coloplast.ca.

Product Code	Page
0089	66
112	44
114	44
116	44
208	43
210	43
212	43
214	43
240	43
305	43
306	43
308	43
310	43
408	44
410	44
412	44
414	44
416	44
418	44
0421	34
0422	34
0423	34
0424	34
450	44
460	44
475	43
0505	65
0506	65
0509	35
608	44
610	44
612	44
614	44
808	44
810	44
812	44
814	44
816	44
818	44
0897	61
0901	61
0925	34
1003	36
1061	58
1063	58
1110	36
1114	44
1116	44
1166	58
1500	36
1644	65
1645	65
1802	62
1814	62
1843	62
1907	34

Product Code	Page
1964	58
1967	58
2041	34
2115	32
2125	32
2151	31
2166	31
2181	31
2320	34
2325	34
2330	34
2340	34
2350	34
2490	29
2501	35
2520	29
2650	34
2655	33
2792	28
2793	28
2801	35
2802	35
2803	35
2807	35
2808	35
2831	27
2832	27
2833	27
2836	27
2837	27
2846	27
2847	27
2881	26
2882	26
2883	26
3210	32
3215	32
3220	32
3410	54
3412	54
3413	54
3416	54
3420	54
3423	54
3430	54
3462	54
3470	55
3472	55
3475	55
3483	54
3485	54
3488	54
3705	56
3710	56
3715	56
3740	56

Product Code	Page
3755	56
3760	56
3765	56
3780	56
3900	58
3906	58
4110	56
4112	56
4115	56
4120	56
4140	56
4142	56
4610	57
4612	57
4615	57
4620	57
4625	57
4630	57
4632	57
4635	57
4639	57
4645	57
5030	45
5030H	45
5035	45
5035H	45
5050	48
5062	48
5063	48
5070	48
5100	46
5105	55
5110	46
5130	46
5135	46
5161	47
5163	47
5167	47
5170	48
5174	48
5200	46
5205	46
5210	46
5212	46
5215	46
5221	45
5225	45
5225H	45
5230	45
5230H	45
5235	45
5235H	45
5290	46
5300	46
5390	46
5400	46

Product Code	Page
5490	46
5580	31
5585	31
5590	46
6105	54
6200	46
6300	46
6400	46
7056	62
7067	63
7068	63
7073	62
7074	62
7075	62
7086	64
7087	64
7090	63
7094	63
7222	63
7223	63
7224	63
7230	62
7231	62
7233	62
7234	62
7251	66
7262	66
7300	61
7301	61
7565	64
7573	64
7583	35
7690	58
7710	64
7712	64
7725	64
7910	65
7912	65
7915	65
8000	47
8009	32
8100	47
8200	47
8205	47
8300	47
8305	47
8400	47
8500	47
9500	35
9505	35
9622	55
9623	55
9625	55
9626	55
9631	55
9632	55

Product Code	Page
9635	55
9641	55
9643	55
10011	20
10012	20
10013	20
10015	21
10016	21
10017	21
10021	20
10022	20
10023	20
10025	21
10026	21
10027	21
10031	20
10032	20
10035	21
10036	21
10041	20
10045	21
10101	22
10102	22
10103	22
10104	23
10105	23
10106	23
10108	22
10111	22
10113	22
10114	22
10116	22
10117	22
10154	21
10155	21
10156	21
10164	21
10165	21
10166	21
10167	21
10451	13
10461	13
10471	13
10472	13
10473	13
10474	13
10475	13
10476	13
10477	13
10478	13
10479	13
10480	13
10481	13
10482	14
10483	14
10484	14

Product Code	Page
10487	14
10489	13
10502	9
10503	9
10504	9
10505	9
10506	9
10512	9
10513	9
10514	9
10515	9
10516	9
10522	9
10523	9
10524	9
10525	9
10531	9
10551	12
10561	12
10562	12
10563	12
10564	12
10565	12
10566	12
10567	12
10568	12
10571	12
10572	12
10573	12
10582	17
10583	17
10584	17
10585	17
10586	17
10870	16
10871	16
10872	16
10873	16
10874	16
10881	16
10883	16
10884	16
10885	16
10904	23
10905	23
10914	23
10915	23
10918	23
11011	20
11012	20
11015	21
11016	21
11021	20
11022	20
11023	20
11024	20

Product Code	Page
11025	21
11026	21
11027	21
11028	21
11031	20
11032	20
11035	21
11036	21
11041	20
11045	21
11114	21
11115	21
11124	21
11125	21
11126	21
11127	21
11135	21
11136	21
11185	21
11186	21
11194	21
11195	21
11196	21
11197	21
11301	22
11302	22
11303	22
11304	23
11305	23
11306	23
11307	22
11308	23
11312	22
11313	22
11314	22
11315	22
11316	22
11401	10
11402	10
11412	10
11413	10
11421	10
11422	10
11425	10
11432	10
11443	10
11461	10
11462	10
11463	10
11471	10
11472	10
11473	10
11474	10
11482	10
11483	10
11484	10

Product Code	Page
11486	10
11487	10
11491	11
11492	11
11496	11
11497	11
11498	11
11499	11
11500	11
11505	23
11514	23
11515	23
11518	23
11612	23
11615	23
11690	23
11801	25
11802	25
11803	25
11804	25
11811	25
11812	25
11814	25
11815	25
11821	23
11822	23
11824	23
11825	23
11827	23
11828	23
11841	22
11842	22
11844	22
11845	22
11846	22
11851	22
11852	22
11854	22
11855	22
11856	22
12000	33
12010	33
12011	33
12014	19
12020	33
12021	33
12030	33
12032	33
12033	33
12035	32
12036	33
12037	32
12038	33
12039	32
12042	33
12045	32

Product Code	Page
12046	33
12047	32
12048	33
12049	32
12050	34
12060	34
12061	34
12070	33
12072	33
12074	33
12076	33
12090	33
12091	33
12092	33
12093	33
12094	33
12095	33
12110	30
12130	30
12140	30
12144	30
12145	30
12146	30
12147	30
12149	30
12170	30
12174	30
12175	30
12176	30
12177	30
12201	13
12202	13
12212	13
12213	13
12222	13
12232	12
12243	12
12262	12
12271	12
12272	12
12274	12
12282	12
12284	12
12291	13
12292	13
12297	13
12298	13
12344	28
12346	28
12354	28
12355	28
12356	28
12374	28
12375	28
12376	28
12384	28

Product Code	Page
12385	28
12386	28
12432	31
12433	31
12471	31
12472	31
12533	29
12595	31
12596	31
12701	26
12702	26
12703	26
12704	26
12705	26
12706	26
12707	26
12708	26
12709	26
12716	26
12719	26
12800	30
12802	30
12805	30
12808	30
12810	30
12820	36
12834	36
12835	36
12836	36
12841	26
12842	26
12844	26
12845	26
12847	26
12848	26
12991	31
12992	31
12993	31
12994	31
12995	31
12996	31
12997	31
13671	14
13672	14
13673	14
13674	14
13675	14
13676	14
13678	17
13679	17
13680	17
13683	16
13684	16
13685	16
13706	29
13840	29

Product Code	Page
13844	29
13845	29
13846	29
13860	17
13870	29
13874	29
13875	29
13876	29
13924	28
13925	28
13964	28
13965	28
13966	28
13974	28
13975	28
13976	28
13984	28
13985	28
13986	28
14010	37
14015	37
14021	37
14031	37
14041	37
14050	37
14060	37
14070	37
14103	29
14106	29
14163	29
14164	29
14165	29
14166	29
14173	29
14174	29
14175	29
14176	29
14196	29
14201	31
14202	31
14224	28
14225	28
14226	28
14227	28
14228	28
14229	28
14232	27
14233	27
14234	27
14235	27
14236	27
14237	27
14238	27
14239	27
14243	27
14246	27

Product Code	Page
14249	27
14251	27
14252	27
14254	27
14255	27
14258	27
14261	26
14262	26
14263	26
14271	26
14272	26
14273	26
14274	26
14275	26
14276	26
14277	26
14278	26
14281	27
14282	27
14283	27
14291	27
14292	27
14293	27
14294	27
14295	27
14296	27
14307	32
14308	32
14404	29
14405	29
14411	29
14412	29
14414	29
14415	29
14422	29
14423	29
14424	29
14441	30
14445	30
14452	30
14453	30
14454	30
14494	28
14498	28
14681	32
14682	32
14691	32
14692	32
14701	31
14702	31
15403	25
15404	25
15405	25
15406	25
15407	25
15408	25

Product Code	Page
15413	25
15414	25
15415	25
15441	25
15450	25
15470	25
15472	25
15473	25
15474	25
15480	25
15482	25
15483	25
15484	25
15521	24
15522	24
15523	24
15524	24
15531	24
15532	24
15533	24
15534	24
15605	24
15606	24
15607	24
15608	24
15633	24
15634	24
15635	24
15695	24
15696	24
15821	24
15981	24
16301	16
16305	16
16306	16
16310	16
16311	16
16312	16
16313	16
16315	16
16322	16
16323	16
16325	16
16326	16
16327	16
16330	16
16331	16
16332	16
16333	16
16335	16
16345	16
16346	16
16347	16
16350	16
16351	16
16352	16

Product Code	Page
16353	16
16471	12
16472	12
16481	12
16483	12
16484	12
16485	12
16486	12
16487	12
16488	12
16491	12
16493	12
16494	12
16705	14
16706	14
16715	14
16716	14
16725	14
16726	14
16727	14
16730	14
16731	14
16732	14
16733	14
16734	14
16735	14
16737	14
16738	14
16739	14
16748	14
16755	14
16756	14
16757	14
16760	14
16761	14
16762	14
16763	14
16764	14
16767	14
16768	14
16770	14
16771	14
16772	14
16773	14
16778	14
16801	17
16805	17
16810	17
16825	17
16826	17
16827	17
16835	17
16836	17
16837	17
16840	17
16841	17

Product Code	Page
16842	17
16855	17
16856	17
16857	17
16865	17
16866	17
16867	17
16870	17
16871	17
16872	17
16901	9
16903	9
16911	9
16913	9
16914	9
16915	9
16921	9
16923	9
16924	9
16931	9
16933	9
16934	9
16941	9
16942	9
16951	9
16953	9
16954	9
16955	9
16961	9
16963	9
16964	9
16971	9
16973	9
16974	9
18151	17
18152	17
18160	17
18161	17
18162	17
18163	17
18170	17
18171	17
18172	17
18173	17
18180	17
18181	17
18182	17
18185	17
18300	15
18301	15
18302	15
18305	15
18310	15
18312	15
18313	15
18314	15

Product Code	Page
18315	15
18320	15
18321	15
18322	15
18323	15
18324	15
18330	15
18331	15
18332	15
18333	15
18340	15
18341	15
18342	15
18343	15
18344	15
18386	13
18389	12
18390	12
18401	18
18402	18
18403	18
18404	18
18405	18
18500	10
18501	10
18510	10
18511	10
18512	10
18513	10
18520	10
18521	10
18522	10
18523	10
18524	10
18530	10
18531	10
18532	10
18600	11
18601	11
18602	11
18603	11
18610	11
18611	11
18612	11
18613	11
18620	11
18621	11
18622	11
18630	11
18631	11
18632	11
18640	11
18641	11

Product Code	Page
18642	11
18650	13
18651	13
18652	13
18655	13
18656	13
18657	13
18660	18
18661	18
18665	18
18666	18
18670	18
18671	18
18675	18
18676	18
18680	18
18681	18
18690	18
18691	18
18700	19
18705	19
18710	20
18711	20
18712	20
18713	20
18720	19
18721	19
18722	19
18723	19
18724	19
18725	19
18730	19
18731	19
18735	19
18747	19
19001	25
19003	25
19010	25
19011	25
19020	25
19021	25
19030	25
19040	22
19041	22
19042	22
19050	22
19051	22
19052	22
21034	47
21054	47
21346	48
21356	48
21365	37

Product Code	Page
22025	45
22028	45
22030	45
22035	45
22040	45
22121	45
22125	45
22130	45
22135	45
25501	47
27408	41
27410	41
27412	41
27414	41
27416	41
27418	41
27490	41
27492	41
27494	41
27496	41
27508	41
27510	41
27512	41
27514	41
27516	41
27706	41
27708	41
27710	41
28110	42
28112	42
28114	42
28422	43
28520	43
28522	43
28524	43
28578	42
28580	42
28582	42
28584	42
28692	42
28810	42
28812	42
28814	42
28916	41
28920	41
28922	41
28924	41
29140	49
29142	49
29143	49
29145	49
29146	49
29147	49

Product Code	Page
29149	49
33110	57
33115	57
33120	57
33146	57
33283	57
33285	57
33400	53
33401	53
33406	53
33434	53
33435	53
33436	53
33437	53
33438	53
33444	53
33445	53
33446	53
33452	53
33453	53
33501	55
33503	55
33506	55
33507	55
33508	55
33530	57
33533	57
33536	57
33537	57
33539	57
33547	57
33548	57
39000*	53
39001*	53
39010*	53
39020	53
39021	53
39022	53
39023	53
39024	53
39025	53
39026	53
39636	54
39637	54
39638	54
39639	54
39640	54
39644	54
39645	54
39650	54
39651	54
39652	54
61761	61

Product Code	Page
61762	61
61768	61
67145	61
383004	48
504311	46
504331	46
504341	46
504351	46
505311	46
505321	46
505331	46
505341	46
505421	46
505431	46
505441	46
505451	46
505711	46
505731	46
505741	46
505751	46
506110	47
506130	47
506140	47
506150	47
506200	47
506220	47
506240	47
506260	47
C5200	46
C5210	46
C5310	46
C5410	46
C6100	46

Coloplast® Care

Support for your ostomy or continence care needs







Website and emails
A dedicated website and regular emails with articles and reliable advice

Self-Assessment Tools
Helping users keep track of their ostomy or continence health through regular self-assessment

Phone support
Advisors ready to answer questions and provide support by phone

Product guidance
Access to product samples on request

Enrolling in Care is easy!

Online at **www.coloplast.ca/careform**

By phone at 1-866-293-6349

Request Coloplast Care enrollment forms from your Coloplast Territory Manager

Velcro® is a registered trademark of Velcro Industries B.V.

Ostomy Care / Continence Care / Wound & Skin Care / Interventional Urology

Coloplast Corp. Oakville, ON L6H 5S9 / 1-888-880-8605
www.coloplast.ca

The Coloplast logo is a registered trademark of Coloplast A/S. © 2022 Coloplast Corp. All rights reserved.

Coloplast